

SENIOR CARE PRODUCTS

summer 2012 www.sproductsmag.com | solutions for independent living & aging in place

Accessible Tubs & Showers FOR THE HME MARKET



also inside:

Products for
Seniors with
Diabetes



FROM THE PUBLISHER OF
HomeCare

MEDICATION MANAGEMENT | HME PROVIDER PROFILE | **TRAVEL WHEELCHAIRS**
HOME MONITORING & PERS | **BEDS & ACCESSORIES** | **MEDTRADE SPRING SCRAPBOOK**



respiratory BREATHING CARE

introduction to pulse oximeters

story by ASHLEY WOOD

How they work and how to select the best products for your clients.

Pulse oximetry is a noninvasive method for monitoring the oxygen saturation levels (SpO_2) in your patients' blood. SpO_2 is the percentage of oxygen that blood carries compared to the maximum it is capable of carrying. A normal rate is 94 to 99 percent in adults. Patients with COPD or other lung diseases generally have blood oxygen levels of 92 percent or lower, causing them to require supplemental oxygen treatment.

WHEN IS A PULSE OXIMETER USED?

In general, pulse oximeters are used to provide insight into the functioning of the patient's respi-

ratory system. It can be used anytime these levels need to be tested or diagnosed.

In the senior care and home care setting, a portable pulse oximeter like the one shown above is most commonly used to treat or diagnose patients with respiratory diseases that cause low levels of oxygen in their blood when they breathe air. A pulse oximeter is an easy, noninvasive way to measure oxygen levels at any time so supplemental oxygen amounts can be adjusted as needed throughout the day.

In addition, different types of pulse oximeters are also commonly used in the hospital setting, often during surgery with anesthesia or in the recovery room. They are also used in sleep labs

as a part of sleep apnea diagnosis. In a nonmedical setting, pilots and athletes commonly require spot checks at high altitudes or during rigorous training.

HOW DO THEY WORK?

Generally, pulse oximeters work by shining two different kinds of light (red and infrared) through a thin part of a patient's body, usually a fingertip or earlobe. Blood absorbs light differently depending on the level of oxygen it contains: oxygenated blood absorbs more infrared light than red light. So by comparing the changes in amounts of red and infrared light received, the instrument can calculate the SpO_2 reading.

As a side benefit, you can also check your pulse, as the increase in the amount of blood with each heartbeat also affects the light.

IS A PRESCRIPTION REQUIRED?

It depends on the device and how it's being used. According to the FDA, if an oximeter is accurate enough to diagnose or treat a disease, a prescription is generally required. If not, a prescription may not be needed. For example, if an athlete or rock climber wanted to monitor their pulse and oxygen levels, they could purchase an OTC model that is not intended for diagnosis and would not require a prescription. To be sure, ask the device manufacturer.

WHAT TYPE OF PULSE OXIMETER DOES MY CLIENT NEED?

There are three common types that fit a variety of

needs. The most common is the portable fingertip pulse oximeter, which is often used at home. It is relatively inexpensive, lightweight and easy to use. However, it doesn't always provide an accurate reading for those with circulatory problems. Also, it is designed for spot checks instead of continuous monitoring, so there are also several other options available.

The handheld pulse oximeter is most commonly used in hospitals and has a clip and wire attached to a handheld monitor. It can be clipped to a finger, earlobe or toe in an emergency situation and can be used for spot checks or continuous monitoring. The monitor can record the patient's information for several hours at a time. Many models also come with a built-in alarm that sounds if the patient's SpO₂ or pulse rate goes outside a designated range. Some also include the ability to connect directly to a printer to prevent

documentation errors.

The final type is a wrist pulse oximeter that has a small fingertip sensor that is attached to a wrist-watch style recording system that continuously monitors the blood oxygen saturation level and pulse rate. It is ideal for monitoring daily activities or for an overnight sleep study because it allows the freedom to perform daily tasks and has a memory that can last for several days at a time.

ABOUT THE AUTHOR: Ashley Wood is the director of operations for Sunset Healthcare Solutions, a national manufacturer and distributor of CPAP, oxygen and other respiratory supplies. She can be reached at ashley@sunsethcs.com or 312-533-2446. For more information on the products mentioned in this article, visit the company's website at www.sunsethcs.com.



SUNSET HEALTHCARE SOLUTIONS COMFORT CANNULA WITH CUSHIONS

Increase patients' comfort and satisfaction with Sunset's Comfort Cannula with built-in soft foam cushions. The cushions protect sensitive ear tissue and are added during the manufacturing process so they are guaranteed to stay on the tube with no assembly required. Simply slide the cushions into place behind the top of the ear for the best fit. All Sunset oxygen nasal cannulas are tested and proven to provide maximum safety and comfort. The tubes are six-channel, kink-free safety tubes to ensure uninterrupted airflow. The nasal prongs are curved and tapered for comfort that feels natural and the material is soft, memory-free and lightweight for long-term use. Ear cushions are available as a separate purchase to add them to your own cannula.

Sunset Healthcare Solutions, 877-578-6738, www.sunsethcs.com

Circle 176 on Product Card



PHILIPS RESPIRONICS ULTRAFILL HOME OXYGEN FILLING SYSTEM

Philips Respironics' UltraFill is an advanced home oxygen system that combines a stationary oxygen concentrator, filling station and standard or optional high-capacity cylinders to meet the needs of a wide range of oxygen patients, including those who are highly active or require continuous flow oxygen. UltraFill's long-lasting 3,000 PSI cylinders are capable of providing patients with greater than 40 percent more oxygen than similar sized 2,000 PSI cylinders. UltraFill also has the needs of home care providers in mind with features intended to save money by reducing deliveries and utilizing existing inventory.

Philips Respironics, 724-387-5000, www.philips.com/respironics

Circle 178 on Product Card

All HCPCS codes subject to change.



RESPCARE ALOHA NASAL PILLOW SYSTEM

RespCare introduces the latest design of its CPAP interface device, the Aloha. The Aloha is a nasal-pillow style patient interface that is lightweight and quiet, and offers a stable and comfortable seal. The Aloha features proprietary Arced-Track technology that allows the angle of the pillow reservoir to be adjusted for a more customized, natural fit. The interface's ball-and-socket elbow adapts to the most active sleepers, accommodating various sleeping positions. The Aloha system comes standard with three different pillow sizes and a headgear with Comfort Wraps. The Aloha contains all of the components required to properly size the patient, eliminating the hassle of opening multiple products to get the right fit.

RespCare, Inc., 800-200-9842, www.respcare.com

Circle 177 on Product Card