

# UltraMax 02

OXYGEN ANALYZER

Instructions for Use

ENGLISH





phone: (800) 748.5355 fax: (801) 973.6090 email: sales@maxtec.com web: www.maxtec.com **(€**-0123

Conforms to:

AAMI STD ES60601-1, ISO STD
80601-2-55, IEC STDS 606011-6, 60601-1-8 & 62366
Certified to:
CSA STD C22.2 No. 60601-1

**NOTE:** The latest edition of this operating manual can be downloaded from our website at **www.maxtec.com** 

**NOTE:** The UltraMax 02 is for use only by trained personnel. Before use, all individuals using the UltraMax 02 should become familiar with the information contained in this Operation Manual. Adherence to these instructions is necessary for safe, effective product performance. Thoroughly read all instructions and labeling provided with this device and any other equipment that will be used.

#### CLASSIFICATION

Protection against electric shock	Internally powered equipment
Protection against water	IPX1
Mode of operation	
Sterilization	
Flammable anaesthetic mixture Not for use in presence of	flammable anaesthetic mixtures
Power specification	1.8-3.2V = = = 32mW10mA

**CAUTION:** Federal law restricts this device to sale by or on the order of a physician or other licensed healthcare practitioner.



#### **Product Disposal Instructions:**

The sensor, batteries, and circuit board are not suitable for regular trash disposal. Return sensor to Maxtec for proper disposal or dispose according to local guidelines. Follow local guidelines for disposal of other components.

#### INDICATIONS FOR USE

The UltraMax 02 Oxygen Analyzer is a tool used to measure oxygen purity, flow and pressure at the outlet of an oxygen concentrator. It is not intended to be used by patients who are prescribed oxygen, nor is it intended to continuously monitor or confirm oxygen delivery to a patient. The UltraMax 02 Oxygen Analyzer is intended to be used in an environment where oxygen concentrators are being serviced or repaired. This includes Hospitals, Nursing Homes, Extended Care Facilities, Patient Homes, and Respiratory Device Service and Repair Centers.

#### WARRANTY

Under normal operating conditions, Maxtec warrants the UltraMax 02 to be free from defects of workmanship or materials for a period of Three (3) years from the date of shipment from Maxtec, provided that the unit is properly operated and maintained in accordance with Maxtec's operating instructions. Based on Maxtec product evaluation, Maxtec's sole obligation under the foregoing warranty is limited to making replacements, repairs, or issuing credit for equipment found to be defective. This warranty extends only to the buyer purchasing the equipment directly from Maxtec or through Maxtec's designated distributors and agents as new equipment.

Routine maintenance items, such as batteries, are excluded from warranty. Maxtec and any other subsidiaries shall not be liable to the purchaser or other persons for incidental or consequential damages or equipment that has been subject to abuse, misuse, mis-application, alteration, negligence or accident.

These warranties are exclusive and in lieu of all other warranties, expressed or implied, including warranty of merchantability and fitness for a particular purpose.

#### PRINCIPLE OF OPERATION

The UltraMax O2 Oxygen Analyzer measures oxygen concentration and flow using ultrasound technology and measures pressure using a piezoresistive silicon pressure sensor.

#### **WARNINGS**

Indicates a potentially hazardous situation, if not avoided, could result in death or serious injury.

- Not for use in an MRI environment.
- Improper use of the UltraMax 02 can cause inaccurate oxygen readings leading to improper treatment and/or patient harm. Follow the procedures outlined in this user manual.
- The UltraMax 02 is for checking oxygen concentrators only.

- **DO NOT** use the UltraMax O2 for continuous oxygen monitoring.
- **DO NOT** use the UltraMax 02 to measure the oxygen concentration of a concentrator when flowing at rates lower than its optimal performance as specified by the concentrator manufacturer; generally 4 LPM or less on concentrators that have a maximum flow of 10 LPM, and 1 LPM or less on concentrators that have a maximum flow of 5 LPM.
- Not for use in anesthesia applications or for measuring oxygen concentration from any sources other than conventional oxygen concentrators.
- Not for use with inhalation agents. Operating the UltraMax 02 in flammable or explosive environments may result in fire or explosion.
- Not suitable for use in the presence of flammable anesthetic mixtures.
- Oxygen rapidly accelerates combustion.
- **DO NOT** smoke while using the UltraMax 02 for checking oxygen concentrators.

Users must become thoroughly familiar with the information contained in this Operation Manual before use. Strict adherence to the operating instructions is necessary for safe, effective product performance. This product will perform only as designed if operated in accordance with the manufacturer's operating instructions.

- Use only genuine Maxtec accessories. Failure to do so may seriously impair the performance
  of the UltraMax 02. Repair or alteration of the UltraMax 02 by anyone other than an authorized Maxtec service representative could cause the product to fail to perform as designed.
- Use of the UltraMax O2 near devices that generate electrical fields may cause erratic readings.
- If the UltraMax 02 is ever exposed to liquids from spills or immersion, immediately remove the batteries and let the device dry completely. When dry, replace the batteries and check for proper operation.
- **DO NOT** autoclave or expose the UltraMax O2 to high temperatures (>60°C).
- **DO NOT** use ethylene oxide sterilization.
- **DO NOT** expose the UltraMax O2 to irradiation, vacuum, steam, or harsh chemicals.
- **DO NOT** expose the UltraMax 02 to pressure greater than 50 psi. Exposure to pressure above 50 psi could cause leaks in the device which may adversely affect performance in flow and pressure readings.

### **CAUTIONS**

Indicates a potentially hazardous situation, if not avoided, could result in minor or moderate injury and property damage.

- Replace the batteries with high quality AA Alkaline or Lithium batteries.
- **DO NOT** use rechargeable batteries.
- When not in use for periods greater than 30 days remove the batteries to protect the UltraMax 02 from potential battery leakage.
- Avoid dropping the UltraMax 02 to prevent damage which may adversely affect its performance. If damage to the device is suspected, perform the calibration verification procedure in Section 2.3 of this operating manual.
- Avoid foreign matter entry into the UltraMax 02.
- **DO NOT** use the UltraMax 02 to check a concentrator with a humidifier in place. Humidity from a humidifier could damage the device.
- **DO NOT** check a concentrator while holding the mode button or the reading will be inaccurate.
- Following storage in extremely hot or cold conditions, allow the gas to flow through the analyzer long enough for the internal sensors to reach the gas stream temperature, or wait for the analyzer to equilibrate to room temperature before use.

#### **Symbol Guide**

The following symbols and safety labels are found on the UltraMax 02:

A	Warning	EC REP	Authorized Representative in the European Commuity
BAT	Low Battery	SN	Serial Number
	Do not throw away. Follow local guidelines for disposal	REF	Catalog Number
ett classified c us us Intertek 9700630	Meets ETL standards	LPM	Liter per minute flow
	Manufacturer	PSI	Pounds per square inch
$\mathbb{Z}^{\mathbb{Z}}$	Date of Manufacture	KPA	Kilopascals
MD	Medical Device	%	Percent
IPX1	Ingress Protection Rating	4	Gas sample inlet
$R_{\!$	Federal law (USA) restricts this device to sale by or on order of a physician.	<b>←</b> □	Gas sample outlet
	Latex free		Direct current
0	On/Off Button	0	DO NOT
4	Mode Button	Â	Caution
	Follow instructions for use		

#### 1.0 SYSTEM OVERVIEW

#### 1.1 Description & Principle of Operation


The UltraMax 02 is an oxygen analyzer designed to check the oxygen concentration, flow and outlet pressure of oxygen concentrators. The UltraMax 02 provides unparalleled performance and reliability from its advanced design that includes the following features and operational benefits:

- Accurate oxygen measurements.
- No in-field calibration required.
- · Convenient ability to measure pressure in PSI or kPa.
- Durable, compact design.
- Large, easy-to-read, liquid crystal display (LCD).
- Shielded, reinforced sample gas inlet port.
- Long battery life with 2 AA batteries.
- Auto-off after 4 minutes.
- Low battery indication.
- Self-diagnostics.
- Easy to clean.

#### 1.2 Indication for Use

The UltraMax 02 Oxygen Analyzer is a tool used to measure oxygen purity, flow and pressure of an oxygen concentrator. The UltraMax 02 Oxygen Analyzer is intended to be used in an environment where oxygen concentrators are being serviced or repaired. This includes Hospitals, Nursing Homes, Extended Care Facilities, Patient Homes, and Respiratory Device Service and Repair Centers.

#### 1.3 Component Identification



- (1) **3 1/2 DIGIT DISPLAY** The LCD provides direct readout of oxygen concentration, gas flow and gas pressure. The LCD also displays error codes as necessary.
- **MODE BUTTON** Switches between measuring the concentration of gas produced by an oxygen concentrator and pure oxygen (for calibration verification).
- (3) **ON/OFF BUTTON** Turns the device on or off.
- (4) **PSI** Indicates the pressure measurement is in units of pounds per square inch.
- (5) **KPA** Indicates the pressure measurement is in units of kilopascals.
- **(6) GAS SAMPLE INLET** Used to receive the gas sample.
- (7) **% SYMBOL** Illuminated next to the concentration measurement.
- (8) LOW BATTERY INDICATOR Indicates the voltage of the batteries is below normal operating levels.
- (9) LPM Illuminated next to the flow measurement. (Not shown when in calibration verification mode).
- GAS SAMPLE OUTLET Used as an outlet for the gas sample and as a trigger for pressure measurement when occluded.
- (11) BATTERY DOOR

**GAS SAMPLE TUBING** — Used to connect to gas sample sources (not shown).

#### 2.0 OPERATING INSTRUCTIONS

#### 2.1 Oxygen, Flow and Pressure Measurement

To check oxygen concentration, flow and pressure of a gas sample from a concentrator:

Connect the gas sample tubing to the gas sample inlet of the UltraMax 02.

- 1. Attach the other end of the gas sample tubing to the oxygen concentrator.
- 2. Initiate the flow of gas to the UltraMax O2 at a rate of 1-10 liters per minute (2 liters per minute is recommended). Ensure the concentrator's output is stable per the concentrator manufacturer's recommendations.
- 3. Turn on the UltraMax 02.
- Allow the oxygen reading to stabilize for approximately 10 seconds before reading the oxygen concentration and flow.
- To check pressure, cover the gas sample outlet with thumb or finger while gas is flowing.
- 6. Wait 5 seconds for the display to read pressure.
- **DO NOT** hold the mode button while checking a concentrator or the reading will be inaccurate.

#### 2.2 Switching Pressure Units of Measure

The UltraMax 02 can measure pressure in PSI or kPa. The UltraMax 02 is factory set to measure in PSI. To switch to kPa:

- Using a #1 Phillips screwdriver loosen the battery door screw and remove the battery door
- 2. Toggle the switch inside the battery compartment.
- 3. Replace the battery door and tighten the battery door screw.

#### 2.3 Calibration Verification Procedure

A calibration verification mode is provided to verify that the UltraMax 02 is functioning properly. To perform the calibration verification:

- 1. Turn on the UltraMax 02.
- 2. Connect a source of pure oxygen (≥99.95%) to the gas sample inlet.
- 3. Flow 2-5 LPM of gas into the UltraMax 02. Ensure that the gas flowing to the UltraMax 02 is at a stable temperature.
- 4. Press and hold the mode button. While holding the mode button, the gas measurement should read between 98.5 and 101.5% oxygen. If the gas measurement is not within this range, call Maxtec Customer Service. Calibration verification mode is indicated by "CAL" and "VER" flashing on screen beneath the gas measurement.

## 3.0 FACTORS INFLUENCING ACCURATE READINGS

#### 3.1 Effects of Temperature

The UltraMax 02 compensates for temperature and will perform within specifications throughout the operating temperature range. However, taking measurements during rapid changes in gas temperature should be avoided.

#### 3.2 Effects of Humidity

The UltraMax O2 has a humidity sensor to detect and compensate for the humidity of gas entering the device. However, high levels (condensing) of humidity can affect the accuracy and reliability of the UltraMax O2. To prevent possible damage:

Avoid usage in environments of greater than 95% relative humidity.

**DO NOT** use this device in a breathing circuit.

**DO NOT** breathe or blow into the UltraMax 02.

#### 3.3 Effects of Other Gases

The UltraMax 02 is designed to measure two different types of gas mixtures:

- Oxygen, nitrogen and argon from oxygen concentrators.
- Pure oxygen during calibration verification mode.

Any other concentrations or combinations of gases will cause the UltraMax 02 to measure oxygen concentration incorrectly.

#### 3.4 Effects of Low Flow

Oxygen concentrators function on the principle of removing nitrogen gas from air, leaving concentrated oxygen and argon at a specific oxygen to argon ratio. This operating principle may be altered when concentrators are set to flow at the low end of their operational range. At low flows they may output a low oxygen concentration, e.g. 85% to 91%, for reasons other than high nitrogen, possibly due to an increase in argon content. The UltraMax 02 requires that the ratio of oxygen to argon remain constant in order to guarantee an accuracy of +/-1.5% oxygen.

**DO NOT** use the UltraMax 02 to measure the oxygen concentration of a concentrator when flowing at rates lower than its optimal performance as specified by the concentrator manufacturer; generally 4 LPM or less on concentrators that have a maximum flow of 10 LPM, and 1 LPM or less on concentrators that have a maximum flow of 5 LPM.

#### 4.0 ERROR CODES

The UltraMax 02 has self diagnostic features built into the software to detect faulty readings outside of normal operating ranges. The codes, descriptions and recommended actions are:

**E01:** Oxygen measurement out of range Hi (≥102.0% calculated by algorithm). Recommended Action: Verify that the UltraMax 02 is being used in the correct mode (Concentrator or Calibration Verification mode). If error code repeats; perform a calibration verification per section 2.3 of this manual. If error code repeats again; contact Maxtec Customer Service.

**E02:** Oxygen measurement out of range Low (≤-2.0% calculated by algorithm). Recommended Action: Verify that the UltraMax 02 is being used in the correct mode (Concentrator or Calibration Verification mode). If error code repeats; perform a calibration verification per section 2.3 of this manual. If error code repeats again; contact Maxtec Customer Service.

**E03:** Device memory corrupt or missing. Recommended Action: Return the UltraMax 02 to the manufacturer for factory repair.

**E04:** Signal reading not stable. Recommended Action: Return the UltraMax 02 to the manufacturer for factory repair.

**E05:** Pressure measurement out of Range Hi (≥50 PSI). Recommended Action: Check the pressure on a known gas source pressure. If error code repeats; contact Maxtec Customer Service.

**E06:** Outside of operating temperature Hi ( $\geq$ 40° C). Recommended action: The UltraMax 02 is too hot, cool the device closer to room temperature before use.

**E07:** Outside of operating temperature Low (≤15° C). Recommended action: The UltraMax 02 is too cold, warm the device closer to room temperature before use.

**E08:** Device self check found error. Recommended Action: Remove and replace the batteries. If error code repeats; return the UltraMax 02 to the manufacturer for factory repair.

#### 5.0 CHANGING THE BATTERIES

Batteries should be changed by service personnel. Use only brand name batteries. Replace with two AA batteries and insert per orientation marked on the device. Batteries should be changed when the DATT icon illuminates. The icon will remain lit until the batteries are changed. If the battery power level is too low the UltraMax O2 will not power on until the batteries are changed.

#### **5.1 Battery Replacement Procedure**

- Using a #1 Phillips screwdriver loosen the battery door screw and remove the battery door
- 2. Remove the batteries.
- 3. Insert new batteries ensuring correct polarity. **DO NOT** use rechargeable batteries.
- 4. Replace the battery door and tighten the battery door screw.
- If the UltraMax 02 does not power on when done verify the batteries are installed correctly and that the batteries are fresh.

#### **6.0 CLEANING AND MAINTENANCE**

Use caution to prevent any fluid from entering the UltraMax 02.

**DO NOT** soak or immerse the UltraMax 02 in fluid.

**DO NOT** autoclave or expose the UltraMax O2 to ethylene oxide sterilization.

#### 6.1 Cleaning

Wipe down the exterior surfaces of the UltraMax 02 with a moist cloth and mild hand or dish soap (pH 6-8).

#### 6.2 Maintenance

Replace the batteries with high quality AA Alkaline or Lithium batteries.

**DO NOT** use rechargeable batteries.

- When not in use for periods greater than 30 days, remove the batteries to protect the UltraMax 02 from potential battery leakage.
- Store the UltraMax 02 between -15°C and 60°C (5°F 140°F)

#### 7.0 SPECIFICATIONS

.... b...

#### Oxygen

Oxygen measurement kange	
(from a concentrator)	20.9 - 96%
Oxygen Measurement Accuracy	±1.5 % of full scale at
,	constant temperature and optimal flow*
Oxygen Measurement Resolution	
Flow	
Flow Measurement Range	0 - 10 LPM

Pressure	
Pressure Measurement Range	
Pressure Measurement Accuracy	±0.5% (PSI), ±0.5% (kPa)
Pressure Measurement Resolution	0.1 (PSI), 0.1 up to 199, 1 from 200 to 344 (kPa)
Response Time	≤17 seconds
Warm-up Time	<1second
Operating Temperature	15°C - 40°C (59°F-104°F)
Storage Temperature	-15°C - 60°C (5°F-140°F)
Pressure	800 - 1000 mBars
Power Requirements	2 AA Alkaline batteries (2 x 1.5 Volts)
Battery Life	≥ 1,100 hours (16,500 read cycles)
Low Battery Indication	"Low Battery" icon displayed on LCD
Dimensions	3.16" x 5.10" x 1.04" (80.3mm x 129.5mm x 26.4mm)
Weight	

#### 8.0 SPARE PARTS AND ACCESSORIES

#### 8.1 Included with Your Unit

PART NUMBER	ITEM
R211M11	Operating Manual and Instructions for Use*
RP46P05	Gas Sample Tubing

#### 8.2 Optional Accessories

PART NUMBER		ITEM	
	R221P15	Soft Cover	

Repair of this equipment must be performed by a qualified service technician experienced in repair of portable hand held medical equipment.

Equipment in need of repair shall be sent to:

Maxtec Customer Service Department 2305 South 1070 West Salt Lake City, UT 84119

(Include RMA number issued by Customer Service)

## 9.0 ELECTROMAGNETIC COMPATIBILITY

The information contained in this section (such as separation distances) is in general specifically written with regard to the UltraMax 02. The numbers provided will not guarantee faultless operation but should provide reasonable assurance of such. This information may not be applicable to other medical electrical equipment; older equipment may be particularly susceptible to interference.

Note: Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in this document and the remainder of the instructions for use this device.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Cables and accessories not specified within the instructions for use are not authorized. Using other cables and/or accessories may adversely impact safety, performance and electromagnetic compatibility (increased emission and decreased immunity).

Care should be taken if the equipment is used adjacent to or stacked with other equipment; if adjacent or stacked use is inevitable, the equipment should be observed to verify normal operation in the configuration in which it will be used.

ELECTROMAGNETIC EMISSIONS		
This equipment is intended for use in the electromagnetic environment specified below.  The user of this equipment should assure that it is used in such an environment.		
EMISSIONS	COMPLIANCE ACCORDING TO	ELECTROMAGNETIC ENVIRONMENT
RF Emissions (CISPR 11)	Group 1	The UltraMax 02 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
CISPR Emissions Classification	Class A	The UltraMax 02 is suitable for use in all establishments other than domestic and
Harmonic Emissions (IEC 61000-3-2)	Class A	those directly connected to the public low- voltage power supply network that supplies
Voltage Fluctuations	Complies	buildings used for domestic purposes.
		NOTE: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radiofrequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

#### **ELECTROMAGNETIC IMMUNITY**

This equipment is intended for use in the electromagnetic environment specified below. The user of this equipment should assure that it is used in such an environment.

IMMUNITY AGAINST	IEC 60601-1-2: (4TH Edition) test level		ELECTROMAGNETIC ENVIRONMENT
	Professional Healthcare Facility Environment	Home Healthcare Environment	
Electrostatic discharge, ESD (IEC 61000-4-2)	Contact discharge: ±8 Air discharge: ±2 kV, ±		Floors should be wood, concrete, or ceramic tile. If floors are covered
Electrical fast transients / bursts (IEC 61000-4-5)	Power supply lines: ±2 Longer input / output		with synthetic material, the relative humidity should be kept at levels to reduce electrostatic
Surges on AC mains lines (IEC 61000-4-5)	Common mode: ±2 kV Diferential mode: ±1 k		charge to suitable levels.  Mains power quality
3 A/m power frequency magnetic field 50/60 Hz (IEC 61000-4-8)	30 A/m 50 Hz or 60 Hz		should be that of a typical commercial or hospital environment.
Voltage dips and short interruptions on AC mains input lines (IEC 61000-4-11)	Dip>95%, 0.5 periods Dip 60%, 5 periods Dip 30%, 25 periods Dip >95%, 5 seconds		Equipment which emits high levels of power line magnetic fields (in excess of 30A/m) should be kept at a distance to reduce the likelihood of interference.
			If user requires continued operation during power mains interruptions, ensure that batteries are installed and charged. Ensure that battery life exceeds longest anticipated power outages or provide an additional uninterruptible power source.

Recommended separation distances between portable and mobile
RF communications equipment and the equipment

RATED MAXIMUM OUTPUT POWER OF	Separation distance according to frequency of transmitters in meters		
TRANSMITTER <b>W</b>	150 kHz to 80 MHz d=1.2/V1] √P	80 MHz to 800 MHz d=1.2/V1] √P	800MHz to 2.5 GHz d=2.3 √P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance

for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propogation is affected by absorption and reflection from structures, objects, and people.

This equipment is intended for use in the electromagnetic environment specified below. The customer or the user of this equipment should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601-1-2: 2014 (4TH Edition) test level		ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
	Professional Healthcare Facility Environment	Home Healthcare Environment	
Conducted RF coupled into lines (IEC 61000-4-6)	3V (0.15 - 80 MHz) 6V (ISM bands)	3V (0.15 - 80 MHz) 6V (ISM & Amateur bands)	Portable and mobile RF communications equipment (including cables) should be used no closer to any part of the recommended separation distance calculated
Radiated RF immunity	3 V/m	10 V/m	from the equation applicable to the frequency of the transmitter as below.
(IEC 61000-4-3)	80 MHz - 2.7 GHz 80% @ 1 KHz AM Modulation	80 MHz - 2.7 GHz 80% @ 1 KHz AM Modulation	Recommended sparation distance: d=1.2 √P d=1.2 √P 80 MHz to 800 MHz d=2.3 √P 800 MHz to 2.7 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b.
			Interference may occur in the vicinity of equipment marked with the following symbol:

The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the equipment.

This page intentionally left blank



2305 South 1070 West Salt Lake City, Utah 84119 (800) 748-5355 www.maxtec.com