Deluxe Nasal Pillow CPAP Mask

REF CM036

Intended Use

The nasal pillow mask is intended to provide an interface for Continuous Positive Airway Pressure (CPAP) or bilevel system. The mask to be used by adult patients (>66lb/30kg) for whom positive airway pressure has been prescribed. The mask is for single-patient re-use in the home environment and multi-patient re-use in the hospital/institutional environment.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Excluding patient population

This mask may not be suitable for patients with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex and hiatal hernia. It should not be used if the patient is unresponsive, obtunded, unable to remove the mask, or uncooperative.

Warnings

This mask is designed for use with CPAP or bilevel systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bilevel system is turned on and operating properly. Do not block or try to seal the exhalation port.

CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the device is turned on and functioning properly, new air from the device pushes the exhaled air out through the attached mask exhalation port. However, when the device is not operating, enough fresh air will not be provided through the mask, and exhaled air may be re-breathed. At low pressures, the air flow through the exhalation port may not be enough to clear all of the exhaled gas (CO₂) from the mask. You may breathe in some of the air that you have exhaled.

If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating.

- When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, nasal interface selection, and the leak rate.

Some users may experience skin redness. If this happens, contact your health care professional.

Cleaning the Mask

The cushion, nasal interface, and tubing should be cleaned daily or as needed.

- Hand wash in warm water with a mild liquid dishwashing detergent. Do not use bleach, alcohol, or cleaning solutions containing alcohol. Do not use cleaners containing conditioners or moisturizers.
- Rinse thoroughly. Pat with a towel and air dry completely before use
- Inspect all parts for damage or wear; replace if necessary.


Cleaning the Headgear

Hand wash the headgear weekly or as needed. Lay the headgear flat or line dry.

Disposal

The nasal pillow mask does not contain any hazardous substances and may be disposed of with your normal household refuse.

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