
SUNSET
HEALTHCARE SOLUTIONS

EN english
ES español
FR français

Instruction Manual

Heated CPAP Tube and Adapter for ResMed AirSense 10

REF **TUB06HS10**
TUB06HS10-ADAPT

Sunset Healthcare Solutions
141 W Jackson Blvd Ste 1950
Chicago IL 60604
www.sunsethcs.com

IFU-00148 REV 1.03
2026-08-19

Introduction

The Heated Tube Adapter for ResMed AirSense 10 is used with the Sunset Heated CPAP Tube, which is designed to prevent condensation of humidified air during CPAP treatment, also known as “rainout.”

When combined with a heated humidifier, the Sunset Heated CPAP Tube and Adapter controls humidification by ensuring that air temperature does not cool prior to inhalation. This supports comfortable and consistent treatment.

Table of Contents

1. Understanding the Adapter and Heated CPAP Tube

1.1 Components....1
1.2 Indications for Use....1
1.3 Warnings and Cautions....2-3

2. Use and Care


2.1 Directions for Use....4-5
2.2 Cleaning....6
2.3 Errors and Troubleshooting....6

3. Additional Information

3.1 Operating Conditions....7
3.2 Storage and Transport Conditions....7
3.3 Replacement....7
3.4 Warranty....7
3.5 Product Specifications....8
3.6 Standards Compliance....8
3.7 EMC Statement....8
3.8 Symbol Glossary....9

1. Understanding the Adapter and Heated CPAP Tube

1.1 Components

1. Machine connector cuff
2. Connector lead
3. Squeeze release tabs
4. Mask-end cuff
5. Flexible tubing
6. Adapter for the ResMed AirSense 10

1.2 Indications for Use

The Sunset Heated CPAP Tube and Adapter is a heated wire breathing tube intended to provide warmed and/or humidified breathing gases before they enter a patient's airway. It is indicated for single-patient reuse in the home and in clinical settings, such as hospitals, institutions, sleep laboratories, and sub-acute care facilities. It may be used with non-invasive ventilation for adult patients.

1.3 ⚠ Warnings and Cautions


Please retain these instructions for use for future reference.

- WARNING: Failure to use a mask or accessory that minimizes rebreathing of carbon dioxide or permits spontaneous breathing can cause asphyxiation.
- WARNING: The user is accountable for the compatibility of the Sunset Heated CPAP Tube and the humidifier and all of the parts and accessories used to connect to the patient before use.
- WARNING: Always connect-disconnect the tube from the humidifier by using the grips on the cuff. Do not pull directly on the tube.
- WARNING: To prevent disconnection of the tubing or tubing system during use, only use accessories in compliance with ISO 5367 or ISO 80601-2-74.
- WARNING: When used routinely, check the tube daily for holes, deformation, or any other form of damage. When damaged, replace the tube immediately.
- WARNING: Operating the Sunset Heated CPAP Tube and Adapter close to other heat sources including lamps, electric blankets, etc., can cause elevated airflow temperatures that can affect the quality of the therapy or cause injury to the patient
- WARNING: The tube and adapter shall be replaced after 90 days of use. Using the Sunset Heated CPAP Tube and Adapter after the maximum 90 days duration of use can affect the quality of therapy or cause injury to the patient.
- WARNING: Do not sterilize the Sunset Heated CPAP Tube or Adapter.
- WARNING: The Sunset Heated CPAP Tube and Adapter is for single patient reuse only. Using the Sunset Heated CPAP Tube and Adapter on multiple patients can lead to contamination and injure the patient.
- WARNING: Do not subject the tube to pressure, crushing or other mechanical stresses, as damage to the tube can affect the quality of the therapy or cause injury to the patient.
- WARNING: Replace damaged (having holes, kinks, tears, exposed heating wires, etc. or is not functional) tube immediately. Do not try to repair the tube yourself. Damaged CPAP tubes can affect quality of the therapy or injure the patient.
- WARNING: Using brushes or objects to clean the inside of the tube can damage the tubing and can affect the quality of the therapy or cause injury to the patient. See cleaning instructions.
- WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is required, this equipment and the other equipment should be observed to verify normal operation.

- WARNING: Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and resulting in improper operation.
- WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Sunset Heated CPAP Tube and Adapter, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- WARNING: MR Unsafe and should remain outside the MRI scanner room.
- WARNING: Beware of electrocution. If liquids are spilled into or onto the device, unplug the device and let the parts dry. Always unplug the device before cleaning and make sure that all parts are dry before plugging it back in.
- CAUTION: Dispose of the Heated CPAP Tube and Adapter in accordance with your local government regulations. Upon permanent disposal, mark the tube and adapter so that it cannot be re-used by another patient.
- CAUTION: When using electrical products, always follow applicable safety precautions.
- CAUTION: The Sunset Heated CPAP Tube and Adapter complies with applicable Radiated RF Emission levels to prevent adverse events to the patients and operator due to electromagnetic disturbances. In order to prevent adverse events to patient and operator due to electromagnetic disturbances do not modify or change The Sunset Heated CPAP Tube and Adapter.
- CAUTION: Patients undergoing CPAP therapy may be susceptible to airway infections. Single patient reuse devices that are not replaced within the prescribed frequency may be a potential source of infection. Proper hand hygiene (e.g. hand washing or hand disinfection) is mandatory before and after replacement of the device.
- CAUTION: Do not use (microwave) ovens, fans, heating devices or any other appliances to dry the tube after washing. See cleaning instructions.

2.1 Directions for Use


The Sunset Heated CPAP Tube features a Mask-end Cuff and an Adapter Connector Cuff.

Before connecting, inspect the tube and adapter for any damage. Discard and replace if damaged.

To Connect

Connect the adapter to the ResMed AirSense 10 air port. Align the adapter with the air outlet port on the CPAP machine. Place the Connector Lead into the mating connector located at the top of the air outlet port. Press until the tabs on each side of the connector click into place.

Connect the tube's Adapter Connector Cuff to the Adapter.


Connect Mask-end Cuff of the tube to the connector on your CPAP mask. For proper placement, please consult instructions included with your mask. Ensure there is airflow through the tubing.

Apply mask and use as directed by your physician or respiratory therapist, and in accordance with the instructions in your CPAP or mask ventilator system user manual.

For extended directions for use, please refer to your ResMed CPAP machine manual.

To Disconnect

Lightly squeeze Release Tabs on side of tube Connector Cuff and carefully disconnect from the adapter. To disconnect the Adapter from the CPAP Machine, lightly squeeze Release Tabs on side of Adapter and carefully remove.

⚠ Caution: To avoid circuit leaks, do not pull or stretch tubing. Always handle tubing by the cuff.

2.2 Cleaning

Clean Sunset Heated CPAP Tube and Adapter before first use and daily during use.

1. Disconnect Heated CPAP Tube from the Adapter and mask. Disconnect the Adapter from the CPAP machine.
2. Thoroughly rinse inside and outside of Heated CPAP Tube and Adapter with cold tap water to remove excess soil.
3. Gently wash Heated CPAP Tube and Adapter by immersing in a solution of warm water and mild unscented liquid dish soap, gently agitating by hand and using a soft bristle brush to remove contaminants.
4. Thoroughly rinse Heated CPAP Tube and Adapter with clean water to remove all soap residue.
5. With connector cuffs directed down, toward the floor, gently tap connector cuffs to remove excess water.
6. Hang Heated CPAP Tube with connector cuffs directed downward, towards the floor, for drying.
7. Allow Heated CPAP Tube and Adapter to air dry thoroughly prior to next use.
8. Visually inspect Heated CPAP Tube and Adapter for cleanliness. Repeat cleaning process if not visually clean.
9. Inspect Heated CPAP Tube for damage or wear. Discard and replace if necessary.

⚠ Caution: Do not use harsh soaps, detergents, solvents, or alcohols to clean the Sunset Heated CPAP Tube and Adapter. Do not attempt to clean Heated CPAP Tube and Adapter in a dishwasher, ultrasonic cleaner or any other cleaning device.

Only reuse the Heated CPAP Tube and Adapter when fully dry.

2.3 Errors and Troubleshooting

Consult your CPAP machine user manual.

3.1 Environmental Operating Conditions

- Ambient Temperature: 41°F – 104°F (5°C – 40°C)
- Relative humidity: 5% – 95%
- Altitude: Air pressure range 700 hPa to 1060 hPa

3.2 Storage and Transport Conditions

Environmental conditions for transport and storage:

This device should be transported and stored in its original packaging and protected from direct sunlight, dust, or any other form of contamination. . Temperature during storage and transport should not exceed 140°F (60°C) or fall below -13°F (-25°C) with Relative humidity can range from 5 % to 95%.

3.3 Replacement

The Heated CPAP Tube and Adapter are intended to be replaced every 3 months, or according to your provider's recommended schedule. If damaged, leaking or malfunctioning, replace your heated tube or adapter immediately.

3.4 Warranty

The Sunset Heated CPAP Tube warranty is 30 days under normal operating conditions. This is based on daily usage of 8 hours.

3.5 Product Specifications

Description	Specification
Sunset Heated CPAP Tube Length	6ft (1.83m)
Internal Diameter	15mm
Flow resistance @60 l min-1	≤ 0.1 kPa
Tubing Compliance Standards	ISO 5367:2014 ISO 80601-2-70:2015 ISO 80601-2-74:2017
Tube Ingress Protection	IP22: Protected from touch by fingers and objects greater than 12 millimeters. Protected from water spray less than 15 degrees from vertical
Compatible CPAP Machine	ResMed Airsense 10
Periodic visual safety inspections that should be performed by the OPERATOR	Daily
HUMIDIFIER classification according to the ISO 80601-2-74	Humidifier classification: 2

3.6 Standards Compliance

These devices are designed to conform to the following standards:

- IEC 60601-1 General Requirements for Basic Safety and Essential Performance of Medical Electrical Equipment
- EN 60601-1-2 Electromagnetic Compatibility
- ISO 80601-2-70 Sleep Apnea Breathing Therapy Equipment
- ISO 80601-2-74 Respiratory humidifying equipment

IEC 60601-1 Classification

- Degree of Protection Against Electric Shock: Type BF Applied Part
- Degree of Protection against Ingress of Water: Drip Proof, IPx2
- Mode of Operation: Continuous

Electrical

- Power Consumption 3.3 Vdc, max. 5A.

3.7 EMC Statement

The Sunset Heated CPAP Tube and Adapter comply with all applicable electromagnetic compatibility requirements (EMC) according to ISO 60601-1-2.

3.8 Symbol Glossary

Symbol	Meaning
	Attention, consult accompanying documents
	Read instruction manual carefully
	Single patient - multiple use
	Catalog item number
	Batch code
	Name and address of manufacturer
	Date of manufacture
	Not made with natural latex rubber
	Not made with phthalates (PHT/DEHP)
	Keep away from direct sunlight
	Keep dry
	Storage and transportation temperature range
	Storage and transportation humidity range
	Do not use if the package is damaged
	MR Unsafe and should remain outside the MRI scanner room
	Type BF applied part

Standard References

ISO 15223-1:2012, Medical Symbols to be used with medical device labels, labeling, and information to be supplied - General requirements

IEC 60601-1, Medical electrical equipment - General requirements for basic safety and essential performance

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