



# SUNSET

HEALTHCARE SOLUTIONS



## DEVILBISS

### DeVilbiss PAP Devices & Portable Oxygen Concentrator

- Less than 1% failure rate and backed by the best-in-industry, 3-year warranty
- Patented auto on/off automatically starts flow of air with patient breath
- Whisper-quiet operation (26 dBA)
- Eight-foot power cord provides greater positioning flexibility – no bulky power brick
- Stacked design creates an exceptionally small footprint
- Reminders alert patient to replace or check parts
- Made in the USA
- Easy and intuitive setup
- Includes attractive travel bag

#### IntelliPAP® With Heated Humidification

Order#

<b>CAPDV51D-HH</b>	Standard CPAP
<b>CAPDV53D-HH</b>	Standard Plus CPAP
<b>CAPDV54D-HH</b>	AutoAdjust CPAP
<b>CAPDV55D-HH</b>	Bilevel S
<b>CAPDV57D-HH</b>	Auto Bilevel

- SmartLink modem available on all models
- Also available without heated humidifier

Visit [www.sunsethcs.com/intellipap](http://www.sunsethcs.com/intellipap) for more info on each device

CALL 877-578-6738  
to place an order today



#### iGo Portable Oxygen Concentrator

Order#

<b>RES306DS-A</b>	With Wheeled Carry Case
<b>RES306D-413</b>	Spare Battery

Visit [www.sunsethcs.com/igo](http://www.sunsethcs.com/igo) for more info on the iGo

- FAA Approved for air travel
- External battery charger
- Battery lasts as long as 5 hours
- Built in oxygen sensing device to ensure accurate oxygen delivery and reduced periodic maintenance schedule
- No maintenance needed for three years
- Integrated PulseDose oxygen conserving technology extends use time, offering increased mobility with improved comfort and increased efficiency
- If the unit is operating on AC power and the power is interrupted, then the iGo will automatically switch to battery operation

©2013 Sunset Healthcare Solutions Inc

SUNSET HEALTHCARE SOLUTIONS - CHICAGO, IL

P 312.997.9980 P 877.578.6738 F 312.997.9985 W [SUNSETHCS.COM](http://SUNSETHCS.COM)