





# DELUXE CPAP MASK SIZE GUIDE

877-578-6738  
www.SunsetHCS.com


LARGE




MEDIUM



NASAL MASK  
CM006




Mouth should be slightly open,  
as if breathing.



FULL FACE MASK  
CM007

**\*IMPORTANT\*** Print without page scaling (print at 100% scale). Turn off any auto scaling or shrink to fit options. Printed area should measure 8 inches by 10.5 inches. Cut inside blue and red guidelines. Place sheet over face.

LARGE



MEDIUM



SMALL

