

SoClean

Compatibilty Chart

Most CPAP machines with a standard hose will work with the SoClean right out of the box. Some CPAP machines may need a specific adapter, particularly units with a heated hose. To see whether you need to have a special adapter to use the SoClean with your CPAP, please refer to our guide below.



| | | N/A = Not Applicable | ✓ = Out of Box | ⊕ = Need Adapter | X = Not Compatible |
|---------------------|---------------------------------|----------------------|-----------------|------------------|--------------------|
| | | PAP | WITH HUMIDIFIER | WITH HEATED HOSE | |
| Aeiomed | | | | | |
| | Everest | ✓ | ✓ | N/A | |
| Apex Medical | | | | | |
| | ICH | ✓ | ✓ | N/A | |
| | XT | ✓ | ✓ | N/A | |
| Carefusion | | | | | |
| | PureSom | ✓ | ✓ | N/A | |
| Devilbiss | | | | | |
| | IntelliPAP | ✓ | ✓ | N/A | |
| Evo | | | | | |
| | ComfortPAP | ✓ | ✓ | N/A | |
| Fisher & Paykel | | | | | |
| | Icon | ⊕ | ⊕ | ⊕ | |
| | SleepStyle 200 Series | ✓ | ✓ | N/A | |
| | SleepStyle 600 Series | ⊕ | ⊕ | ⊕ | |
| Philips Respironics | | | | | |
| | Dreamstation | ✓ | ✓ | ⊕ | |
| | PR System One REMstar 60 Series | ✓ | ✓ | ⊕ | |
| | PR System One REMstar | ✓ | ✓ | ⊕ | |
| | SleepEasy | ✓ | ✓ | N/A | |
| | SleepEasy II | ✓ | ✓ | N/A | |
| | M-Series V 2 | ✓ | ✓ | N/A | |
| | M-Series V 3 | ✓ | ✓ | N/A | |
| | Remstar | ✓ | ✓ | N/A | |
| PMI Probasic | | | | | |
| | Zzz-PAP | ✓ | ✓ | N/A | |
| Puritan Bennett | | | | | |
| | GoodKnight 420 | ✓ | ✓ | N/A | |
| | Sandman | ✓ | ✓ | N/A | |
| RESmart | | | | | |
| | 600 | ✓ | ✓ | N/A | |
| | 700 | ✓ | ✓ | N/A | |
| ResMed | | | | | |
| | S8 | ✓ | ✓ | N/A | |
| | S9 | ⊕ | ⊕ | ⊕ | |
| | AirSense I0 | ⊕ | ⊕ | ⊕ | |
| Transcend | | | | | |
| | Transcend II | ✓ | (HME) X | N/A | |