



# SUNSET

HEALTHCARE SOLUTIONS

**ADVANCED**

## FINGER PULSE OXIMETER

Order# RES5101



- FDA approved
- 2-Year Warranty
- 6 display modes
- Dual color OLED display
- SpO2, PR, pulse bar and waveform
- Adult and pediatric use
- Light and compact
- Easy and convenient to operate and carry
- Ergonomic design
- Low power consumption (more than 30 hours of operation on 2 AAA batteries - included)

### SpO2

Measurement range  
70-99%

Resolution  
1%

Measurement accuracy  
80%-99%;  $\pm 2\%$   
70%-79 %;  $\pm 3\%$   
0%-69%;unspecified

### Pulse Rate

Measurement range  
30-235 bpm

Resolution  
1bpm

Measurement accuracy  
30- 99;  $\pm 2$  bpm  
101- 235;  $\pm 2\%$

### Environmental

Operation Temperature  
41°F - 104°F

Storage Temperature  
-4°F - 131°F

Operation Humidity  
RH less than or equal to 80%

Storage Humidity  
RH less than or equal to 93%

### Dimensions

2.28in (W) x 1.25in (H) x 1.33in(D)

### Weight

1oz (without batteries)

Device works normally when mixed noise produced by  
BIO-TEK INDEX Pulse Oximeter tester.

CALL 877-578-6738

© 2012 Sunset Healthcare Solutions

SUNSET HEALTHCARE SOLUTIONS - CHICAGO, IL

P 312.997.9980 P 877.578.6738 F 312.997.9985 W SUNSETHCS.COM