



# SUNSET

HEALTHCARE SOLUTIONS

**IMPROVED**



## BUBBLE HUMIDIFIER

Order# **RES050**  
6psi  
20/case

- Disposable, durable plastic
- Unique diffuser designed to provide patient with comfortable, therapeutic humidity and minimal noise
- Defined maximum and minimum water level lines
- Recessed nipple reduces possibility of breaking off outlet connection
- Holds 500 mL
- Preset audible alarm at 6 psi pressure relief
- Allows for use of longer lengths of oxygen tubing without the audible alarm sounding

Visit [www.SunsetHCS.com](http://www.SunsetHCS.com) or call 877-578-6738 to see our full catalog or place an order.