

## **O2 TUBING & CANNULAS**



## 6-channel safety tubing is kink-free



- All tubes have 6 channels for added strength and durability
- The material is soft, memory free and lightweight for longterm use.
- Green tubing is always the same price as clear —
  add a "G" to the end of the order# to purchase green

## Order#

**RES3007** 7ft Oxygen Tube – 50/case

RES3015 15ft Oxygen Tube - 50/case

RES3025 25ft Oxygen Tube – 25/case

RES3040 40ft Oxygen Tube - 20/case

RES3050 50ft Oxygen Tube - 20/case

Visit www.SunsetHCS.com to see our full line of tubes cannulas

© 2011 Sunset Healthcare Solutions Inc

- Sunset's oxygen nasal cannulas are tested and proven to provide maximum safety and comfort.
- All tubes are kink-free to ensure uninterrupted airflow.
- The cannulas are contoured to match the shape of the patient's face to help keep it in place.
- The nasal prongs are curved and tapered for comfort that feels natural.
- The material is soft, memory-free and lightweight for long-term use.

## Order#

RESI 100 Adult Cannula - 50/case

RES I 104 Adult Cannula with 4ft Tube - 50/case

RES 1 107 Adult Cannula with 7ft Tube— 50/case

RESIII5 Adult Cannula with 15ft Tube - 50/case

RESI125 Adult Cannula with 25ft Tube - 25/case

Call 877-578-6738 to place an order today