

SUNSET

HEALTHCARE SOLUTIONS

OXYGEN REGULATORS



CGA 870 Style

Order# RES042
0 to 8 L/Min Barbed Outlet

Order# RES043
0 to 15 L/Min Barbed Outlet

CGA 540 Style

Order# RES044
0 to 8 L/Min DISS Outlet

Order# RES045
0 to 15 L/Min DISS Outlet

Click reg with mini type oxygen regulator inlet
Adjustable flow 0-8 Lpm, 0-15 Lpm
Heavy duty brass core
Durable anodized aluminum body
Input pressure 2200psi
Output pressure 50psi

© 2012 Sunset Healthcare Solutions