

# CPAP Batteries

- No additional adapters needed
- Outlet functions like a wall socket
- Universal 110-240V recharging conversion
- 1-Year Warranty

## Zopec Explore

- Run any electronics with power requirement of 120 watts or under
- AutoSwitch and Bypass Technology for uninterrupted sleep when power outages occur
- Unlimited linking of multiple batteries together for longer use before recharging



Zopec Explore 5700  
Item# **CAP1024-2157**

Zopec Explore 8200  
Item# **CAP1024-2182**



Zopec Power Outlet 2  
Item# **CAP1024-21ZDP02**

- Includes 45W USB-C



Portable Outlet  
Item# **CAP1008**

Item#	<b>CAP1008</b>	<b>CAP1024-21ZDP02</b>	<b>CAP1024-2157</b>	<b>CAP1024-2182</b>
<b>Name</b>	Portable Outlet	Power Outlet2	Explore 5700	Explore 8200
<b>Manufacturer</b>	Portable Outlet	Zopec Medical	Zopec Medical	Zopec Medical
<b>OEM #</b>	155W	21ZDP02	215700C	218200C
<b>Dimensions (L x W x H)</b>	7.0" x 5.0" x 1.5"	7.0" x 5.0" x 1.5"	8.5" x 6.0" x 1.0"	7.75" x 5.5" x 1.75"
<b>Volume (cubic inches)</b>	53	53	51	74.5
<b>Weight</b>	2.0 lb	2.0 lb	2.5 lb	3.5 lb
<b>Capacity</b>	42,000 mAh	40,000 mAh	54,000mAh	80,000 mAh
<b>Watt Hours (Wh)</b>	155Wh	146Wh	160Wh	288Wh
<b>AC Output</b>	100W	100W	120W	120W
<b>UPS Function (Autoswitch)</b>	No	No	Yes	Yes
<b>USB Type C</b>	No	PD 45W Smart Hub	15W (5V/3A)	15W (5V/3A)
<b>USB Type A</b>	5V/2.1A	5V/2.1A	5V/3.1A	5V/3.1A
<b>Battery Cell</b>	LG	LG	Panasonic	Panasonic
<b>Capacity Indicator</b>	4 lights	4 lights	Digital Display	Digital Display
<b>Unlimited Linking Technology</b>	No	No	Yes	Yes
<b>Recharge Time</b>	Less than 4 Hours	Less than 4 Hours	Less than 4 Hours	Less than 5 Hours
<b>Car Charger Included</b>	No	No	Yes (45W)	Yes (45W)
<b>Optional Solar Panel Charger</b>	No	Yes	Yes	Yes
<b>Air Travel</b>	Yes	Yes	Yes	No
<b>Charge Cycles</b>	Over 500 Cycles	Over 500 Cycles	Over 500 Cycles	Over 500 Cycles
<b>Warranty</b>	1 Year	1 Year	1 Year	1 Year