



**CPAP BATTERY**

*Estimated Machine Run Times*

Machine	Pressure	Hours (Estimated)	Machine	Pressure	Hours (Estimated)
Respironics M Series and REMstar	6	20	ResMed S8 and S9	6	17
	7	18		7	15
	8	17		8	14
	9	16		9	13
	10	14		10	12
	11	13		11	11
	12	11		12	10
	13	10		13	9
	14	9		14	9
	15	8		15	8
PR System One	6	24			
	7	19			
	8	16			
	9	14			
	10	12			
	11	11			
	12	11			
	13	10			
	14	10			
	15	9			

©2014 Sunset Healthcare Solutions Inc