




Acute Care Product Catalog



New!

See Anestand
on pages 4-5



POLICIES

Ordering Options

Three Easy Ways To Order

Phone

Toll free (877) 578-6738

Our customer service professionals are here to assist you M-F 8:00 AM – 5:00 PM Central Time.

Fax

312-997-9985 — Accepted 24 hours a day.

Email

sales@sunsethcs.com

Shipping

Drop shipping available — ask your sales rep about sending your order to multiple locations.

UPS

Whenever possible, orders ship via UPS. If another carrier is desired, please contact your sales rep.

All orders are packed in the smallest boxes available and combined with other products on the order to reduce shipping costs.

All products ship prepaid and are added to the invoice.

All products ship from Bloomingdale, IL or Oxnard, CA. Please refer to UPS.com for transit times.

Damaged Merchandise

Call us within 10 days. Keep original cartons for UPS inspection.

We proudly partner with



Payment Methods

Payment Portal and ACH

Visit www.sunsethcs.com/paymentportal to pay online with ACH or credit card.

Net 30-Day Payment Terms

For qualified companies. A 1.5% monthly finance charge will be assessed on all past due invoices.

COD

For non-qualified companies.

Credit Card

We accept all major credit cards. A processing fee may apply. After 10 days of invoice date, a surcharge of 2.5% will be applied.



Customer Service

Returns

Call us for a return authorization number. There is no restocking charge on orders returned within 20 days of invoice. Orders returned after 20 days may be assessed a restocking fee of 15%. No returns after 90 days. Toll free (833) 693-0490

Central Distribution Center

Sunset Healthcare Solutions
279 Madsen Dr Ste 101
Bloomingdale, IL 60108

Western Distribution Center

Sunset Healthcare Solutions
2450 Statham Blvd
Oxnard, CA 93033

Problems

Please call us immediately if you have any complaints (it is helpful to ask for the customer service representative with whom you placed your order).

Or email us at: customerservice@sunsethcs.com

Questions?

Our experienced customer service staff will be glad to provide you with information. Just give us a call.

Pricing

All prices are in U.S. Dollars.

Taxes are added where applicable.

We make every effort to hold pricing throughout the year. However, prices are subject to change without notice.

Contact us for more information on fulfillment and compliance solutions that best fit your needs.

CONTENT

Anestand

Anesthesia workstation

Page

4-5

PPE

Masks and Face Shields

6

Isolation Gowns

6

Sanitizers

7

Non-Contact Thermometers

7

Respiratory

Ventilator Circuits

8

Ventilator Connectors, Adapters, and Filters

9

Nebulizers

10

Aerosol

11

Medication Swallowing Aid

11

Suction

12

Tracheostomy

13

Oxygen

Firesafe

14

Masks

15

Pulse Oximeters

15

Cannulas

16

Tubing

17

Connectors

18

Bacteria Filters

19

Accessories

19

CPAP/NIV

Masks and Replacement Cushions

Page

20-22

Tubing

23

Tubing Connectors and Adapters

23

Chinstraps

24

Cleaners

25

Comfort Care

25



Customized anesthesia workstation for every anesthetizing location

Item#

ANE77006

Anestand Anesthesia Stand, Complete Tray and Clamp

ANE77001 Replacement Tray

ANE77005 Replacement Clamp

Benefits to Providers


- Keep medications near IV port
- Ergonomical work space
- Makes wiping IV port easy
- Extra set of hands when help not around
- Helps make RSI safe
- Creates valuable space for non-OR anesthesia

Value to Hospital

- Helps with regulatory compliance
- Sharps management
- Could help reduce HAIs
- Demonstrates support for anesthesia team
- Reusable/single fixed cost

1. The recessed tray holds and organizes supplies for induction of anesthesia and most anesthetic procedures
2. Clips keep IV tubing securely in place, so it's easy to "scrub the hub"
3. Corner grips for suction tubing and more
4. The adjustable support post lets you position it just where you need it
5. The anesthesia circuit clips keep your circuit fixed and firmly in place
6. Patented clamp attaches to OR bed, IV pole, or inpatient bed

**Anestand can help health systems save lives,
lower costs and improve compliance.
All at the same time.**



The first and only anesthesia stand to safely and securely keep all supplies where they're needed can also be one of your best weapons in combatting infections.

Anestand Facts

- A small investment could help in compliance, save time and quickly pay for itself.
- By having supplies and medications organized and easily available with IV port ready, Anestand aids in patient safety.
- The bottom line: Anestand could help prevent costly infections and save money by helping anesthesia providers save OR time.

Of all the many safety and security features Anestand offers, perhaps the single greatest is how it helps health systems to lower HAIs.

- Having IV port ready for site alcohol scrub makes disinfection (and compliance) easier — Anestand keeps your IV port ready for scrubbing.
- By keeping medications immediately available and infusion ports clearly defined Anestand provides additional safety.
- By keeping supplies immediately available, Anestand helps make RSI (rapid sequence inductions) safe.
- During COVID, Anestand assists anesthesia providers as an extra pair of hands, while other healthcare personnel are being asked to keep distance from the airway.

Anestand can provide value in safety, compliance, and convenience. As soon as you use it, you'll wonder how you ever managed without it.

- CDC, JCAHO, SHEA, and other organizations require cleaning the port with every injection. Anestand helps your facility and providers comply with these regulations.
- Patients with an HAI stay 19 additional days, or \$11.2 million in HAI treatment and/or lost revenue. For a small fixed cost your facility could get a large return on investment.
- Keeping supplies organized on a tray may prevent dropped items and reduce replacement costs. Anestand makes it far easier to keep supplies organized and right at hand.



PPE

KN95 Mask

Item#
PPE1001KU
Box of 30 Masks

- 5 layers



Face Shield

Item#
PPE3001FS-F
Disposable Face Shield
with Elastic Strap
Each


- High-density foam forehead rest
- Clear, anti-fog vision



Disposable Level 2 Surgical Face Mask

Item#
PPE1002SM
Box of 50 Masks

- Medical classification and 510(k)
- 98% bacterial filtration efficiency
- Rated for moderate risk levels of aerosols, spray, and/or fluids



Disposable 3-Ply Face Mask

Item#
PPE1002PS
Adult • 6 7/8 x 3 3/4 in
Box of 50 Masks

PPE1003PPB
Pediatric • 5 5/8 x 3 3/4 in
Box of 10 Masks

- Hypoallergenic
- Comfortable
- Breathable
- Fiberglass-free
- 3 layers of polypropylene (a meltblown layer between two soft non-woven layers)



Disposable Tearaway Gown

Item#
PPE5002
25/pack

- Level 3 ASTM F3352
- Individually sealed
- Made in the USA



Disposable Gown

Item#
PPE5003L
100/case

- Meets the standard performance requirement of AAMI PB70 Level 1 for isolation gowns, including the seam of the neckline and cuff.
- 20 bags of 5 gowns per case



PPE

Alcohol Wipes

Item#
CAP1003A
Each (Pack of 50 Wipes)

- 75% Isopropyl Alcohol



Hand Sanitizer

Item#
CAP1022-16
16.9oz Pump Bottle
Each

- Kills 99.9% of Germs
- 70% Ethyl Alcohol



Mini Pistol Grip Non-Contact Thermometer

Item#
PPE4001-J
Infrared
Each

- Forehead mode for all ages
- Object measurement mode
- HD LCD display
- Battery: AAAx2
- Fever warning: 6 short double beeps
- Memory storage: 20 sets of measureable data
- Fahrenheit and Celsius
- Automatic shutdown: 10 seconds
- Professional accuracy (plus or minus 0.6 degrees F when measuring forehead temperatures from 89.6 through 107.9 degrees F)



CleanSmart Disinfectant

Item#
CAP1019-C
23oz Spray Bottle
6/case

- CleanSmart is EPA listed as a broad-spectrum and hospital disinfectant.
- EPA Registration Number 89896-2
- CleanSmart is approved for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment.
- CleanSmart is included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2, the virus that causes COVID-19.
- Active ingredient Hypochlorous Acid is safe and breaks down to saline solution.



CleanSmart Alcohol-Free Hand Cleanser

Item#
CAP1021-HC8
8oz Spray Bottle
6/pack

CAP1021-HC2
2oz Spray Bottle
9/pack

- Same active ingredient as CleanSmart cleaner, in a hand cleanser.



Mini Wand Style Non-Contact Thermometer

Item#
PPE4002-A
Infrared
Each

- One second reading
- One click operation
- Fahrenheit and Celsius
- ¾ to 2-inch measuring distance



Visit www.sunsethcs.com

SUNSET
HEALTHCARE SOLUTIONS

VENTILATOR CIRCUITS

Passive Heated Breathing Circuit

Item# **VC1114**

For use with
MR850 Humidifier
10/case



Passive Breathing Circuit

Item# **VC1100**

With fixed exhalation port
For use with MR850
Humidifier
10/case




Breathing Circuits


For use with MR850 Humidifier • 22mm cuffs
30/case

Item#	Mode	Capacity	Size	Tube Diameter
VC1001	Active	Non-heated	Adult	22mm
VC1001P	Active	Non-heated	Pediatric	15mm
VC1002	Active	Heated	Adult	22mm
VC1002P	Active	Heated	Pediatric	15mm
VC1003	Passive	Heated	Adult	22mm
VC1003P	Passive	Heated	Pediatric	15mm
VC1004	Passive	Non-heated	Adult	22mm
VC1004P	Passive	Non-heated	Pediatric	15mm


VC1001
VC1001P




VC1002
VC1002P



VC1003
VC1003P



VC1004
VC1004P



VENTILATOR FILTERS

Bacterial/Viral HME Filter — Adult

Item# **BF104**

15 mm ID x 15 mm ID/22 mm OD

Each

HCPCS# A7507



- Bacterial efficiency 99.9999+%
- Viral efficiency 99.99+%
- Gas sampling port
- 35.5mg/L @ 500ml VT

Bacteria Filter

Item# **BF2200**

22mm Bacteria
Disk Filter

5/pack



VENTILATOR CONNECTORS AND ADAPTERS

Valved T-Adapter — Adult

Item# **VEN001**

22F-18F-22M/15F

Each



Swivel Connector Silencer

Item# **VEN002**

For Ventilator Circuits
and CPAP

Each



Fixed Exhalation Port

Item# **VEN003**

22mm OD x 15mm ID

10/pack



Circuit Connectors

Item# **RES013**

15mm ID x 22mm OD

10/pack

Item# **RES014**

Aerosol Connector

22mm OD x 22mm OD with Lip

10/pack



VENTILATOR CONNECTORS, ADAPTERS, AND ANESTHESIA

NEBULIZERS

Disposable Nebulizer Kit with T-Piece


Item# **RES091**

Disposable
Single Patient
50/case

HCPCS# A7003

Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Reusable Nebulizer Kit with T-Piece

Item# **RES093**

Reusable
Single Patient
10/case

HCPCS# A7005

Includes:

- Jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

Disposable Nebulizer Kit with Mask

Item# **RES090**

Adult Size • **50/case**

Item# **RES092**

Pediatric Size • **50/case**

HCPCS# A7003



Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble

Compressor Nebulizer

Item#

NEB100

- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size:
respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM



AEROSOL

Aerosol Mask

Item# **RES2110**
Adult • 50/case

Item# **RES2210**
Pediatric • 50/case

HCPCS# A7015



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

Aerosol Tube

Item# **TUB100-6**
Corrugated Aerosol Tubing
100ft in 6in segments



Drainage Bag with T-Adapter



Item# **RES070**
750cc • 50/case

- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

Drainage Bag with Y-Adapter



Item# **RES071**
750cc • 50/case

SperaChamber

Item# **RES095**
Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler
Each

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



MEDICATION SWALLOWING AID

Phazix

Item#
RES6109
2.5oz Tube
6/case



- A swallowing gel designed to help those with difficulty swallowing tablets, capsules, and powders
- A clear, moderately thick, slippery gel that is placed on a spoon with the medication and consumed in a single swallow
- Lubricates the mouth and throat to smoothly carry medication to the stomach
- Has an all natural, pleasant vanilla flavor and aroma to mask any unpleasant taste or smell of medications
- Does not impact absorption rates and breaks down quickly in the stomach to maintain the therapeutic effect of the medication

Visit www.sunsethcs.com



SUCTION

Suction Kit

Item# **RES026**

RES025 Suction Tubing 1/4in x 6ft
 RES024 Suction Tubing Connector 3.5in
 BF3812 Bacteria Filter
 RES023 Suction Canister

HCPCS#
 CANISTER - A7000
 TUBING - A7002



Suction Canister

Item#

RES023

800cc • Single Use

10/case

HCPCS# A7000



- Float shut off valve for protection against overflow and premature shut-off
- Leak-free seal to minimize vacuum loss

Suction Tubing

1/4in

Length	Item#	Qty
6ft	RES025	10/case



Suction Tube Connector

10/pack

Size	Item#
3.5in	RES024



Rubber Suction Elbow Connector

Item#

RES028

6/pack



Suction Machine Inline Filters

Item# **BF1312**

1/3in to 1/2in
1/pack



Item# **BF3812**

3/8in to 1/2in
10/pack



Item# **BF516TS**

5/16in Barbed with Threaded NPT
10/pack



TRACHEOSTOMY

Tracheostomy Mask

Item#
RES2130
Adult • With O₂ Connector
50/case

RES2230
Pediatric • With O₂ Connector
50/case

HCPCS# A7525



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing

Tracheostomy Swivel Elbow Connectors

15mm ID x 15mm OD
10/pack



Heat Moisture Exchangers



RES027
T-HME



RES027OP
T-HME with
Oxygen Port



RES027IL
Inline HME —
Large Volume



RES027NE
Inline HME —
Neonatal

TV Range (mL)	250-1000	250-1000	15-50	15-50
Moisture Output (mg H ₂ O/L)	27.3, Vt = 0.5L	26, Vt = 0.5L	31, Vt = 0.1L	30, Vt = 0.2L
Weight (g)	5	6.5	10	5
Resistance (cm H ₂ O)	0.7@60 LPM	0.4@30 LPM	0.4@10 LPM	0.5@10 LPM
Dead Space (mL)	12.5	12	13	2.5
Connection	15F	15mmF, 5.5mm Oxygen port	15M, 22M/15F	15M, 15F

Flex HME

Item#
RES027FLEX
40/case
Foam
22mm Male
15mm Female



TV Range (mL)	300-1500
Moisture Output (mg H ₂ O/L)	33.8 mg/L at Vt = 0.5L
Weight (g)	17
Resistance (cm H ₂ O)	0.50@30 LPM, 1.34@60 LPM
Dead Space (mL)	90

Adapter for T-HME

Item#
RES027A
Each



OXYGEN PRODUCTS

FIRESAFE

Firesafe Thermal Shut-Off Valve

Item#
RES010
Each


- Bidirectional



- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire

FireSafe Tubing Connector

Item#
RES010TC 100/case



MASKS

Simple O₂ Mask



- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

With 7ft Safety Tubing

Item# **RES2100**
Adult Size
50/case

Item# **RES2200V**
Pediatric Size
50/case

Item# **RES2300V**
Infant Size
50/case



HCPCS# A4620

Partial Rebreather O₂ Mask



- 5-15 LPM
- High concentration

With 7ft Safety Tubing & Reservoir Bag

Item# **RES2101V**
Adult Size
50/case

Item# **RES2201V**
Pediatric Size
50/case



Item# **RES2301V**
Infant Size
50/case

HCPCS# A4620

Non-Rebreather O₂ Mask



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

With 7ft Safety Tubing & Reservoir Bag

Item# **RES2102**
Adult Size
50/case

Item# **RES2202V**
Pediatric Size
50/case



HCPCS# A4620

PULSE OXIMETERS

Advanced Finger Pulse Oximeter

Item# **RES5101**
Sunset
Each

Item# **RES5101V2**
Generic
100/case



- 6 display modes
- Dual color OLED display
- SpO₂, PR, pulse bar and waveform

Both Pulse Oximeters

- 2-year warranty
- Ergonomic design
- Adult and pediatric use
- Low power consumption — automatically powers off
- Easy and convenient to operate and carry
- More than 30 hours of operation on 2 AAA batteries — included
- Light and compact

Standard Finger Pulse Oximeter

Item# **RES5100**
Sunset
Each



- 1 display mode
- Digital LED display
- SpO₂, PR, and pulse bar

OXYGEN CANNULAS

Cannula with Tubing

Adult Cannula with Tubing

Length	Item#	Qty
4ft	RES1104	50/case
7ft	RES1107	50/case
15ft	RES1115	50/case
25ft	RES1125	25/case

HCPCS# A4615



- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort

Item#

RES1207V

Pediatric Cannula
with 7ft Supply Tube
50/case



RES1307V

Infant Cannula
with 7ft Supply Tube
50/case



Cannula without Tubing

Item# **RES11000**

Adult Cannula without
Supply Tubing
50/case

HCPCS# A4615


HCPCS# A4615



Dual-Lumen (Demand) Cannula

Length	Item#	Qty
4ft	RES4144V	25/case
7ft	RES4147V	25/case

HCPCS# A4615




- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

Secure Cannula Stickers

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



OXYGEN TUBING

Oxygen Supply Tube

Length	Item#	Qty
7ft	RES3007	50/case
15ft	RES3015	50/case
25ft	RES3025	25/case
40ft	RES3040	20/case
50ft	RES3050	20/case

HCPCS# A4616



Oxygen Supply Tube — Green

Length	Item#	Qty
7ft	RES3007G	50/case
15ft	RES3015G	50/case
25ft	RES3025G	25/case
40ft	RES3040G	20/case
50ft	RES3050G	20/case

HCPCS# A4616



- Crush resistant with low memory to prevent kinking
- Great durability and comfort level
- Latex-free
- 1-6 LPM
- 6 channels

OXYGEN CONNECTORS AND ACCESSORIES

Christmas Tree Connector

Item# **RES001**
Green without Swivel
10/pack



Item# **RES002**
Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Y-Connector

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack



Tubing Reducer Connector

Item# **RES011**
Fits 3/16in, 1/4in,
5/16in, 3/8in, and 7/16in
ID Tubing
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

Oxygen Tubing Connector — Non-Swivel

Item# **RES012**
Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Oxygen Tubing Connector — Swivel

Item# **RES018**
10/pack



- Enhances patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

Oxygen Tubing Connector — Swivel, Female

Item#
RES018KIT-1
Each



Brass Washer with Rubber Ring

Item#
RES035
25/pack




OXYGEN BACTERIA FILTERS AND ACCESSORIES

Final Bacteria Filter




Item# **BF516**
10/pack
 5/16in Barbed
 on Both Sides

Final Bacteria Filter — 1/4in to 3/8in



Item# **BF1438**
10/pack
 1/4in to 3/8in
 Stepped Barb
 on Both Sides

Final Bacteria Filter — 1/4in to 1/2in



Item# **BF1412**
10/pack
 1/4in to 1/2in Stepped
 Barb on Both Sides

Final Bacteria Filter — 3/16in



Item# **BF316**
10/pack
 3/16in Barbed
 on Both Sides

Oxygen Liter Meter

Item# **RES034**
 0-10 Liter/Minute
Each



- Oxygen pen
- Works with tanks and concentrators

Moisture Therapy

Item# **RES6104L**
 1oz Tube • **18/case**

- Apply before beginning oxygen therapy and as needed during treatment
- Moisturizes the nose, lip, and over the ear where friction may cause discomfort from cannula
- Petroleum-free with aloe vera and vitamins A and E



SUNSET CPAP MASKS

Sunset Deluxe Full Face CPAP Masks

Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	CM005S
Medium	CM005M
Large	CM005L

Unobtrusive design offers a clear line of sight



CM SIZE01

Deluxe Sizing Guide

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031

Mask with Adjustable Forehead T, Headgear, and Replaceable Cushion

Size	Item#
Small	CM105S
Medium	CM105M
Large	CM105L



Replacement Cushions

Size	Item#
Small	CUCM005S
Medium	CUCM005M
Large	CUCM005L

- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

Sunset Deluxe Nasal CPAP Masks

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM006S
Medium	CM006M
Large	CM006L

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



HG003 Headgear - Fits Most
HG003XL Headgear - XL

CM SIZE01

Deluxe Sizing Guide

HCPCS#
 MASK - A7034
 HEADGEAR - A7035
 CUSHION - A7032

Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	CM106S
Medium	CM106M
Large	CM106L



Replacement Cushions

Size	Item#
Small	CUCM006S
Medium	CUCM006M
Large	CUCM006L

SUNSET CPAP MASKS

Mask Series	Mask Type	Replaceable Cushion	Forehead Support	Flex Tube Included	Headgear Included
CM005	Full Face	✓		✓	✓
CM006	Nasal	✓	✓	✓	✓
CM106	Nasal	✓	✓ (Adjustable)		✓
CM007	Full Face	✓			✓
CM008	Nasal	✓			✓
CM110	Nasal	✓			✓
CM036	Nasal Pillow	✓		✓	✓

Sunset Deluxe Nasal Pillow CPAP Mask

Item#

CM036

Mask with Headgear and Small, Medium, and Large Pillows

HCPCS#

MASK - A7034

HEADGEAR - A7035

PILLOWS - A7033

- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design



Sunset Classic Full Face CPAP Mask



Mask with Headgear

Size	Item#
Small	CM007S
Medium	CM007M
Large	CM007L
Extra Large	CM007XL



1/2 THE WEIGHT OF LEADING MASKS

- Conventional forehead pad removed to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

Sunset Classic Nasal CPAP Mask



Mask with Headgear

Size	Item#
Medium	CM008M
Large	CM008L

HCPCS#

MASK - A7034

HEADGEAR - A7035

- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

SLEEPNET CPAP MASKS


MiniMe Pediatric Nasal Mask

Mask with Headgear

Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214
Large	CM004L	60215



HCPCS#
MASK - A7034
HEADGEAR - A7035



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

MiniMe 2 Pediatric Nasal Mask

Vented Mask with Headgear

Size	Item#
Small-Medium	CM027S
Extra Small	CM027XS



Non-Vented Mask with Headgear

Size	Item#
Small-Medium	CM027SNV
Extra Small	CM027XSNV

HCPCS#
MASK - A7034
HEADGEAR - A7035



iQ Nasal Mask

Item#
CM001
iQ Nasal Mask
with Headgear
OEM# TMS-889

HCPCS#
MASK - A7034
HEADGEAR - A7035



- One size fits most
- Flexible shell can be shaped to conform to facial contours, creating comfort and optimal seal
- Small footprint reduces skin irritation and is great for active sleepers or those with facial hair

CPAP/NIV TUBING, CONNECTORS, AND ADAPTERS

Durable CPAP Tubing — Gray


Item#

TUB06

6ft with 22mm Cuffs

HCPCS# A7037

- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free



Slim CPAP Tube

Item#

TUB006SS

6ft x 15mm with 22mm cuffs

HCPCS# A7037

Highly durable

Strain relief cuff decreases separation and tearing

Virtually no memory

Designed to be the most flexible tube available

Minimal weight

One of the lightest tubes on the market at 2.6 ounces

Great value

Affordable and sustainable



Circuit Connectors

Item# **RES013**

22mm OD x
15mm ID/
22mm OD
10/pack



Item# **RES014**

22mm OD x
22mm OD
With Lip
10/pack



Oxygen Bleed-In Adapter

Item# **RES019**

22mm with
6mm Stem
10/pack



- Barbed rotating stem
- Used for oxygen bleed-in

Humidifier Tubing Adapter

Item# **RES016**

22mm ID x
22mm OD
with Lip
10/pack



Swivel Connector Silencer

Item# **VEN002**

For CPAP and
Ventilator Circuits
Each



CHINSTRAPS

All chinstraps are latex-free • Hand wash only • Lay flat to dry


Universal Neoprene Chinstrap



Item# **CS001**
 Universal Chinstrap
 HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby
 OEM# TMS-09 • HCPCS# A7036

- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Adjustable Neoprene Chinstrap



Size	Item#
6in	CS003S
7in	CS003R
8in	CS003M
9in	CS003L
10in	CS003XL

Compare to Puresom Adjustable Ruby
 OEM# TMS-09 ADJ
 HCPCS# A7036

- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head.
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

Deluxe Chinstrap



Size	Item#
26in	CS004S
28in	CS004L

Compare to Respirationics Deluxe
 OEM# 302425
 HCPCS# A7036

- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

Comfort Chinstrap



Item# **CS005**
 Compare to ResMed Sullivan
 OEM# 16015
 HCPCS# A7036

- One size fits most
- Easy to use—attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

Premium Chinstrap



Item# **CS006**
 Compare to Respirationics Premium
 OEM# 1012911
 HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

CLEANERS

CPAP Cleaning Wipes — Flowpack

Item#
CAP1003S

12/case
Flowpack of 64 wipes

- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks



CPAP Cleaning Wipes — Single Packs

Item#
CAP1005S

Each (Box of 30 single-pack wipes)

- Convenient individually wrapped wipes are great for patient samples or travel use



CleanSmart CPAP Disinfectant Spray

Item#
CAP1018-L 16oz **6/case**

CAP1018-S 2oz **9/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Zero chemical residue with active ingredient Hypochlorous Acid
- EPA Registration Number 89896-2
- Approved for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



Moisture Therapy

Item# **RES6104L**
1oz Tube • **18/case**

- Prevent chafing, irritation, dryness, or cracking from using a CPAP mask or other delivery device
- Apply to the face at the beginning of treatment
- Petroleum-free with Aloe Vera and Vitamins A & E



Visit www.sunsethcs.com



Our goal is simple:
To help your facility thrive and be
your acute care partner

We're just easier to work with

- Easy account setup
- One stop shop for respiratory disposables
- A real person answers the phone
- Free samples
- Ability to source new products quickly
- Custom kitting



The Chicago headquarters team



Phone (877) 578-6738
Fax (312) 997-9985
Email sales@sunsethcs.com
Visit www.sunsethcs.com

Corporate Headquarters
Sunset Healthcare Solutions
180 N Michigan Ave Ste 2000
Chicago, IL 60601