



● Portable Suction Machine AC and battery power

Item#
SU100DC

Includes internal rechargeable 12-volt lithium battery, which can be used for up to 1 hour of operation per charge.

● Suction Machine AC power

Item#
SU100AC

Both machines include

- Machine
- AC Adapter and Power Cord
- Carrying Bag
- 800cc Suction Canister
- 1.5in Suction Tube Connector
- 10in Suction Tube Connector
- 6ft Suction Tubing
- Hydrophobic Filter

Both machines feature

- Effective removal of respiratory fluids
- Durability
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty
- HCPCS# E0600

Specifications

Specialty	Small/Lightweight
Carry Case	Yes
Dimensions (L x H x W in)	12.83" x 6.37" x 8.12"
Weight (lbs.)	SU100DC: 5.29 lbs / SU100AC: 4.08 lbs
Vacuum Range (mm Hg)	50-560
LPM	25
Sound Level (dBA)	60 dBA
Canister Type w/ purchase of machine	Disposable
Uses Suction Canister Size (cc)	800 /1000/1200
Machine Warranty	2 years
Can be operated inside carry case?	Yes

Accessories and replacement parts

RES024S-KIT	Replacement Filter Kit (connectors and filter)
RES026S-SS	Replacement Suction Kit (canister, tubing, connectors, and filter)
SU100DC-Adapt	Replacement Power Adapter (without cord)
SU100DC-Cord	Replacement Power Cord (without adapter)
SU100DC-BAG	Carrying/Storage Bag

