



NEW!

*Laser-etched
private label*

Talk to your rep for details.



**2
YEAR
WARRANTY**

CGA 870 Style

Item# RES042

Adult 0-8 L/Min Barbed Outlet
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

Item# RES043

Adult 0-15 L/Min Barbed Outlet
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

Item# RES046

Pediatric 0-4 L/Min Barbed Outlet
0, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4 LPM

CGA 540 Style

Item# RES044

Adult 0-8 L/Min DISS Outlet
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

Item# RES045 /RES045L (Long Stem)

Adult 0-15 L/Min DISS Outlet
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

Item# RES047

Pediatric 0-4 L/Min DISS Outlet
0, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2, 3, 4 LPM

- Click setting wheel
- Heavy duty brass core
- Durable anodized aluminum body
- Input pressure 2200psi
- Output pressure 50psi
- HCPCS# E1353