

CPAP Batteries

- No additional adapters needed
- Outlet functions like a wall socket
- Universal 110-240V recharging conversion
- 1-Year Warranty

Zopec Explore

- Run any electronics with power requirement of 120 watts or under
- AutoSwitch and Bypass Technology for uninterrupted sleep when power outages occur
- Unlimited linking of multiple batteries together for longer use before recharging



Zopec Explore 5700
Item# **CAP1024-2157**

Zopec Explore 8200
Item# **CAP1024-2182**



Zopec Power Outlet 2
Item# **CAP1024-21ZDP02**

- Includes 45W USB-C



Portable Outlet
Item# **CAP1008**

Item#	CAP1008	CAP1024-21ZDP02	CAP1024-2157	CAP1024-2182
Name	Portable Outlet	Power Outlet2	Explore 5700	Explore 8200
Manufacturer	Portable Outlet	Zopec Medical	Zopec Medical	Zopec Medical
OEM #	155W	21ZDP02	215700C	218200C
Dimensions (L x W x H)	7.0" x 5.0" x 1.5"	7.0" x 5.0" x 1.5"	8.5" x 6.0" x 1.0"	7.75" x 5.5" x 1.75"
Volume (cubic inches)	53	53	51	74.5
Weight	2.0 lb	2.0 lb	2.5 lb	3.5 lb
Capacity	42,000 mAh	40,000 mAh	54,000mAh	80,000 mAh
Watt Hours (Wh)	155Wh	146Wh	160Wh	288Wh
AC Output	100W	100W	120W	120W
UPS Function (Autoswitch)	No	No	Yes	Yes
USB Type C	No	PD 45W Smart Hub	15W (5V/3A)	15W (5V/3A)
USB Type A	5V/2.1A	5V/2.1A	5V/3.1A	5V/3.1A
Battery Cell	LG	LG	Panasonic	Panasonic
Capacity Indicator	4 lights	4 lights	Digital Display	Digital Display
Unlimited Linking Technology	No	No	Yes	Yes
Recharge Time	Less than 4 Hours	Less than 4 Hours	Less than 4 Hours	Less than 5 Hours
Car Charger Included	No	No	Yes (45W)	Yes (45W)
Optional Solar Panel Charger	No	No	Yes	Yes
Air Travel	Yes	Yes	Yes	No
Charge Cycles	Over 500 Cycles	Over 500 Cycles	Over 500 Cycles	Over 500 Cycles
Warranty	1 Year	1 Year	1 Year	1 Year