


# PreserveTech™ 360° Swivel Bath Chair

Item#  
**BSE001**



Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"-22", and the padded armrests can be raised out of the way.

Drive PreserveTech™ products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

\*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

## Product Specs

Product Dimensions	18" (W) x 16" (L) x 15"-21" (H)
Product Weight	8.80 lbs.
Color	Grey
Weight Capacity	350 lbs.
Warranty	Limited Lifetime



## Box Specs

Retail Box Dimensions	7.09" (W) x 20.08" (L) x 16.14" (H)
Retail Box Weight	10.30 lbs.