

Bubble Humidifier

Item#
RES050
20/case



- Multi-threaded lid ensures a secure fit every-time, preventing leaks and service calls
- Compatible with oxygen flow up to 15 LPM
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem: it eliminates shredded aluminum diffuser so often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle has a 500 ml water capacity
- Bottle fill line is very easy to see
- There are finger grips at the top of the bottle, making the lid easier to screw on/off for the patient