

Drive DeVilbiss 10 Liter Oxygen Concentrator

Item# **RES1025DS**

HCPCS: **E1390**

- Built upon the reliable features of the 525 concentrator, the 1025 Oxygen Concentrator is designed to provide optimal oxygen delivery for patients with varying flow requirements
- Assembled in the USA and built in the same shell as the 5 Liter, the 10 Liter Concentrator is also compatible with cylinder transfill systems
- With a weight of 42 lbs. and an identical footprint to the 525 unit, the 1025 is substantially smaller and lighter than competitive 10 Liter concentrators
- Capable of delivering 87% to 96% of O₂ purity from 2 to 10 LPM, the 10 Liter Concentrator is an ideal investment for your patients and business alike
- Exclusive DeVilbiss OSD® (Oxygen Sensing Device), standard on every unit, monitors the oxygen produced and ensures patient safety and reliability for longer service intervals



Sold separately, the RES039-RACKD Transfill Caddy is designed with a steel frame to simply connect the 18 x 14 inch shelf, designed to support up to 40 lbs, to the DeVilbiss 10L Concentrator.

SPECIFICATIONS	Description
Delivery Rate	2.0 to 10.0 LPM
Outlet Pressure	20.0 PSI ± 1.0 PSI
Oxygen Percentage	2-10 LPM, 87-96%
Power Consumption	120VAC, 60Hz
Product Dimensions	24.5" (H) x 13.5" (W) x 12" (D)
Product Weight	42 lbs
Warranty	3 Years