



Mini Mesh Nebulizer with Rechargeable Battery

Item#
NEB401



Small particle size
MMAD 2.2µm*



Convenient
Small and Portable



Efficient
Low Residual Liquid



Quiet
≤50dB (at one meter)




Modern
Rechargeable Lithium Battery

Included

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual

Power Input	5VDC Min1.5A
Rechargeable Lithium Battery	3.7VDC 1800mAh
Frequency	110KHz±10%
Nebulization Rate	≥0.2mL/min
Sound Level	≤50dB (at one meter)
Medicine Liquid Temperature	≤140°F



*According to the Standard EN 13544-1:2007+A1:2009
Respiratory therapy equipment - Part 1: Nebulizing
systems and their components