

# BASKETBALL

## Pediatric Compressor Nebulizer

Item#

### NEB300BASKETBALL

Comes with Disposable Neb Kit, Pediatric Aerosol Mask, & Filters  
HCPCS# E0570 + A7003 + A7015 + A7013

### NEB300BASKETBALL-R

Comes with Disposable and Reusable Neb Kits,  
Pediatric Aerosol Mask, & Filters  
HCPCS# E0570 + A7003 + A7005 + A7015 + A7013

Neb kits include nebulizer bottle, t-piece, mouthpiece,  
6in corrugated reservoir tube, and 7ft supply tube.



Great value

**FAST  
TREATMENT  
TIMES**



Comes in a fun retail box



## Specifications

<b>Power</b>	120V AC 60 Hz
<b>Power consumption</b>	200VA (75 W Max)
<b>Compressor flow</b>	6 L
<b>Particle size</b>	0.5 to 6 microns
<b>MMAD</b>	1.663 um
<b>Average mobilization rate</b>	≥0.15 ml/min
<b>Remain amount</b>	0.15 mL
<b>Noise level</b>	Below 63 dBA
<b>Compressor pressure range</b>	30 psig or greater / 206 kPa
<b>Operating temperature range</b>	50°F to 104°F (10°C to 40°C)
<b>Operation humidity</b>	10% to 95% RH
<b>Operating pressure range</b>	8 ~ 13 psi / 55kPa ~ 90kPa
<b>Storage temperature range</b>	-4°F to 158°F (-10°C to 70°C)
<b>Storage humidity range</b>	10% to 95% RH
<b>Dimension (L x W x H)</b>	6.3 x 4.6 x 7.9 in (161 x 117 x 200mm)
<b>Product weight</b>	3.3 lbs (1.5 kg)