

Splash Defense™ Transfer Bench

Item#
BSE002



Splash Defense™ Transfer Bench is the first of its kind to allow the shower curtain liner to remain inside the bath tub while a user is bathing. This unique product enhancement is essential to help prevent water from spilling out onto the bathroom floor, creating a danger zone for users and caretakers. The Splash Defense™ Transfer Bench provides security for the caretaker, while giving users back their independence and allowing them to bathe safely and privately again.

REVERSIBLE TO ACCOMMODATE ANY BATHROOM

DURABLE ALUMINUM FRAME SUPPORTS UP TO 400 LB.



SEAT HEIGHT ADJUSTMENTS
17.5-22.5"

PATENTED SHOWER CURTAIN GUARD

Product Specs

Outside Leg Width/Depth	33.50" (W) x 23.00" (D)
Seat Width / Depth	28.00" (W) x 16.00" (D)
Product Weight	12.70 lb.
Color	White
Weight Capacity	400 lb.
Warranty	Limited Lifetime*

*Refer to instruction guide or website for warranty details.



Box Specs

Retail Box Dimensions	9.30" (W) x 21.10" (L) x 17.50" (H)
Retail Box Weight	17 lb.