

Handheld Compressor Nebulizer

1
YEAR
WARRANTY




Item#

NEB200

Comes with compressor, medication cup, AC power adapter, air hose, elbow, mouthpiece, and mask.

Replacement Parts

NEB200CUP	Medication Cup
NEB200MP	Mouthpiece
NEB200RE	Elbow



Nebulization Rate ≥ 0.25 ml/min
Medication Capacity 6 ml (maximum) / 2 ml (minimum)

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No filter required

Model Number: NEB200	IP Classification: IP22
Dimension & Weight: 115 x 85 x 207 mm; 205 g	Mode of Operation: 30 mins on/ 30 mins off
Power Source: DC 6V or 4 AA batteries	Nebulization Rate: ≥ 0.25 ml/min
Power Consumption: 6 W	Medication Capacity: 6 ml (maximum) / 2 ml (minimum)
Compressor Pressure Range: 90 to 138 kPa	MMAD: Approximately 3 microns
Operating Pressure Range: 35 to 69 kPa	Aerosol Output: 1.8 ± 0.05 (ml) ^{**}
Maximum Flow Rate: ≥ 3 LPM(l/min)	Aerosol Output Rate: 0.32 ± 0.003 (ml/min) ^{**}
Noise Level: 60 dBA	
Input (AC adapter): 100-240V / 50/60Hz	

Operating conditions: 41°F to 104°F (5°C to 40°C), 15% to 93% relative humidity, 700 hPa to 1060 hPa
Storage conditions: -13°F to 158°F (-25°C to 70°C), 10% to 95% relative humidity