

# Mini Mesh Nebulizer

Item#  
**NEB400**



## Features

- Fast treatment times
- Quiet
- Small particle size
- Attainable for more patients
- Runs on USB or 2 AA batteries

## Specifications

<b>Vibrating Frequency</b>	Approx. 117kHz
<b>Power Source</b>	2 AA Batteries or USB 5V
<b>Power Consumption</b>	<1.5W
<b>Nebulization Rate</b>	≥0.25 ml/min
<b>Particle Size</b>	Approx. 5µm (MMAD)
<b>Capacity of Medicine Cup</b>	8ml
<b>Dimensions</b>	Approx. L78 x W41 x H73 mm
<b>Weight</b>	Approx. 2.6 oz (excluding batteries)
<b>Operation Temperature Range</b>	50–100°F/10–40°C
<b>Operating Humidity Range</b>	30%–85% RH
<b>Storage Temperature Range</b>	-68–158°F/-20–70°C
<b>Storage Humidity Range</b>	20%–75% RH

## Included

- Main unit
- Medication cup
- Adult mask
- Pediatric mask
- Mouthpiece
- USB cord
- Carry bag
- User manual

