

# SleepWeaver<sup>®</sup> Élan

Soft cloth CPAP mask kit  
by Circadiance



SleepWeaver's loose fit provides ultimate comfort, a perfect seal, and eliminates strap marks

**All SleepWeaver masks should be fitted loosely for the best seal and comfort**

### Kit Includes

- Mask
- Regular Cushion
- Large Cushion
- Headgear
- Feather Weight Tube
- Overhead Tube Tether Strap
- Quick Start Illustration Guide

### Features

- Replaceable cushion design
- Machine washable
- Three cushion sizes included: Small, Regular and Large
- Detachable 90° elbow with 360° rotation for over the head hose placement
- Soft cloth material inflates like a balloon to create seal
- Clear line of site makes it ideal for wearing glasses, reading and watching TV
- Latex free
- 30-day money back trial

Item#  
**CM022B-KIT**  
Blue

Item#  
**CM022T-KIT**  
Tan