

# Finger Pulse Oximeter

Item#  
**RES5102**



## Product Specs

### DISPLAY RANGE

Functional Oxygen Saturation (SpO2)	0-99%
Pulse Rate (PR)	25-250 bpm
Perfusion Index (PI)	0.1-20%

### ACCURACY

SpO2 Declared Accuracy Range	70%~100%
SpO2 Accuracy	±2%
PR Declared Accuracy range	25-250 bpm
PR Accuracy	3 bpm

### MEASUREMENT RESOLUTION

Oxygen Saturation (%SpO2)	1%
Pulse Rate (PR)	1 bpm

### POWER REQUIREMENTS

Specification of alkaline batteries	2 AAA(LR03)
Operating current	25-50mA
Operating mode	Spot checking
Degree of protection against hazards of explosion	IP22

### PHYSICAL SPECIFICATION

Product dimensions	57 x 30 x 31 mm
Product net weight	28g
Gift box dimensions	9 x 6 x 4 cm
Package	100pcs/ctn
Carton box size	47 x 32 x 20 cm
GW(with battery)	5.5kgs
GW(without battery)	7.8kgs

## Features

- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual



SCAN TO LEARN MORE