




Product Catalog



vizient[®]
Awarded Supplier

Contract

MS1097
MS1049
MS1087



Area Business Director
Jackson Hanft
C: 971-416-6695
jackson@sunsethcs.com

Senior Area Business Director
Bryan Eckard
C: 215-375-5765
beckard@sunsethcs.com


Director, Strategic Accounts
Lisa Iannazzone
C: 636-222-3134
lisa@sunsethcs.com

Sales Representative
Egle Sventnickaite
O: 872-258-1654
esventnickaite@sunsethcs.com

Sales Representative
Owen Ruth
O: 312-533-2451
oruth@sunsethcs.com

Director, Sales Operations
Mark Lipson
C: 248-981-9669
O: 312-533-2465
mlipson@sunsethcs.com

Sales Representative
Christian Straussner
O: 872-282-0746
C: 847-525-6340
cstraussner@sunsethcs.com



Corporate Office

Phone 872-275-3210

Fax 312-997-9985

180 N Michigan Ave Ste 2000

Central Returns/Distribution Center

279 Madsen Dr Ste 101
Bloomington, IL 60108

West Returns/Distribution Center

CONTENT

Contract MS1097

MDD

Page

Nebulizer Kits	4
Nebulizer Compressor	5
Inhaler Chamber	5
Peak Flow Meter	5
Valved T-Adapter	5

Contract MS1049

NIV

Page

Masks	6-8
Mask Elbows	9
Connectors/Ports	9
Headgear	9

Contract MS1087

Respiratory Disposables

Page

Breathing Circuits	10
Bacterial/Viral HME Filter	10
Connectors and Adapters	10
Heat Moisture Exchangers	11
Tracheostomy Masks	11
Berman and Guedel Airways	12
Spirometers	13
Aerosol Mask, Tubing, Valved T-Adapter	13
Oxygen Cannula	14-15
Oxygen Masks	16
Oxygen Tubing	17


Disposable Nebulizer Kits with Mouthpiece

50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- T-piece
- Mouthpiece
- Corrugated tube
- 7ft supply tube



REF	Size	Accessories
RES091L	Adult	All
RES091T	Adult	No corrugated tube

- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Disposable Nebulizer Kits with Mask

50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- Mask
- 7ft supply tube



REF	Size	Accessories
RES090NM	Adult	No Mask
RES092	Pediatric	All

- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble

Medicine Cup

Item#

RES090A

50/case

Nebulizer Compressor

Item# **NEB100**
Each



Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Peak Flow Meter

Item# **RES052**
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**
Pediatric Mouthpiece
for Peak Flow Meter
10/case

Valved T-Adapter

Item# **VEN001**
22F-18F-22M/15F
10/pack



SperaChamber

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Item# **RES195-10**
10/case

- Inhalation valve
- Universal adapter
- Spiral vents

RES195-50
50/case

- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask
10/case

REF	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)
RES195ML	Large (Adult)

Tracheostomy SperaChamber

Item# **RES195T**
10/case

Full Face Mask

With Standard Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041S
Medium	CM041M
Large	CM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEE
Medium	CM041MEE
Large	CM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV



Standard Elbow, Non-Vented

Size	Item#
Small	VCM041S
Medium	VCM041M
Large	VCM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	VCM041SEE
Medium	VCM041MEE
Large	VCM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	VCM041SEEV
Medium	VCM041MEEV
Large	VCM041LEEV



- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031

With Cap Style Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041SC
Medium	CM041MC
Large	CM041LC

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEEC
Medium	CM041MEEC
Large	CM041LEEC

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEVC
Medium	CM041MEEVC
Large	CM041LEEVC

Nasal Mask

With Standard Headgear

Size	Item#
Small	CM008SV
Medium	CM008MV
Large	CM008LV



Size	Item#
Small	VCM008S
Medium	VCM008M
Large	VCM008L



- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

Total Face Mask

Standard Elbow

Size	Item#
Small	CM040S
Medium	CM040M
Large	CM040L

EE Elbow

Size	Item#
Medium	CM040MEE
Large	CM040LEE



Veraseal 2 Disposable Gel Mask with Headgear



Standard Elbow, Non-Vented



Size	Item#
Small	CM50933
Medium	CM50934
Large	CM50935
Extra Large	CM50936



Size	Item#
Small	VCM053SZ
Medium	VCM053MZ
Large	VCM053LZ
Extra Large	VCM053XLZ



Single Limb Elbow, Non-Vented, AAV



Size	Item#
Small	CM50937
Medium	CM50938
Large	CM50939
Extra Large	CM50940



Size	Item#
Small	VCM054SZ
Medium	VCM054MZ
Large	VCM054LZ
Extra Large	VCM054XLZ



Single Limb Elbow, Vented, AAV



Size	Item#
Small	CM50827
Medium	CM50828
Large	CM50829
Extra Large	CM50830



Size	Item#
Small	VCM052SZ
Medium	VCM052MZ
Large	VCM052LZ
Extra Large	VCM052XLZ

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.

- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.

- The active headgear connector moves as the patient moves for a reliable seal.

- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

Innova

Nasal Mask with Headgear

Vented

Size	Item#
Small/Medium	CM50162
Medium/Large	CM50163

Non-Vented

Size	Item#
Small/Medium	CM50594
Medium/Large	CM50595



Innova® CPAP masks feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

MiniMe 2 Pediatric Nasal Mask



Vented

Size	Item #
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented

Size	Item #
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261



Vented

Size	Item #
Extra Small	VCM055XSZ
Small	VCM055SZ
Medium	VCM055MZ
Large	VCM055LZ

Non-Vented

Size	Item #
Extra Small	VCM056XSZ
Small	VCM056SZ
Medium	VCM056MZ
Large	VCM056LZ



- Maximum comfort, minimum weight. Patented AIR^ogel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA

Ray Pediatric Masks Comparable to Nonny masks

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL

Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031



Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	CM060-KITS
L&XL Cushions	CM060-KITL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

HCPCS#
 MASK - A7034
 HEADGEAR - A7035
 CUSHION - A7032



Standard Elbow

NIV Contract MS1049

Standard Elbow, Non-Vented

Item# **CAP009**
10/case



Single Limb Elbow, Non-Vented, AAV

Item# **CAP008**
10/case




Single Limb Elbow, Vented, AAV

Item# **CAP010**
10/case



Swivel Connector Silencer

Item# **VEN002**
10/case



Single Use Exhalation Port

Item# **CAP011**
10/case



Standard 4-Point Headgear

Item# **HG003**



Cap Style Headgear

Item# **HG002**



RESPIRATORY DISPOSABLES Contract MS1087



Passive Heated Breathing Circuit

Item#

VC1114Z

For use with MR850 Humidifier

10/case



Non-Heated Dual Limb Breathing Circuits

10/case

REF	Size	Humidifier Limb	Water Trap
VC2001PZ	Pediatric	No	No
VC2001Z	Adult	No	No
VC2002PZ	Pediatric	Yes	No
VC2002Z	Adult	Yes	No
VC2003Z	Adult	Yes	Yes

Bacterial/Viral HME Filter

Item#

BF2200Z

5/pack



Xmas Tree Connector

Item# **RES001-10**

10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Connector

Item# **RES012-10**

Fits 3/16in and 1/4in ID Tubing

10/pack



- Barbed ends ensure secure, leak-free connections

Circuit Connectors

Item# **RES014-10**

22mm OD x
22mm OD
With Lip

10/pack



Humidifier Tubing Adapter

Item# **RES016**

22mm ID x
22mm OD
with Lip

10/pack



RESPIRATORY DISPOSABLES Contract MS1087

Flex HME

Item# **RES027FLEX-40**
40/case

TV Range (mL)	300-1500
Moisture Output (mg H ₂ O/L)	33.8, Vt = 0.5L
Weight (g)	17
Resistance (cm H ₂ O)	0.5@30 LPM, 1.34@60 LPM
Dead Space (mL)	90
Connection (mm)	22M, 15F



Heat Moisture Exchangers

30/case



RES027-30
T-HME



RES0270P-30
T-HME with
Oxygen Port



RES0271L-30
Inline HME —
Large Volume



RES0271LOHME
Inline HME —
Large Volume
with Oxygen Port




RES027NE-10
Inline HME —
Neonatal

TV Range (mL)	250-1000	250-1000	15-50	50-250	15-50
Moisture Output (mg H ₂ O/L)	27.3, Vt = 0.5L	26, Vt = 0.5L	31, Vt = 0.1L	31.2, Vt = 0.1L	30, Vt = 0.2L
Weight (g)	5	6.5	10	10	5
Resistance (cm H ₂ O)	0.7@60 LPM	0.4@30 LPM	0.4@10 LPM	0.4@10 LPM	0.5@10 LPM
Dead Space (mL)	12.5	12	13	13	2.5
Connection (mm)	15F	15mmF 5.5mm Oxygen port	15M, 22M/15F	15M, 22M/15F 5.5mm Oxygen port	15M, 15F

Adapter for T-HME

Item# **RES027AZ**
Each



Tracheostomy Mask

50/case

Item# **RES2130Z**
Adult Regular Size

RES2130L
Adult Large Size

RES2230Z
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing

RESPIRATORY DISPOSABLES Contract MS1087

Berman Oral Airway

50/pack

REF	Size	Color
AOA40BC-50	40mm	Pink
AOA50BC-50	50mm	Blue
AOA60BC-50	60mm	Black
AOA70BC-50	70mm	White
AOA80BC-50	80mm	Green
AOA90BC-50	90mm	Yellow
AOA100BC-50	100mm	Red
AOA110BC-50	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

Oropharyngeal Guedel Airway

10/pack

REF	Size	Color
AOP500GC	50mm	Blue
AOP601GC	60mm	Black
AOP702GC	70mm	White
AOP803GC	80mm	Green
AOP904GC	90mm	Yellow
AOP1005GC	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

RESPIRATORY DISPOSABLES Contract MS1087

Volumetric Incentive Spirometer

12/case

Item# **RES055**
2500ml

RES056
5000ml

RES055P
Pediatric



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.

Three Ball Spirometer

12/case

Item# **RES057**

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

Aerosol Mask

50/case

Item# **RES2110Z**
Adult Size

Item# **RES2110F**
Adult Size
w/ Flex Tube


Item# **RES2210Z**
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

Aerosol Tubing

Item# **TUB100-6**
Corrugated Aerosol Tubing
100ft in 6in segments



Valved T-Adapter, Aerosol, 22mm ID x 22mm OD

Item# **NEB9001T**
30/case



RESPIRATORY DISPOSABLES Contract MS1087

Adult Oxygen Cannula

REF	Tubing Length	Connector	Qty
RES1100Z	None	Standard	50
RES1104Z	4 ft	Standard	50
RES1107U	7 ft	Universal	50
RES1107Z	7 ft	Standard	50
RES11110	10 ft	Standard	50
RES1114N	14 ft	Standard	50
RES1114U	14 ft	Universal	50
RES1115U	15 ft	Universal	50
RES1115Z	15 ft	Standard	50
RES1125Z	25 ft	Standard	25
RES1135	35 ft	Standard	20
RES1150U	50 ft	Universal	20
RES1150Z	50 ft	Standard	20



Soft Cannula


REF	Tubing Length	Connector	Qty
RES1107SU	7 ft	Universal	50
RES1107SZ	7 ft	Standard	50
RES1114S	14 ft	Standard	50
RES1114SU	14 ft	Universal	50
RES1125S	25 ft	Standard	25
RES1125SU	25 ft	Universal	25

Flared Cannula

REF	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

Adult Cannula with Ear Cushions

REF	Tubing Length	Connector	Touch	Tip	QTY
RES1100ECZ	None	Standard	Standard	Standard	50
RES1104EC	None	Standard	Standard	Standard	50
RES1107ECU	7 ft	Universal	Standard	Standard	50
RES1107ECZ	7 ft	Standard	Standard	Standard	50
RES1114EC	14 ft	Standard	Standard	Standard	50
RES1125EC	25 ft	Standard	Standard	Standard	25



RESPIRATORY DISPOSABLES Contract MS1087

Pediatric, Infant, and Neonate Cannula

50/case

REF	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No


Humidoflo High Flow Nasal Cannula

10/case

REF	Size
HFN001AS	Small
HFN001AM	Medium
HFN001AL	Large

- Includes standard ISO connector and F&P adaptor so it can be used with almost any machine on the market
- Kink-free prongs and flexible, low-noise tubing
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)
- Can be used with F&P 950, 850, or Airvo system

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O₂ concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.



RESPIRATORY DISPOSABLES Contract MS1087

Medium Concentration Oxygen Mask

50/case

REF	Size	Tubing Length	Connector
RES2100NT	Adult	None	Standard
RES2100U	Adult	7 ft	Universal
RES2100Z	Adult	7 ft	Standard
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203I	Infant	None	Standard
RES2203I7	Infant	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

High Concentration Oxygen Mask

50/case

REF	Size	Tubing Length
RES2202A	Adult	None
RES2202A7	Adult	7 ft
RES2202I	Infant	None
RES2202I7	Infant	7 ft
RES2202P	Pediatric	None
RES2202P7	Pediatric	7 ft

Non-Rebreather Oxygen Mask

50/case

REF	Size	Tubing Length	Connector	Safety Vent
RES2102I	Infant	7 ft	Standard	No
RES2102IM	Infant	None	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No
RES2102S	Adult	7 ft	Standard	Yes
RES2102SM	Adult	None	Standard	Yes
RES2102U	Adult	7 ft	Universal	No
RES2102Z	Adult	7 ft	Standard	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



RESPIRATORY DISPOSABLES Contract MS1087

Oxygen Supply Tubing

REF	Length	Connector	Color	Bore	Qty
RES3007SB	7 ft	Standard	None	Smooth	50
RES3007U	7 ft	Universal	None	Channeled	50
RES3007Z	7 ft	Standard	None	Channeled	50
RES3015GZ	15 ft	Standard	Green	Channeled	50
RES3015Z	15 ft	Standard	None	Channeled	50
RES3025GZ	25 ft	Standard	Green	Channeled	25
RES3025U	25 ft	Universal	None	Channeled	25
RES3025Z	25 ft	Standard	None	Channeled	25
RES3030	30 ft	Standard	None	Channeled	20
RES3030G	30 ft	Standard	Green	Channeled	20
RES3040GZ	40 ft	Standard	Green	Channeled	20
RES3050GZ	50 ft	Standard	Green	Channeled	20
RES3050Z	50 ft	Standard	None	Channeled	20



Universal Connector

