



Sol

Full Face CPAP Mask



Easy pop cushions
Cushions attach or detach from mask in an easy pop; no need to manually align pieces




Quiet vent
Honeycomb design reduces airflow erosion intensity for gentle and quiet flow



Quick release connector
Rapidly and easily separates or connects mask and the tube



Soft and light headgear
Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment



Full Face CPAP Mask with Headgear and Cushion

CM205S	Small
CM205M	Medium
CM205L	Large

Mask without Headgear

CM205SNH	Small
CM205MNH	Medium
CM205LNH	Large


Replacement Cushions

CUCM205S	Small
CUCM205M	Medium
CUCM205L	Large

Replacement Headgear

HG205

Sol comes in attractive retail packaging



SCAN TO LEARN MORE



Sunset Healthcare Solution's new Sol Full Face CPAP Mask offers premium features without the premium price tag. Maximize reimbursement without sacrificing quality or compliance and help patients wake up to a brighter day with Sol!