


# Ray Pediatric CPAP/NIV Masks

Comparable to Nonny masks

Scan to see online



Available size guide

- Easily adjustable forehead T-piece
- Small, light and comfortable one piece molded silicone mask
- Convenient quick release clips
- Vented for CPAP, BiPAP, and NIV
- Dual 360 degree swivel at the mask attachment and at the standard 22mm tubing attachment point for more freedom

## Full Face

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

### Replacement Headgear

Size	Item#
Small	<b>HG061S</b>
Large	<b>HG061L</b>

## Nasal

Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

- Comes with either Small & Medium cushions or Large & XL cushions to help fit the patient more easily

### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>



### Replacement Headgear

Size	Item#
Small	<b>HG060S</b>
Large	<b>HG060L</b>