



# Acute Care Product Catalog



*Sales Representative*

**Egle Sventnickaite**

O: 872-258-1654

esventnickaite@sunsethcs.com

*Sales Representative*

**Christian Straussner**

O: 872-282-0746

C: 847-525-6340

cstraussner@sunsethcs.com

*Sales Representative*

**Owen Ruth**

O: 312-533-2451

oruth@sunsethcs.com

*Director, Sales Operations*

**Mark Lipson**

O: 312-533-2465

mlipson@sunsethcs.com

*Sales and Marketing Manager*

**Jackson Hanft**

C: 971-416-6695

jackson@sunsethcs.com

*Director, Strategic Accounts*

**Lisa Iannazzone** RN, BSN

C: 636-222-3134

lisa@sunsethcs.com

*RRT Clinical Specialist*


**Carlos Mejia** MS-QSM, RRT, RPGST, LSSB

cmejia@sunsethcs.com

*Clinical Specialist*

**Marion Sachs** RN, MSN, CNOR, CVAHP

msachs@sunsethcs.com



### Corporate Office

Phone 872-275-3210

Fax 312-997-9985

180 N Michigan Ave Ste 2000  
Chicago, IL 60601

### Central Returns/Distribution Center

279 Madsen Dr Ste 101  
Bloomington, IL 60108

### West Returns/Distribution Center

6750 W Van Buren St  
Phoenix, AZ 85043

# CONTENT

## CPAP

### Page

Masks and Replacement Parts	4-12
Sunset Healthcare Solutions	4-6
Sleepnet	7-9
Circadiance	10-11
InnoMed	12
Airway Management	12
Sleep Apnea Appliances	13
Tubing	14-15
Filters	16-24
Respironics	16-17
ResMed	18-19
Fisher & Paykel	20
3B	20
AEIOMed	21
Breas/GE	21
CareFusion/Viasys	22
Covidien Nellcor Puritan Bennett	22
DeVilbiss	23
Evo Medical	23
Healthydyne	24
Invacare	24
ProBasics	24
Humidifier Chambers	25
Chin straps	26-27
Headgear	27
Cash Items	28-35
Cleaners and Accessories	28-29
Battery	30
Mask Liners and Soft Headgear Covers	32
Self Care	33
Connectors, Adapters, and Inline Filter	33
Pillows	34
Storage	34
Pressure Meter	35
Tubing Wraps	35

## Oxygen

### Page

Pulse Oximeter	36
Firesafe	37
Concentrators	38
Batteries	38
Masks	39
Cannula	40-41
Tubing	41
Regulators and Conserver	42
Cylinder Carriers	44-45
Connectors and Adapters	46
Accessories	46-48
Cleaners	48
Bacteria and HEPA Filters	49
Concentrator Cabinet Filters	50-55
Caire/AirSep/SeQual	50
CHAD	51
Covidien Nellcor Puritan Bennett	51
DeVilbiss/Drive	52
Inogen	53
Invacare	53
Jumao, Lifestyle, Attendant	54
Nidek	54
3B	54
Respironics/Healthydyne	55

## Other

### Page

Aerosol	57
Tabletop Nebulizers	58-59
Handheld Nebulizers	59
Nebulizer Kits	60
Nebulizer and Anesthesia Filters	61
Suction	62-65
Tracheostomy	66
Ventilator Filters, Masks, and Connectors	67
Ventilator Circuits	68-69
StrongArm Comfort Cane	70
Bath Safety	71
PPE	72-74
Mini TENS	74
Maternity	75

## Innovative Technology

### Page

GuideBlade	76
PeDIA	77
Anestand	78-79

# SUNSET HEALTHCARE SOLUTIONS

## Sol Full Face CPAP Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM205S</b>
Medium	<b>CM205M</b>
Large	<b>CM205L</b>

**NEW**

Ask us about our  
30-day mask guarantee

### Replacement Cushions

Size	Item#
Small	<b>CUCM205S</b>
Medium	<b>CUCM205M</b>
Large	<b>CUCM205L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



### A premium mask without the premium price

#### Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

#### Quiet vent


Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

#### Quick release connector

Rapidly and easily separates or connects mask and the tube

#### Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment



Sol comes in attractive retail packaging



## Ray Pediatric CPAP/NIV Masks

### Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031



**NEW**

### Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032



• Comparable to Nonny masks



# SUNSET HEALTHCARE SOLUTIONS

## Sunset Deluxe Full Face CPAP Masks

### Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	<b>CM005S</b>
Medium	<b>CM005M</b>
Large	<b>CM005L</b>

Unobtrusive design offers a clear line of sight


**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPSC#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031

### Mask without Headgear

Size	Item#
Small	<b>CM005SNH</b>
Medium	<b>CM005MNH</b>
Large	<b>CM005LNH</b>



### Replacement Cushions

Size	Item#
Small	<b>CUCM005S</b>
Medium	<b>CUCM005M</b>
Large	<b>CUCM005L</b>

- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

### Mask with Adjustable Forehead T, Headgear, and Replaceable Cushion

Size	Item#
Small	<b>CM105S</b>
Medium	<b>CM105M</b>
Large	<b>CM105L</b>



## Sunset Deluxe Nasal CPAP Masks

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM006S</b>
Medium	<b>CM006M</b>
Large	<b>CM006L</b>

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPSC#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032

### Mask without Headgear

Size	Item#
Small	<b>CM006SNH</b>
Medium	<b>CM006MNH</b>
Large	<b>CM006LNH</b>



### Replacement Cushions

Size	Item#
Small	<b>CUCM006S</b>
Medium	<b>CUCM006M</b>
Large	<b>CUCM006L</b>

### Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	<b>CM106S</b>
Medium	<b>CM106M</b>
Large	<b>CM106L</b>



# SUNSET HEALTHCARE SOLUTIONS

## Sunset Deluxe Nasal Pillow CPAP Mask

Item#

**CM036**

Mask with Headgear and Small, Medium, and Large Pillows

**HG036**

Headgear Only

HCPCS#

MASK - A7034

HEADGEAR - A7035

PILLOWS - A7033



### Replacement Pillows

Size	Item#
Small	<b>CM036SP</b>
Medium	<b>CM036MP</b>
Large	<b>CM036LP</b>



### Mask without Headgear

**CM036NH**

- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

## Sunset Clearsight Deluxe Nasal Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM110S</b>
Medium	<b>CM110M</b>
Large	<b>CM110L</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

### Replacement Cushions

Size	Item#
Small	<b>CUCM110S</b>
Medium	<b>CUCM110M</b>
Large	<b>CUCM110L</b>

HCPCS#

MASK - A7034

HEADGEAR - A7035



- Clearsight design provides unobstructed view

Mask Series	Mask Type	Replaceable Cushion	Forehead Support	Flex Tube Included	Headgear Included
<b>CM005</b>	Full Face	✓		✓	✓
<b>CM105</b>	Full Face	✓	✓ (Adjustable)		✓
<b>CM006</b>	Nasal	✓	✓	✓	✓
<b>CM106</b>	Nasal	✓	✓ (Adjustable)		✓
<b>CM007</b>	Full Face				✓
<b>CM008</b>	Nasal				✓
<b>CM110</b>	Nasal	✓			✓
<b>CM036</b>	Nasal Pillow	✓		✓	✓

## Sunset Classic Full Face CPAP Mask

### Mask with Headgear

Size	Item#
Small	<b>CM007S</b>
Medium	<b>CM007M</b>
Large	<b>CM007L</b>
Extra Large	<b>CM007XL</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

**CMSIZE02**

Classic Sizing Guide



HCPCS#  
MASK - A7030  
HEADGEAR - A7035



**1/2 THE WEIGHT  
OF LEADING MASKS**

- Conventional forehead pad removed
- to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

## Sunset Classic Nasal CPAP Mask

### Mask with Headgear

Size	Item#
Medium	<b>CM008M</b>
Large	<b>CM008L</b>

**HG003** Headgear - Fits Most

**HG003XL** Headgear - XL

**CMSIZE02**

Classic Sizing Guide



HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

# SLEEPNET

## Phantom®2 and iQ®2

### Phantom®2

Item#  
**CM050**

Replacement Cushion  
**CUCM050**



**NEW**

### iQ®2

Item#  
**CM049**

Replacement Cushion  
**CUCM049**



Replacement Frame/Tubing Assembly  
**CM051**

Ask us about our  
30-day mask guarantee

- Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel® allow the mask to seal with minimal tension. By design AIR°gel® is twice as soft as competing gel products.
- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.
- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS

# SLEEPNET

## Innova

### Nasal Mask with Headgear

Size	Item#
Small/Medium	<b>CM034SM</b>
Medium/Large	<b>CM034ML</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

### Full Face Mask with Headgear

Size	Item#
Small	<b>CM035S</b>
Medium	<b>CM035M</b>
Large	<b>CM035L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035



Innova® CPAP masks feature patented AIR<sup>o</sup>gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

## MiniMe Pediatric Nasal Mask

### Mask with Headgear

Headgear Size	Item#	OEM#
Small	<b>CM004S</b>	60213
Medium	<b>CM004M</b>	60214
Large	<b>CM004L</b>	60215




### Replacement Headgear

Size	Item#	OEM#
Small	<b>HG010S</b>	55241
Medium	<b>HG010M</b>	55242
Large	<b>HG010L</b>	55243

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Non-vented version available



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

## MiniMe 2 Pediatric Nasal Mask

### Mask with Headgear

Size	Item#	OEM#
Extra Small	<b>CM027XS</b>	60216
Small	<b>CM027S</b>	60217
Medium	<b>CM027M</b>	60218
Large	<b>CM027L</b>	60221



### Replacement Headgear

Size	Item#
Extra Small/Small	<b>HG027S</b>
Medium/Large	<b>HG027L</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

- Non-vented version available





## Ascend

### Nasal Mask with Headgear and All Cushions

Item# **CM037N**

#### Replacement Cushions

Size	Item#
Small	<b>CM037NS</b>
Medium	<b>CM037NCM</b>
Large	<b>CM037NCL</b>

### Nasal Mask with Headgear and 1 Cushion

Size	Item#
Small	<b>CM037NS</b>
Medium	<b>CM037NM</b>
Large	<b>CM037NL</b>

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032

### Full Face Mask with Headgear and All Cushions

Item# **CM037F**

#### Replacement Cushions


Size	Item#
Small	<b>CM037CS</b>
Medium	<b>CM037CM</b>
Large	<b>CM037CL</b>

### Full Face Mask with Headgear and 1 Cushion

Size	Item#
Small	<b>CM037S</b>
Medium	<b>CM037M</b>
Large	<b>CM037L</b>

HCPCS#  
MASK - A7030  
HEADGEAR - A7035  
CUSHION - A7031

- One frame with three sizes of replaceable cushions: S, M, L
- Unique flexible spacebar design provides superior stability, holding a seal without contacting the skin
- Sleepnet's patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation while offering maximum comfort with minimum weight
- Lightweight, breathable quick-release headgear slides off easily with one hand
- 360-degree rotating elbow for optimal tube positioning and easier mobility
- Active Headgear Connector moves as the patient moves for a reliable seal
- Made in the USA



## iQ Nasal Mask

Item# **CM001**

iQ Nasal Mask with Headgear  
OEM# TMS-889

**HG005**

3-Point Headgear Only  
OEM# TMS-89-55046

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

**HG009**

Holey Cap Style Headgear Only



## Mojo 2 Full Face Mask

### Mask with Headgear

Size	Item#	OEM#
Small	<b>CM50817</b>	50817
Medium	<b>CM50818</b>	50818
Large	<b>CM50819</b>	50819
Extra Large	<b>CM50820</b>	50820

- Maximum comfort, minimum weight. Patented AIR<sup>o</sup>gel minimizes pressure points and skin irritation.
- Forehead adjustment screw removes pressure from the bridge of the nose for personalized fit.
- The shapeable shell of each mask can be easily molded to fit individual needs. Bend it. Mold it. Shape it.
- Lightweight, breathable headgear disconnects from mask easily with one hand thanks to the quick-release magnetic headgear clips.
- Oxygen port allows easy delivery of supplemental oxygen.
- 360° rotating elbow for optimal tube positioning and easier mobility.



## Mojo Full Face Mask

### Mask with Headgear

Size	Item#	OEM#
Small	<b>CM003SH</b>	TMS-835
Medium	<b>CM003MH</b>	TMS-836
Large	<b>CM003LH</b>	TMS-837

### Replacement Headgear

**HG003**  
OEM# TMS-89

HCPCS#  
MASK - A7030  
HEADGEAR - A7035



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

# CIRCADIANCE

## SleepWeaver® Advance Nasal Mask



- One size fits most (Small Tan available for smaller faces)
- Includes Advance headgear that can be adjusted to fit a broader range of patients



### Mask with Adjustable Headgear

Mask Color	Item#	OEM#
Blue	<b>CM017B</b>	100274
Camo	<b>CM017C</b>	100281
Leopard	<b>CM017L</b>	100285
Pink	<b>CM017P</b>	100277
Tan	<b>CM017T</b>	100332

**CM017T-S**  
Small Tan  
OEM# 101224

HCPCS#  
MASK - A7034  
HEADGEAR - A7035

**HG017**  
Headgear Only  
OEM# 100435

**CM017TNH**  
Tan — No Headgear  
OEM# 100336

- Easy to store and pack for travel

## SleepWeaver® Élan Nasal Mask

### Mask Kit

Mask Color	Item#	OEM#
Tan	<b>CM022T-KIT</b>	100563
Blue	<b>CM022B-KIT</b>	100782

### Kits Include:

- Mask
- Small Cushion
- Regular Cushion
- Large Cushion
- Headgear
- Feather Weight Tube®
- Overhead Tube Tether Strap
- Quick Start Illustration Guide
- Instructional DVD

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
TUBE - A7037

**TUB0018F**  
Replacement  
Feather Weight Tube®  
OEM# 100477

**HG022**  
Headgear Only  
OEM# 100630

### Mask with Headgear only

Size and Color	Item#	OEM#
Small Blue	<b>CM022B-S</b>	100728
Regular Blue	<b>CM022B-R</b>	100519
Large Blue	<b>CM022B-L</b>	100482
Small Tan	<b>CM022T-S</b>	100727
Regular Tan	<b>CM022T-R</b>	100518
Large Tan	<b>CM022T-L</b>	100481

## All SleepWeaver Masks

- Made of soft cloth
- Allow sleeping on side, stomach, and even wearing glasses
- Easy to store and pack for travel
- Hand or machine washable gentle cycle, mild detergent, hang to dry

### Replacement Cushions

Size and Color	Item#	OEM#
Small Blue	<b>CUCM022B-S</b>	100743
Regular Blue	<b>CUCM022B-R</b>	100654
Large Blue	<b>CUCM022B-L</b>	100664
Small Tan	<b>CUCM022T-S</b>	100742
Regular Tan	<b>CUCM022T-R</b>	100653
Large Tan	<b>CUCM022T-L</b>	100663
Small Leopard	<b>CUCM022L-S</b>	100746
Reg. Leopard	<b>CUCM022L-R</b>	100657
Large Leopard	<b>CUCM022L-L</b>	100667



- Detachable 90° elbow with 360° rotation allows for over-the-head hose placement

## SleepWeaver® Anew™ Full Face Mask

### Blue Mask with Headgear

Size	Item#	OEM#
Small	<b>CM026B-S</b>	100955
Regular	<b>CM026B-R</b>	100960
Large	<b>CM026B-L</b>	100965

### Headgear

Size	Item#	OEM#
Small	<b>HGCM026-S</b>	101145
Reg/Lg	<b>HGCM026-R</b>	101084

### TUB0018F

Feather Weight Tube®  
OEM# 100477

HCPCS#  
MASK - A7027  
HEADGEAR - A7035  
TUBE - A7037



- Soft material wicks away moisture and leaves no pressure marks on the face
- Loose fit provides ultimate comfort and eliminates strap marks
- Open line of sight accommodates eyeglasses and is great for people who like to fall asleep while reading or watching TV
- 90° elbow swivel connector rotates 360° allowing users to sleep in more natural positions; easy to remove and clean

## Pediatric SleepWeaver® CPAP Mask

Item#  
**CM101325**

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.



## SleepWeaver® 3D®

Item#  
**CM027**

**NEW**

### Blue Mask with Headgear

Headgear	Item#
With Large Headgear	<b>CM027-LB</b>
With Regular Headgear	<b>CM027-RB</b>
Nasal Mask (no headgear)	<b>CM027</b>

### 100/packs

Headgear	Item#
With Large Headgear	<b>CM027-L-100</b>
With Regular Headgear	<b>CM027-R-100</b>



- New wing-shaped stay routes the headgear under the ears, ensuring comfort and a great fit.
- Permanent 360° swivel provides reliable performance in any position, including over the headboard.
- Breathable cloth design eliminates fogging that is common in hard plastic masks. Nothing breathes like cloth, keeping skin dry and comfortable.
- For those who want to read or watch TV while in bed, SleepWeaver 3D can be worn with glasses.
- Users prefer a soft cloth against the face over rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad.



# INNOMED

## Aloha® Nasal Pillows System

Item#  
**CM014**

All Sizes Kit  
OEM# ALO100

### Replacement Headgear

**HG014**  
Aloha Headgear  
OEM# ALO130

### Replacement Pillows

Size	Item#	OEM#
Small	<b>CM014SP</b>	ALO121
Medium	<b>CM014MP</b>	ALO123
Large	<b>CM014LP</b>	ALO125



### Frame

Size	Item#	OEM#
Small	<b>CM014FS</b>	ALO131
Medium	<b>CM014FM</b>	ALO133
Large	<b>CM014FL</b>	ALO135

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
PILLOWS - A7033

- One of the highest consumer-rated CPAP masks on the market
- Limited headgear reduces points of contact on the face and prevents interference with the eyes
- Removable micro-fleece headgear wraps are included for comfort
- Pillow angle can be adjusted for a more customized and natural fit

## AIRWAY MANAGEMENT

### Tap Pap

### No headgear needed!

Item#  
**CM024NH**

Mask Kit with Small, Medium, and Large Pillows  
OEM# PAP-NP1-000

### Replacement Pillows

Size	Item#	OEM#
Small	<b>CUCM024S</b>	PAP-NP1-101
Medium	<b>CUCM024M</b>	PAP-NP1-102
Large	<b>CUCM024L</b>	PAP-NP1-103

**HG024**  
Tap Pap Headgear  
OEM# PAP-NP1-300

**CM024MP**  
Replacement Mouthpiece  
OEM# PAP-NP1-201

HCPCS#  
MASK - A7034  
PILLOWS - A7033  
HEADGEAR - A7035



- Mouthpiece molds to the upper teeth to provide a custom and secure fit even if the bottom teeth are not in contact with the mouthpiece
- Unique nut and bolt joint allows the nasal pillows to be set to the desired angle
- Provides a stable fit regardless of the sleep position
- Patients report not being able to feel the mask on their face
- 3 sizes of nasal pillows included (Small, Medium, and Large)
- 90-day compliance guarantee for new users of the mask, 30-day guarantee for experienced users



# AIRWAY MANAGEMENT

## myTAP



Item# **CM039**

**Precision-fit, immediate, simple and affordable solution for snoring and mild to moderate sleep apnea.**

**NEW**

- Offers immediate treatment
- Can be fitted in 15 minutes in a sleep lab, physician's office or HME
- Clinically proven in over 40 independent peer reviewed studies to treat snoring and obstructive sleep apnea (OSA).
- Patients control and remold trays with Precision-Fit, enabling overall comfort, compliance and success. 90% of patients achieve ideal fit upon first try. No bite registration needed.
- Thin, sleek trays create a comfortable, low-profile appliance
- Simple design makes it an affordable and convenient option
- Includes an AM Aligner to help offset potential bite changes
- Includes a mouth shield accessory to prevent dry mouth or excess saliva. A secondary benefit is the promotion of nasal breathing, the best natural way to breathe.

## AIRAVANT

### Bongo Rx EPAP

Item# **CAPBNG500**

Bongo Rx all sizes starter kit

- All four sizes: small, medium, large, and extra large
- Travel case
- Drying stand
- Headgear

**NEW**



#### Annual Replenishment Pack

Includes 4 pairs, travel case, drying stand, and headgear

Size	Item#
Small	<b>CAPBNG502</b>
Medium	<b>CAPBNG503</b>
Large	<b>CAPBNG504</b>
Extra Large	<b>CAPBNG505</b>



The Bongo Rx is a simple device with soft nasal seals and is available in various sizes (S, M, L, and XL). The Bongo Rx is an EPAP (expiratory positive airway pressure) device and works as follows:

- During inhalation, small valves open to allow normal breathing through the device.
- During exhalation, the valves close and direct exhaled air through small vents to generate EPAP, which keeps the airway open.

- No machine, tube, electricity, or batteries needed
- Use at home, traveling, camping... use anytime, anywhere
- Small enough to fit in a shirt pocket
- Reusable and easy to clean with simply soap and water
- Designed, molded, and assembled in the USA
- Clinically tested to demonstrate safety and effectiveness
- Supplemental lifestyle option for your compliant patients
- Patient friendly option for your non-compliant patients

Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS

# CPAP TUBING

## Durable CPAP Tubing — Gray



Item#

**TUB06**

6ft

19mm with 22mm cuffs

HCPCS# A7037

- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free


### All Lengths

<b>TUB0018</b>	18in
<b>TUB002</b>	2ft
<b>TUB004</b>	4ft
<b>TUB06</b>	6ft
<b>TUB006.5</b>	6ft 6in — Compatible with ResMed tube wraps
<b>TUB007</b>	7ft
<b>TUB008</b>	8ft
<b>TUB009</b>	9ft
<b>TUB010</b>	10ft

Compatible with Sunset's new quick release connector!

Item#

**CAP1016**



## Slim CPAP Tube

Item#

**TUB006SS**

6ft x 15mm with  
22mm cuffs

HCPCS# A7037



### Highly durable

Strain relief cuff decreases separation and tearing

### Virtually no memory

Designed to be the most flexible tube available

### Minimal weight

One of the lightest tubes on the market at 2.6 ounces


### Great value

Affordable and sustainable

## Pressure Port Tubing

Item# **TUB006PP**

6ft with Pressure Sensor Lines  
22mm Cuffs • HCPCS# A7037



- Built-in sensor line for use with the Puritan Bennett: GoodKnight 420S, 420E and 425 and KnightStar 330

# HEATED CPAP TUBING

## Heated Tube for Respironics

Item#


**TUB06HPR**

Heated Tube for  
Philips Respironics  
Dreamstation Series  
and System One

OEM# HT15

HPCPS# A4604

- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off



**TUBHC900ICON208**

Heated Tube for Fisher & Paykel ICON

OEM# 900ICON208

## Hybernite Heated CPAP Tube 6ft

Item#

**TUB006ROC**

*MOQ may apply*

**NEW**

- Copper wire embedded in the tubing wall generates uniform heat and maintains the air temperature inside the tube to prevent rainout.
- Stand-alone kit can be used in the sleep lab or at home with any CPAP or BiPAP machine with humidifier.
- Includes heated circuit tube and independent power supply.

Compatible with machines that use standard CPAP Tubing.  
Not compatible with machines that use slim 15mm tubing  
such as newer Respironics and ResMed machines.



**TUB006HBT**

Hybernite Replacement Heated CPAP Tube 6ft

## DreamStation 2

**NEW**



Item# **CF1012**  
Disposable Filter  
**6/pack**

OEM# 1142831  
HCPCS# A7038  
**CF1012-1** 1/pack  
**CF1012-2** 2/pack  
**CF1012-3** 3/pack  
**CF1012-50** 50/pack

Item# **CF1012F**  
Reusable Filter  
**2/pack**  
OEM# 1142687  
HCPCS# A7039  
**CF1012F-1** 1/pack

## DreamStation



Item# **CF1010**  
Disposable Filter  
**6/pack**

OEM# 1122519  
HCPCS# A7038  
**CF1010-1** 1/pack  
**CF1010-2** 2/pack

Item# **CF1010F**  
Reusable Filter  
**2/pack**  
OEM# 112446 (1-pack)  
HCPCS# A7039  
**CF1010F-1** 1/pack  
**CF1010F-6** 6/pack

## DreamStation Go Filters



Item# **CF101130D**  
30-Day Disposable Filter  
**6/pack**

OEM# 1133744  
HCPCS# A7038

Item# **CF101190D**  
90-Day Reusable Filter  
**2/pack**

OEM# 1133743  
HCPCS# A7039

## PR System One



Item# **CF1009NT**  
Ultrafine Filter - No Tab  
**6/pack** • 1 3/4" x 7/8"

OEM# 1029331  
HCPCS# A7038  
**CF1009NT-1** 1/pack  
**CF1009NT-2** 2/pack

Item# **CF1009**  
Ultrafine Filter - Tabbed  
**6/pack** • 1 3/4" x 7/8"

HCPCS# A7038  
**CF1009-1** 1/pack  
**CF1009-2** 2/pack

Item# **CF1009F**  
Foam Filter  
**2/pack** • 1 5/8" x 7/8"

OEM# 1029330  
HCPCS# A7039  
**CF1009F-1** 1/pack

## M-Series



Item# **CF1007**  
Ultrafine Filter - Tabbed  
**6/pack** • 1 3/4" x 7/8"

HCPCS# A7038  
**CF1007-1** 1/pack  
**CF1007-2** 2/pack

Item# **CF1007F**  
Foam Filter  
**2/pack** • 1 5/8" x 7/8"

HCPCS# A7039  
**CF1007F-1** 1/pack



## REMstar Pro, REMstar Plus



**Item# CF1006**  
Ultrafine Filter  
**2/pack** • 3 1/2" x 1 1/2"  
OEM# 1005945  
HCPCS# A7038

**Item# CF1006F**  
Foam Filter  
**2/pack** • 3 5/8" x 1 5/8"  
OEM# 1005964  
HCPCS# A7039

## BiPAP Pro, Plus, Duet LX, Synchrony



**Item# CF1005**  
Ultrafine Filter  
**6/pack** • 4 3/4" x 3 1/4"  
OEM# 622219  
HCPCS# A7038

**Item# CF1005F**  
Foam Filter  
**2/pack** • 4 3/4" x 2 3/8"  
OEM# 622220  
HCPCS# A7039

## Aria, Virtuoso, Duet, Quartet



**Item# CF1004**  
Ultrafine Filter  
**6/pack** • 3 1/2" x 2 3/16"  
OEM# 532311  
HCPCS# A7038

**Item# CF1004F**  
Foam Filter  
**2/pack** • 3 1/4" x 2 1/8"  
OEM# 532247  
HCPCS# A7039


## REMstar



**Item# CF1003**  
Ultrafine Filter  
**6/pack** • 5 1/2" x 4 3/4"  
OEM# 362522  
HCPCS# A7038

**Item# CF1003F**  
Foam Filter  
**1/pack** • 5 3/4" x 5 1/4"  
OEM# 362521  
HCPCS# A7039

## Solo, Aria LX, Virtuoso LX, REMstar LX



**Item# CF1002**  
Ultrafine Filter  
**6/pack** • 3 3/8" x 2 3/4"  
OEM# 622017  
HCPCS# A7038

**Item# CF1002F**  
Foam Filter  
**2/pack** • 3 5/16" x 2 5/16"  
OEM# 622018  
HCPCS# A7039

## SleepEasy/BiPAP



**Item# CF1001**  
Ultrafine Filter  
**6/pack** • 7 3/4" x 4 7/8"  
OEM# 302064  
HCPCS# A7038

*All filters can be ordered in any pack quantity.  
Contact your sales rep for details.*

## AirSense 11



Item#  
**CF2111**

Disposable Filter  
**3/pack**

HCPCS# A7038

**CF2111-1** 1/pack OEM# 39300

**CF2111-2** 2/pack

**NEW**

Item#

**CF2111HA-3**

Hypoallergenic Filter  
**3/pack**

1 1/4" x 1 3/16"

HCPCS# A7038

**CF2111HA-1** 1/pack

**CF2111HA-2** 2/pack

## AirMini



Item#

**CF2108**

Disposable Filter

**2/pack** • 1 1/4" x 1 3/16"

OEM# 38815

HCPCS# A7038

## Adapt SV



Item#

**CF2106**

Disposable Filter

**3/pack** • 1 1/4" x 1 1/2"

OEM# 26926

HCPCS# A7038

**CF2106-50** 50/pack

Item#

**CF2106HA**

Hypoallergenic Filter

**3/pack** • 1 1/4" x 1 1/2"

OEM# 26927

HCPCS# A7038

**CF2106HA-50** 50/pack

## S9 and AirSense 10 Series



Item#

**CF2107**

Disposable Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107-1** 1/pack OEM# 36850

**CF2107-2** 2/pack OEM# 36851

**CF2107-50** 50/pack OEM# 36853

Item#

**TUB06S9COMBO**

S9 Filter and CPAP Tube Kit

6 - CF2107

1 - 6ft CPAP Tube

**KIT**

Item#

**CF2107HA**

Hypoallergenic Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107HA-1** 1/pack OEM# 36855

**CF2107HA-2** 2/pack OEM# 36856

**CF2107HA-50** 50/pack

OEM# 36858

## Auto Set T



Item# **CF2104**

Disposable Filter


**3/pack** • 4 1/4" x 3 1/4"

OEM# 17937

HCPCS# A7038

**CF2104-50** 50/pack OEM# 17938

## S8 Series



Item#  
**CF2105**  
Disposable Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
OEM# 33916  
HCPCS# A7038

**CF2105-1** 1/pack OEM# 33915  
**CF2105-2** 2/pack OEM# 33973  
**CF2105-50** 50/pack OEM# 33918



Item#  
**TUB06S8COMBO**  
S8 Filter and CPAP Tube Kit  
6 - CF2105  
1 - 6ft CPAP Tube



Item#  
**CF2105HA**  
Hypoallergenic Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>8</sub>"  
OEM# 33920  
HCPCS# A7038  
  
**CF2105HA-1** 1/pack OEM# 33919  
**CF2105HA-2** 2/pack OEM# 33975  
**CF2105HA-50** 50/pack OEM# 33922

## S7 Series



Item#  
**CF2103**  
Disposable Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 30907  
HCPCS# A7038

**CF2103-1** 1/pack OEM# 30919  
**CF2103-2** 2/pack  
**CF2103-50** 50/pack OEM# 30908



Item#  
**CF2103HA**  
Hypoallergenic Filter  
**3/pack** • 1<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
HCPCS# A7038

**CF2103HA-1** 1/pack OEM# 30933  
**CF2103HA-50** 50/pack

## S6 Series



Item#  
**CF2102**  
Disposable Filter  
**3/pack** • 5<sup>3</sup>/<sub>4</sub>" x 5<sup>5</sup>/<sub>8</sub>"  
OEM# 21935  
HCPCS# A7038

**CF2102-50** 50/pack  
OEM# 21936

# FISHER & PAYKEL

## SleepStyle Series



Item#  
**CF4004-2**  
Disposable Filter  
**2/pack**  
HCPCS# A7038



## ICON Series



Item#  
**CF4003**  
Disposable Filter  
**4/pack**  
2 1/16" x 5/8"  
HCPCS# A7038




**CF4003-1** 1/pack  
OEM# 900ICON503

**CF4003-6** 6/pack  
OEM# 900ICON212

## HC230 Series, HC600 Series



Item#  
**CF4002**  
Disposable Filter  
**4/pack** • 2 3/4" x 7/8"  
OEM# 900HC240  
HCPCS# A7038  
**CF4002-2** 2/pack



## HC200, HC220, HC221



Item#  
**CF4001**  
Disposable Filter  
**4/pack** • 5" x 3/4"  
OEM# 900HC222  
HCPCS# A7038  
**CF4001-2** 2/pack



## 3B

### Luna G3

**NEW**



Item# **CF8006-6**  
Disposable Filter  
**6/pack**  
HCPCS# A7038  
**CF8006-1** 1/pack  
**CF8006-2** 2/pack



Item# **CF8006F-1**  
Reusable Filter  
**1/pack**  
HCPCS# A7039

### Luna and Luna II



Item#  
**CF8005**  
Disposable Filter  
**12/pack** • 2" x 1 13/16"  
HCPCS# A7038




Item#  
**CF8005F-1**  
Reusable Filter  
**1/pack** • 2" x 1 13/16"  
HCPCS# A7039  
**CF8005F-2** 2/pack

## RESmart



Item#  
**CF8004**  
Disposable Filter  
**3/pack** • 3 5/16" x 1 5/8"  
HCPCS# A7038





# AEIOMed

## Everest Aura, Everest III and Integrated



Item# **CF8001**  
Foam Filter  
**2/pack** • 2" x 1<sup>9</sup>/<sub>16</sub>"  
OEM# 100446  
HCPCS# A7039

# BEYOND

## ResPlus



Item# **CF9005F-2**  
**NEW** Reusable Filter  
**2/pack**  
HCPCS# A7039  
**CF9005F-1** 1/pack

# BREAS/GE

## iSleep 20 Series



Item# **CF7003**  
Ultrafine Filter  
**2/pack** • 2<sup>1</sup>/<sub>4</sub>" x 7<sup>7</sup>/<sub>8</sub>"  
OEM# 004153  
HCPCS# A7038

Item# **CF7003F**  
Foam Filter  
**Each** • 2<sup>3</sup>/<sub>4</sub>" x 1"  
OEM# 004154  
HCPCS# A7039

## PV10



Item# **CF7002**  
Ultrafine Filter  
**2/pack** • 2<sup>3</sup>/<sub>4</sub>" x 7<sup>7</sup>/<sub>8</sub>"  
OEM# 2B01975  
HCPCS# A7038

Item# **CF7002F**  
Foam Filter  
**Each** • 2<sup>5</sup>/<sub>8</sub>" x 7<sup>7</sup>/<sub>8</sub>"  
OEM# 2B01972  
HCPCS# A7039

## PV1



Item# **CF7001**  
Ultrafine Filter  
**2/pack** • 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
HCPCS# A7038

Item# **CF7001F**  
Foam Filter  
**3/pack** • 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
HCPCS# A7039

# CAREFUSION/VIASYS

## PureSom



Item#  
**CF9003**  
Ultrafine Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7511  
HCPCS# A7038



Item#  
**CF9003F**  
Foam Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7510  
HCPCS# A7039

## Orion and Pegasus




Item#  
**CF8002**  
Ultrafine Filter  
**6/pack** • 2<sup>7</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"  
OEM# 774688  
HCPCS# A7038



Item#  
**CF8002F**  
Foam Filter  
**Each** • 2<sup>7</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>4</sub>"  
OEM# 774689  
HCPCS# A7039

# COVIDIEN NELLCOR PURITAN BENNETT

## GoodKnight 314/318



Item#  
**CF6005**  
Ultrafine Filter  
**2/pack** • 6<sup>1</sup>/<sub>8</sub>" x 9<sup>9</sup>/<sub>16</sub>"  
OEM# S-126924-00  
HCPCS# A7038



Item#  
**CF6005F**  
Foam Filter  
**2/pack** • 6<sup>1</sup>/<sub>4</sub>" x 1"  
OEM# S-290609-00  
HCPCS# A7039

## GoodKnight 420/425



Item#  
**CF6003**  
Ultrafine Filter  
**2/pack** • 1<sup>3</sup>/<sub>16</sub>" x 3<sup>3</sup>/<sub>16</sub>"  
OEM# M-413950-04  
HCPCS# A7038



Item#  
**CF6003F**  
Foam Filter  
**Each** • 1<sup>3</sup>/<sub>16</sub>" x 3<sup>3</sup>/<sub>16</sub>"  
OEM# M-413950-02  
HCPCS# A7039

## GoodKnight 418 A/G/P/S




Item#  
**CF6002**  
Ultrafine Filter  
**4/pack** • 5<sup>3</sup>/<sub>16</sub>" x 3<sup>3</sup>/<sub>4</sub>"  
OEM# Y-413560-01  
HCPCS# A7038



Item#  
**CF6002F**  
Foam Filter  
**6/pack** • 5<sup>1</sup>/<sub>4</sub>" x 3<sup>3</sup>/<sub>4</sub>"  
OEM# Y-413560-02  
HCPCS# A7039

## GoodKnight 418



Item#  
**CF6001**  
Ultrafine Filter  
**4/pack** • 2<sup>1</sup>/<sub>4</sub>" x 1<sup>2</sup>/<sub>5</sub>"  
OEM# S-645030-00  
HCPCS# A7038



Item#  
**CF6001F**  
Foam Filter  
**6/pack** • 2<sup>7</sup>/<sub>16</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# S-645029-00  
HCPCS# A7039

## IntelliPAP



Item#  
**CF5004**  
Ultrafine Filter  
**4/pack** • 1<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
OEM# DV51D-603  
HCPCS# A7038

Item#  
**CF5004F**  
Foam Filter  
**4/pack** • 1<sup>5</sup>/<sub>8</sub>" x 1<sup>3</sup>/<sub>16</sub>"  
OEM# DV51D-602  
HCPCS# A7039

## Horizon LT



Item#  
**CF5003**  
Ultrafine Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-603  
HCPCS# A7038

Item#  
**CF5003F**  
Foam Filter  
**4/pack** • 4<sup>1</sup>/<sub>4</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
OEM# 8000D-602  
HCPCS# A7039

# Evo MEDICAL

## RemRest



Item#  
**CF8003F**  
Foam Filter  
**3/pack** • 2<sup>3</sup>/<sub>16</sub>" x 1<sup>1</sup>/<sub>16</sub>"  
OEM# 21117  
HCPCS# A7039

## Comfort PAP



Item#  
**CF9003**  
Ultrafine Filter  
**3/pack** • 1<sup>3</sup>/<sub>16</sub>" x 1<sup>7</sup>/<sub>16</sub>"  
OEM# 7511  
HCPCS# A7038

Item#  
**CF9003F**  
Foam Filter  
**3/pack** • 1<sup>3</sup>/<sub>8</sub>" x 1<sup>1</sup>/<sub>4</sub>"  
OEM# 7510  
HCPCS# A7039

Enhance your branding!

## FREE PRIVATE LABELING

Sunset's free service allows you to put your logo and contact information on our products.

Increase brand awareness and allow patients to easily contact you.

Work with our marketing team to customize the layout with your logo, phone number, website address, and reorder information.

Available on all CPAP tubes and filters.  
Minimum orders apply.  
Contact your sales rep for details.




# HEALTHDYNE

## Quest



Item# **CF3003**  
Ultrafine Filter  
with Curves  
**5/pack** • 4 1/2" x 4 3/8"  
OEM# 7302  
HCPCS# A7038



Item# **CF3003F**  
Foam Filter  
with Curves  
**5/pack** • 4 1/2" x 4 1/4"  
OEM# 7070  
HCPCS# A7039

## Tranquility MPR



Item# **CF3002**  
Ultrafine Filter  
Rectangle/Curve  
**5/pack** • 4 5/16" x 1 7/8"  
OEM# 7270  
HCPCS# A7038



Item# **CF3002F**  
Foam Filter  
**5/pack**  
5" x 2 1/4" x 2 1/8"  
OEM# 7301  
HCPCS# A7039

## Tranquility



Item# **CF3001**  
Ultrafine Filter  
**5/pack** • 5 1/8" x 1 7/8"  
OEM# 7170  
HCPCS# A7038

# RESVENT

## iBreeze



HCPCS# A7038

Ultrafine Filter  
**CF9901-1** 1/pack  
**CF9901-2** 2/pack  
**CF9901-6** 6/pack  
**CF9901-50** 50/pack

# INVACARE

## Polaris LT



Item# **CF9002**  
Ultrafine Filter  
**3/pack** • 5" x 2"  
OEM# 1082428  
HCPCS# A7038



Item# **CF9002F**  
Foam Filter  
**3/pack** • 5" x 2"  
OEM# 1082995  
HCPCS# A7039

## Polaris EX




Item# **CF9001**  
Ultrafine Filter  
**3/pack** • 2 3/16" x 1"  
HCPCS# A7038



Item# **CF9001F**  
Foam Filter  
**3/pack** • 2 1/4" x 1 1/16"  
OEM# 1124211  
HCPCS# A7039

## Polaris TR



Item# **CF8001**  
Foam Filter  
**2/pack** • 2" x 1 9/16"  
OEM# 100446  
HCPCS# A7039

# PROBASICS

## Zzz-PAP



Item# **CF9003**  
Ultrafine Filter  
**3/pack** • 1 3/16" x 1 7/16"  
OEM# 7511  
HCPCS# A7038



Item# **CF9003F**  
Foam Filter  
**3/pack** • 1 3/16" x 1 7/16"  
OEM# 7510  
HCPCS# A7039



# FISHER & PAYKEL

## Humidifier Chamber for ICON Series



Item# **HUMHC900ICON200**

Replacement Chamber for the Fisher & Paykel ICON Series

Each OEM# 900ICON200 • HCPCS# A7046

Compatible with:

- ICON Premo
- ICON Novo
- ICON Auto



## Humidifier Chambers for HC200 Series



Item# **HUMHC355**

Extended Life Humidifier Chamber for the Fisher & Paykel HC200 Series

Each OEM# HC 355



Item# **HUMHC385**

Replacement Chamber for the Fisher & Paykel SleepStyle HC200 Series

Each OEM# HC 385S

HCPCS# A7046



Compatible with:

• HC200	• HC231	• HC236
with: • HC201	• HC232	• HC238
• HC220	• HC233	• HC242
• HC221	• HC234	• HC244
		• HC254

## Humidifier Chambers for HC100 Series



Item# **HUMHC300**

Extended Life Humidifier Chamber for the Fisher & Paykel HC100 Series

Each OEM# HC 300



Item# **HUMHC325S**

Replacement Chamber for the Fisher & Paykel HC100 Series

Each OEM# HC 325S

HCPCS# A7046

Compatible with:

- HC100
- HC105
- HC150



## Humidifier Chamber for HC600 Series



Item# **HUMHC365**

Replacement Chamber for the Fisher & Paykel SleepStyle 600 Series

Each OEM# HC 365S • HCPCS# A7046

Compatible with:

- HC604
- HC608



# RESPIRONICS

## Humidifier Chamber for PR System One



Item# **HUM1136502**

Humidifier Chamber for PR System One 60 Series

Each HCPCS# A7046

## Humidifier Chamber for DreamStation



Item# **HUM1122520**

Humidifier Chamber for DreamStation

Each HCPCS# A7046

**NEW**

# DEVILBISS

## Humidifier Chamber for IntelliPAP



Item# **HUMDV5C**

Humidifier Chamber for DeVilbiss IntelliPAP Series

Each HCPCS# A7046

# SUNSET CHINSTRAPS

All chinstraps are latex-free • Hand wash only • Lay flat to dry

## Universal Neoprene Chinstrap



Item#

**CS001**

Universal Chinstrap

HCPSC# A7036



- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

## Basic Chinstrap



Item# **CS002**

Compare to  
Respironics Basic

OEM# 302175

HCPSC# A7036



- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

## Adjustable Neoprene Chinstrap



Size

Item#

6in

**CS003S**

7in

**CS003R**

8in

**CS003M**

9in

**CS003L**

10in

**CS003XL**

Compare to Puresom Ad-  
justable Ruby

OEM# TMS-09 ADJ

HCPSC# A7036

For fitting, measure along one side of the head, as indicated by the dashed line above.

- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head.
- Great for hard-to-fit patients

## Deluxe Chinstrap



Size

Item#

26in

**CS004S**

28in

**CS004L**

Compare to  
Respironics Deluxe

OEM# 302425

HCPSC# A7036



- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

## Comfort Chinstrap



Item# **CS005**

Compare to  
ResMed Sullivan

OEM# 16015

HCPSC# A7036



- One size fits most
- Easy to use —Attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

## Premium Chinstrap



Item# **CS006**

Compare to  
Respironics Premium

OEM# 1012911

HCPSC# A7036




- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

# SUNSET CHINSTRAPS AND HEADGEAR

All chinstraps and headgear are latex-free • Hand wash only • Lay flat to dry

## Ruby Style Chinstrap



Size	Item#
21.5-22.5in	<b>CS007S</b>
22-23in	<b>CS007M</b>
23-24in	<b>CS007L</b>

Compare to Puresom Ruby  
OEM# TMS-09 • HCPCS# A7036

- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

## Halo Style Chinstrap



Item# **CS025**  
Compare to BreatheWear Halo Chinstrap  
OEM# H3407  
HCPCS# A7036

## Ruby Style Adjustable Chinstrap



Size	Item#
Regular	<b>CS007-ADJ</b>
Extra Large	<b>CS007-ADJXL</b>

HCPCS# A7036

- Fully adjustable
- Super soft materials for extra comfort



## Universal Mesh Headgear



Item# **HG001**  
Compare to Respironics SimpleStrap Headgear  
OEM# 302328  
HCPCS# A7035

- One size fits most
- Can be used with most 4-point nasal or full face masks



## Universal Headgear



Size	Length	Item#
Regular	22.5in	<b>HG003</b>
Extra Large	24.5in	<b>HG003XL</b>

HCPCS# A7035

- Can be used with most 4-point nasal or full face masks



# CLEANERS AND ACCESSORIES

## Sunset Cleaning Wipes — Flowpack

Item#  
**CAP1003S**

**12/case**  
Flowpack of 64 wipes

**CAP1003S-80**

**12/case**  
Flowpack of 80 wipes

- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks and equipment



## Private Label CPAP Wipes

Order the flowpack wipes with your logo and contact info on a lid sticker, with a minimum quantity commitment. Talk to your sales rep for more info.



## Citrus II Cleaning Wipes

Item# **CAP1003**  
Canister of 62 wipes  
**12/case**



## Cleaning Wipes — Single Packs

**CAP1005**  
Citrus II Box of 12 single-pack wipes  
**12/case**

- Convenient individually wrapped wipes are great for patient samples or travel use



## CleanSmart CPAP Disinfectant Spray

Item#  
**CAP1018-L**  
16oz • **6/case**

**CAP1018-S**  
2oz • **9/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Active ingredient Hypochlorous Acid does not leave a chemical residue
- EPA Registration Number 89896-2
- Formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



## Citrus II Mask Cleaning Spray

Item# **CAP1001**  
1.5oz • **24/case**

Item# **CAP1002**  
8oz • **12/case**

- All-natural cleanser for daily mask cleaning
- Available in 1.5oz travel size or 8oz non-aerosol spray cans






# CLEANERS AND ACCESSORIES

## CPAP Tube Cleaning Brush


Size	Item#
Standard 6ft	<b>CAP1013</b>
Slim 6ft	<b>CAP1013SS</b>

- Brushes all the way through a 6ft CPAP tube without leaving debris in the middle
- Stainless steel wire and soft nylon bristles



## Filter Replacement Kit for SoClean 2

Item# **CAP1007S-KIT**  
Each



- With filter, check valve, retail packaging and instruction manual
- Replace every 6 months

## Cleaner Accessories

<b>ZCD100KIT</b>	Filter and Check Valve Kit
<b>CAP1016S</b>	Quick Release Connector
<b>ZCD100A-DS</b>	Heated Tube Adapter for DreamStation/PR System One
<b>ZCD100A-A10</b>	Heated Tube Adapter Lid for ResMed AirSense 10
<b>ZCD100A-S9</b>	Heated Tube Adapter Lid for ResMed S9
<b>ZCD100TC</b>	Replacement Tube Connector (used when cleaning machine)
<b>ZCD100A-P</b>	Replacement Ozone Plug Adapter (used when not cleaning machine)
<b>ZCD100IT</b>	Replacement Injection Tube with Check Valve
<b>ZCD100F</b>	Replacement Filter
<b>ZCD100PA</b>	Replacement Power Adapter

# CPAP BATTERY

## Portable Outlet UPS CPAP Battery

Item#

**CAP10082**

Includes battery, wall charger, and travel bag

- Autoswitch technology: When power goes out, your CPAP stays on
- Three-prong AC outlet
- Two USB ports
- Digital display
- Convenient travel bag
- Meets TSA/FAA battery guidelines
- 1-year warranty
- Rated for over 500 charge cycles
- Recharges in less than 4 hours
- May be used for other devices powered by AC or USB
- AC and USB ports may be used simultaneously

### Specs

- 2.875 lbs
- 5.375 in wide x 9 in long x 1.5 in tall
- Rated Capacity: 159Wh/14.4Ah
- Rated Voltage: 14.5V
- Charge Input: 16.8V/2A
- USB Output: 5V/4.2A
- AC Output: 110v AC 60hz/100w
- AC Adapter Input: 100V-240V/50-60Hz



## CPAP HCPCS AND REPLACEMENT SCHEDULE

HCPCS	Item	Recommended replacement schedule
A4604	Heated Circuit Tubing	1 per 3 months
A7027	Combination Oral/Nasal Mask	1 per 3 months
A7028	Replacement Oral Cushion for Combination Oral/Nasal Mask	2 per month
A7029	Replacement Nasal Pillows for Combination Oral/Nasal Mask — Pair	2 per month
A7030	Full Face Mask	1 per 3 months
A7031	Replacement Full Face Mask Cushion or Liner	1 per month
A7032	Replacement Nasal Cushion or Liner	2 per month
A7033	Replacement Pillows — Pair	2 per month
A7034	Nasal or Pillows Mask	1 per 3 months
A7035	Headgear	1 per 6 months
A7036	Chinstrap	1 per 6 months
A7037	Tubing	1 per 3 months
A7038	Filter — Disposable	2 per month
A7039	Filter — Nondisposable	1 per 6 months
A7044	Oral Interface	
A7045	Replacement Exhalation Port	
A7046	Replacement Water Chamber	1 per 6 months
A9279	Monitoring Feature/Device	
E0470	BiLevel Device without Backup Rate	
E0471	BiLevel Device with Backup Rate	
E0561	Humidifier – Non-Heated	
E0562	Humidifier – Heated	
E0601	CPAP Device	

*Due to the rapid changes in reimbursement, this information is for reference only. Sunset Healthcare Solutions cannot guarantee its comprehensiveness or accuracy. Sunset urges you to contact a reimbursement specialist for up-to-date information and advice.*

# CASH ITEMS


## Snugz CPAP Mask Liners

Item#  
**CAP7001**

Machine-Washable Liners  
for Most Full Face Masks  
**2/pack**

**CAP7002**

Machine-Washable Liners  
for Most Nasal Masks  
**2/pack**



- 2 per pack, 90-day supply
- Machine washable
- Fits MOST nasal or full face CPAP masks (Not recommended for ResMed Mirage Quattro or Respironics Amara View and Dreamwear)
- Blend of spandex, moisture-wicking polyester and lycra fibers
- Latex-free

Snugz is a machine washable, one-size-fits-most liner for CPAP masks that is designed to help reduce noisy leaks, prevent redness and irritation, and feel great against the skin. The liner helps to grip the skin and keep the mask in place around the nose and mouth, reducing the need for a chinstrap. The improved seal is especially beneficial for high-pressure/high-flow PAP patients and can reduce nightly sleep disturbances and AHI.


## Soft Headgear Strap Covers

Item#  
**CAP2004**

**NEW**



- Reduces strap marks
- No bulky and uncomfortable velcro
- Reduces strap irritation on sensitive skin
- Thick padding for extra comfort
- No need for clips – just slip wraps around headgear straps
- Easy and helpful for patients with arthritis in hands or fingers
- Can be machine washed and dried
- Wraps securely grip headgear and slide to adjust
- One size fits all
- Made from 230g Coral Velvet



## RemZzzs CPAP Mask Liners

Item#

**CAP3001(S,M,L)**

Liners for Full Face ResMed  
& Respirationics

OEM#	K1-FS	SMALL
	K2-FM	MEDIUM
	K3-FL	LARGE

**CAP3002(S,M,L)**

Liners for Full Face Fisher &  
Paykel, Hans Rudolph, and  
Sunset HCS

OEM#	K4-FS	SMALL
	K5-FM	MEDIUM
	K6-FL	LARGE

**CAP3003(S,M,L)**

Liners for Nasal  
Respironics

OEM#	K7-NS	SMALL
	K8-NM	MEDIUM
	K9-NL	LARGE

**CAP3004(S,M,L)**

Liners for Nasal ResMed  
and Fisher & Paykel

OEM#	K10-NS	SMALL
	K11-NM	MEDIUM
	K12-NL	LARGE

**CAP3005L**

Liners for Respironics  
FitLife Mask

OEM#	6A-FLK	SMALL AND LARGE
	6B-FLK	EXTRA LARGE

**CAP3006**

Liners for Nasal Pillows

OEM# 7A-NPK

**CAP3007**

Liners for Respironics  
DreamWear and Amara View

OEM# 6C-AVK

### 30-Day Supply

Boxes sold individually.

Full face and nasal mask RemZzzs include 30 facial liners and 2 forehead liners.

- Increases compliance by adding a comfortable and absorbent pad to CPAP mask
- Reduces air leaks and protects skin from red marks





# CASH ITEMS

## CPAP Moisture Therapy

Item# **CAP6100L**  
1oz Tube • 18/cas

- Formulated to prevent chafing, irritation, dryness or cracking from using a CPAP mask or other delivery device
- Apply to the face at the beginning of treatment and as needed
- Petroleum-free with aloe vera, vitamins A & E, and argan oil



## Pur-Sleep CPAP Aromatherapy

Item# **CAP5001**  
VaporClear Aromatherapy Starter Kit  
OEM# VAPOR

### Kit Includes:

- 1 aluminum diffusion stand
- 4 reusable diffusion pads
- 1.5ml bottle of Clear (peppermint, lime)
- 1.5ml bottle of Calm (eucalyptus, lavender, rosemary)



### 30ml Aromatic Refills

Scent/OEM#	Item#
Black Cherry	<b>CAP5004-BCH</b>
Calm	<b>CAP5004-CLM</b>
Clear	<b>CAP5004-CLR</b>
Dark Chocolate	<b>CAP5004-DCH</b>
Deep	<b>CAP5004-DEE</b>
French Vanilla	<b>CAP5004-FRV</b>
Fresh	<b>CAP5004-FSH</b>
Peace	<b>CAP5004-PCE</b>
Spice	<b>CAP5004-SPC</b>
Strawberry	<b>CAP5004-STR</b>

Item# **CAP5002**  
Reusable Diffusion Pads  
10/pack • OEM# DP10

Item# **CAP5003**  
Replacement Aluminum Diffusion Tray  
Each • OEM# TRAY

- Aromatherapy products specifically for use with CPAP
- Set up the diffuser right below the CPAP filter intake, add a drop of essential oil to the pad and the fragrance is drawn into the CPAP machine



## H2Doze Distilled Water for CPAP


Item# **CAPH20-6**  
16.9oz Bottles  
6/pack

- TSA cleared bottle when accompanied with CPAP machine
- Convenient 16.9 oz bottle perfect for travel
- Enough water for multiple nights of use
- Distilled water designed to protect CPAP equipment
- Allows patient to breathe healthy, safe air throughout the night



## Quick Tube Release Connector

Item# **CAP1016**  
Each



## CPAP Final Bacteria Filter

Item# **BF2200**  
22mm Bacteria Disk Filter  
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

## Circuit Connectors

Item# **RES013**  
22mm OD x 15mm ID/  
22mm OD  
10/pack



Item# **RES014**  
22mm OD x 22mm OD  
With Lip  
10/pack



## Humidifier Tubing Adapter

Item# **RES016**  
22mm ID x 22mm OD  
with Lip  
10/pack



## Oxygen Bleed-In Adapter

Item# **RES019** 10/pack  
22mm with 6mm  
Rotating Stem

**RES019-CAP** 10/pack  
With Stem Cap



## Swivel Connector Silencer

Item# **VEN002**  
For CPAP and Ventilator Circuits



# CASH ITEMS

## CPAP Pillows

Item# **CAP4003**

CPAPMax Pillow

Each OEM# 15-501R-1-920



**CPAPMax Replacement Cover**

**CAP4003RC**

Each OEM# 1-630-101-1-933

**Replacement Pillow Case**

**CAP4003C**


Each OEM# 1-626-900-1-02



**Cool 3D Mesh Air-Flow**  
Reduces heat and perspiration



**Fiberfill**  
Soft microfiber pillow top with hypoallergenic fiberfill



**Charcoal Infused High Resiliency Foam**

Charcoal neutralizes odors for freshness every night

### BEFORE

Traditional pillows interfere with mask performance and can cause a stiff neck, shifting, leaks, and facial discomfort.



### AFTER

The CPAPMax Pillow offers superb mask accommodation, excellent support, and unparalleled comfort.



## CPAP Bedside Storage Bag

Item# **CAP1015**



- Well-made and stylish
- Designed to hold CPAP supplies such as a mask and tube
- Weighted lid allows bag to hang off top of nightstand or tabletop



## CPAP Holder

Item# **CAP1010**



- Sits between mattress and boxspring
- Saves space by allowing machine to rest in the area between bed and nightstand
- CPAP machine is kept close to patient's head, allowing more freedom to move throughout the night
- Keeps machine below mask, to help prevent rainout



# CASH ITEMS

## Handheld CPAP Pressure Meter

Item#  
**CAP003**

- Measures flow 0-18 cm/H2O
- Connects directly to the CPAP tubing without an adapter



## Comfort Tubing Wrap

Item#  
**CAP2002**  
6ft Nylon  
Tubing Wrap  
**Each**  
Compare to  
ResMed  
Tubing Wrap



- 6ft nylon wrap with fleece lining maximizes heat retention and reduces noise
- Zipper closure for ease of use

## Deluxe Tubing Wrap

Item#  
**CAP2001**  
6ft Velour  
Tubing Wrap  
**Each**  
Compare to  
Respironics  
Tubing Wrap



- 6ft plush velour material adds comfort while preventing rainout
- Zipper closure for ease of use

**Both Tubing Wraps come in attractive retail packaging**



Item#  
**CAP6100L**  
1oz Tube  
**18/case**

- Prevent chafing, irritation, dryness or cracking from using a CPAP or oxygen delivery device
- Apply to the face at the beginning of treatment
- Petroleum-free with Aloe Vera, Vitamins A & E, and Argan oil



Item#  
**RES6104L**  
1oz Tube  
**18/case**



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)



**Oxygen****Page**

Pulse Oximeter	40
Firesafe	41
Concentrators	42
Batteries	42
Masks	43
Cannula	44-45
Tubing	45
Regulators and Conserver	46
Cylinder Carriers	48-49
Connectors and Adapters	50
Accessories	50-52
Cleaners	52
Bacteria and HEPA Filters	53
Concentrator Cabinet Filters	54-59
Caire/AirSep/SeQual	54
CHAD	55
Covidien Nellcor Puritan Bennett	55
DeVilbiss/Drive	56
Inogen	57
Invacare	57
Jumao, Lifestyle, Attendant	58
Nidek	58
3B	58
Respironics/Healthdyne	59

**OXYGEN  
PRODUCTS****Finger Pulse Oximeter**Item# **RES5102**

- Oversized 1.3" OLED display face is the largest of any finger pulse oximeter
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual

**NEW**



# FIRESAFE

## Firesafe Thermal Shut-Off Valve

Item#  
**RES010**  
Each



Bidirectional



**NEW**  
**IMPROVED**

- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years

- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

### Fatalities per 100,000 patients



\* 1 death was reported among 73% of the patient population (2013-17)

## FireSafe Kits

Item#  
**RES010KIT**  
2 Firesafe Valves  
1 Tubing Connector

Item#  
**RES010CANNULAKIT**  
2 Firesafe Valves  
1 Tubing Connector  
1 Adult Soft Cannula



## FireSafe Tubing Connector

Item#  
**RES010TC** 100/case  
**RES010TC-10** 10/case



## FireSafe Nozzle

Item#  
**RES010FND-20**  
20/pack



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS

FIRESAFE

# OXYGEN CONCENTRATORS

## Drive DeVilbiss Oxygen Concentrators

Item#

**RES525DS**

**5 Liter**

Item#

**RES1025DS**

**10 Liter**

- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter and 10 Liter Oxygen Concentrators were designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device



## Drive DeVilbiss iG02 Portable Oxygen Concentrator

Item#

**RES125D**

**iG02**

Item#

**RES125D-XB**

**iG02 with Extra Battery**

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market



# OXYGEN BATTERIES

## Zopec OxygenPro 4XL & 8XL

Item#

**CAP1024-4XL**

**4XL**

Item#

**CAP1024-8XL**

**8XL**

**NEW**

- Provide ~4 or ~8 hours backup at 5 L/M continuous flow for Oxygen Concentrators.
- Relax with powerful new backup for oxygen concentrators—and more.
- Batteries fully charge in 10 hours.



# OXYGEN MASKS

## Simple O<sub>2</sub> Mask

With 7ft Safety Tubing

Item# **RES2100**

Adult Size

50/case

Item# **RES2200V**

Pediatric Size

50/case



HCPCS# A4620

- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## Partial Rebreather O<sub>2</sub> Mask

With 7ft  
Safety Tubing  
& Reservoir Bag

- 5-15 LPM
- High concentration

Item# **RES2101V**

Adult Size

50/case

HCPCS# A4620



## Non-Rebreather O<sub>2</sub> Mask

With 7ft Safety Tubing  
& Reservoir Bag

- 10+ LPM
- Low-resistance  
check valves prevent  
rebreathing and allow  
exhaled gas to escape

Item# **RES2102**

Adult Size

50/case

Item# **RES2202V**

Pediatric Size

50/case



HCPCS# A4620



# OXYGEN CANNULA

## Cannula with Tubing

### Adult Cannula with Tubing

Length	Material	Item#	Qty
4ft	Standard	<b>RES1104</b>	<b>50/case</b>
4ft	Soft	<b>RES1104S</b>	<b>50/case</b>
7ft	Standard	<b>RES1107</b>	<b>50/case</b>
7ft	Soft	<b>RES1107S</b>	<b>50/case</b>
15ft	Standard	<b>RES1115</b>	<b>50/case</b>
25ft	Standard	<b>RES1125</b>	<b>25/case</b>
25ft	Soft	<b>RES1125S</b>	<b>25/case</b>

HCPCS# A4615



- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort

Item#

**RES1207V**

Pediatric Cannula  
with 7ft Supply Tube  
**50/case**



**RES1307V**

Infant Cannula  
with 7ft Supply Tube  
**50/case**



## High Flow Soft Cannula

Length	Item#	Qty
4ft	<b>RES1104SHF</b>	<b>25/case</b>
7ft	<b>RES1107SHF</b>	<b>25/case</b>
15ft	<b>RES1115SHF</b>	<b>25/case</b>
25ft	<b>RES1125SHF</b>	<b>25/case</b>

HCPCS# A4615



- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

## Cannula without Tubing


### Adult Cannula without Tubing

Port	Material	Item#	Qty
Male	Standard	<b>RES1100</b>	<b>50/case</b>
Male	Soft	<b>RES1100S</b>	<b>50/case</b>
Female	Soft	<b>RES1100FS</b>	<b>50/case</b>

HCPCS# A4615



Male port



Female port

## Soft Cannula with Universal Connector

Item# **RES1107SUC**

Adult Soft Cannula with  
Universal Connector  
and 7ft Supply Tube  
**50/case**

HCPCS# A4615




**NEW**

## Dual-Lumen (Demand) Cannula

Length	Item#	Qty
4ft	<b>RES4144V</b>	<b>25/case</b>
7ft	<b>RES4147V</b>	<b>25/case</b>

HCPCS# A4615



- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

## Tabbed Cannula

Item# **RES1007V**

Adult Tabbed Cannula  
with 7ft Supply Tube  
**50/case**

HCPCS# A4615






# OXYGEN CANNULA AND TUBING

## Humidoflo High Flow Nasal Cannula

10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>

**NEW**




- Includes standard ISO connector and F&P adaptor so it can be used with almost any machine on the market
- Kink-free prongs and flexible, low-noise tubing
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)
- Can be used with F&P 950, 850, or Airvo system

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O2 concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

## Secure Cannula Stickers

Item# **RES080**  
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



## Comfort Cannula with Ear Cushions

Item# **RES1107EC**  
7ft Cannula with  
Ear Cushions  
**25/case**  
HCPCS# A4615



- 1-6 LPM
- Soft foam cushions add comfort and increase patient satisfaction
- Cushions are added during the manufacturing process so they stay in place and no assembly is required

Item# **RES020**  
Cushions Only  
**50 Pairs/case**

- Add cushions to existing cannula



## OXYGEN TUBING

### Oxygen Supply Tube



Length	Item#	Qty
7ft	<b>RES3007</b>	<b>50/case</b>
15ft	<b>RES3015</b>	<b>50/case</b>
25ft	<b>RES3025</b>	<b>25/case</b>
40ft	<b>RES3040</b>	<b>20/case</b>
50ft	<b>RES3050</b>	<b>20/case</b>

HCPCS# A4616

- Crush resistant with low memory to prevent kinking
- Great durability and comfort level
- Latex-free
- 1-6 LPM
- 6 channels

### Oxygen Supply Tube — Green



Length	Item#	Qty
7ft	<b>RES3007G</b>	<b>50/case</b>
15ft	<b>RES3015G</b>	<b>50/case</b>
25ft	<b>RES3025G</b>	<b>25/case</b>
40ft	<b>RES3040G</b>	<b>20/case</b>
50ft	<b>RES3050G</b>	<b>20/case</b>

HCPCS# A4616

Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS

# OXYGEN REGULATORS AND CONSERVER

## Mini Oxygen Regulator CGA 870 Style



"T" Handle for easy mounting — no tools needed

Barbed hose outlet for easy connection to the oxygen tubing

For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.

### Item# RES042

Adult • Barbed  
0-8 LPM • **Each**  
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

### Item# RES043

Adult • Barbed  
0-15 LPM • **Each**  
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

### Item# RES046

Pediatric • Barbed  
0-4 LPM • **Each**  
0,  $\frac{1}{32}$ ,  $\frac{1}{16}$ ,  $\frac{1}{8}$ ,  $\frac{1}{4}$ ,  $\frac{1}{2}$ ,  $\frac{3}{4}$ , 1, 1  $\frac{1}{2}$ , 2, 3, 4 LPM



HCPCS# E1353

## Mini Oxygen Regulator CGA 540 Style



Threaded DISS outlet

For use with M60 and larger cylinders: M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.

### Item# RES044

Adult • DISS  
0-8 LPM • **Each**  
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

### Item# RES045

Adult • DISS  
0-15 LPM • **Each**  
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

### Item# RES047

Pediatric • DISS  
0-4 LPM • **Each**  
0,  $\frac{1}{32}$ ,  $\frac{1}{16}$ ,  $\frac{1}{8}$ ,  $\frac{1}{4}$ ,  $\frac{1}{2}$ ,  $\frac{3}{4}$ , 1, 1  $\frac{1}{2}$ , 2, 3, 4 LPM



HCPCS# E1353

### Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty



**Private label laser etching now available!**

Ask your rep for details.

## Oxygen Regulator CGA 540 Style Long Stem



### Item# RES045L

Adult • DISS  
0-15 LPM • **Each**  
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

## Bonsai Pneumatic Conserver

### Item# RES049

Single Lumen • Barbed Outlet  
**Each**




- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5, and 6
- Two continuous flow (CF) settings preset @ 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses breath/delivers oxygen in first half of inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy to operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag
- Limited 2-year warranty

## Sample Oxygen Cylinder Specifications\*

Use the handy chart below to determine which accessories to purchase for your oxygen cylinder.

\* This information is intended to be used as a guide. Dimensions may vary by manufacturer. Sunset Healthcare Solutions does not sell oxygen tanks.



Older Name		A	B			C	D	JD	E		
Newer Name	M2	M4	M6	ML6	M7	M9	M15	M22	M24	M60	M/MM/M122
Diameter (in.)	2.5	3.2	3.2	4.3	4.3	4.3	4.3	5.3	4.3	7.3	8
Height (in.)	5.3	8.5	11.5	7.6	9.1	11	16.5	16.5	25.5	23	36
Empty Weight (lb)	0.7	1.6	2.2	2.8	3.3	3.7	5.3	8	7.9	22.3	39.5
Capacity (L)	42	113	164	165	198	255	425	640	680	1738	3455
Transport Method	CSBM6, CSBM6BP	CSBM6, CSBM6BP	CSBM6, CSBM6BP, CSBM6Flex, CSBM6-Carrier	CSBML6M7, CSBML6-Hardshell	CSBML6M7	CSBM9, CSBM6BP, CSBM6Flex, RES038B	CSBD, CSBDDBP, CSBD1Roller, CSBD2Roller, RES038B, RES060		CSBE, CSBD1Roller, CSBD2Roller, RES038B, RES060	Not Portable	Not Portable
Regulator Type	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES044, RES045	RES044, RES045

## OXYGEN HCPCS

A4615	Cannula, nasal
A4616	Oxygen tubing, per foot
A4620	Oxygen mask, variable concentration
E0424	Stationary compressed gaseous oxygen system, rental
E0431	Portable gaseous oxygen system, rental
E0443	Portable oxygen contents, gaseous (for use only with portable gaseous systems when no stationary gas or liquid system is used), 1 month supply = 1 unit
E0445	Oximeter device for measuring blood oxygen levels non-invasively
E0560	Humidifier, durable for supplemental humidification during IPPB treatment or oxygen delivery
E1353	Oxygen supply regulator
E1354	Oxygen supply wheeled cart
E1355	Oxygen supply stand/rack
E1390	Oxygen concentrator, single delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1391	Oxygen concentrator, dual delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1392	Portable oxygen concentrator, rental
K0738	Portable gaseous oxygen system, rental; home compressor used to fill portable oxygen cylinders; includes portable containers, regulators, flow meter, humidifier, cannula or mask and tubing

# CYLINDER BAGS

## Comfort Shoulder Bags

Cylinder	Item#
ML6/M7	<b>CSBML6M7</b>
C/M9	<b>CSBM9</b>
B/M6	<b>CSBM6</b>
D	<b>CSBD</b>



**CASH ITEM**

- Comfort flex shoulder strap
- Double zipper for easy access
- Machine washable



CSBML6M7



CSBM9



CSBM6



CSBD

## Flex 02 Bag

Item#  
**CSBM6Flex**  
Reconfigurable Bag  
for M6 and M9  
Cylinders  
**Each**

- Adjusts to 5 carrying modes



**CASH ITEM**



## Backpack

Cylinder	Item#
B/M6/C/M9	<b>CSBM6BP</b>
D	<b>CSBDBP</b>



**CASH ITEM**



- Fully adjustable and comfortable backpack straps
- Machine washable
- Converts to either backpack or shoulder bag

## Wheelchair/Scooter Carrier

Item# **CSBE**  
For E Cylinders  
**Each**



**CASH ITEM**



- Multi-part stabilizer straps for easy chair attachment
- Machine washable



# CYLINDER CARRIERS

## Single Tank Rolling Bag/Backpack

Item#  
**CSBD1Roller**  
For D/E Cylinders  
Each



**CASH ITEM**



## Double Tank Rolling/Shoulder Bag

Item#  
**CSBD2Roller**  
For D/E Cylinders  
Each



**CASH ITEM**

- May also be used with shoulder strap



## ML6 Hardshell Tote

Item#  
**CSBML6-Hardshell**  
For Homefill Cylinders  
Each



**CASH ITEM**



## Cylinder Rack

Item# **RES038B**  
Each  
For D, E, or M9 Cylinders —  
Lightweight, Durable  
Plastic Composite Rack



## Oxygen Cylinder Carrier

Item#  
**CSBM6-Carrier**  
Each  
Soft Polyester Tote  
for M6 Cylinders



- Tote handles
- Padded shoulder strap
- Fastening strap for mobility
- Holds up to 8 M6 tanks

## D/E Cart

Item#  
**RES060**  
4/case

- Accommodates one D or E tank
- Sturdy yet lightweight (4.18 LBS welded metal construction)
- Adjustable height handle with PVC grip
- Large EVA Foam wheels have just the right amount of resilience for good traction and maneuverability
- Elegant, grey satin finish
- Bottom fork rest with rubber caps



# OXYGEN CONNECTORS AND ACCESSORIES

## Christmas Tree Connector

Item# **RES001**  
Green without Swivel  
**10/pack**

Item# **RES002**  
Green with Swivel  
**10/pack**



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Accessory Kit

Item# **RES1000**  
RES002 Christmas Tree Connector - Swivel  
RES012 Tubing Connector - Non-Swivel  
RES017 Tubing Y-Connector  
RES018 Tubing Connector - Swivel  
RES032 Wrench  
RES035 Brass Washer

Each



**CUSTOMIZABLE KIT**



- Order all common oxygen accessories in one easy-to-handle kit
- Contents can be customized to fit your needs

## Oxygen Tubing Y-Connector

Item# **RES017**  
Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)  
**10/pack**



- Connects 2 oxygen sources for higher flows

## Oxygen Tubing Connector — Non-Swivel

Item# **RES012**  
Fits 3/16in and 1/4in ID Tubing  
**10/pack**



- Barbed ends ensure secure, leak-free connections

## Tubing Reducer Connector

Item# **RES011**  
Fits 3/16in, 1/4in, 5/16in, 3/8in, and 7/16in ID Tubing  
**10/pack**



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Oxygen Tubing Connector — Swivel

Item# **RES018**  
Male/Male **10/pack**

**RES018-25** 25/pack

**RES018KIT-1**  
Female/Male **Each**



- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

## Oxygen In Use Signs

Item# **RES003**  
6in x 9in With String  
**50/pack**



- Double-sided (English and Spanish)

## Washers with Rubber Ring

Item# **RES035**  
Brass  
**25/pack**

Item# **RES036**  
Aluminum  
**25/pack**



# OXYGEN ACCESSORIES

## O<sub>2</sub> Starter Kit

### Item# **RES02KIT**

RES1107 Cannula with 7ft Tube x2  
RES3015 15ft Supply Tube x2  
RES002 Christmas Tree Connector — Swivel  
RES003 Oxygen In Use Sign  
RES012 Tubing Connector — Non-Swivel  
RES017 Tubing Y-Connector  
RES018 Tubing Connector — Swivel  
RES032 Wrench  
RES035 Brass Washer



**CUSTOMIZABLE KIT**

## Disposable Bubble Humidifier 500 ml

### Item#

### **RES050-50**

Up to 15 LPM • 6 psi • 500 ml

**50/case**

HCP# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

## Wrenches

### Item# **RES031**

With Chain **10/pack**

### Item# **RES032**

Without Chain **10/pack**

### Item# **RES031M**

Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators **5/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

## Humidifier Connectors

### Item# **HUM008**


L-Connector 10in • **10/pk**

### Item# **HUM008L**

L-Connector 21in • **10/pk**

### Item# **HUM008S**

Straight Connector 9in • **10/pk**



## Water Traps

### Item# **RES037**

Push On End Caps • **Each**

### Item# **RES037T**

Threaded End Caps • **Each**



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)



# ACCESSORIES AND CLEANERS

## Homefill Tubing Connector


Item# **RES039**  
10/pack

- Compatible with Invacare, Drive, and Respiration homefill machines
- Durable, coiled tubing




## Homefill Racks

Item# **RES039-RACKB**  
For Invacare Perfecto2  
4/case



**RES039-RACKW**  
For Invacare Platinum  
4/case

**RES039-RACKD**  
For DeVilbiss 525/1025 Series  
4/case



Easily attaches onto concentrator and is designed to conveniently hold Homefill and transfill oxygen systems.

## Oxygen Liter Meter

Item# **RES034**  
0-10 Liter/Minute  
Each



- Oxygen pen
- Works with tanks and concentrators

## RoEzIt Dermal Care

Item# **RES6104L**  
1oz Tube • 18/cas



**CASH ITEM**

- Apply before beginning oxygen therapy and as needed during treatment
- Moisturizes the nose, lip, and over the ear where friction may cause discomfort from cannula
- Petroleum-free with aloe vera and vitamins A and E

## Germicidal Deodorizing Cleaner

Item# **RES6100**  
12/case  
22oz spray bottles



- Heavy duty cleaner contains no bleach, alcohol, or phenols
- An effective bactericide, virucide, and fungicide
- Formulated to disinfect hard, non-porous surfaces

## Germicidal Deodorizing Cleaner — Gallon

Item# **RES6102**  
1 Gallon



- No-scrub, hospital strength germicidal cleaner
- Cleans, deodorizes, and disinfects surfaces, equipment and non-critical instruments



# OXYGEN BACTERIA AND HEPA FILTERS

## 6in HEPA Capsule with Removable Sticker

Item# **BFH6S**  
Threaded  
Each



## Intake Bacteria Filter

Item# **BF2222**  
22mm fittings  
5/pack




## EverFlo Intake Filter

Item# **BF1001**  
Each



## Extended Life Bacteria Filter with Muffler

Item# **OF6011**  
Each



## Extended Life Bacteria Filter

Item# **BF5000**  
Nonthreaded  
Each



## Extended Life Intake Bacteria Filter

Item# **OF1025D**  
Each



## Extended Life Intake Bacteria Filter

Item# **BF1103**  
Each



## Intake Bacteria Filter

Item# **BF7000**  
Each  
For Jumao, Attendant,  
and Lifestyle 5L




## Final Bacteria Filter




Item# **BF516**  
10/pack  
1/4in Barbed  
on Both Sides

## Final Bacteria Filter — 1/4in to 3/8in



Item# **BF1438**  
10/pack  
1/4in to 3/8in  
Stepped Barb  
on Both Sides

## Final Bacteria Filter — 1/4in to 1/2in



Item# **BF1412**  
10/pack  
1/4in to 1/2in  
Stepped  
Barb on Both Sides

## Final Bacteria Filter — 3/16in



Item# **BF316**  
10/pack  
3/16 Stepped Barb and  
1/4in Straight Barb  
on Both Sides

# CAIRE/AIRSEP/SEQUAL


## CAIRE Companion 5



Item# **BF1103**  
Extended Life Intake  
Bacteria Filter  
**Each**  
OEM# 2185

**NEW**

## AirSep Newlife, QuietLife, Intensity



Item# **OF8002**  
Foam Cabinet Filter  
**10/pack**  
5 1/4" x 4 7/8" x 7/8"  
OEM# F1-50-08-05



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# F1-30-03



Item# **OF8001**  
Felt Filter  
**10/pack**  
1 7/8" OD x 3/4" ID x 1"  
OEM# M1 161-1



**KIT**

Item# **OF8002SP**  
Newlife Kit  
2 - OF8002  
1 - BF1438  
1 - OF8001  
Each

## SeQual Integra



Item# **OF11005**  
Foam Cabinet Filter  
**10/pack**  
7 5/8" x 3 5/8" x 3/4"  
OEM# 2185



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**  
OEM# 2033

Item# **OF11005SP**  
SeQual Kit  
2 - OF11005  
1 - BF5000  
1 - BF516  
Each




**KIT**



Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
**10/pack**

## SeQual Eclipse



Item# **OF11006**  
Foam Cabinet Filter  
**10/pack**  
2 1/2" x 2 1/2" x 3/4"

# CHAD

## Total O<sub>2</sub>



Item# **OF11007**  
Foam Cabinet Filter  
**10/pack**  
43/8" x 4" x 3/4"



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**



**KIT**

Item#  
**OF11007SP**  
CHAD Kit  
2 - OF11007  
1 - BF2222  
1 - OF3003  
Each



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
**10/pack**  
25/8" OD x 3/4" ID x 23/8"  
OEM# H613

# COVIDIEN NELLCOR PURITAN BENNETT

## 590 and 592A



Item# **OF6001**  
Foam Cabinet Filter  
**10/pack**  
63/4" x 41/4" x 3/8"  
OEM# 492672



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**



**KIT**

Item#  
**OF6001SP**  
590 and 592 Kit  
2 - OF6001  
1 - BF5000  
1 - BF316  
Each



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**  
OEM# 492141

## SolAiris III, 303 Series



Item#  
**OF5004SP**  
SolAiris III Kit  
2 - OF5004  
1 - BF1438  
1 - BF5000  
Each



KIT



Item# **OF5004**  
Foam Cabinet Filter  
**10/pack**  
3 3/4" x 3 3/4" x 1/2"  
OEM# 30302605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**


## 1025 Series



Item# **OF1025D**  
Extended Life Bacteria  
Filter  
**Each**



Item# **OF5004**  
Foam Cabinet Filter  
**10/pack**  
3 3/4" x 3 3/4" x 1/2"  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651

## SolAiris V, 505, 515 Series



Item#  
**OF5005SP**  
SolAiris V Kit  
2 - OF5005  
1 - BF6011  
1 - BF1438  
Each



KIT




Item# **OF5005**  
Foam Cabinet Filter  
**10/pack**  
5 1/2" x 3 5/8" x 5/8"  
OEM# 50502694



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# MC29D-651



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**



Item# **OF6011**  
Extended Life  
Bacteria Filter  
22mm fittings  
**Each**  
Use if SolAiris was  
manufactured  
before Nov '03


## 525 DS Series



Item# **OF5004**  
Foam Cabinet Filter  
**10/pack**  
3 3/4" x 3 3/4" x 1/2"  
OEM# 303DZ-605



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



Item# **OF6011**  
Extended Life Bacteria  
Filter — 22mm fittings  
**Each**


## All MC Units



Item# **OF5001**  
Foam Cabinet Filter  
**10/pack**  
6 1/8" x 3 1/2" x 1/4"  
OEM# 444-514



Item# **BF5000**  
Extended Life  
Bacteria Filter • **Each**



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# MC29D-651



KIT

Item#  
**OF5001SP**  
MC Kit  
2 - OF5001  
1 - BF5000  
1 - BF1438  
Each



# INOGEN

## One



Item# **OF11008**  
Foam Cabinet Filter  
**10/pack**  
2 1/4" x 3/4" x 1/2"  
OEM# RP-100

# INVACARE

## Mobilaire 3 and 5



Item# **OF9001**  
Foam Cabinet Filter  
**10/pack**  
6 3/4" x 3 3/8" x 1/2"  
OEM# 2000489

Item# **OF9001SP**  
Mobilaire Kit  
2 - OF9001  
1 - BFH6S  
1 - BF516  
Each



KIT



Item# **BFH6S**  
HEPA Filter  
6in Capsule  
**Each**  
OEM# 2000806



Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
**10/pack**  
OEM# P0750X030

## Platinum 5 and 10



Item# **OF9002**  
Foam Cabinet Filter  
**10/pack**  
9" x 27/8" x 1/2"  
OEM# 1107412

Item# **OF9002SP**  
Platinum Kit  
2 - OF9002  
1 - BFH6S  
1 - BF516  
Each



KIT



Item# **BFH6S**  
HEPA Filter 6in  
Capsule with Sticker  
**Each**  
OEM# 2000806



Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
**10/pack**  
OEM# P0750X030



Item# **HUM008S**  
Straight Humidifier  
Connector with  
9in Tube  
**10/pk**

## Perfecto2



Item# **OF9004**  
Foam Cabinet Filter  
**10/pack**  
6 7/8" x 3 3/4" x 1/2"  
OEM# 1143492

Item# **OF9004SP**  
Perfecto2 Kit  
2 - OF9004  
1 - BFH6S  
1 - BF516  
Each



KIT




Item# **BFH6S**  
HEPA Filter 6in Capsule  
**Each**  
OEM# 2000806



Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
**10/pack**  
OEM# P0750X030

## HomeFill



Item# **OF9001**  
Foam Cabinet Filter  
**10/pack**  
6 3/4" x 3 3/8" x 1/2"

# JUMAO, LIFESTYLE, ATTENDANT

## Jumao, Lifestyle, Attendant



Item# **OF7000**  
Foam Cabinet Filter  
**10/pack**

**NEW**



Item# **BF7000**  
Bacteria Filter  
**Each**

## NIDEK

### Nuvo



Item# **OF11003SP**  
Nuvo Kit  
2 - OF11003  
1 - BF1438  
1 - BF5000  
**Each**



**KIT**



Item# **OF11003**  
Foam Cabinet Filter  
**10/pack**  
9 1/2" x 4" x 1/4"  
OEM# 9250-1025



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in Barbed  
**10/pack**  
OEM# 7631-1053



Item# **BF5000**  
Extended Life  
Bacteria Filter  
**Each**  
OEM# 9250-1028

### Mark IV, V, V+




Item# **OF11004SP**  
Mark Kit  
2 - OF11004  
1 - BF1438



**KIT**



Item# **OF11004**  
Foam Cabinet Filter  
**10/pack**  
10" x 4" x 1/2"  
OEM# 9007-6G37



Item# **BF1438**  
Final Bacteria Filter  
1/4in to 3/8in  
Barbed  
**10/pack**  
OEM# 7631-1053

## 3B

### Stratus 5



Item# **OF6030**  
Foam Cabinet Filter  
**10/pack**  
OEM# STR1030



Item# **BF6020**  
HEPA Filter  
**10/pack**  
OEM# STR1020



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and 1/4in Straight Barb  
**10/pack**  
OEM# H621

# RESPIRONICS/HEALTHDYNE

## EverGo



Item#  
**OF1005SP**  
EverGo Kit  
2 - OF1005  
1 - BF316  
Each

KIT



Item# **OF1005**  
Foam Cabinet Filter  
**10/pack**

2<sup>3</sup>/<sub>16</sub>" x 11<sup>1</sup>/<sub>16</sub>" x 3<sup>3</sup>/<sub>8</sub>"  
OEM# 900-150



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**

OEM# H621

## Millennium



Item#  
**OF1001SP**  
Millennium Kit  
2 - OF1001  
1 - BF2222  
1 - OF3003  
1 - BF316  
Each



KIT



Item# **OF1001**  
Foam Cabinet Filter  
**10/pack**

8<sup>1</sup>/<sub>4</sub>" x 4" x 3<sup>3</sup>/<sub>8</sub>"  
(3/16" notches)  
OEM# H619



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**

OEM# H620



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
**10/pack**

2<sup>5</sup>/<sub>8</sub>" OD x 3<sup>3</sup>/<sub>4</sub>" ID x 2<sup>3</sup>/<sub>8</sub>"  
OEM# H618



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb  
and 1/4in Straight  
Barb  
**10/pack**

OEM# H621



Item# **OF1002**  
Felt Filter  
**6/pack**

NEW



Item# **OF2000**  
Sieve Bed Filter  
4.125in  
**5/pack**

NEW



Item# **OF2001**  
Sieve Bed Filter  
2.875in  
**5/pack**

## Alliance



Item#  
**OF3002SP**  
Alliance Kit

2 - OF3002  
1 - BF2222  
1 - OF3003  
1 - BF316  
Each



KIT



Item# **OF3002**  
Foam Cabinet Filter  
**10/pack**

8<sup>1</sup>/<sub>4</sub>" x 3<sup>3</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>2</sub>"  
OEM# 100-0078-10



Item# **BF2222**  
Intake Bacteria Filter  
22mm fittings  
**5/pack**

OEM# 100-6631-00



Item# **OF3003**  
Foam Intake Filter  
for use with BF2222  
**10/pack**

2<sup>5</sup>/<sub>8</sub>" OD x 3<sup>3</sup>/<sub>4</sub>" ID x 2<sup>3</sup>/<sub>8</sub>"  
OEM# H618



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**

OEM# H621

## EverFlo



Item#  
**OF1003SP**  
EverFlo Kit

1 - BF1001  
1 - BF316  
Each



KIT



Item# **OF1006SP**  
Cabinet Filter Kit  
(No Cover)  
**2 Filters, 1 Casing**  
OEM# 1055625



Item# **OF1006**  
Cabinet Filter  
**10/pk**  
5<sup>11</sup>/<sub>16</sub>" x 4<sup>1</sup>/<sub>2</sub>" x 7<sup>1</sup>/<sub>16</sub>"



Item# **OF1006-CASE**  
Foam Cabinet Filter  
Casing  
**10/pack**



Item# **BF1001**  
EverFlo Intake Filter  
**Each**  
OEM# 1038831



Item# **BF316**  
Final Bacteria Filter  
3/16 Stepped Barb and  
1/4in Straight Barb  
**10/pack**  
OEM# H621



Item# **OF2000**  
Sieve Bed Filter 4.125in  
**5/pack**

**Other****Page**

Aerosol	61
Tabletop Nebulizers	62-63
Handheld Nebulizers	63
Nebulizer Kits	64
Nebulizer and Anesthesia Filters	65
Suction	66-69
Tracheostomy	70
Ventilator Filters, Masks, and Connectors	71
Ventilator Circuits	72-73
StrongArm Comfort Cane	74
Bath Safety	75
PPE	76-78
Maternity	79

**OTHER****OTHER HCPCS**

A6402	Gauze - sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A4483	Heat and Moisture Exchanger — Disposable
A4626	Tracheostomy Cleaning Brush
A4628	Oropharyngeal Suction Catheter (Yankauer)
A4629	Tracheostomy Care Kit - for established tracheostomy
A7000	Suction Canister — Disposable
A7002	Tubing for Suction Pump
A7003	Nebulizer Set - Small Volume—Disposable
A7005	Nebulizer Set - Small Volume—Nondisposable
A7013	Nebulizer Compressor Filter— Disposable
A7015	Aerosol Mask used with Nebulizer
A7525	Tracheostomy Mask
E0570	Nebulizer with Compressor
E0600	Respiratory Suction Pump - home model, portable or stationary, electronic
E2601	General Use Wheelchair Seat Cushion - width less than 22 inches, any depth



# AEROSOL

## Aerosol Mask

Item# **RES2110**  
Adult Size  
**50/case**

Item# **RES2210**  
Pediatric Size  
**50/case**

HCPCS# A7015



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



## Peak Flow Meter

Item# **RES052**  
**10/case**

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



**NEW**

Item# **RES052MPD**  
Pediatric Mouthpiece  
for Peak Flow Meter  
**10/case**

## Drainage Bag with T-Adapter

Item# **RES070**  
750cc • **50/case**



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

## Drainage Bag with Y-Adapter

Item# **RES071**  
750cc • **50/case**



## SperaChamber

Anti-Static Valved Inhaler Chamber  
for Metered Dose Inhaler

Chamber Only

Item#  
**RES195-10 10/pack**  
**RES195-50 50/pack**

Chamber with Mask  
**10/pack**

Item#  
**RES195MS Small**  
**RES195MM Medium**

**RES195ML Large (Adult)**

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free
- Dog mask appeals to kids



# TABLETOP COMPRESSOR NEBULIZERS

## Compressor Nebulizers

Item#

### NEB100

- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

### NEB100-R

- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



**NEB100-BAG**  
Carry Bag

## Sunny the Bear

### NEB300BEAR

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

### NEB300BEAR-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

## Gordon the Fish

### NEB300FISH

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

### NEB300FISH-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)

**NEW**



## Basketball Pediatric Compressor Nebulizer

Item#

### NEB300BASKETBALL-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)

**NEW**



Sunset's Basketball Nebulizer brings some play to treatment time! The Basketball Nebulizer is packaged in a vivid and fun retail box that appeals to kids.



# HANDHELD NEBULIZERS

## Handheld Mesh Nebulizer

Item#  
**NEB402**

**With internal rechargeable battery**

**Included**  
Medication cup  
USB cord  
Mask  
Mouthpiece  
Portable bag  
Manual

**NEW**



**Particle Size** MMAD < 5µm  
**Nebulization Rate** 0.25ml/min-0.9ml/min  
**Medication Cup Capacity** 10ml (Max)  
**Device Weight** Approx. 1/4 lb  
**Device Size** Approx. 1 1/2 x 2 x 3 7/8 inches

## Mini Mesh Nebulizer

Item#  
**NEB400**

**Included**  
USB cord  
Medicine cup  
Mouthpiece  
Adult mask  
Pediatric mask  
Carrying bag

- Runs on USB or 2 AA batteries
- Quiet
- Fast treatment times
- Small particle size
- Attainable for more patients

### Replacement Parts

<b>NEB400CUP</b>	Medication Cup
<b>NEB400MP</b>	Mouthpiece
<b>NEB400MASKP</b>	Pediatric Mask
<b>NEB400MASKA</b>	Adult Mask



## Handheld Compressor Nebulizer

Item#  
**NEB200**

**Included**  
Compressor  
Medication cup  
AC power adapter  
Air hose  
Elbow  
Mouthpiece  
Mask  
4 AA batteries



- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required


Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

# NEBULIZER KITS

## Disposable Nebulizer Kit with T-Piece

Item# **RES091**  
Disposable  
Single Patient  
**50/case**  
HCPCS# A7003

- Includes:
- 6cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

## Reusable Nebulizer Kit with T-Piece

Item# **RES093**  
Reusable — Single Patient  
**10/case**  
HCPCS# A7005

- Includes:
- 10cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

## Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**  
Disposable  
Single Patient  
**50/case**  
HCPCS# A7003



- Includes:
- 6cc jet nebulizer bottle
  - Angled mouthpiece
  - 7ft supply tube



## Disposable Nebulizer Kit with Mask

Item# **RES090**  
Adult Size  
**50/case**

Item# **RES092**  
Pediatric Size  
**50/case**

HCPCS# A7003



### Includes

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

Disposable single patient use for long term delivery procedures

6cc bottle is easy to assemble and disassemble



# NEBULIZER FILTERS

## Plastic Intake Filter

Item# **BF015**  
**10/pack**  
 HCPCS# A7013



## Plastic Intake Filter with Head

Item# **BF016**  
**10/pack**  
 HCPCS# A7013



## Plastic Intake Filter with Head

Item# **BF017**  
**5/pack**  
 HCPCS# A7013



## Foam Filter

Item# **BF018**  
**10/pack**  
 HCPCS# A7013



## Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**  
**10/pack**  
 HCPCS# A7013

- Mini
- MicroElite
- MiniElite



## Nebulizer Filter Compatibility

### Caire

Medi-Mist	BF018
-----------	-------

### DeVilbiss

PulmoNeb	BF017
PulmoNeb LT	BF017
Sunrise Compact	BF018
Pulmo-Aide	BF017
Pulmo-Aide Compact	BF017
PulmoMate	BF017

### Evo

AeroMax	BF018
SportsNeb	BF018
Super SportsNeb	BF018

### Invacare

Passport	BF018
Envoy	BF018
Envoy Jr.	BF018
Select	BF018

### John Bunn

Nebulite EV	BF018
LT2	BF018
LX	BF018

### MedLine

SportsMist	BF018
------------	-------

### Respironics

Inspiration 323	BF015
Inspiration 626 & 929	BF016
MisterNeb	BF016
Inspiration Elite	BF016
MiniElite	BF018
InnoSpire	BF026S

### Salter

AIRE Plus	BF018
-----------	-------

### Schuco

2000	BF018
3000	BF018

### Sunset

NEB100 Series	BF018
NEB300 Series	BF018

# SUCTION

## Suction Machine

Item#

**SU100DC**

With Internal  
Rechargeable  
Battery

**SU100AC**

Without Battery



### Comes with

Machine  
Carry Bag  
800cc Suction Canister  
AC Adapter & Power Cord  
6ft Suction Tubing  
3.5in Suction Tube Connector  
4.5in Suction Tube Connector  
BF1438S Bacteria Filter

- Effective removal of respiratory fluids
- Durable yet lightweight
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty

### Accessories and replacement parts

**SU100DC-PA** DC Mobile Power Adapter (for car power outlet)

**SU100DC-Adapt** Replacement Power Adapter (without cord)

**SU100DC-Cord** Replacement Power Cord (without adapter)

**SU100DC-BAG** Carrying/Storage Bag

**SU100-PA** Car Charger (DC Mobile Power Adapter)

**SU100DC-BATTERY** Internal rechargeable battery replacement

## Suction Kit for SU100DC

Item# **RES026S-SS**

RES023AW Suction Canister  
RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438S Bacteria Filter  
RES025 Suction Tubing 1/4in x 6ft

HCPCS#  
CANISTER - A7000  
TUBING - A7002



## Suction Filter Component Kit for SU100DC

Item# **RES024-MKIT**


RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438S Bacteria Filter



## Suction Filter Component Kit

Item# **RES024XL-3KIT**

RES024XL Suction Tubing Connector 16in  
BF3812 Filter  
RES028 Elbow Connector



## Suction Canister

Item#

**RES023AW**

800cc • Single Use

**48/case**

HCPCS# A7000

Item#

**RES025AW1200**

1200cc • Single Use

**48/case**



- Float valve shutoff

## Suction Kit

Item# **RES026S-E**

RES023AW Suction Canister  
RES024M Suction Tubing Connector 4.5in  
RES025 Suction Tubing 6ft  
RES028 Elbow Connector  
BF3812 Bacteria Filter

HCPCS#  
CANISTER - A7000  
TUBING - A7002



# SUCTION

## Suction Machine Inline Filters

Item# **BF1438S**  
1/4in to 3/8in  
10/pack



Item# **BF3812**  
3/8in to 1/2in  
10/pack



Item# **BF1312**  
1/3in to 1/2in  
1/pack



Item# **BF516TS**  
5/16in Barbed with  
Threaded NPT  
10/pack



## Suction Tubing

### 1/4in

Length	Item#	Qty
18in	<b>RES025S</b>	<b>50/case</b>
6ft	<b>RES025</b>	<b>10/case</b>

### 3/16in

Length	Item#	Qty
18in	<b>RES025-31618</b>	<b>100/case</b>
20in	<b>RES025-31620</b>	<b>100/case</b>
6ft	<b>RES025-31672</b>	<b>50/case</b>

HCPCS# A7002



## Suction Tube Connectors

### 10/pack

Size	Item#
1.5in	<b>RES024S</b>
3.5in	<b>RES024</b>
4.5in	<b>RES024M</b>
6in	<b>RES024L</b>
10in	<b>RES024L10</b>
16in	<b>RES024XL</b>

**NEW**



## Rubber Suction Elbow Connector Drive DeVilbiss Style

Item# **RES028**  
6/pack



## Suction Elbow Connector Precision Medical Style

Item# **RES033**  
6/pack



## Bulb Tip Yankauers

Item# **TRAA001S**  
Vented • 50/case

Item# **TRAA002S**  
Non-Vented • 50/case

HCPCS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber

## CLOSED SUCTION

NEW

## Standard

## 24 Hour Adult

## Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

## Tracheostomy 305mm

VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr
VEN46011400M	14Fr, MDI Port
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

## 72 Hour Adult

## Endotracheal 540mm

VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

## Tracheostomy 305mm

VEN56061000A	10Fr
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port



## Neonatal/Pediatric


## With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000600	305mm, 6Fr
VEN66000700	305mm, 7Fr
VEN66000800	305mm, 8Fr
VEN66001000	405mm, 10Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000M	405mm, 10Fr, MDI Port

## Deluxe

## Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single handed operation design, which allows you to keep your eye on the patient and do no harm. One handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.



## Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.

**MDI Port** provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

**Irrigation Port** offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

**Isolation Valve** seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

**Black Ring Indicator** ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



# CLOSED SUCTION

## Deluxe

### 24 Hour Adult

#### Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN160160	16Fr
VEN160160M	16Fr, MDI Port

#### With Wedge + Adapter

Item#	Description
VEN162100A	10Fr
VEN162100AM	10Fr, MDI Port
VEN162120	12Fr
VEN162120M	12Fr, MDI Port
VEN162140	14Fr
VEN162140M	14Fr, MDI Port
VEN162160	16Fr
VEN162160M	16Fr, MDI Port

#### Tracheostomy 305mm

VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN161120	12Fr
VEN161140	14Fr
VEN161160	16Fr
VEN161160M	16Fr, MDI Port

#### With Wedge + Adapter

VEN163100A	10Fr
VEN163100AM	10Fr, MDI Port
VEN163120	12Fr
VEN163140	14Fr
VEN163160	16Fr

### 72 Hour Adult

#### Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN365160	16Fr
VEN365160M	16Fr, MDI Port

#### With Wedge + Adapter

Item#	Description
VEN367100A	10Fr
VEN367100AM	10Fr, MDI Port
VEN367120	12Fr
VEN367120M	12Fr, MDI Port
VEN367140	14Fr
VEN367140M	14Fr, MDI Port
VEN367160	16Fr
VEN367160M	16Fr, MDI Port

#### Tracheostomy 305mm

VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN366160	16Fr
VEN366160M	16Fr, MDI Port

#### With Wedge + Adapter

VEN368100A	10Fr
VEN368100AM	10Fr, MDI Port
VEN368120	12Fr
VEN368120M	12Fr, MDI Port
VEN368140	14Fr
VEN368140M	14Fr, MDI Port
VEN368160	16Fr
VEN368160M	16Fr, MDI Port

### Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port

# TRACHEOSTOMY

## Tracheostomy Mask

Item#

**RES2130**

Adult Size

**50/case**

**RES2230**

Pediatric Size

**50/case**

HCPCS# A7525



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## HMEs

### Heat Moisture Exchangers



**RES027FLEX**

Flex HME —  
Foam



**RES027**

T-HME



**RES027OP**

T-HME with  
Oxygen Port



**RES027IL**

Inline HME —  
Large Volume



**RES027NE**

Inline HME —  
Neonatal

	<b>RES027FLEX</b> Flex HME — Foam	<b>RES027</b> T-HME	<b>RES027OP</b> T-HME with Oxygen Port	<b>RES027IL</b> Inline HME — Large Volume	<b>RES027NE</b> Inline HME — Neonatal
TV Range (mL)	300-1500	250-1000	250-1000	50-250	15-50
Moisture Output (mg H <sub>2</sub> O/L)	33 mg/L	30, Vt = 0.1L		31, Vt = 0.1L	30 mg/L
Weight (g)	17	4.5	7	10	4.5
Resistance (cm H <sub>2</sub> O)	0.5 @ 30 LPM 1.34 @ 60 LPM	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM
Dead Space (mL)	75	12.5	12	13	2.5
Connection (mm)	22M,15F	15F	15F, 5.5mm Oxygen port	15M,22M/15F	15M,15F

### Tracheostomy Ties with Fastener Tabs

**20/case**


One-piece

- Easy on

Item#

**RES029** Adult

**RES029P** Pediatric




Two-piece

- Easy off

Item#

**RES030** Adult

**RES030P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton



## Tracheostomy Swivel Elbow Connectors

15mm ID x 15mm OD

**10/pack**



Item# **RES021**

With suction port

Item# **RES022**



## Adapter for T-HME

Item# **RES027A**

**Each**



# VENTILATOR FILTERS, MASKS, AND CONNECTORS

## NIV Mask with Headgear — Non-vented

Each

Item#

**CM038M** Medium

**CM038L** Large

- Disposable, single patient use
- PC frame with PVC air cushion seal
- Low dead space for low CO2 buildup
- Clear mask allows a good view of mouth and nasal area
- Headgear with Velcro fasteners provide maximum versatility and movement
- Dual 360 degree swivels at the mask attachment and at the standard 22mm tubing attachment point



## Bacteria Filter

Item#

**BF2200**

22mm Bacteria

Disk Filter

**5/pack**




## Respironics Trilogy Ventilator Filter

Item#

**CF1009F**

**2/pack**

OEM# 1029330



## Air Inlet Vent Filter for Philips Respironics V60

Item#

**VF002**

For use with Philips Respironics V60

**5/pack**



## Inlet Vent Filter for Pulmonetic LTV Series

Item#

**VF001**

For use with Pulmonetic LTV Series

**5/pack**



## Valved T-Adapter — Adult

Item#

**VEN001**

22F-18F-22M/15F

**Each**



## Ventilator Felt Filter

Item#

**VF006**

For use with NPB LP6 & LP10

**10/pack**



## Swivel Connector Silencer

Item#

**VEN002**

For Ventilator Circuits and CPAP

**Each**



## Circuit Connectors

Item#

**RES013**

22mm OD x

15mm ID/22mm OD

**10/pack**

**RES014**

22mm OD x 22mm OD

**10/pack**



## Fixed Exhalation Port

Item#

**VEN003**

For Ventilator Circuits and CPAP

22 mm OD x 15 mm ID

**Each**

**NEW**



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS



# VENTILATOR CIRCUITS WITHOUT HUMIDIFIER CHAMBER

## Passive Heated Breathing Circuit

Item#  
**VC1114**

For use with MR850 Humidifier  
**10/case**

OEM# RT1114



## Passive Non-Heated Breathing Circuits

For use with Respiration Trilogy

Item#  
**VC1000**

With swivel exhalation port  
**Each**

OEM# 1073228 & 1073236



**VC1100**

With fixed exhalation port  
**Each**

OEM# 1069467 & 1069468

## Cough Assist Kits

Item#  
**CAKIT001**  
**Each**

OEM# 1090835

- 1 RES013 Circuit Connector  
22mm OD, 15mm ID x  
22mm OD
- 1 BF2200 Bacteria/Viral  
Filter
- 1 TUB06 Circuit Tube
- 
- 

**CAKIT002**  
**Each**

- 1 RES013 Circuit Connector  
22mm OD, 15mm ID x  
22mm OD
- 1 BF2200 Bacteria/Viral  
Filter
- 1 TUB06 Circuit Tube
- 1 RES019 Oxygen  
Enrichment Adapter





# VENTILATOR CIRCUITS WITHOUT HUMIDIFIER CHAMBER


## Non-Heated Breathing Circuits

For use with Respironics Trilogy  
22mm cuffs  
**30/case**


**NEW**

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
<b>VC1001</b>	1073223 & 1073232	Active	No	Adult	22mm
<b>VC1001P</b>	1073220 & 1073229	Active	Yes	Pediatric	15mm
<b>VC1004</b>	1069443 & 1069444	Passive	Yes	Adult	22mm
<b>VC1004P</b>	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#  
**VC1001**  
Active Adult  
OEM#  
1073223  
1073232




Item#  
**VC1001P**  
Active Pediatric  
OEM#  
1073220  
1073229




To use active circuits with Respironics Trilogy, simply cut ends off pressure and exhalation lines.



Item#  
**VC1004**  
Passive Adult  
OEM#  
1069443  
1069444



Item#  
**VC1004P**  
Passive Pediatric  
OEM#  
1073224  
1073233



# STRONGARM

## StrongArm® Comfort Cane

NEW

- Ergonomic handgrip supports your forearm and relieves pressure from your wrist
- Offset bend centers grip over tip for maximum stability when cane is at the user's side
- Lightweight & strong
- 12 height adjustments
- Comfort Fit Handle
- Supports up to 500 lbs
- For left & right arm
- Self-standing base: Larger 'footprint' greatly improves traction and overall stability
- Safely walk on grass, gravel, sand, carpet, hard surfaces, and stairs



10/case

**BM001** Black

**BM001-BL** Blue

**BM001-TI** Titanium

One of each color pack

**BM001-3MC**

STRONGARM COMFORT CANE



# BATH SAFETY

## PreserveTech™ 360° Swivel Bath Chair

Item#  
**BSE001**

Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"-22", and the padded armrests can be raised out of the way.



## Splash Defense™ Transfer Bench

Item#  
**BSE002**

Splash Defense™ Transfer Bench is the first of its kind to allow the shower curtain liner to remain inside the bath tub while a user is bathing. This unique product enhancement is essential to help prevent water from spilling out onto the bathroom floor, creating a danger zone for users and caretakers. The Splash Defense™ Transfer Bench provides security for the caretaker, while giving users back their independence and allowing them to bathe safely and privately again.



## PreserveTech™ 360° Swivel Bath Chair

Item#  
**BSE003**

Introducing the PreserveTech™ Universal Raised Toilet Seat. With our first ever, truly universal raised toilet seat users can set three different seat heights (2, 4 or 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space. This item also includes a lid with a recessed grip for effortless opening.





# PPE

## Mini Pistol Grip Non-Contact Thermometer

Item#  
**PPE4001-J**  
Infrared  
Each



- Forehead mode for all ages
- Object measurement mode
- HD LCD display
- Battery: AAAx2
- Fever warning: 6 short double beeps
- Memory storage: 20 sets of measurable data
- Fahrenheit and Celsius
- Automatic shutdown: 10 seconds
- Professional accuracy (plus or minus 0.6 degrees F when measuring forehead temperatures from 89.6 through 107.9 degrees F)

## Simpli Germicidal Disposable Wipes

Item#  
**RES6108W**  
Canister of 160 wipes

**NEW**

- Bactericidal, Tuberculocidal, Fungicidal, Virucidal\*
- Premoistened wipes containing a quaternary/alcohol-based solution.
- Made in the USA.
- Kills \*SARS-CoV-2, Covid 19, Mycobacterium bovis – BCG (Tuberculosis—TB) in 1 Minute
- Kills Candida Auris in 2 Minutes



## Mini Wand Style Non-Contact Thermometer

Item#  
**PPE4002-A**  
Infrared  
Each



- One second reading
- One click operation
- Fahrenheit and Celsius
- ¾ to 2-inch measuring distance

## CleanSmart Disinfectant

Item#  
**CAP1019-C**  
23oz Spray Bottle  
6/case

- CleanSmart is EPA listed as a broad-spectrum and hospital disinfectant
- EPA Registration Number 89896-2
- CleanSmart is formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment
- CleanSmart is included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2, the virus that causes COVID-19
- Active ingredient Hypochlorous Acid does not leave a chemical residue



## Nitrile Gloves

Ask us about our gloves and item numbers, which may vary based on available supply.



## CleanSmart Alcohol-Free Hand Cleanser

Item#  
**CAP1021-HC8**  
8oz Spray Bottle  
6/pack

**CAP1021-HC2**  
2oz Spray Bottle  
9/pack

- Same active ingredient as CleanSmart cleaner, in a hand cleanser.





# PPE

## Disposable 3-Ply Face Mask

Item#

**PPE1002PS**

Adult • 6 7/8 x 3 3/4 in

**Box of 50 Masks**

- Hypoallergenic
- Comfortable
- Breathable
- Fiberglass-free
- 3 layers of polypropylene (a meltblown layer between two soft non-woven layers)



## Disposable Face Shields with Elastic Strap

Item#

**PPE3001FS-F**

Adult

**Each**

- High-density foam forehead rest
- Clear, anti-fog vision



**PPE3002PFS-M**

Pediatric

Monster Design

**Each**

**PPE3002PFS-A**

Pediatric

Animal Design

**Each**



## Disposable Level 2 Surgical Face Mask

Item#

**PPE1002SM**

**Box of 50 Masks**

- Medical classification and 510(k)
- 98% bacterial filtration efficiency
- Rated for moderate risk levels of aerosols, spray, and/or fluids



## Disposable Tearaway Gown

Item#

**PPE5002**

**25/pack**


- Level 3 ASTM F3352
- Individually sealed
- Made in the USA



## BOLT Clean Air System

**NEW**

Item#  
**BOL1001**  
Each



BOLT's technology is a recommended method to reduce exposures to SARS-CoV-2, the virus that causes COVID-19.

Included:

- 1 BOLT Louver System
- 1 UV-C Germicidal Light Bulb
- All necessary mounting hardware
- Operating instructions

BOLT is the 360° commercial unit that can be used safely while people are in the room, allowing for a continuous and comprehensive germicidal UVC disinfection that kills 99.9% of viruses. BOLT uses ultraviolet light to create an eradication zone, destroying viruses, bacteria, and other pathogens at a DNA level, and disabling their ability to multiply and cause disease.

BOLT safely eradicates viruses with shielded UVC light. Its state-of-the-art engineered louver system protects people from direct UVC light rays. This means BOLT provides real time air sanitation and protection for customers, clients, employees, and students while they work, shop, learn, and dine indoors. BOLT's patent pending design comes from the team's decades of niche experience with plastics.

# MATERNITY

## Twist Double Electric Breast Pump

Item#  
**MAT2300T1**

**NEW**

The Motif Twist is ideal for moms looking for portable convenience. The Twist works at home or on-the-go, with an AC adapter or battery power option. Users simply twist the knob to select different suction and massage settings. The blue LED light timer tracks pumping sessions. The Twist allows moms to customize their pumping experience to fit their comfort level.

- Super lightweight design. Weighs less than ½ a pound.
- Single or double pumping session.
- Closed system. Backflow prevention keeps the motor and tubing clean, preventing contamination.
- Dual power options. Use it stationary or on battery power (4 AA batteries).
- Quiet motor for discreet pumping.
- Adjustable suction pattern designed for maximum letdown.



## Motif BiliTouch™ Phototherapy Blanket

Item#  
**MAT2300EAB**

**NEW**

The Motif BiliTouch™ phototherapy blanket is a biliblanket designed by healthcare professionals to treat newborn infants diagnosed with jaundice (hyperbilirubinemia). This portable, convenient phototherapy treatment device allows your patients to continue to hold, feed, and care for their baby while reducing high bilirubin levels at home.

- Convenient, ultra-portable phototherapy treatment with battery-powered option
- Patient-centered, lightweight design facilitates bonding with newborn infant
- LED light source, no fiber optic cables for easy positioning



Visit [www.sunsethcs.com/acuteacare](http://www.sunsethcs.com/acuteacare)

**SUNSET**  
HEALTHCARE SOLUTIONS



Precise, on the wire incisions for guidewire assisted catheter insertion

## GuideBlade—Precision Incisions, Every Time!

### Evaluation of the efficacy and safety of a novel dermatotomy device for central venous cannulation

Worasak Keeyapaj, MD., Albert T. Cheung, MD.

Department of Anesthesiology, Perioperative and Pain Medicine, Stanford University, Stanford, California

Outcome	Number (%)
<b>Primary outcome:</b> Successful lines insertion without additional instrument	188 (100%)
Wire damage	0 (0%)
Operator injury	0 (0%)



### About the GuideBlade™

The GuideBlade™ is the only disposable scalpel on the market engineered specifically to streamline over the wire vascular access placement. The GuideBlade™ allows clinicians to perfectly align the dermatotomy incision made during vascular access placement with the in situ guidewire. This alignment prevents residual skin tags and eliminates the need for dermatotomy revisions.

### GuideBlade™ Clinical Benefits

- Less insertion site bleeding.
- Maintenance of insertion site skin integrity.
- Decreased risk of catheter bacterial colonization.
- Decreased risk of wire related complications (wire fracture and wire kinking).
- Decreased scarring/improved cosmetic result.

### GuideBlade™ In Practice

- Appropriate for use with guidewires up to and including 0.035" diameter.
- Useful for placement in a wide range of venous access devices including central venous catheters, PICC lines, venous introducer sheaths, hemodialysis catheters, and percutaneous ECMO cannulas.
- Particularly useful in procedures placed while wearing protective faceshields where condensation and glare can obscure visualization of the guidewire insertion site.
- Particularly useful in obese or bleeding patients where wire entry site can be obscured by excess soft tissue or entry site bleeding.
- Increased benefit when used by trainees and less experienced practitioners.

### The GuideBlade™ Value Proposition

- Decreases the need to replace kit components (like guidewires, dilators or catheters) that are easily damaged by residual skin tags often seen with standard disposable scalpel use.
- Decreases insertion site bleeding, thereby decreases the need for frequent dressing changes.
- Maintains skin integrity and thereby decreases risk of catheter colonization and associated infection.

Misaligned skin incisions make passage of dilators and catheters difficult and can lead to guidewire kinking and fracture. Revision of incisions may increase the risk of insertion site bleeding, scarring and catheter associated bloodstream infection.





## Stop “brutane” mask induction and turn panic into play

For children, the process of undergoing surgery is already unsettling, and the process of mask-based induction is often the most traumatizing part. Often kids will kick and scream rather, believing they’re being trapped by the gas or suffocated by the mask. The experiences create negative emotional scarring that can last weeks, months or even years.

With the introduction of the revolutionary Pediatric Device for Induction of Anesthesia device, or “PeDIA” for short, that stress and anxiety can be replaced with something colorful and fun.

As the world’s first mask-free system and method for inhalation induction, the PeDIA provides a solution to this problem that anesthesia providers have been researching for decades. It replaces the mask with a colorful balloon, which actually allows the child to put themselves to sleep rather than being forced to submit to an unfamiliar mask.

PeDIA resolves the problems of traditional mask induction while fostering a sense of trust and patient satisfaction.

### Item#

<b>ANEP001B</b>	Blue
<b>ANEP001Y</b>	Yellow
<b>ANEP001P</b>	Purple

Instructional videos for children available in English and Spanish





## Customized anesthesia workstation for every anesthetizing location

Item#

**ANE77006**

Anestand Anesthesia Stand, Complete Tray and Clamp

**ANE77001** Replacement Tray

**ANE77005** Replacement Clamp

### Benefits to Providers



- Keep medications near IV port
- Ergonomical work space
- Makes wiping IV port easy
- Extra set of hands when help not around
- Helps make RSI safe
- Creates valuable space for non-OR anesthesia

### Value to Hospital

- Helps with regulatory compliance
- Sharps management
- Could help reduce HAIs
- Demonstrates support for anesthesia team
- Reusable/single fixed cost

1. The recessed tray holds and organizes supplies for induction of anesthesia and most anesthetic procedures
2. Clips keep IV tubing securely in place, so it's easy to "scrub the hub"
3. Corner grips for suction tubing and more
4. The adjustable support post lets you position it just where you need it
5. The anesthesia circuit clips keep your circuit fixed and firmly in place
6. Patented clamp attaches to OR bed, IV pole, or inpatient bed

**Anestand can help health systems save lives,  
lower costs and improve compliance.  
All at the same time.**





Scan to learn more  
←

**The first and only anesthesia stand to safely and securely keep all supplies where they're needed can also be one of your best weapons in combatting infections.**

### **Anestand Facts**

- A small investment could help in compliance, save time and quickly pay for itself.
- By having supplies and medications organized and easily available with IV port ready, Anestand aids in patient safety.
- The bottom line: Anestand could help prevent costly infections and save money by helping anesthesia providers save OR time.

**Of all the many safety and security features Anestand offers, perhaps the single greatest is how it helps health systems to lower HAIs.**

- Having IV port ready for site alcohol scrub makes disinfection (and compliance) easier — Anestand keeps your IV port ready for scrubbing.
- By keeping medications immediately available and infusion ports clearly defined Anestand provides additional safety.
- By keeping supplies immediately available, Anestand helps make RSI (rapid sequence inductions) safe.
- During COVID, Anestand assists anesthesia providers as an extra pair of hands, while other healthcare personnel are being asked to keep distance from the airway.

**Anestand can provide value in safety, compliance, and convenience.  
As soon as you use it, you'll wonder how you ever managed without it.**

- CDC, JCAHO, SHEA, and other organizations require cleaning the port with every injection. Anestand helps your facility and providers comply with these regulations.
- Patients with an HAI stay 19 additional days, or \$11.2 million in HAI treatment and/or lost revenue. For a small fixed cost your facility could get a large return on investment.
- Keeping supplies organized on a tray may prevent dropped items and reduce replacement costs. Anestand makes it far easier to keep supplies organized and right at hand.



