

Product Catalog

Contract PP-NS-1571
**Closed Suction
Catheters**



Contract PP-NS-1753
Medication Delivery



Contract PP-NS-1741
Respiratory Disposables



Contract PP-NS-1762
Resuscitation Devices



Contract PP-OR-2217
Laryngeal Masks





Sales Representative

Egle Sventnickaite

O: 872-258-1654

esventnickaite@sunsethcs.com

Sales Representative

Christian Straussner

O: 872-282-0746

C: 847-525-6340

cstraussner@sunsethcs.com

Sales Representative

Owen Ruth

O: 312-533-2451

oruth@sunsethcs.com

Director, Sales Operations

Mark Lipson

O: 312-533-2465

mlipson@sunsethcs.com

Sales and Marketing Manager

Jackson Hanft

C: 971-416-6695

jackson@sunsethcs.com

Director, Strategic Accounts

Lisa Iannazzone RN, BSN

C: 636-222-3134

lisa@sunsethcs.com

RRT Clinical Specialist


Carlos Mejia MS-QSM, RRT, RPGST, LSSB

cmejia@sunsethcs.com

Clinical Specialist

Marion Sachs RN, MSN, CNOR, CVAHP

msachs@sunsethcs.com



Corporate Office

Phone 872-275-3210

Fax 312-997-9985

180 N Michigan Ave Ste 2000
Chicago, IL 60601

Central Returns/Distribution Center

279 Madsen Dr Ste 101
Bloomington, IL 60108

West Returns/Distribution Center

6750 W Van Buren St
Phoenix, AZ 85043

CONTENT

Contract PP-NS-1571

Closed Suction Catheters **Page**

Standard Switch	4
Deluxe Switch	5-6

Contract PP-NS-1753

Medication Delivery **Page**

Nebulizers	7
Nebulizer Kits	8
Aerosol Masks	8
Nebulizer Filters	9
Inhaler Chamber	9
Peak Flow Meter	9
Adapters	9

Contract PP-NS-1741

Respiratory Disposables **Page**

Suction Catheter Kits	10
Suction Catheters	11
Trach Ties	11
Connectors	11
HMEs	12
Trach Masks	12
Spirometers	13
Drainage Bags	13
Oxygen Filters	14
Oxygen Connectors	15
Firesafe	16
Oxygen Cannula and Accessories	17-18
Oxygen Masks	19
Oxygen Supply Tubing	20
Oxygen Accessories	21
Vent Circuits and Filters	22

Contract PP-NS-1762

Resuscitation Devices **Page**

Resuscitators	23
Manometers	24
PEEP Valves	24
Reservoir Bags	25
Pocket Resuscitator Set	25
Breathing Masks	25

Contract PP-OR-2217


Laryngeal Masks **Page**

Silicone Disposable	26
Silicone Reusable	26
Flexible Silicone Disposable	26
Flexible Silicone Reusable	26
PVC Disposable	27
Flexible PVC Disposable	27

Standard Switch

Sterile and
latex-free

ISO color coded for
instant recognition of
catheter sizes.



Once the airway is clear, withdraw the catheter completely into the plastic sheath. A black mark on the catheter indicates full withdrawal.

Smooth tubing surface and atraumatic rounded catheter tip for minimal tracheal wall injury.

With the unique isolation valve, 72 hour closed suction system features safe and simple access for mechanically ventilated patients for 72 hours of continuous use.

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr
VEN46011400M	14Fr, MDI Port
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN56061000A	10Fr
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port

Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000600	305mm, 6Fr
VEN66000700	305mm, 7Fr
VEN66000800	305mm, 8Fr
VEN66001000	405mm, 10Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000M	405mm, 10Fr, MDI Port

Deluxe Switch (One-handed operation)

Sterile and latex-free

Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

Catheter Tip


Atraumatic rounded catheter tip for minimal tracheal wall injury.

Isolation Valve

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



Deluxe Switch (One-handed operation)

Sterile and latex-free

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN160160	16Fr
VEN160160M	16Fr, MDI Port

With Wedge + Adapter

Item#	Description
VEN162100A	10Fr
VEN162100AM	10Fr, MDI Port
VEN162120	12Fr
VEN162120M	12Fr, MDI Port
VEN162140	14Fr
VEN162140M	14Fr, MDI Port
VEN162160	16Fr
VEN162160M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN161120	12Fr
VEN161140	14Fr
VEN161160	16Fr
VEN161160M	16Fr, MDI Port

With Wedge + Adapter

Item#	Description
VEN163100A	10Fr
VEN163100AM	10Fr, MDI Port
VEN163120	12Fr
VEN163140	14Fr
VEN163160	16Fr

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN365160	16Fr
VEN365160M	16Fr, MDI Port

With Wedge + Adapter

Item#	Description
VEN367100A	10Fr
VEN367100AM	10Fr, MDI Port
VEN367120	12Fr
VEN367120M	12Fr, MDI Port
VEN367140	14Fr
VEN367140M	14Fr, MDI Port
VEN367160	16Fr
VEN367160M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN366160	16Fr
VEN366160M	16Fr, MDI Port

With Wedge + Adapter

Item#	Description
VEN368100A	10Fr
VEN368100AM	10Fr, MDI Port
VEN368120	12Fr
VEN368120M	12Fr, MDI Port
VEN368140	14Fr
VEN368140M	14Fr, MDI Port
VEN368160	16Fr
VEN368160M	16Fr, MDI Port

Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port

MEDICATION DELIVERY Contract PP-NS-1753

Compressor Nebulizers

Item#
NEB100
Each

Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

NEB100-R

- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Sunny the Bear

NEB300BEAR

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)



NEB100-BAG
Carry Bag

Gordon the Fish

NEB300FISH

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)



Handheld Compressor Nebulizer

Item#
NEB200
Each

Included

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (approximately 3 microns)
- No replacement filter required



Handheld Mesh Nebulizer

Item#
NEB401
Each

With internal rechargeable battery

Included

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual

- Small particle size: MMAD 2.2µm*
- Efficient: Low Residual Liquid
- Convenient: Small and Portable
- Quiet: ≤50dB (at one meter)
- Modern: Rechargeable Lithium Battery



MEDICATION DELIVERY Contract PP-NS-1753

Disposable Nebulizer Kit

Item# **RES090**

Adult Size

50/case

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

Item# **RES090A**

No Accessories

50/case

Item# **RES090NM**

- 6cc jet nebulizer bottle
- 7ft supply tube

50/case

Item# **RES092**

Pediatric Size

50/case

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



Disposable Nebulizer Kit with T-Piece

Item# **RES091**

Disposable
Single Patient

50/case

Item# **RES091NA**

Disposable
Single Patient,
No Aerosol

50/case

HCPCS# A7003



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

Reusable Nebulizer Kit with T-Piece

Item# **RES093**

Reusable — Single Patient

10/case

HCPCS# A7005

Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**

Disposable
Single Patient

50/case

HCPCS# A7003

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



Aerosol Mask

50/case

Item# **RES2110**

Adult Size

RES2210

Pediatric Size

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



MEDICATION DELIVERY Contract PP-NS-1753

Plastic Intake Filter with Head

Item# **BF017**

5/pack

HCPCS# A7013

Compatible with DeVilbiss Pulmo Series



Foam Filter

Item# **BF018**

10/pack

HCPCS# A7013

Compatible with Caire Medi-Mist, DeVilbiss Sunrise Compact, Evo AeroMax/SportsNeb/Super SportsNeb, Invacare Passport/Envoy/Envoy Jr./Select, John Bunn Nebulite EV/LT2/LX, MedLine SportsMist, Respironics MiniElite, Salter AIRE Plus, Schuco 2000/3000, Sunset NEB100/NEB300 Series



Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**

10/pack

HCPCS# A7013

- Mini
- MicroElite
- MiniElite



SperaChamber

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item# **RES195-10**
10/case

- Inhalation valve
- Universal adapter
- Spiral vents

RES195-50
50/case

- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask
10/case

Item #	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)
RES195ML	Large (Adult)

Peak Flow Meter

Item# **RES052**
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**
Pediatric Mouthpiece for Peak Flow Meter
10/case

Valved T-Adapter, Aerosol, 22mm ID x 22mm OD

Item# **NEB9001T**
10/case



Valved T-Adapter — Adult

Item# **VEN001**
Each
22F-18F-22M/15F



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Suction Catheter Kits

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	TRATK305
6fr	Y	DeLee	Pediatric	TRATK306
8fr	Y	DeLee	Pediatric	TRATK308
8fr	Y	Whistle	Adult	TRATK108
8fr	Chimney	Whistle	Adult	TRATK108CV
10fr	Y	Whistle	Adult	TRATK110
10fr	Chimney	Whistle	Adult	TRATK110CV
12fr	Y	Whistle	Adult	TRATK112
12fr	Chimney	Whistle	Adult	TRATK112CV
14fr	Y	Whistle	Adult	TRATK114
14fr	Chimney	Whistle	Adult	TRATK114CV
16fr	Y	Whistle	Adult	TRATK116
16fr	Chimney	Whistle	Adult	TRATK116CV

- Sterile
- Includes Catheter and Gloves



Suction Catheter Kits with Cup

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	TRATK405
6fr	Y	DeLee	Pediatric	TRATK406
8fr	Y	DeLee	Pediatric	TRATK408
8fr	Y	Whistle	Adult	TRATK208
8fr	Chimney	Whistle	Adult	TRATK208CV
10fr	Y	Whistle	Adult	TRATK210
10fr	Chimney	Whistle	Adult	TRATK210CV
12fr	Y	Whistle	Adult	TRATK212
12fr	Chimney	Whistle	Adult	TRATK212CV
14fr	Y	Whistle	Adult	TRATK214
14fr	Chimney	Whistle	Adult	TRATK214CV
16fr	Y	Whistle	Adult	TRATK216
16fr	Chimney	Whistle	Adult	TRATK216CV

- Sterile
- Includes Catheter and Gloves



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Suction Catheters

Size	Valve Type	Tip	Intended Use	Packaged	Item#	Qty
5fr	Y	DeLee	Pediatric	Coiled	TRASC305C	50
5fr	Y	DeLee	Pediatric	Straight	TRASC305S	100
6fr	Y	DeLee	Pediatric	Coiled	TRASC306C	50
6fr	Y	DeLee	Pediatric	Straight	TRASC306S	100
8fr	Y	DeLee	Pediatric	Coiled	TRASC308C	50
8fr	Y	DeLee	Pediatric	Straight	TRASC308S	100
10fr	Y	Whistle	Adult	Coiled	TRASC110C	50
10fr	Y	Whistle	Adult	Straight	TRASC110S	100
12fr	Y	Whistle	Adult	Coiled	TRASC112C	50
12fr	Y	Whistle	Adult	Straight	TRASC112S	100
14fr	Y	Whistle	Adult	Coiled	TRASC114C	50
14fr	Y	Whistle	Adult	Straight	TRASC114S	100
16fr	Y	Whistle	Adult	Coiled	TRASC116C	50
16fr	Y	Whistle	Adult	Straight	TRASC116S	100




- Sterile
- Single patient use
- Latex-free

Tracheostomy Ties with Fastener Tabs

20/case

One-piece
• Easy on



Item#
RES029 Adult
RES029P Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

Two-piece
• Easy off

Item#
RES030 Adult
RES030P Pediatric



Tracheostomy Swivel Elbow Connectors

10/pack
15mm ID x 15mm OD

Item# **RES021**
With suction port



Visit SunsetHCS.com/AcuteCare

RESPIRATORY DISPOSABLES Contract PP-NS-1741

Flex HME

Item# **RES027FLEX-40**
40/case

TV Range (mL)	300-1500
Moisture Output (mg H ₂ O/L)	33 mg/L
Weight (g)	17
Resistance (cm H ₂ O)	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (mL)	75
Connection (mm)	22M,15F



Heat Moisture Exchangers

30/case



RES027-30
T-HME



RES027OP-30
T-HME with
Oxygen Port



RES027IL-30
Inline HME —
Large Volume



RES027NE-10
Inline HME —
Neonatal

	RES027-30 T-HME	RES027OP-30 T-HME with Oxygen Port	RES027IL-30 Inline HME — Large Volume	RES027NE-10 Inline HME — Neonatal
TV Range (mL)	250-1000	250-1000	50-250	15-50
Moisture Output (mg H ₂ O/L)	30, Vt = 0.1L	29.50	31, Vt = 0.1L	30
Weight (g)	4.5	7	10	4.5
Resistance (cm H ₂ O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM
Dead Space (mL)	12.5	12	13	2.5
Connection (mm)	15F	15F 5.5mm Oxygen port	15M,22M/15F	15M,15F

Tracheostomy Mask with O₂ Connector

50/case

Item# **RES2130**
Adult Regular Size

RES2130L
Adult Large Size

RES2230
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing

Adapter for T-HME

Item# **RES027A**
Each

RES027A-30
30/pack



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Volumetric Incentive Spirometer

12/case

Item# **RES055**
2500ml

RES056
5000ml

RES055P
Pediatric



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design

Three Ball Spirometer

10/case

Item# **RES057-10**

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

Drainage Bag with T-Adapter

Item# **RES070**
50/case
750cc



Drainage Bag with Y-Adapter

Item# **RES071**
50/case
750cc




- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

RESPIRATORY DISPOSABLES Contract PP-NS-1741

EverFlo Intake Filter

Item# **BF1001**
Each



Final Bacteria Filter — 3/16in

Item# **BF316**
10/pack
3/16 Stepped Barb and
1/4in Straight Barb
on Both Sides




Final Bacteria Filter

Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack
OEM# P0750X030



CPAP Bacteria Filter

Item# **OF6011**
Extended Life
Bacteria Filter
22mm fittings
Each
Use if SolAiris was manufactured
before Nov '03



Foam Cabinet Filter

Item# **OF9001**
Foam Cabinet Filter
10/pack
6 3/4" X 3 3/8" X 1/2"
OEM# 2000489



Oxygen Filter Compatibility

Covidien	
590 and 592A	BF316

SeQual	
Integra	BF516

DeVilbiss/Drive	
SolAiris V, 505, 515 Series	OF6011
525 DS Series	OF6011

Invacare	
Mobilaire 3 and 5	OF9001
Mobilaire 3 and 5	BF516
PerfectO ₂	BF516
Platinum 5 and 10	BF516
HomeFill	OF9001

3B	
Stratus 5	BF316

Respironics/Healthdyne	
EverGo	BF316
Millennium	BF316
Alliance	BF316
EverFlo	BF1001
EverFlo	BF316

RESPIRATORY DISPOSABLES Contract PP-NS-1741

Christmas Tree Connector

Item# **RES001**
Green without Swivel
10/pack

RES001-30 30/pack



Item# **RES002**
Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Connector — Non-Swivel

Item# **RES012**
Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Tubing Reducer Connector

Item# **RES011**
Fits 3/16in, 1/4in,
5/16in, 3/8in, and 7/16in
ID Tubing
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

Humidifier Tubing Adapter

Item# **RES016**
22mm ID x
22mm OD
with Lip
10/pack



Circuit Connectors

Item# **RES013**
22mm OD x
15mm ID/
22mm OD
10/pack



Item# **RES014**
22mm OD x
22mm OD
With Lip
10/pack



Oxygen Tubing Connector — Swivel

Item# **RES018**
Male/Male **10/pack**

RES018-25 25/pack



Oxygen Tubing Y-Connector

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack

- Connects two oxygen sources for higher flows



Oxygen Bleed-In Adapter

Item# **RES019** **10/pack**
22mm with 6mm
Rotating Stem



RESPIRATORY DISPOSABLES Contract PP-NS-1741

FireSafe Nozzle

Item#
RES010FND
10/pack




Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients



* 1 death was reported among 73% of the patient population (2013-17)

- A thermal fuse designed to stop the flow of oxygen in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Adult Oxygen Cannula

Item #	Tubing Length	Connector	Qty
RES1100	None	Standard	50
RES1104	4 ft	Standard	50
RES1107	7 ft	Standard	50
RES1107U	7 ft	Universal	50
RES1110	10 ft	Standard	50
RES1114N	14 ft	Standard	50
RES1114U	14 ft	Universal	50
RES1115	15 ft	Standard	50
RES1115U	15 ft	Universal	50
RES1125	25 ft	Standard	25
RES1150U	50 ft	Universal	20



Soft Cannula

Item #	Tubing Length	Connector	Qty
RES1100S	None	Standard	50
RES1107S	7 ft	Standard	50
RES1107SU	7 ft	Universal	50
RES1114S	14 ft	Standard	50
RES1114SU	14 ft	Universal	50
RES1115S	15 ft	Standard	50
RES1125S	25 ft	Standard	25



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Pediatric, Infant, and Neonate Cannula

50/case

Item #	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No

Adult Cannula with Ear Cushions

Item #	Tubing Length	Connector	Touch	Tip	QTY
RES1100EC	None	Standard	Standard	Standard	50
RES1107EC	7 ft	Standard	Standard	Standard	50
RES1107ECU	7 ft	Universal	Standard	Standard	50
RES1114EC	14 ft	Standard	Standard	Standard	50
RES1115EC	15 ft	Standard	Standard	Standard	50
RES1125EC	25 ft	Standard	Standard	Standard	25

Ear Cushions for Cannula

Item# **RES020**
Cushions Only
50 Pairs/case



• Add cushions to existing cannula



Secure Cannula Stickers

Item# **RES080**
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



RESPIRATORY DISPOSABLES Contract PP-NS-1741


Simple O₂ Mask

With 7ft Safety Tubing

Item# **RES2100**
Adult Size
50/case

HCPCS# A4620

- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



High Concentration Oxygen Mask

50/case

Item #	Size	Tubing Length
RES2202A	Adult	None
RES2202A7	Adult	7 ft
RES2202I7	Infant	7 ft

Medium Concentration Oxygen Mask

50/case

Item #	Size	Tubing Length	Connector
RES2100NT	Adult	None	Standard
RES2100U	Adult	7 ft	Universal
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203A	Adult	None	Standard
RES2203A7	Adult	7 ft	Standard
RES2203I7	Infant	7 ft	Standard
RES2203P7	Pediatric	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather Oxygen Mask

50/case

Item #	Size	Tubing Length	Connector	Safety Vent
RES2102	Adult	7 ft	Standard	No
RES2102I	Infant	7 ft	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No
RES2102S	Adult	7 ft	Standard	Yes
RES2102SM	Adult	None	Standard	Yes
RES2102U	Adult	7 ft	Universal	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Oxygen Supply Tubing

Item #	Length	Connector	Color	Qty
RES3007	7 ft	Standard	None	50
RES3007G	7 ft	Standard	Green	50
RES3007U	7 ft	Universal	None	50
RES3015	15 ft	Standard	None	50
RES3025	25 ft	Standard	None	25
RES3025G	25 ft	Standard	Green	25
RES3025U	25 ft	Universal	None	25
RES3040G	40 ft	Standard	Green	20
RES3050G	50 ft	Standard	Green	20

- 6-channel
- PHT/DEHP-free



Universal Connector



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Standard Finger Pulse Oximeter

Item# **RES5100**
Each

- 1 display mode
- Digital LED display
- SpO2, PR, and pulse bar



Advanced Finger Pulse Oximeter

Item# **RES5101**
Each

- 6 display modes
- Dual color OLED display
- SpO2, PR, pulse bar and waveform



Wrenches

Item# **RES031**
With Chain **10/pack**



Water Traps

Item# **RES037**
Push On End Caps • **Each**

- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position



Oxygen Liter Meter

Item# **RES034**
0-10 Liter/Minute
Each

- Oxygen pen
- Works with tanks and concentrators



Humidifier Connectors

Item# **HUM008**
L-Connector 10in • **10/pk**

Item# **HUM008L**
L-Connector 21in • **10/pk**


Item# **HUM008S**
Straight Connector 9in • **10/pk**



RESPIRATORY DISPOSABLES Contract PP-NS-1741

Aerosol Tubing

Item# **TUB100-6B**
Corrugated Aerosol Tubing
100ft in 6in segments
Each



Passive Non-Heated Breathing Circuits

For use with Respironics Trilogy


Item# **VC1100**
With fixed exhalation port
Each
OEM# 1069467 & 1069468



Non-Heated Breathing Circuits

For use with Respironics Trilogy
22mm cuffs
30/case

Item# **VC1001**
Active Adult,
No Water Trap
OEM#
1073223
1073232




Single Use Exhalation Port

Item# **CAP011**
10/case



Passive Heated Breathing Circuit

Item# **VC1114Z**
For use with MR850 Humidifier
10/case



Bacterial/Viral HME Filter — Adult

Item# **BF104**
15 mm ID x 15 mm ID/
22 mm OD
Each
HCPCS# A7507



Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item# **BF2200**
22mm Bacteria
Disk Filter
5/pack



Attach between the machine and tube for additional protection against bacteria and allergens

Filter Specifications

Bacterial Filtration Efficiency: (BFE) 99.99+%
Viral Filtration Efficiency: (VFE) 99.9+%
NaCl Particulate Filtration Efficiency: 96+%
Flow Resistance: Approximately 8.8 mmH2O at 30 LPM
Internal Volume: 66.3 mL without fittings / 78.3 mL with fittings
Connections:
• Patient End: 22 mm OD / 15 mm ID
• Machine End: 22 mm ID

RESUSCITATION DEVICES Contract PP-NS-1762

Adult Resuscitator

12/case

Item #	Pop-Off	Mask
RES8100	No	None
RES8100H	Yes	None
RES8101	Yes	#4
RES8102	Yes	#5



Child Resuscitator

12/case

Item #	Pop-Off	Mask
RES8110H	Yes	None
RES8111	Yes	#2
RES8112	Yes	#3



Infant Resuscitator

12/case

Item #	Pop-Off	Mask
RES8120H	Yes	None
RES8121	Yes	#0
RES8122	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

RESUSCITATION DEVICES Contract PP-NS-1762

MANOMETERS AND PEEP VALVES

Disposable Pressure Manometer

25/case
Item# **RES8201**



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

Disposable Pressure Manometer, L-Shaped

25/case
Item# **RES8201L**



Reusable PEEP Valve, Orange

10/case
2-10cm H2O

Item# **RES8203R**
15mm ID

Item# **RES8203RL**
30mm ID



Reusable PEEP Valve, Blue

10/case
5-20cm H2O

Item# **RES8204B**
15mm ID

Item# **RES8204BL**
30mm ID



- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

Disposable PEEP Valve, Orange

10/case
2.5-10cm H2O

Item# **RES8205R**
15mm ID

Item# **RES8205RL**
30mm ID



Disposable PEEP Valve, Blue

10/case
5-20cm H2O

Item# **RES8205B**
15mm ID

Item# **RES8205BL**
30mm ID



- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

RESUSCITATION DEVICES Contract PP-NS-1762

Reservoir Bag, PVC

25/case

Item #	Size
RES8300P	0.5L
RES8301P	1L
RES8302P	2L



Reservoir Bag, Silicone

25/case

Item #	Size
RES8300S	0.5L
RES8301S	1L
RES8302S	2L



Suitable for reusable resuscitator.

CPR Pocket Resuscitator Set

12/case

Item# **RES8200**

- High transparency for better visualization
- Non-reusable, single use only
- 100% Medical level PVC material
- Latex-free components
- Economical choice
- Blow molding



Breathing Mask

20/case

Item #	Size	Hooks
SBC01N	Neonatal (#0)	No
SBC01I	Infant (#1)	No
SBC01P	Pediatric (#2)	No
SBC01S	Adult Small (#3)	No
SBC01M	Adult Medium (#4)	No
SBC01L	Adult Large (#5)	No
SBC01XL	Adult Extra Large (#6)	No
SBC02N	Neonatal (#0)	Yes
SBC02I	Infant (#1)	Yes
SBC02P	Pediatric (#2)	Yes
SBC02S	Adult Small (#3)	Yes
SBC02M	Adult Medium (#4)	Yes
SBC02L	Adult Large (#5)	Yes
SBC02XL	Adult Extra Large (#6)	Yes



LARYNGEAL MASKS Contract PP-OR-2217

Silicone Disposable Laryngeal Masks

Item #	Size
ALM1001	#1
ALM1015	#1.5
ALM1002	#2
ALM1025	#2.5
ALM1003	#3
ALM1004	#4
ALM1005	#5



- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use

Silicone Reusable Laryngeal Masks


Item #	Size
ALM1001R	#1
ALM1015R	#1.5
ALM1002R	#2
ALM1025R	#2.5
ALM1003R	#3
ALM1004R	#4
ALM1005R	#5



- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch

Flexible Silicone Disposable Laryngeal Masks

Item #	Size
ALM3001	#1
ALM3015	#1.5
ALM3002	#2
ALM3025	#2.5
ALM3003	#3
ALM3004	#4
ALM3005	#5



- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness

Flexible Silicone Reusable Laryngeal Masks

Item #	Size
ALM3001R	#1
ALM3015R	#1.5
ALM3002R	#2
ALM3025R	#2.5
ALM3003R	#3
ALM3004R	#4
ALM3005R	#5




LARYNGEAL MASKS Contract PP-OR-2217

PVC Disposable Laryngeal Masks

Item #	Size
ALM2001	#1
ALM2015	#1.5
ALM2002	#2
ALM2025	#2.5
ALM2003	#3
ALM2004	#4
ALM2005	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



Flexible PVC Disposable Laryngeal Masks

Item #	Size
ALM4001	#1
ALM4015	#1.5
ALM4002	#2
ALM4025	#2.5
ALM4003	#3
ALM4004	#4
ALM4005	#5



Size	Patient Weight and Age	
#1	Up to 5 kg	Pediatric
#1.5	5-10 kg	Pediatric
#2	10-20 kg	Pediatric
#2.5	20-30 kg	Pediatric
#3	30-50 kg	Pediatric
#4	50-70 kg	Adult
#5	70-100 kg	Adult

