




# Product Catalog



**vizient**<sup>®</sup>  
Awarded Supplier

Contract

MS1097  
MS1049  
MS1087



Sales Representative  
**Egle Sventnickaite**  
O: 872-258-1654  
esventnickaite@sunsethcs.com

Director, Sales Operations  
**Mark Lipson**  
O: 312-533-2465  
mlipson@sunsethcs.com

Director, Strategic Accounts  
**Lisa Iannazzone** RN, BSN  
C: 636-222-3134  
lisa@sunsethcs.com

Sales Representative  
**Christian Straussner**  
O: 872-282-0746  
C: 847-525-6340  
cstrausner@sunsethcs.com


Sales and Marketing Manager  
**Jackson Hanft**  
C: 971-416-6695  
jackson@sunsethcs.com

RRT Clinical Specialist  
**Carlos Mejia** MS-QSM, RRT, RPGST, LSSB  
cmejia@sunsethcs.com

Sales Representative  
**Owen Ruth**  
O: 312-533-2451  
oruth@sunsethcs.com

Clinical Specialist  
**Marion Sachs** RN, MSN, CNOR, CVAHP  
msachs@sunsethcs.com

Contract Manager, Strategic Accounts  
**Gina Rosenberg**  
C: 914-837-5215  
grosenberg@sunsethcs.com



**Corporate Office**

Phone 872-275-3210

Fax 312-997-9985

180 N Michigan Ave Ste 2000  
Chicago, IL 60601

**Central Returns/Distribution Center**

279 Madsen Dr Ste 101  
Bloomington, IL 60108

**West Returns/Distribution Center**

6750 W Van Buren St  
Phoenix, AZ 85043

# CONTENT

Contract MS1097

## **MDD**

**Page**

Nebulizer Kits	4
Nebulizer Compressor	5
Inhaler Chamber	5
Peak Flow Meter	5
Valved T-Adapter	5

Contract MS1049

## **NIV**

**Page**

Masks	6-8
Mask Elbows	9
Connectors/Ports	9
Headgear	9

Contract MS1087

## **Respiratory Disposables**

**Page**

Breathing Circuits	10
Bacterial/Viral HME Filter	10
Connectors and Adapters	10
Heat Moisture Exchangers	11
Tracheostomy Masks	11
Berman and Guedel Airways	12
Spirometers	13
Aerosol Mask, Tubing, Valved T-Adapter	13
Oxygen Cannula	14-15
Oxygen Masks	16
Oxygen Tubing	17


## Disposable Nebulizer Kits with Mouthpiece

50/case

All include 6cc jet nebulizer medicine cup

**Accessory options**

- T-piece
- Mouthpiece
- Corrugated tube
- 7ft supply tube



REF	Size	Accessories
<b>RES091L</b>	Adult	All
<b>RES091T</b>	Adult	No corrugated tube

- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal


## Disposable Nebulizer Kits with Mask

50/case

All include 6cc jet nebulizer medicine cup

**Accessory options**

- Mask
- 7ft supply tube



REF	Size	Accessories
<b>RES090NM</b>	Adult	No Mask
<b>RES092</b>	Pediatric	All

- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble

### Medicine Cup

Item#

**RES090A**

50/case

## Nebulizer Compressor

Item# **NEB100**  
Each



### Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

### Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

## Peak Flow Meter

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece  
for Peak Flow Meter  
10/case

## Valved T-Adapter

Item# **VEN001**  
22F-18F-22M/15F  
10/pack



## SperaChamber

Anti-Static Valved Inhaler Chamber  
for Metered Dose Inhaler

Item# **RES195-10**  
10/case

Item# **RES195-50**  
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask  
10/case

REF	Size
<b>RES195MS</b>	Small (Infant)
<b>RES195MM</b>	Medium (Pediatric)
<b>RES195ML</b>	Large (Adult)



Adult Mask

**Full Face Mask**

**With Standard Headgear**

**Standard Elbow, Non-Vented**

Size	Item#
Small	<b>CM041S</b>
Medium	<b>CM041M</b>
Large	<b>CM041L</b>

**Single Limb Elbow, Non-Vented, AAV**

Size	Item#
Small	<b>CM041SEE</b>
Medium	<b>CM041MEE</b>
Large	<b>CM041LEE</b>

**Single Limb Elbow, Vented, AAV**

Size	Item#
Small	<b>CM041SEEV</b>
Medium	<b>CM041MEEV</b>
Large	<b>CM041LEEV</b>



**Standard Elbow, Non-Vented**

Size	Item#
Small	<b>VCM041S</b>
Medium	<b>VCM041M</b>
Large	<b>VCM041L</b>

**Single Limb Elbow, Non-Vented, AAV**

Size	Item#
Small	<b>VCM041SEE</b>
Medium	<b>VCM041MEE</b>
Large	<b>VCM041LEE</b>

**Single Limb Elbow, Vented, AAV**

Size	Item#
Small	<b>VCM041SEEV</b>
Medium	<b>VCM041MEEV</b>
Large	<b>VCM041LEEV</b>



- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

HCPCS#  
 MASK - A7030  
 HEADGEAR - A7035  
 CUSHION - A7031

**With Cap Style Headgear**

**Standard Elbow, Non-Vented**

Size	Item#
Small	<b>CM041SC</b>
Medium	<b>CM041MC</b>
Large	<b>CM041LC</b>

**Single Limb Elbow, Non-Vented, AAV**

Size	Item#
Small	<b>CM041SEEC</b>
Medium	<b>CM041MEEC</b>
Large	<b>CM041LEEC</b>

**Single Limb Elbow, Vented, AAV**

Size	Item#
Small	<b>CM041SEEVC</b>
Medium	<b>CM041MEEVC</b>
Large	<b>CM041LEEVC</b>

**Nasal Mask**

**With Standard Headgear**

Size	Item#
Small	<b>CM008SV</b>
Medium	<b>CM008MV</b>
Large	<b>CM008LV</b>



Size	Item#
Small	<b>VCM008S</b>
Medium	<b>VCM008M</b>
Large	<b>VCM008L</b>



- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

**Total Face Mask**

**Standard Elbow**

Size	Item#
Small	<b>CM040S</b>
Medium	<b>CM040M</b>
Large	<b>CM040L</b>

**EE Elbow**

Size	Item#
Medium	<b>CM040MEE</b>
Large	<b>CM040LEE</b>



**Veraseal 2 Disposable Gel Mask with Headgear**



**Standard Elbow, Non-Vented**

**Sleepnet**

Size	Item#
Small	<b>CM50933</b>
Medium	<b>CM50934</b>
Large	<b>CM50935</b>
Extra Large	<b>CM50936</b>

**novaplus**

Size	Item#
Small	<b>VCM053S</b>
Medium	<b>VCM053M</b>
Large	<b>VCM053L</b>
Extra Large	<b>VCM053XL</b>

**Single Limb Elbow, Non-Vented, AAV**

**Sleepnet**

Size	Item#
Small	<b>CM50937</b>
Medium	<b>CM50938</b>
Large	<b>CM50939</b>
Extra Large	<b>CM50940</b>

**novaplus**

Size	Item#
Small	<b>VCM054S</b>
Medium	<b>VCM054M</b>
Large	<b>VCM054L</b>
Extra Large	<b>VCM054XL</b>

**Single Limb Elbow, Vented, AAV**

**Sleepnet**

Size	Item#
Small	<b>CM50827</b>
Medium	<b>CM50828</b>
Large	<b>CM50829</b>
Extra Large	<b>CM50830</b>

**novaplus**

Size	Item#
Small	<b>VCM052S</b>
Medium	<b>VCM052M</b>
Large	<b>VCM052L</b>
Extra Large	<b>VCM052XL</b>

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.
- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.
- The active headgear connector moves as the patient moves for a reliable seal.
- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

**PEDIATRIC NIV MASKS**

**Innova**

**Nasal Mask with Headgear**

**Vented**

Size	Item#
Small/Medium	<b>CM50180</b>
Medium/Large	<b>CM50181</b>

**Non-Vented**

Size	Item#
Small/Medium	<b>CM50594</b>
Medium/Large	<b>CM50595</b>

**novaplus**

Size	Item#
Small/Medium	<b>VCM034SM</b>
Medium/Large	<b>VCM034ML</b>

**novaplus**

Size	Item#
Small/Medium	<b>VCM034NVSM</b>
Medium/Large	<b>VCM034NVML</b>

Innova® CPAP masks feature patented AIR°gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA



## MiniMe 2 Pediatric Nasal Mask



### Vented

Size	Item #
Extra Small	<b>CM60216</b>
Small	<b>CM60217</b>
Medium	<b>CM60218</b>
Large	<b>CM60221</b>

### Non-Vented

Size	Item #
Extra Small	<b>CM60256</b>
Small	<b>CM60257</b>
Medium	<b>CM60258</b>
Large	<b>CM60261</b>



### Vented

Size	Item #
Extra Small	<b>VCM055XS</b>
Small	<b>VCM055S</b>
Medium	<b>VCM055M</b>
Large	<b>VCM055L</b>

### Non-Vented

Size	Item #
Extra Small	<b>VCM056XS</b>
Small	<b>VCM056S</b>
Medium	<b>VCM056M</b>
Large	<b>VCM056L</b>



- Maximum comfort, minimum weight. Patented AIR<sup>o</sup>gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA

## Ray Pediatric Masks Comparable to Nonny masks

### Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

HCPCS#  
 MASK - A7030  
 HEADGEAR - A7035  
 CUSHION - A7031



### Nasal Mask with Headgear, Tube, and Replaceable Cushions

Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>

HCPCS#  
 MASK - A7034  
 HEADGEAR - A7035  
 CUSHION - A7032





# Standard Elbow

**NIV** Contract MS1049

## Standard Elbow, Non-Vented

Item# **CAP009**  
10/case



## Single Limb Elbow, Non-Vented, AAV

Item# **CAP008**  
10/case



## Single Limb Elbow, Vented, AAV

Item# **CAP010**  
10/case



## Swivel Connector Silencer

Item# **VEN002**  
10/case



## Single Use Exhalation Port

Item# **CAP011**  
10/case



## Standard 4-Point Headgear

Item# **HG003**



## Cap Style Headgear

Item# **HG002**



# RESPIRATORY DISPOSABLES Contract MS1087

BREATHING CIRCUITS, HME FILTER, CONNECTORS AND ADAPTERS



## Passive Heated Breathing Circuit

Item#

**VC1114Z**

For use with MR850 Humidifier

10/case



## Non-Heated Dual Limb Breathing Circuits

10/case

REF	Size	Humidifier Limb	Water Trap
<b>VC2001PZ</b>	Pediatric	No	No
<b>VC2001Z</b>	Adult	No	No
<b>VC2002PZ</b>	Pediatric	Yes	No
<b>VC2002Z</b>	Adult	Yes	No
<b>VC2003Z</b>	Adult	Yes	Yes

## Bacterial/Viral HME Filter

Item#

**BF2200Z**

5/pack



## Xmas Tree Connector

Item# **RES001-10**

10/pack




- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Tubing Connector

Item# **RES012-10**

Fits 3/16in and 1/4in ID Tubing

10/pack



- Barbed ends ensure secure, leak-free connections

## Circuit Connectors

Item# **RES014-10**

22mm OD x  
22mm OD  
With Lip  
10/pack



## Humidifier Tubing Adapter

Item# **RES016**

22mm ID x  
22mm OD  
with Lip  
10/pack



# RESPIRATORY DISPOSABLES Contract MS1087

## Flex HME






Item# **RES027FLEX-40**  
40/case

TV Range (mL)	300-1500
Moisture Output (mg H2O/L)	33 mg/L
Weight (g)	17
Resistance (cm H2O)	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (mL)	75
Connection (mm)	22M,15F



## Heat Moisture Exchangers

30/case

	 <b>RES027-30</b> T-HME	 <b>RES0270P-30</b> T-HME with Oxygen Port	 <b>RES0271L-30</b> Inline HME — Large Volume	 <b>RES0271LOHME</b> Inline HME — Large Volume with Oxygen Port	 <b>RES027NE-10</b> Inline HME — Neonatal
TV Range (mL)	250-1000	250-1000	50-250	50-250	15-50
Moisture Output (mg H2O/L)	30, Vt = 0.1L	29.50	31, Vt = 0.1L	31, Vt = 0.1L	30
Weight (g)	4.5	7	10	10	4.5
Resistance (cm H2O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM
Dead Space (mL)	12.5	12	13	13	2.5
Connection (mm)	15F	15F 5.5mm Oxygen port	15M,22M/15F	15M,22M/15F 5.5mm Oxygen port	15M,15F

## Adapter for T-HME

Item# **RES027AZ**  
Each



## Tracheostomy Mask

50/case

Item# **RES2130Z**  
Adult Regular Size

**RES2130L**  
Adult Large Size

**RES2230Z**  
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

# RESPIRATORY DISPOSABLES Contract MS1087

## Berman Oral Airway

50/pack

REF	Size	Color
<b>AOA40BC-50</b>	40mm	Pink
<b>AOA50BC-50</b>	50mm	Blue
<b>AOA60BC-50</b>	60mm	Black
<b>AOA70BC-50</b>	70mm	White
<b>AOA80BC-50</b>	80mm	Green
<b>AOA90BC-50</b>	90mm	Yellow
<b>AOA100BC-50</b>	100mm	Red
<b>AOA110BC-50</b>	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Oropharyngeal Guedel Airway

10/pack

REF	Size	Color
<b>AOP500GC</b>	50mm	Blue
<b>AOP601GC</b>	60mm	Black
<b>AOP702GC</b>	70mm	White
<b>AOP803GC</b>	80mm	Green
<b>AOP904GC</b>	90mm	Yellow
<b>AOP1005GC</b>	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

# RESPIRATORY DISPOSABLES Contract MS1087

## Volumetric Incentive Spirometer

12/case

Item# **RES055**  
2500ml

**RES056**  
5000ml

**RES055P**  
Pediatric



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.

## Three Ball Spirometer

12/case

Item# **RES057**

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

## Aerosol Mask

50/case

Item# **RES2110Z**  
Adult Size

Item# **RES2110F**  
Adult Size  
w/ Flex Tube


Item# **RES2210Z**  
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

## Aerosol Tubing

Item# **TUB100-6**  
Corrugated Aerosol Tubing  
100ft in 6in segments



## Valved T-Adapter, Aerosol, 22mm ID x 22mm OD

Item# **NEB9001T**  
30/case



# RESPIRATORY DISPOSABLES Contract MS1087

## Adult Oxygen Cannula

REF	Tubing Length	Connector	Qty
<b>RES1100Z</b>	None	Standard	50
<b>RES1104Z</b>	4 ft	Standard	50
<b>RES1107U</b>	7 ft	Universal	50
<b>RES1107Z</b>	7 ft	Standard	50
<b>RES11110</b>	10 ft	Standard	50
<b>RES1114N</b>	14 ft	Standard	50
<b>RES1114U</b>	14 ft	Universal	50
<b>RES1115U</b>	15 ft	Universal	50
<b>RES1115Z</b>	15 ft	Standard	50
<b>RES1125Z</b>	25 ft	Standard	25
<b>RES1135</b>	35 ft	Standard	20
<b>RES1150U</b>	50 ft	Universal	20
<b>RES1150Z</b>	50 ft	Standard	20



## Soft Cannula

REF	Tubing Length	Connector	Qty
<b>RES1107SU</b>	7 ft	Universal	50
<b>RES1107SZ</b>	7 ft	Standard	50
<b>RES1114S</b>	14 ft	Standard	50
<b>RES1114SU</b>	14 ft	Universal	50
<b>RES1125S</b>	25 ft	Standard	25
<b>RES1125SU</b>	25 ft	Universal	25

## Flared Cannula

REF	Tubing Length	Connector	Qty
<b>RES1100F</b>	None	Standard	50
<b>RES1107F</b>	7 ft	Standard	50
<b>RES1107FU</b>	7 ft	Universal	50
<b>RES1114F</b>	14 ft	Standard	50
<b>RES1125F</b>	25 ft	Standard	25
<b>RES1150F</b>	50 ft	Standard	20

## Adult Cannula with Ear Cushions

REF	Tubing Length	Connector	Touch	Tip	QTY
<b>RES1100ECZ</b>	None	Standard	Standard	Standard	50
<b>RES1104EC</b>	None	Standard	Standard	Standard	50
<b>RES1107ECU</b>	7 ft	Universal	Standard	Standard	50
<b>RES1107ECZ</b>	7 ft	Standard	Standard	Standard	50
<b>RES1114EC</b>	14 ft	Standard	Standard	Standard	50
<b>RES1125EC</b>	25 ft	Standard	Standard	Standard	25



# RESPIRATORY DISPOSABLES Contract MS1087

## Pediatric, Infant, and Neonate Cannula

50/case

REF	Size	Tubing Length	Connector	Touch	Ear Cushions
<b>RES1100P</b>	Pediatric	None	Standard	Standard	No
<b>RES1107I</b>	Infant	7 ft	Standard	Standard	No
<b>RES1107NE</b>	Neonate	7 ft	Standard	Standard	No
<b>RES1107P</b>	Pediatric	7 ft	Standard	Standard	No
<b>RES1107PEC</b>	Pediatric	7 ft	Standard	Standard	Yes
<b>RES1107PECU</b>	Pediatric	7 ft	Universal	Standard	Yes
<b>RES1107SI</b>	Infant	7 ft	Standard	Soft	No
<b>RES1107SNE</b>	Neonate	7 ft	Standard	Soft	No
<b>RES1107SP</b>	Pediatric	7 ft	Standard	Soft	No
<b>RES1107SPU</b>	Pediatric	7 ft	Universal	Soft	No
<b>RES1114SP</b>	Pediatric	14 ft	Standard	Soft	No
<b>RES1114SPU</b>	Pediatric	14 ft	Universal	Soft	No

## Humidoflo High Flow Nasal Cannula


10/case

REF	Size
<b>HFN001AS</b>	Small
<b>HFN001AM</b>	Medium
<b>HFN001AL</b>	Large



- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 - 60 L/min and O<sub>2</sub> concentrations range from 21 - 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.



# RESPIRATORY DISPOSABLES Contract MS1087

## Medium Concentration Oxygen Mask

50/case

REF	Size	Tubing Length	Connector
<b>RES2100NT</b>	Adult	None	Standard
<b>RES2100U</b>	Adult	7 ft	Universal
<b>RES2100Z</b>	Adult	7 ft	Standard
<b>RES2203P</b>	Pediatric	None	Standard
<b>RES2100P</b>	Pediatric	7 ft	Standard
<b>RES2203I</b>	Infant	None	Standard
<b>RES2203I7</b>	Infant	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

## High Concentration Oxygen Mask

50/case

REF	Size	Tubing Length
<b>RES2202A</b>	Adult	None
<b>RES2202A7</b>	Adult	7 ft
<b>RES2202I</b>	Infant	None
<b>RES2202I7</b>	Infant	7 ft
<b>RES2202P</b>	Pediatric	None
<b>RES2202P7</b>	Pediatric	7 ft

## Non-Rebreather Oxygen Mask

50/case

REF	Size	Tubing Length	Connector	Safety Vent
<b>RES2102I</b>	Infant	7 ft	Standard	No
<b>RES2102IM</b>	Infant	None	Standard	No
<b>RES2102P</b>	Pediatric	7 ft	Standard	No
<b>RES2102PS</b>	Pediatric	7 ft	Standard	Yes
<b>RES2102PSM</b>	Pediatric	None	Standard	Yes
<b>RES2102PSU</b>	Pediatric	7ft	Universal	Yes
<b>RES2102PU</b>	Pediatric	7 ft	Universal	No
<b>RES2102S</b>	Adult	7 ft	Standard	Yes
<b>RES2102SM</b>	Adult	None	Standard	Yes
<b>RES2102U</b>	Adult	7 ft	Universal	No
<b>RES2102Z</b>	Adult	7 ft	Standard	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape





# RESPIRATORY DISPOSABLES Contract MS1087

## Oxygen Supply Tubing

REF	Length	Connector	Color	Bore	Qty
<b>RES3007SB</b>	7 ft	Standard	None	Smooth	50
<b>RES3007U</b>	7 ft	Universal	None	Channeled	50
<b>RES3007Z</b>	7 ft	Standard	None	Channeled	50
<b>RES3015GZ</b>	15 ft	Standard	Green	Channeled	50
<b>RES3015Z</b>	15 ft	Standard	None	Channeled	50
<b>RES3025GZ</b>	25 ft	Standard	Green	Channeled	25
<b>RES3025U</b>	25 ft	Universal	None	Channeled	25
<b>RES3025Z</b>	25 ft	Standard	None	Channeled	25
<b>RES3030</b>	30 ft	Standard	None	Channeled	20
<b>RES3030G</b>	30 ft	Standard	Green	Channeled	20
<b>RES3040GZ</b>	40 ft	Standard	Green	Channeled	20
<b>RES3050GZ</b>	50 ft	Standard	Green	Channeled	20
<b>RES3050Z</b>	50 ft	Standard	None	Channeled	20



**Universal Connector**



