

SperaChamber

**Anti-static valved inhaler chamber
for metered dose inhaler**

SperaChamber Only

Item#

RES195-10 10/pack

RES195-50 50/pack

SperaChamber with Dog Mask — Pediatric (3-6 Years)

Item#

RES195MS Small

RES195MM Medium



- Latex Free
- DEHP Free
- Universal adapter with spiral vents
- Mask made from soft and comfortable material which easily contours to the patient's face to minimize medication loss and reduce waste
- Mask's exhalation safety valve releases CO₂ when patient breathes out
- Chamber features one-way valve that opens to allow inhalation of aerosol and closes to prevent patient from exhaling into the chamber
- Drug delivery is further enhanced by the antistatic quality of the chamber, a feature that helps prevent medication from clinging to the chamber wall

Specifications

Chamber Volume 160cm³

Weight 65 grams

Mouthpiece Diameter 22mm

Dimensions 53mm x 150mm

Storage Temperature 14°–113°F (-10°–45°C)

Storage Relative Humidity 30%–80%

Single Patient Use Yes



Adult size mask
also available

RES195ML

