

## Experience the Transcend Micro™ Difference

The smallest, lightest CPAP on the market redefines convenience with its exceptional portability, user-friendly interface and compatibility with any CPAP mask\*. Leave your bulky CPAP at home and enjoy your new favorite travel companion.

|                                        | Transcend Micro™<br>Auto CPAP | ResMed AirMini™<br>Portable CPAP | ResMed AirSense™<br>10 AutoSet™ CPAP<br>Card-to-Cloud | ResMed AirSense™<br>11 AutoSet™ CPAP |
|----------------------------------------|-------------------------------|----------------------------------|-------------------------------------------------------|--------------------------------------|
| <b>Dimensions</b>                      | <b>3.6"x3.6"x2.4"</b>         | 5.4"x3.3"x2.1"                   | 10.04"x4.57"x5.91"                                    | 10.21"x3.72"x5.45"                   |
| <b>Weight</b>                          | <b>0.48 lbs</b>               | 0.66 lbs                         | 2.75 lbs                                              | 2.5 lbs                              |
| <b>Footprint</b>                       | <b>12.96 sq inches</b>        | 17.82 sq inches                  | 45.88 sq inches                                       | 37.98 sq inches                      |
| <b>Volume</b>                          | <b>31.1 cubic inches</b>      | 37.4 cubic inches                | 276.1 cubic inches                                    | 200.5 cubic inches                   |
| <b>Sound Volume/Level</b>              | <b>27 dBA</b>                 | 29 dBA                           | 26.6 dBA                                              | 27 dBA                               |
| <b>Mask Compatibility*</b>             | <b>Yes</b>                    | No                               | Yes                                                   | Yes                                  |
| <b>Portable Battery Option</b>         | <b>Yes</b>                    | No                               | No                                                    | No                                   |
| <b>Adapts to User Breathing Cycle</b>  | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Sleep Tracking App</b>              | <b>MySleepDash™</b>           | AirMini™ by ResMed               | myAir™ (SD Card Reader Required)                      | myAir™                               |
| <b>Bluetooth Connectivity</b>          | <b>Yes</b>                    | Yes                              | No                                                    | Yes                                  |
| <b>Waterless Humidification</b>        | <b>Yes</b>                    | Yes                              | No (Fillable Water Tub)                               | No (Fillable Water Tub)              |
| <b>Muffler</b>                         | <b>Yes</b>                    | No                               | No                                                    | No                                   |
| <b>Dry Cycle</b>                       | <b>Yes</b>                    | No                               | No                                                    | No                                   |
| <b>Warranty</b>                        | <b>2 years</b>                | 2 years                          | 2 years                                               | 2 years                              |
| <b>Leak Compensation</b>               | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Ramp Time</b>                       | <b>0-45 minutes</b>           | 0-45 minutes                     | 0-45 minutes                                          | 0-45 minutes                         |
| <b>Auto Ramp</b>                       | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Expiratory Pressure Relief</b>      | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Operating Pressure Range</b>        | <b>4-20 cm H2O</b>            | 4-20 cm H2O                      | 4-20 cm H2O                                           | 4-20 cm H2O                          |
| <b>Max Operating Altitude</b>          | <b>8,000 feet</b>             | 8,500 feet                       | 8,500 feet                                            | 9,870 feet                           |
| <b>FAA-Certified for In-Flight Use</b> | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Travel Case Accessory</b>           | <b>Yes</b>                    | Yes                              | Yes                                                   | Yes                                  |
| <b>Operation Mode</b>                  | <b>CPAP/APAP</b>              | CPAP/APAP                        | CPAP/APAP                                             | CPAP/APAP                            |
| <b>Power</b>                           | <b>100-240V, 50-60Hz</b>      | 100-240V, 50-60Hz                | 100-240V, 50-60Hz                                     | 100-240V, 50-60Hz                    |
| <b>Release Date</b>                    | <b>July 2022</b>              | May 2017                         | August 2014                                           | October 2021                         |

\* Compatible with any CPAP mask with a standard 22mm connection