



Portable Suction Machine AC and battery power

Item#
SU100DC

Includes internal rechargeable 12-volt lithium battery, which can be used for up to 1 hour of operation per charge.

Suction Machine AC power

Item#
SU100AC

Both machines include

- Machine
- Carry Bag
- 800cc Suction Canister
- AC Adapter & Power Cord
- 6ft Suction Tubing
- 3.5in Suction Tube Connector
- 4.5in Suction Tube Connector
- Bacteria Filter

Both machines feature

- Effective removal of respiratory fluids
- Durability
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty*
- HCPCS# E0600

*The SU100DC internal rechargeable battery is covered by a 90-day warranty

Specifications

Specialty	Small/Lightweight
Carry Case	Yes
Dimensions (L x H x W in)	12.83" x 6.37" x 8.12"
Weight (lbs.)	SU100DC: 5.29 lbs / SU100AC: 4.08 lbs
Vacuum Range (mm Hg)	50-560
LPM	25
Sound Level (dBA)	60 dBA
Canister Type w/ purchase of machine	Disposable
Uses Suction Canister Size (cc)	800 /1000/1200
Machine Warranty	2 years
Can be operated inside carry case?	Yes

Accessories and replacement parts

RES024-MKIT	Replacement Filter Kit (connectors and filter)
RES026S-SS	Replacement Suction Kit (canister, tubing, connectors, and filter)
SU100DC-Adapt	Replacement Power Adapter (without cord)
SU100DC-Cord	Replacement Power Cord (without adapter)
SU100DC-BAG	Carrying/Storage Bag
SU100DC-BATTERY	Replacement internal rechargeable battery for SU100DC



SCAN TO LEARN MORE



Suction Kit for SU100DC

Item# **RES026S-SS**

RES023AW Suction Canister
RES024 Suction Tubing Connector 3.5in
RES024M Suction Tubing Connector 4.5in
BF1438 Bacteria Filter
RES025 Suction Tubing 1/4in x 6ft

HCPCS#
CANISTER - A7000
TUBING - A7002



Suction Canister

Item# **RES023AW**

800cc • Single Use
48/case
HCPCS# A7000

Item# **RES025AW1200**

1200cc • Single Use
48/case


Float valve shutoff



Suction Filter Component Kit

Item# **RES024XL-3KIT**

RES024XL Suction Tubing Connector 16in
BF3812 Filter
RES028 Elbow Connector



Suction Filter Component Kit for SU100DC

Item# **RES024-MKIT**

RES024 Suction Tubing Connector 3.5in
RES024M Suction Tubing Connector 4.5in
BF1438S Bacteria Filter



Suction Kit

Item# **RES026S-E**

RES023AW Suction Canister
RES024M Suction Tubing Connector 4.5in
RES025 Suction Tubing 6ft
RES028 Elbow Connector
BF3812 Bacteria Filter

HCPCS#
CANISTER - A7000
TUBING - A7002



To order or for more information

P 312.997.9980 E sales@sunsethcs.com W SUNSETHCS.COM

