

# Product Catalog



Contract PP-NS-1571  
**Closed Suction  
Catheters**



Contract PP-NS-1741  
**Respiratory Disposables**

Contract PP-NS-1762  
**Resuscitation Devices**



Contract PP-NS-1753  
**Medication Delivery**



Contract PP-OR-2217  
**Laryngeal Masks**



*Sales Representatives*



**Owen Ruth**  
312-533-2451  
oruth@sunsethcs.com



**Christian Straussner**  
872-282-0746  
cstraussner@sunsethcs.com




**Cindy Geighes**  
312-533-2478  
cgeighes@sunsethcs.com



**Egle Sventnickaite**  
872-258-1654  
esventnickaite@sunsethcs.com



**Hunter Baudino**  
262-788-1855  
hbaudino@sunsethcs.com



*Director, Sales Operations*  
**Mark Lipson**  
312-533-2465  
mlipson@sunsethcs.com

*Contract Manager, Strategic Accounts*  
**Gina Rosenberg**  
914-837-5215  
grosenberg@sunsethcs.com

*RRT Clinical Specialist*  
**Carlos Mejia** MS-QSM, RRT, RPGST, LSSB  
cmejia@sunsethcs.com

*Director, Strategic Accounts*  
**Lisa Iannazzone** RN, BSN  
636-222-3134  
lisa@sunsethcs.com

*Sales and Marketing Manager*  
**Jackson Hanft**  
971-416-6695  
jackson@sunsethcs.com

*Registered RT*  
**Krystle Bodine** BS, RRT  
kbodine@sunsethcs.com

**Corporate Office**

Phone 872-275-3210

Fax 312-997-9985

180 N Michigan Ave Ste 2000  
Chicago, IL 60601

**Central Returns/Distribution Center**

279 Madsen Dr Ste 101  
Bloomington, IL 60108

**West Returns/Distribution Center**

6750 W Van Buren St  
Phoenix, AZ 85043

# CONTENT

## Contract PP-NS-1741

### Respiratory Disposables Page

Lynx Venturi Oxygen Delivery	4
Spirometers	5
Drainage Bags	5
Oxygen Filters	5
Vent Circuits and Filters	6
Aerosol Tubing	6
Exhalation Ports	6
Oxygen Cannula and Accessories	7-8
Oxygen Masks	9-10
Air Entrainment Masks	10
Oxygen Supply Tubing	11
Firesafe	11
Oxygen Connectors and Exhalation Ports	12
Oxygen Accessories	13
Berman And Guedel Airways	14
HMEs	15
Trach Masks and Adapters	15
Suction Catheter Kits	16
Suction Catheters	17
Trach Ties	17
Connectors	17

## Contract PP-NS-1753

### Medication Delivery Page

Inhaler Chamber	18
Compressor Nebulizers	18
Peak Flow Meter	18
Adapters	18
Nebulizer Kits	19
Aerosol Masks	19

## Contract PP-NS-1571

### Closed Suction Catheters Page

Standard Switch	20
Deluxe Switch	21-22

## Contract PP-NS-1762

### Resuscitation Devices Page

Resuscitators	23
Manometers	24
PEEP Valves	24
Reservoir Bags	25
Pocket Resuscitator Set	25
Breathing Masks	25

## Contract PP-OR-2217

### Laryngeal Masks Page

Silicone Disposable	26
Silicone Reusable	26
Flexible Silicone Disposable	26
Flexible Silicone Reusable	26
PVC Disposable	27
Flexible PVC Disposable	27

# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Lynx Venturi Oxygen Delivery


### Low Concentration

Item# **RES2400**

Each

Adapter Only **RES240LC**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO<sub>2</sub> delivery from 24-60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO<sub>2</sub> washout



### Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

### Do no harm

Quicker de-escalation of O<sub>2</sub> levels and prevention of excessive CO<sub>2</sub> in the patient.

### Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

### Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

### High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO<sub>2</sub> in the EMS ambulance or ED to prevent CO<sub>2</sub> buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub>) allows early intervention and better outcomes
- Excellent CO<sub>2</sub> washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Volumetric Incentive Spirometer

12/case

Item# **RES055**  
2500ml

**RES056**  
5000ml

**RES055P**  
Pediatric



SCAN FOR VIDEO



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design

## Drainage Bag with Y-Adapter

Item# **RES071**  
750cc  
50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

## Drainage Bag with T-Adapter

Item# **RES070**  
750cc  
50/case



## Three Ball Spirometer

Item# **RES057-10**  
10/case



- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient

- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

## Final Bacteria Filter

Item# **BF516**  
Final Bacteria Filter  
1/4in Barbed  
10/pack  
OEM# P0750X030



## Final Bacteria Filter — 3/16in

Item# **BF316**  
3/16 Stepped Barb and  
1/4in Straight Barb  
on Both Sides  
10/pack



SPIROMETERS, DRAINAGE BAGS, AND BACTERIA FILTERS


# RESPIRATORY DISPOSABLES Contract PP-NS-1741

VENT CIRCUITS, AEROSOL TUBING, EXHALATION PORTS, AND FILTERS

## Non-Heated Breathing Circuits

For use with Respironics Trilogy  
22mm cuffs

Item#  
**VC1001**  
Active Adult,  
No Water Trap  
**30/case**  
OEM#  
1073223  
1073232



## Passive Non-Heated Breathing Circuits

For use with Respironics Trilogy

Item#  
**VC1100**  
With fixed exhalation port  
**Each**  
OEM# 1069467 & 1069468




## Aerosol Tubing

Item#  
**TUB100-6**  
Corrugated Aerosol Tubing  
100ft in 6in segments  
**Each**




Item#  
**TUB100-6B**  
Blue Corrugated Aerosol Tubing  
100ft in 6in segments  
**Each**



## Passive Heated Breathing Circuit

Item#  
**VC1114**  
For use with  
MR850 Humidifier  
**10/case**



## Single Use Exhalation Port

Item# **CAP011**  
**10/case**



## Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item#  
**BF2200**  
22mm Bacteria  
Disk Filter  
**5/pack**



Attach between the machine and tube for additional protection against bacteria and allergens

### Filter Specifications

Bacterial Filtration Efficiency: (BFE) 99.99+%  
Viral Filtration Efficiency: (VFE) 99.9+%  
NaCl Particulate Filtration Efficiency: 96+%  
Flow Resistance: Approximately 8.8 mmH2O at 30 LPM  
Internal Volume: 66.3 mL without fittings / 78.3 mL with fittings  
Connections:  
• Patient End: 22 mm OD / 15 mm ID  
• Machine End: 22 mm ID

# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Humidoflo High Flow Nasal Cannula

10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>



- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respirationics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerated due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O<sub>2</sub> concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

## Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty
<b>RES1100</b>	None	Standard	50
<b>RES1104</b>	4 ft	Standard	50
<b>RES1107</b>	7 ft	Standard	50
<b>RES1107U</b>	7 ft	Universal	50
<b>RES1110</b>	10 ft	Standard	50
<b>RES1114N</b>	14 ft	Standard	50
<b>RES1114U</b>	14 ft	Universal	50
<b>RES1115</b>	15 ft	Standard	50
<b>RES1115U</b>	15 ft	Universal	50
<b>RES1125</b>	25 ft	Standard	25
<b>RES1150U</b>	50 ft	Universal	20



## Soft Cannula

Item#	Tubing Length	Connector	Qty
<b>RES1100S</b>	None	Standard	50
<b>RES1107S</b>	7 ft	Standard	50
<b>RES1107SU</b>	7 ft	Universal	50
<b>RES1114S</b>	14 ft	Standard	50
<b>RES1114SU</b>	14 ft	Universal	50
<b>RES1115S</b>	15 ft	Standard	50
<b>RES1125S</b>	25 ft	Standard	25



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
<b>RES1107I</b>	Infant	7 ft	Standard	Standard	No
<b>RES1107P</b>	Pediatric	7 ft	Standard	Standard	No
<b>RES1107PEC</b>	Pediatric	7 ft	Standard	Standard	Yes
<b>RES1107PECU</b>	Pediatric	7 ft	Universal	Standard	Yes
<b>RES1107SI</b>	Infant	7 ft	Standard	Soft	No
<b>RES1107SNE</b>	Neonate	7 ft	Standard	Soft	No
<b>RES1107SP</b>	Pediatric	7 ft	Standard	Soft	No
<b>RES1107SPU</b>	Pediatric	7 ft	Universal	Soft	No

## Adult Cannula with Ear Cushions

Item#	Tubing Length	Connector	Touch	Tip	QTY
<b>RES1100EC</b>	None	Standard	Standard	Standard	50
<b>RES1107EC</b>	7 ft	Standard	Standard	Standard	50
<b>RES1107ECU</b>	7 ft	Universal	Standard	Standard	50
<b>RES1114EC</b>	14 ft	Standard	Standard	Standard	50
<b>RES1115EC</b>	15 ft	Standard	Standard	Standard	50
<b>RES1125EC</b>	25 ft	Standard	Standard	Standard	25

## Ear Cushions for Cannula

Item# **RES020**  
Cushions Only  
**50 Pairs/case**



- Add cushions to existing cannula



## Secure Cannula Stickers

Item# **RES080**  
**25 pairs/case**

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs





# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Simple O<sub>2</sub> Mask

With 7ft Safety Tubing

Item# **RES2100**

Adult Size

50/case

HCPCS# A4620

- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector
<b>RES2100NT</b>	Adult	None	Standard
<b>RES2100U</b>	Adult	7 ft	Universal
<b>RES2203P</b>	Pediatric	None	Standard
<b>RES2100P</b>	Pediatric	7 ft	Standard
<b>RES2203A</b>	Adult	None	Standard
<b>RES2203A7</b>	Adult	7 ft	Standard
<b>RES2203I7</b>	Infant	7 ft	Standard
<b>RES2203P7</b>	Pediatric	7 ft	Standard

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length
<b>RES2202A</b>	Adult	None
<b>RES2202A7</b>	Adult	7 ft
<b>RES2202I7</b>	Infant	7 ft



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent
<b>RES2102</b>	Adult	7 ft	Standard	No
<b>RES2102I</b>	Infant	7 ft	Standard	No
<b>RES2102P</b>	Pediatric	7 ft	Standard	No
<b>RES2102PS</b>	Pediatric	7 ft	Standard	Yes
<b>RES2102PSM</b>	Pediatric	None	Standard	Yes
<b>RES2102PSU</b>	Pediatric	7ft	Universal	Yes
<b>RES2102PU</b>	Pediatric	7 ft	Universal	No
<b>RES2102S</b>	Adult	7 ft	Standard	Yes
<b>RES2102SM</b>	Adult	None	Standard	Yes
<b>RES2102U</b>	Adult	7 ft	Universal	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



## Adult Air Entrainment Mask

50/case

Item#	Size	Kit	Tubing Length	Connector	Flex Tube
<b>RES2400A</b>	Adult	Mask Only	7 ft	Standard	Yes
<b>RES2400AK</b>	Adult	Yes	7 ft	Standard	Yes
<b>RES2400AKU</b>	Adult	Yes	7 ft	Universal	Yes
<b>RES2400AU</b>	Adult	Mask Only	7 ft	Universal	Yes
<b>RES2400P</b>	Pediatric	Mask Only	7 ft	Standard	Yes
<b>RES2400PK</b>	Pediatric	Yes	7 ft	Standard	Yes
<b>RES2400PKU</b>	Pediatric	Yes	7 ft	Universal	Yes
<b>RES2400PU</b>	Pediatric	Mask Only	7 ft	Universal	Yes



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Oxygen Supply Tubing

Item#	Length	Connector	Color	Qty
<b>RES3007</b>	7 ft	Standard	None	50
<b>RES3007G</b>	7 ft	Standard	Green	50
<b>RES3007U</b>	7 ft	Universal	None	50
<b>RES3015</b>	15 ft	Standard	None	50
<b>RES3025</b>	25 ft	Standard	None	25
<b>RES3025G</b>	25 ft	Standard	Green	25
<b>RES3025U</b>	25 ft	Universal	None	25
<b>RES3040G</b>	40 ft	Standard	Green	20
<b>RES3050G</b>	50 ft	Standard	Green	20


- 6-channel
- PHT/DEHP-free

### Universal Connector



## FireSafe Nozzle

Item#  
**RES010FND**  
10/pack



- A thermal fuse designed to stop the flow of oxygen in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire

# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Christmas Tree Connector

Item# **RES001**  
Green without Swivel  
**10/pack**



**RES001-30** 30/pack

Item# **RES002**  
Green with Swivel  
**10/pack**



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Tubing Connector — Non-Swivel

Item# **RES012**  
Fits 3/16in and 1/4in ID Tubing  
**10/pack**



- Barbed ends ensure secure, leak-free connections

## Tubing Reducer Connector

Item# **RES011**  
Fits 3/16in, 1/4in, 5/16in, 3/8in, and 7/16in ID Tubing  
**10/pack**



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Humidifier Tubing Adapter

Item# **RES016**  
22mm ID x 22mm OD with Lip  
**10/pack**



## Circuit Connectors

Item# **RES013**  
22mm OD x 15mm ID/ 22mm OD  
**10/pack**



Item# **RES014**  
22mm OD x 22mm OD With Lip  
**10/pack**



## Oxygen Tubing Connector — Swivel

Item# **RES018**  
Male/Male  
**10/pack**



**RES018-25** 25/pack

## Oxygen Tubing Y-Connector

Item# **RES017**  
Fits 3/16in and 1/4in ID Tubing (All 3 sides are the same)  
**10/pack**



- Connects two oxygen sources for higher flows

## Oxygen Bleed-In Adapter

Item# **RES019**  
22mm with 6mm Rotating Stem  
**10/pack**



## Single Use Exhalation Port

Item# **CAP011**  
**10/case**



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Standard Finger Pulse Oximeter

Item# **RES5100**  
Each

- 1 display mode
- Digital LED display
- SpO2, PR, and pulse bar



## Advanced Finger Pulse Oximeter

Item# **RES5101**  
Each

- 6 display modes
- Dual color OLED display
- SpO2, PR, pulse bar and waveform



## Wrenches

Item# **RES031**  
With Chain  
10/pack



## Water Traps

Item# **RES037**  
Push On End Caps  
Each

- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position



## Oxygen Liter Meter

Item# **RES034**  
0-10 Liter/Minute  
Each

- Oxygen pen
- Works with tanks and concentrators



## Humidifier Connectors

Item# **HUM008**  
L-Connector 10in  
10/pk

**HUM008L**  
L-Connector 21in  
10/pk

**HUM008S**  
Straight Connector 9in  
10/pk



# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Color Coded Dual Channel Berman Oral Airway

50/pack

Item#	Size	Color
<b>A0A40BC-50</b>	40mm	Pink
<b>A0A50BC-50</b>	50mm	Blue
<b>A0A60BC-50</b>	60mm	Black
<b>A0A70BC-50</b>	70mm	White
<b>A0A80BC-50</b>	80mm	Green
<b>A0A90BC-50</b>	90mm	Yellow
<b>A0A100BC-50</b>	100mm	Red
<b>A0A110BC-50</b>	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Oropharyngeal Guedel Airway

10/pack

Item#	Size	Color
<b>A0P500GC</b>	50mm	Blue
<b>A0P601GC</b>	60mm	Black
<b>A0P702GC</b>	70mm	White
<b>A0P803GC</b>	80mm	Green
<b>A0P904GC</b>	90mm	Yellow
<b>A0P1005GC</b>	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Flex HME

Item# **RES027FLEX**  
40/case



<b>TV Range (mL)</b>	300-1500
<b>Moisture Output (mg H2O/L)</b>	33 mg/L
<b>Weight (g)</b>	17
<b>Resistance (cm H2O)</b>	0.5 @ 30 LPM 1.34 @ 60 LPM
<b>Dead Space (mL)</b>	75
<b>Connection (mm)</b>	22M,15F

## Heat Moisture Exchangers

30/case



Item#	<b>RES027</b> T-HME	<b>RES027OP</b> T-HME with Oxygen Port	<b>RES027IL</b> Inline HME — Large Volume	<b>RES027NE</b> Inline HME — Neonatal	<b>RES027HV2</b> Foam HME — Adult/Child
<b>TV Range (mL)</b>	250-1000	250-1000	50-250	15-50	250-1000
<b>Moisture Output (mg H2O/L)</b>	30, Vt = 0.1L	29.50	31, Vt = 0.1L	30	29
<b>Weight (g)</b>	4.5	7	10	4.5	21
<b>Resistance (cm H2O)</b>	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
<b>Dead Space (mL)</b>	12.5	12	13	2.5	29
<b>Connection (mm)</b>	15F	15F 5.5mm Oxygen port	15M,22M/15F	15M,15F	22F/15M, 22M/15F

## Tracheostomy Mask with O<sub>2</sub> Connector

50/case

Item#  
**RES2130**  
Adult Regular Size

**RES2130L**  
Adult Large Size

**RES2230**  
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## Adapter for T-HME

Item#  
**RES027A**  
Each

**RES027A-30**  
30/pack



HMEs, TRACH MASKS, AND ADAPTERS

# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Suction Catheter Kits

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK305</b>
6fr	Y	DeLee	Pediatric	<b>TRATK306</b>
8fr	Y	DeLee	Pediatric	<b>TRATK308</b>
8fr	Y	Whistle	Adult	<b>TRATK108</b>
8fr	Chimney	Whistle	Adult	<b>TRATK108CV</b>
10fr	Y	Whistle	Adult	<b>TRATK110</b>
10fr	Chimney	Whistle	Adult	<b>TRATK110CV</b>
12fr	Y	Whistle	Adult	<b>TRATK112</b>
12fr	Chimney	Whistle	Adult	<b>TRATK112CV</b>
14fr	Y	Whistle	Adult	<b>TRATK114</b>
14fr	Chimney	Whistle	Adult	<b>TRATK114CV</b>
16fr	Y	Whistle	Adult	<b>TRATK116</b>
16fr	Chimney	Whistle	Adult	<b>TRATK116CV</b>

- Sterile
- Includes Catheter and Gloves



## Suction Catheter Kits with Cup

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK405</b>
6fr	Y	DeLee	Pediatric	<b>TRATK406</b>
8fr	Y	DeLee	Pediatric	<b>TRATK408</b>
8fr	Y	Whistle	Adult	<b>TRATK208</b>
8fr	Chimney	Whistle	Adult	<b>TRATK208CV</b>
10fr	Y	Whistle	Adult	<b>TRATK210</b>
10fr	Chimney	Whistle	Adult	<b>TRATK210CV</b>
12fr	Y	Whistle	Adult	<b>TRATK212</b>
12fr	Chimney	Whistle	Adult	<b>TRATK212CV</b>
14fr	Y	Whistle	Adult	<b>TRATK214</b>
14fr	Chimney	Whistle	Adult	<b>TRATK214CV</b>
16fr	Y	Whistle	Adult	<b>TRATK216</b>
16fr	Chimney	Whistle	Adult	<b>TRATK216CV</b>

- Sterile
- Includes Catheter and Gloves





# RESPIRATORY DISPOSABLES Contract PP-NS-1741

## Suction Catheters

Size	Valve Type	Tip	Intended Use	Packaged	Item#	Qty
5fr	Y	DeLee	Pediatric	Coiled	<b>TRASC305C</b>	50
5fr	Y	DeLee	Pediatric	Straight	<b>TRASC305S</b>	100
6fr	Y	DeLee	Pediatric	Coiled	<b>TRASC306C</b>	50
6fr	Y	DeLee	Pediatric	Straight	<b>TRASC306S</b>	100
8fr	Y	DeLee	Pediatric	Coiled	<b>TRASC308C</b>	50
8fr	Y	DeLee	Pediatric	Straight	<b>TRASC308S</b>	100
10fr	Y	Whistle	Adult	Coiled	<b>TRASC110C</b>	50
10fr	Y	Whistle	Adult	Straight	<b>TRASC110S</b>	100
12fr	Y	Whistle	Adult	Coiled	<b>TRASC112C</b>	50
12fr	Y	Whistle	Adult	Straight	<b>TRASC112S</b>	100
14fr	Y	Whistle	Adult	Coiled	<b>TRASC114C</b>	50
14fr	Y	Whistle	Adult	Straight	<b>TRASC114S</b>	100
16fr	Y	Whistle	Adult	Coiled	<b>TRASC116C</b>	50
16fr	Y	Whistle	Adult	Straight	<b>TRASC116S</b>	100



- Sterile
- Single patient use
- Latex-free

## Tracheostomy Ties with Fastener Tabs

**20/case**


One-piece

- Easy on

Item#

**RES029** Adult

**RES029P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton


Two-piece

- Easy off

Item#

**RES030** Adult

**RES030P** Pediatric



## Tracheostomy Swivel Elbow Connectors

Item# **RES021**

With suction port  
15mm ID x 15mm OD  
**10/pack**



Visit [SunsetHCS.com/AcuteCare](http://SunsetHCS.com/AcuteCare)

# MEDICATION DELIVERY Contract PP-NS-1753

## SperaChamber

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item# **RES195-10**  
10/case

- Inhalation valve
- Universal adapter
- Spiral vents

Item# **RES195-50**  
50/case

- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask  
10/case

Item#	Size
<b>RES195MS</b>	Small (Infant)
<b>RES195MM</b>	Medium (Pediatric)
<b>RES195ML</b>	Large (Adult)



Adult Mask

## Compressor Nebulizer

Item# **NEB100**  
Each

- Includes**
- Compressor
  - Disposable T-Piece Kit
  - Foam Filters



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

## Peak Flow Meter

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**  
Pediatric Mouthpiece for Peak Flow Meter  
10/case

## Valved T-Adapter, Aerosol, 22mm ID x 22mm OD

Item# **NEB9001T**  
10/case



## Valved T-Adapter — Adult

Item# **VEN001**  
22F-18F-22M/15F  
Each



# MEDICATION DELIVERY Contract PP-NS-1753

## Disposable Nebulizer Kit

Item#

**RES090**

Adult Size

**50/case**

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

**RES090A**

No Accessories

**50/case**

**RES090NM**

- 6cc jet nebulizer bottle
- 7ft supply tube

**50/case**

**RES092**

Pediatric Size

**50/case**

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



## Reusable Nebulizer Kit with T-Piece

Item# **RES093**

Reusable — Single Patient

**10/case**

HCPCS# A7005

Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube
- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°



## Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**

Disposable  
Single Patient

**50/case**

HCPCS# A7003

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



## Disposable Nebulizer Kit with T-Piece

Item#

**RES091**

Disposable  
Single Patient

**50/case**

**RES091NA**

Disposable  
Single Patient,  
No Aerosol

**50/case**

HCPCS# A7003



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

**Includes:**

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube

## Aerosol Mask

**50/case**

Item#

**RES2110**

Adult Size

**RES2210**

Pediatric Size


- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



# Standard Switch

Sterile and latex-free

ISO color coded for instant recognition of catheter sizes.



Once the airway is clear, withdraw the catheter completely into the plastic sheath. A black mark on the catheter indicates full withdrawal.

Smooth tubing surface and atraumatic rounded catheter tip for minimal tracheal wall injury.

With the unique isolation valve, 72 hour closed suction system features safe and simple access for mechanically ventilated patients for 72 hours of continuous use.

## 24 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr
VEN46011400M	14Fr, MDI Port
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

## 72 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN56061000A	10Fr
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port

## Neonatal/Pediatric

### With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000600	305mm, 6Fr
VEN66000700	305mm, 7Fr
VEN66000800	305mm, 8Fr
VEN66001000	405mm, 10Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000M	405mm, 10Fr, MDI Port



# Deluxe Switch (One-handed operation)

Sterile and latex-free

## Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

## Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

## MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

## Catheter Tip


Atraumatic rounded catheter tip for minimal tracheal wall injury.

## Isolation Valve

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

## Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



# Deluxe Switch (One-handed operation)

Sterile and latex-free

## 24 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN160160	16Fr
VEN160160M	16Fr, MDI Port

### With Wedge + Adapter

Item#	Description
VEN162100A	10Fr
VEN162100AM	10Fr, MDI Port
VEN162120	12Fr
VEN162120M	12Fr, MDI Port
VEN162140	14Fr
VEN162140M	14Fr, MDI Port
VEN162160	16Fr
VEN162160M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN161120	12Fr
VEN161140	14Fr
VEN161160	16Fr
VEN161160M	16Fr, MDI Port

### With Wedge + Adapter

Item#	Description
VEN163100A	10Fr
VEN163100AM	10Fr, MDI Port
VEN163120	12Fr
VEN163140	14Fr
VEN163160	16Fr

## 72 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN365160	16Fr
VEN365160M	16Fr, MDI Port

### With Wedge + Adapter

Item#	Description
VEN367100A	10Fr
VEN367100AM	10Fr, MDI Port
VEN367120	12Fr
VEN367120M	12Fr, MDI Port
VEN367140	14Fr
VEN367140M	14Fr, MDI Port
VEN367160	16Fr
VEN367160M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN366160	16Fr
VEN366160M	16Fr, MDI Port

### With Wedge + Adapter

Item#	Description
VEN368100A	10Fr
VEN368100AM	10Fr, MDI Port
VEN368120	12Fr
VEN368120M	12Fr, MDI Port
VEN368140	14Fr
VEN368140M	14Fr, MDI Port
VEN368160	16Fr
VEN368160M	16Fr, MDI Port

## Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port



## Adult Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8100</b>	No	None
<b>RES8100H</b>	Yes	None
<b>RES8101</b>	Yes	#4
<b>RES8102</b>	Yes	#5



## Child Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8110H</b>	Yes	None
<b>RES8111</b>	Yes	#2
<b>RES8112</b>	Yes	#3



## Infant Resuscitator

12/case

Item#	Pop-Off	Mask
<b>RES8120H</b>	Yes	None
<b>RES8121</b>	Yes	#0
<b>RES8122</b>	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

# RESUSCITATION DEVICES Contract PP-NS-1762

## Disposable Pressure Manometer

Item# **RES8201**  
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

## Disposable Pressure Manometer, L-Shaped

Item# **RES8201L**  
25/case



## Reusable PEEP Valve, Orange

2-10cm H2O  
10/case

Item#  
**RES8203R**  
15mm ID



**RES8203RL**  
30mm ID

## Reusable PEEP Valve, Blue

5-20cm H2O  
10/case

Item#  
**RES8204B**  
15mm ID



**RES8204BL**  
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

## Disposable PEEP Valve, Orange

2.5-10cm H2O  
10/case

Item#  
**RES8205R**  
15mm ID



**RES8205RL**  
30mm ID

## Disposable PEEP Valve, Blue

5-20cm H2O  
10/case

Item#  
**RES8205B**  
15mm ID



**RES8205BL**  
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel



# RESUSCITATION DEVICES Contract PP-NS-1762

## Reservoir Bag, PVC

25/case

Item#	Size
<b>RES8300P</b>	0.5L
<b>RES8301P</b>	1L
<b>RES8302P</b>	2L



## Reservoir Bag, Silicone

25/case

Item#	Size
<b>RES8300S</b>	0.5L
<b>RES8301S</b>	1L
<b>RES8302S</b>	2L



Suitable for reusable resuscitator.

## CPR Pocket Resuscitator Set

12/case

Item# **RES8200**

- High transparency for better visualization
- Non-reusable, single use only
- 100% Medical level PVC material
- Latex-free components
- Economical choice
- Blow molding



## Breathing Mask

20/case

Item#	Size	Hooks
<b>SBC01N</b>	Neonatal (#0)	No
<b>SBC01I</b>	Infant (#1)	No
<b>SBC01P</b>	Pediatric (#2)	No
<b>SBC01S</b>	Adult Small (#3)	No
<b>SBC01M</b>	Adult Medium (#4)	No
<b>SBC01L</b>	Adult Large (#5)	No
<b>SBC01XL</b>	Adult Extra Large (#6)	No
<b>SBC02N</b>	Neonatal (#0)	Yes
<b>SBC02I</b>	Infant (#1)	Yes
<b>SBC02P</b>	Pediatric (#2)	Yes
<b>SBC02S</b>	Adult Small (#3)	Yes
<b>SBC02M</b>	Adult Medium (#4)	Yes
<b>SBC02L</b>	Adult Large (#5)	Yes
<b>SBC02XL</b>	Adult Extra Large (#6)	Yes



# LARYNGEAL MASKS Contract PP-OR-2217

## Silicone Disposable Laryngeal Masks

Item#	Size
<b>ALM1001</b>	#1
<b>ALM1015</b>	#1.5
<b>ALM1002</b>	#2
<b>ALM1025</b>	#2.5
<b>ALM1003</b>	#3
<b>ALM1004</b>	#4
<b>ALM1005</b>	#5

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



## Silicone Reusable Laryngeal Masks

Item#	Size
<b>ALM1001R</b>	#1
<b>ALM1015R</b>	#1.5
<b>ALM1002R</b>	#2
<b>ALM1025R</b>	#2.5
<b>ALM1003R</b>	#3
<b>ALM1004R</b>	#4
<b>ALM1005R</b>	#5


- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



## Flexible Silicone Disposable Laryngeal Masks

Item#	Size
<b>ALM3001</b>	#1
<b>ALM3015</b>	#1.5
<b>ALM3002</b>	#2
<b>ALM3025</b>	#2.5
<b>ALM3003</b>	#3
<b>ALM3004</b>	#4
<b>ALM3005</b>	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



## Flexible Silicone Reusable Laryngeal Masks

Item#	Size
<b>ALM3001R</b>	#1
<b>ALM3015R</b>	#1.5
<b>ALM3002R</b>	#2
<b>ALM3025R</b>	#2.5
<b>ALM3003R</b>	#3
<b>ALM3004R</b>	#4
<b>ALM3005R</b>	#5



# LARYNGEAL MASKS Contract PP-OR-2217

## PVC Disposable Laryngeal Masks

Item#	Size
<b>ALM2001</b>	#1
<b>ALM2015</b>	#1.5
<b>ALM2002</b>	#2
<b>ALM2025</b>	#2.5
<b>ALM2003</b>	#3
<b>ALM2004</b>	#4
<b>ALM2005</b>	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



## Flexible PVC Disposable Laryngeal Masks

Item#	Size
<b>ALM4001</b>	#1
<b>ALM4015</b>	#1.5
<b>ALM4002</b>	#2
<b>ALM4025</b>	#2.5
<b>ALM4003</b>	#3
<b>ALM4004</b>	#4
<b>ALM4005</b>	#5



Size	Patient Weight and Age	
#1	Up to 5 kg	Pediatric
#1.5	5-10 kg	Pediatric
#2	10-20 kg	Pediatric
#2.5	20-30 kg	Pediatric
#3	30-50 kg	Pediatric
#4	50-70 kg	Adult
#5	70-100 kg	Adult

