



Micro™ Power Options

Item#
CAPSOM505001

Transcend Micro
PowerAway®
Battery



Extended Run Time


- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability

Item#
CAPSOM503056

Transcend Micro
Solar Battery
Charger



The durable, military-grade Transcend Portable Solar Battery Charger can recharge P4 and P8 batteries with help from the sun. It works even in low-light conditions. No matter what type of outdoor adventure you're embarking on next, remember to take this portable battery charger with you wherever you go.

*Not compatible with Transcend 365 miniCPAP devices, P10 Battery or Transcend 365 miniCPAP power supply.

Item#
CAPSOM503055

Transcend Micro
Plug Adapters
(AU/EU/UK)



Travel with ease and confidence with the Transcend Plug Adapters. Designed for use with the Transcend Micro CPAP machine and its PSA5 power supply, these adapters ensure you can continue your sleep therapy wherever your journey takes you.

- **Universal Compatibility:** The Transcend Plug Adapters are designed to work with the Transcend PSA5 power supply that comes with the Transcend Micro.
- **International Use:** These adapters are designed for international use, fitting outlets in the EU, AU and UK. They allow you to stay powered up wherever your travels take you.
- **All-in-One Solution:** This package combines multiple plugs into one compact unit, eliminating the need to worry about compatibility.
- **Package Includes:** Each package contains three Transcend Plug Adapters, offering great value.

Item#
CAPSOM505010

Transcend Micro
Power Supply



Replacement power supply set for Transcend Micro. This unit is a direct replacement for the original accessory, and it can quickly and easily charge the powerful P8 battery. This replacement set comes with the AC power brick and US wall plug.