


Low Concentration Oxygen Delivery

Item# **RES2400** Adapter Only Item# **RES240**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO_2 delivery from 24–60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO_2 washout


Six O_2 ports deliver concentrations of:



High Concentration Oxygen Delivery

Item# **RES2401** Adapter Only Item# **RES241**

- Allows rapid de-escalation of FiO_2 in the EMS ambulance or ED to prevent CO_2 buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO_2/FiO_2) allows early intervention and better outcomes
- Excellent CO_2 washout



Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

Do no harm

Quicker de-escalation of O_2 levels and prevention of excessive CO_2 in the patient.

Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.