








Heat/Moisture Exchangers (HMEs)



	 RES027 T-HME	 RES027OP T-HME with Oxygen Port	 RES027IL Inline HME — Large Volume	 RES027UOP HME with Oxygen Port	 RES027NE Inline HME — Neonatal	 RES027FLEX Flex HME
TV Range (ml)	250-1000	250-1000	50-250	50-1000	15-50	300-1500
Moisture Output (mg H ₂ O/L)	30, Vt = 0.1L	29.50	31, Vt = 0.1L	13.2 mg/L	30	33
Weight (g)	4.5	7	10	5	4.5	17
Resistance (cm H ₂ O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.5 @ 10 LPM	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (ml)	12.5	12	13	17	2.5	75
Connection (mm)	15F	15F, 5.5mm Oxygen port	15M, 22M/15F	15F	15M, 15F	22M, 15F