


Lynx Venturi

High Concentration Oxygen Delivery

Item# **RES2401** Adapter Only Item# **RES241HC**

- Allows rapid de-escalation of FiO_2 in the EMS ambulance or ED to prevent CO_2 buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia ($\text{PaO}_2/\text{FiO}_2$) allows early intervention and better outcomes
- Excellent CO_2 washout



Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Unprecedented Simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

Do No Harm

Quicker de-escalation of O_2 levels and prevention of excessive CO_2 in the patient.

Unmatched precision

Provides specific levels of oxygen with less than 3% variance. Comparable products vary by as much as 30%.

Unparalleled customer service

Sunset's industry renowned customer service combines big capabilities with highly customized, in-house sales support.