



# ASC

## Product Catalog



**Contracted Supplier**



**vizient.**  
Awarded Supplier





Sales Representatives



Owen Ruth
312-533-2451
oruth@sunsethcs.com



Christian Straussner
872-282-0746
cstraussner@sunsethcs.com




Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com



General Manager
Mark Lipson
312-533-2465
mlipson@sunsethcs.com

Contract Manager, Strategic Accounts
Gina Rosenberg
914-837-5215
grosenberg@sunsethcs.com

RRT Clinical Specialist
Carlos Mejia MS-QSM, RRT, RPGST, LSSB
cmejia@sunsethcs.com

Director, Strategic Accounts
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com

Sales and Marketing Manager
Jackson Hanft
971-416-6695
jackson@sunsethcs.com

RRT Clinical Specialist
Krystle Bodine BS, RRT
kbodine@sunsethcs.com

Corporate Office

Phone 872-275-3210
Fax 312-997-9985
141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix, AZ 85043

# CONTENT

## Medication Delivery

### Page

Nebulizer Kits	4
Aerosol Masks	5
Inhaler Chambers	5

## Respiratory Disposables

### Page

Oxygen Cannula	6-8
Non-Rebreather & Partial Rebreather Masks	9
Suction Catheter Kits	10
Suction Catheters	11
Suction Canisters and Yankauers	11

## Resuscitation Devices

### Page

Resuscitators	12
Reservoir Bags	13
Breathing Masks	13

# MEDICATION DELIVERY

## Disposable Nebulizer Kit

Item# **RES090** V P

Adult Size

**50/case**

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube



Item# **RES090A** V P

No Accessories

**50/case**

Item# **RES090NM** P

- 6cc jet nebulizer bottle
- 7ft supply tube

**50/case**

Item# **RES092** V P

Pediatric Size

**50/case**

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube

## Disposable Nebulizer Kit with T-Piece V P

Item# **RES091**

Disposable  
Single Patient  
**50/case**

Item# **RES091NA**

Disposable  
Single Patient,  
No Aerosol  
**50/case**

HCPCS# A7003



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

- Includes:**
- 6cc jet nebulizer bottle
  - T-piece with mouthpiece and corrugated reservoir tube
  - 7ft supply tube

## Reusable Nebulizer Kit with T-Piece P

Item# **RES093**

Reusable — Single Patient

**10/case**

HCPCS# A7005



Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube
- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

## Disposable Nebulizer Kit with Angled Mouthpiece V P

Item# **RES096**

Disposable  
Single Patient  
**50/case**

HCPCS# A7003



Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube

# MEDICATION DELIVERY

**V** Contract MS1097

**P** Contract PP-NS-1753

## Aerosol Mask

50/case

Item#

**RES2110** **V** **P**

Adult Size

**RES2110F** **P**

Adult Size w/ Flex Tube

**RES2210** **V** **P**

Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

## SperaChamber **V** **P**

Anti-Static Valved Inhaler Chamber for Metered Dose Inhaler

Item#

**RES195-10**

10/case

**RES195-50**

50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



## SperaChamber with Mask **V** **P**

10/case

Item#

Size

**RES195MS** Small (Infant)

**RES195MM** Medium (Pediatric)

**RES195ML** Large (Adult)



Adult Mask

AEROSOL MASK AND INHALER CHAMBER

# RESPIRATORY DISPOSABLES

## Humidoflo High Flow Nasal Cannula V P

10/case

Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>




- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respirationics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerated due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O2 concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

## Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty	<span>V</span> <span>P</span>
<b>RES1100</b>	None	Standard	50	<span>V</span> <span>P</span>
<b>RES1100N</b>	None	Standard	50	<span>V</span>
<b>RES1104</b>	4 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1107</b>	7 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1107N</b>	7 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1107NU</b>	7 ft	Universal	50	<span>V</span>
<b>RES1107U</b>	7 ft	Universal	50	<span>V</span> <span>P</span>
<b>RES1110</b>	10 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1114N</b>	14 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1114U</b>	14 ft	Universal	50	<span>V</span> <span>P</span>
<b>RES1115</b>	15 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1115U</b>	15 ft	Universal	50	<span>V</span> <span>P</span>
<b>RES1125</b>	25 ft	Standard	25	<span>V</span> <span>P</span>
<b>RES1135</b>	35 ft	Standard	20	<span>V</span>
<b>RES1150</b>	50 ft	Standard	20	<span>V</span> <span>P</span>
<b>RES1150U</b>	50 ft	Universal	20	<span>V</span> <span>P</span>



## Soft Cannula

Item#	Tubing Length	Connector	Qty	<span>V</span> <span>P</span>
<b>RES1100S</b>	None	Standard	50	<span>V</span> <span>P</span>
<b>RES1107S</b>	7 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1107SU</b>	7 ft	Universal	50	<span>V</span> <span>P</span>
<b>RES1114S</b>	14 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1114SU</b>	14 ft	Universal	50	<span>V</span> <span>P</span>
<b>RES1115S</b>	15 ft	Standard	50	<span>V</span> <span>P</span>
<b>RES1125S</b>	25 ft	Standard	25	<span>V</span> <span>P</span>
<b>RES1125SU</b>	25 ft	Universal	25	<span>V</span>



# RESPIRATORY DISPOSABLES

V Contract MS1087

P Contract PP-NS-1741

## Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions	
RES1100P	Pediatric	None	Standard	Standard	No	V
RES1107I	Infant	7 ft	Standard	Standard	No	V P
RES1107NE	Neonate	7 ft	Standard	Standard	No	V
RES1107P	Pediatric	7 ft	Standard	Standard	No	V P
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes	V P
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes	V P
RES1107SI	Infant	7 ft	Standard	Soft	No	V P
RES1107SNE	Neonate	7 ft	Standard	Soft	No	V P
RES1107SP	Pediatric	7 ft	Standard	Soft	No	V P
RES1107SPU	Pediatric	7 ft	Universal	Soft	No	V P
RES1114P	Pediatric	14 ft	Standard	Standard	No	V
RES1114PU	Pediatric	14 ft	Universal	Standard	No	V
RES1114SP	Pediatric	14 ft	Standard	Soft	No	V
RES1114SPU	Pediatric	14 ft	Universal	Soft	No	V

## Adult Cannula with Ear Cushions

Item#	Tubing Length	Connector	Qty	
RES1100EC	None	Standard	50	V P
RES1104EC	4 ft	Standard	50	V
RES1107EC	7 ft	Standard	50	V P
RES1107ECU	7 ft	Universal	50	V P
RES1114EC	14 ft	Standard	50	V P
RES1115EC	15 ft	Standard	50	V P
RES1125EC	25 ft	Standard	25	V P

## Flared Cannula V

Item#	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1114FU	14 ft	Universal	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

## High Flow Soft Cannula V

25/case

Item#	Tubing Length
RES1104SHF	4ft
RES1107SHF	7ft
RES1115SHF	15ft
RES1125SHF	25ft



- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

Visit [SunsetHCS.com/AcuteCare](http://SunsetHCS.com/AcuteCare)




# RESPIRATORY DISPOSABLES

## 7 ft Cannula With CO<sub>2</sub> Sampling Line

50/case

Item#	Size	Connector	Luer	Filter
<b>RES8007</b>	Adult	Standard	Male	No
<b>RES8007F</b>	Adult	Standard	Male	Yes
<b>RES8007I</b>	Infant	Standard	Male	No
<b>RES8007IF</b>	Infant	Standard	Male	Yes
<b>RES8007IU</b>	Infant	Universal	Male	No
<b>RES8007P</b>	Pediatric	Standard	Male	No
<b>RES8007PF</b>	Pediatric	Standard	Male	Yes
<b>RES8007PU</b>	Pediatric	Universal	Male	No
<b>RES8007U</b>	Adult	Universal	Male	No
<b>RES8107</b>	Adult	Standard	Female	No
<b>RES8107F</b>	Adult	Standard	Female	Yes
<b>RES8107I</b>	Infant	Standard	Female	No
<b>RES8107IF</b>	Infant	Standard	Female	Yes
<b>RES8107P</b>	Pediatric	Standard	Female	No
<b>RES8107PF</b>	Pediatric	Standard	Female	Yes



# RESPIRATORY DISPOSABLES

**V** Contract MS1087

**P** Contract PP-NS-1741

NON-REBREATHER MASKS AND PARTIAL REBREATHER MASKS

## Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent	
<b>RES2102</b>	Adult	7 ft	Standard	No	<b>V</b> <b>P</b>
<b>RES2102I</b>	Infant	7 ft	Standard	No	<b>V</b> <b>P</b>
<b>RES2102IM</b>	Infant	7 ft	Standard	No	<b>V</b>
<b>RES2102P</b>	Pediatric	7 ft	Standard	No	<b>V</b> <b>P</b>
<b>RES2102PS</b>	Pediatric	7 ft	Standard	Yes	<b>V</b> <b>P</b>
<b>RES2102PSM</b>	Pediatric	None	Standard	Yes	<b>V</b> <b>P</b>
<b>RES2102PSU</b>	Pediatric	7ft	Universal	Yes	<b>V</b> <b>P</b>
<b>RES2102PU</b>	Pediatric	7 ft	Universal	No	<b>V</b> <b>P</b>
<b>RES2102S</b>	Adult	7 ft	Standard	Yes	<b>V</b> <b>P</b>
<b>RES2102SM</b>	Adult	None	Standard	Yes	<b>V</b> <b>P</b>
<b>RES2102U</b>	Adult	7 ft	Universal	No	<b>V</b> <b>P</b>

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



## High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	
<b>RES2202A</b>	Adult	None	<b>V</b> <b>P</b>
<b>RES2202A7</b>	Adult	7 ft	<b>V</b>
<b>RES2202I</b>	Infant	None	<b>V</b> <b>P</b>
<b>RES2202I7</b>	Infant	7 ft	<b>V</b> <b>P</b>
<b>RES2202P</b>	Pediatric	None	<b>V</b>
<b>RES2202P7</b>	Pediatric	7 ft	<b>V</b>

## Simple O<sub>2</sub> Mask **V** **P**

With 7ft Safety Tubing  
Elongated

Item# **RES2100**

Adult Size

50/case

- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



Visit [SunsetHCS.com/AcuteCare](http://SunsetHCS.com/AcuteCare)



# RESPIRATORY DISPOSABLES

SUCTION CATHETER KITS

## Suction Catheter Kits P

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK305</b>
6fr	Y	DeLee	Pediatric	<b>TRATK306</b>
8fr	Y	DeLee	Pediatric	<b>TRATK308</b>
8fr	Y	Whistle	Adult	<b>TRATK108</b>
8fr	Chimney	Whistle	Adult	<b>TRATK108CV</b>
10fr	Y	Whistle	Adult	<b>TRATK110</b>
10fr	Chimney	Whistle	Adult	<b>TRATK110CV</b>
12fr	Y	Whistle	Adult	<b>TRATK112</b>
12fr	Chimney	Whistle	Adult	<b>TRATK112CV</b>
14fr	Y	Whistle	Adult	<b>TRATK114</b>
14fr	Chimney	Whistle	Adult	<b>TRATK114CV</b>
16fr	Y	Whistle	Adult	<b>TRATK116</b>
16fr	Chimney	Whistle	Adult	<b>TRATK116CV</b>

- Sterile
- Includes Catheter and Gloves



## Suction Catheter Kits with Cup P

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK405</b>
6fr	Y	DeLee	Pediatric	<b>TRATK406</b>
8fr	Y	DeLee	Pediatric	<b>TRATK408</b>
8fr	Y	Whistle	Adult	<b>TRATK208</b>
8fr	Chimney	Whistle	Adult	<b>TRATK208CV</b>
10fr	Y	Whistle	Adult	<b>TRATK210</b>
10fr	Chimney	Whistle	Adult	<b>TRATK210CV</b>
12fr	Y	Whistle	Adult	<b>TRATK212</b>
12fr	Chimney	Whistle	Adult	<b>TRATK212CV</b>
14fr	Y	Whistle	Adult	<b>TRATK214</b>
14fr	Chimney	Whistle	Adult	<b>TRATK214CV</b>
16fr	Y	Whistle	Adult	<b>TRATK216</b>
16fr	Chimney	Whistle	Adult	<b>TRATK216CV</b>

- Sterile
- Includes Catheter and Gloves



# RESPIRATORY DISPOSABLES

## Suction Catheters **P**

Size	Valve Type	Tip	Intended Use	Packaged	Item#	Qty
5fr	Y	DeLee	Pediatric	Coiled	<b>TRASC305C</b>	50
5fr	Y	DeLee	Pediatric	Straight	<b>TRASC305S</b>	100
6fr	Y	DeLee	Pediatric	Coiled	<b>TRASC306C</b>	50
6fr	Y	DeLee	Pediatric	Straight	<b>TRASC306S</b>	100
8fr	Y	DeLee	Pediatric	Coiled	<b>TRASC308C</b>	50
8fr	Y	DeLee	Pediatric	Straight	<b>TRASC308S</b>	100
10fr	Y	Whistle	Adult	Coiled	<b>TRASC110C</b>	50
10fr	Y	Whistle	Adult	Straight	<b>TRASC110S</b>	100
12fr	Y	Whistle	Adult	Coiled	<b>TRASC112C</b>	50
12fr	Y	Whistle	Adult	Straight	<b>TRASC112S</b>	100
14fr	Y	Whistle	Adult	Coiled	<b>TRASC114C</b>	50
14fr	Y	Whistle	Adult	Straight	<b>TRASC114S</b>	100
16fr	Y	Whistle	Adult	Coiled	<b>TRASC116C</b>	50
16fr	Y	Whistle	Adult	Straight	<b>TRASC116S</b>	100



- Sterile
- Single patient use
- Latex-free

## Suction Canister

Item# **RES023AW**  
 800cc • Single Use  
**48/case**  
 HCPCS# A7000  
 • Float valve shutoff



## Bulb Tip Yankauers

Item# **TRAA001S**  
 Vented • **50/case**  
 Item# **TRAA002S**  
 Non-Vented • **50/case**  
 HCPCS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber

SUCTION CATHETERS, CANISTERS AND YANKAUSERS

# RESUSCITATION DEVICES

## Adult Resuscitator **P**

12/case

Item#	Pop-Off	Mask
<b>RES8100</b>	No	None
<b>RES8100H</b>	Yes	None
<b>RES8101</b>	Yes	#4
<b>RES8102</b>	Yes	#5



## Child Resuscitator **P**

12/case

Item#	Pop-Off	Mask
<b>RES8110H</b>	Yes	None
<b>RES8111</b>	Yes	#2
<b>RES8112</b>	Yes	#3



## Infant Resuscitator **P**

12/case

Item#	Pop-Off	Mask
<b>RES8120H</b>	Yes	None
<b>RES8121</b>	Yes	#0
<b>RES8122</b>	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

# RESUSCITATION DEVICES

## Reservoir Bag, PVC **P**

25/case

Item#	Size
<b>RES8300P</b>	0.5L
<b>RES8301P</b>	1L
<b>RES8302P</b>	2L



## Reservoir Bag, Silicone **P**

25/case

Item#	Size
<b>RES8300S</b>	0.5L
<b>RES8301S</b>	1L
<b>RES8302S</b>	2L



Suitable for reusable resuscitator.

## Breathing Mask **P**

20/case

Item#	Size	Hooks
<b>SBC01N</b>	Neonatal (#0)	No
<b>SBC01I</b>	Infant (#1)	No
<b>SBC01P</b>	Pediatric (#2)	No
<b>SBC01S</b>	Adult Small (#3)	No
<b>SBC01M</b>	Adult Medium (#4)	No
<b>SBC01L</b>	Adult Large (#5)	No
<b>SBC01XL</b>	Adult Extra Large (#6)	No
<b>SBC02N</b>	Neonatal (#0)	Yes
<b>SBC02I</b>	Infant (#1)	Yes
<b>SBC02P</b>	Pediatric (#2)	Yes
<b>SBC02S</b>	Adult Small (#3)	Yes
<b>SBC02M</b>	Adult Medium (#4)	Yes
<b>SBC02L</b>	Adult Large (#5)	Yes
<b>SBC02XL</b>	Adult Extra Large (#6)	Yes



RESERVOIR BAGS AND BREATHING MASKS

