



Product Catalog

Contract MS1049

NIV



Contract MS1087

Respiratory Disposables



Contract MS1097

MDD



vizient[®]
Awarded Supplier



Sales Representatives



Owen Ruth
312-533-2451
oruth@sunsethcs.com



Christian Straussner
872-282-0746
cstraussner@sunsethcs.com




Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com



 *General Manager*
Mark Lipson
312-533-2465
mlipson@sunsethcs.com

 *Contract Manager, Strategic Accounts*
Gina Rosenberg
914-837-5215
grosenberg@sunsethcs.com

 *RRT Clinical Specialist*
Carlos Mejia MS-QSM, RRT, RPGST, LSSB
cmejia@sunsethcs.com

 *Director, Strategic Accounts*
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com

 *Sales and Marketing Manager*
Jackson Hanft
971-416-6695
jackson@sunsethcs.com

 *RRT Clinical Specialist*
Krystle Bodine BS, RRT
kbodine@sunsethcs.com

Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix, AZ 85043

CONTENT

Contract MS1049

NIV

Page

Adult NIV Masks	4-6
Pediatric NIV Masks	6-7
CPAP Tubing	7
Connectors	8
Headgear	8

Contract MS1087

Respiratory Disposables

Page

Lynx Venturi Oxygen Delivery	9
Filters	10
Connectors	10
Breathing Circuits	11-12
HME Filters	12
Connectors and Adapters	12
HMEs and HME Adapters	13
Tracheostomy Masks	13
Trach Ties	14
Berman and Guedel Airways	14
Pulse Oximeters	15
Firesafe	16
Spirometers	17
Aerosol Masks and Tubing	17
Oxygen Cannula and Cushions	18-20
Oxygen Masks	21
Oxygen Supply Tubing	22
Air Entrainment Masks	22

Contract MS1097

MDD

Page

Nebulizer Kits	23
Nebulizer Compressor	24
Peak Flow Meter	24
Valved T-Adapter	24
Inhaler Chamber	24

Full Face Mask

With Standard Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041S
Medium	CM041M
Large	CM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEE
Medium	CM041MEE
Large	CM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

novaplus⁺

Standard Elbow, Non-Vented

Size	Item#
Small	VCM041S
Medium	VCM041M
Large	VCM041L

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	VCM041SEE
Medium	VCM041MEE
Large	VCM041LEE

Single Limb Elbow, Vented, AAV

Size	Item#
Small	VCM041SEEV
Medium	VCM041MEEV
Large	VCM041LEEV



- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

With Cap Style Headgear

Standard Elbow, Non-Vented

Size	Item#
Small	CM041SC
Medium	CM041MC
Large	CM041LC

Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM041SEEC
Medium	CM041MEEC
Large	CM041LEEC

Single Limb Elbow, Vented, AAV

Size	Item#
Small	CM041SEEV
Medium	CM041MEEV
Large	CM041LEEV

Nasal Mask

With Standard Headgear

Size	Item#
Small	CM008S
Medium	CM008M
Large	CM008L

novaplus⁺

Size	Item#
Small	VCM008S
Medium	VCM008M
Large	VCM008L



Total Face Mask

Standard Elbow

Size	Item#
Small	CM040S
Medium	CM040M
Large	CM040L

EE Elbow

Size	Item#
Small	CM040SEE
Medium	CM040MEE
Large	CM040LEE



- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy


Sol Full Face CPAP Mask

Mask with Headgear and Replaceable Cushion

10/case

Size	Item#
Small	CM205S
Medium	CM205M
Large	CM205L

Easy Pop Cushions



A premium mask without the premium price

Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

Quick release connector

Rapidly and easily separates or connects mask and the tube

Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

Ray Pediatric Masks Comparable to Nonny masks

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL

Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

Headgear Only

Size	Item#
Small	HG061S
Large	HG061L



Nasal Mask with Headgear, Tube, and Replaceable Cushions


Comes with	Item#
S&M Cushions	CM060-KITS
L&XL Cushions	CM060-KITL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

Headgear Only

Size	Item#
Small	HG060S
Large	HG060L



Available Size Guide

Innova

Nasal Mask with Headgear

Vented

Size	Item#
Small/Medium	CM50180
Medium/Large	CM50181



Size	Item#
Small/Medium	VCM034SM
Medium/Large	VCM034ML

Non-Vented

Size	Item#
Small/Medium	CM50594
Medium/Large	CM50595



Size	Item#
Small/Medium	VCM034NVSM
Medium/Large	VCM034NVML

Innova® CPAP masks feature patented AIR®gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA



MiniMe 2 Pediatric Nasal Mask



Vented

Size	Item#
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented

Size	Item#
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261



Vented

Size	Item#
Extra Small	VCM055XSZ
Small	VCM055SZ
Medium	VCM055MZ
Large	VCM055LZ

Non-Vented

Size	Item#
Extra Small	VCM056XSZ
Small	VCM056SZ
Medium	VCM056MZ
Large	VCM056LZ

- Maximum comfort, minimum weight. Patented AIR®gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA



Veraseal 2 Disposable Gel Mask with Headgear



**Standard Elbow,
Non-Vented**



Size	Item#
Small	CM50933
Medium	CM50934
Large	CM50935
Extra Large	CM50936



Size	Item#
Small	VCM053S
Medium	VCM053M
Large	VCM053L
Extra Large	VCM053XL



**Single Limb Elbow,
Non-Vented, AAV**



Size	Item#
Small	CM50937
Medium	CM50938
Large	CM50939
Extra Large	CM50940



Size	Item#
Small	VCM054S
Medium	VCM054M
Large	VCM054L
Extra Large	VCM054XL



**Single Limb Elbow,
Vented, AAV**



Size	Item#
Small	CM50827
Medium	CM50828
Large	CM50829
Extra Large	CM50830



Size	Item#
Small	VCM052S
Medium	VCM052M
Large	VCM052L
Extra Large	VCM052XL

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.
- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.
- The active headgear connector moves as the patient moves for a reliable seal.
- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

NIV Tubing



Item#
TUB06 6ft
TUB010 10ft
19mm with 22mm cuffs
10/case

- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

Standard Elbow, Non-Vented

Item# **CAP009**
10/case




Single Limb Elbow, Non-Vented, AAV

Item# **CAP008**
10/case



Swivel Connector Silencer

Item# **VEN002**
10/case



Single Use Exhalation Port

Item# **CAP011**
10/case



Standard 4-Point Headgear

Item#
HG003 Standard
HG003XL Extra Large



Cap Style Headgear

Item# **HG002**



RESPIRATORY DISPOSABLES Contract MS1087

Lynx Venturi Oxygen Delivery


Low Concentration

Item# **RES2400**

Each

Adapter Only **RES240LC**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO₂ delivery from 24–60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO₂ washout



Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

Do no harm

Quicker de-escalation of O₂ levels and prevention of excessive CO₂ in the patient.

Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO₂ in the EMS ambulance or ED to prevent CO₂ buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO₂/FiO₂) allows early intervention and better outcomes
- Excellent CO₂ washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

RESPIRATORY DISPOSABLES Contract MS1087

Final Bacteria Filter

Item# **BF516**
1/4in Barbed
on Both Sides
10/pack



Final Bacteria Filter — 3/16in

Item# **BF316**
3/16 Stepped Barb and
1/4in Straight Barb
on Both Sides
10/pack



Tubing Reducer Connector

Item# **RES011**
Fits 3/16in, 1/4in,
5/16in, 3/8in, and 7/16in
ID Tubing
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

Circuit Connectors

Item# **RES013**
22mm OD x 15mm ID/22mm OD
10/pack



Humidifier Connectors

Item# **HUM008**
L-Connector 10in
10/pk

Item# **HUM008L**
L-Connector 21in
10/pk

Item# **HUM008S**
Straight Connector 9in
10/pk



Oxygen Tubing Y-Connector

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack



- Connects 2 oxygen sources for higher flows

Oxygen Tubing Connector — Swivel

Item# **RES018**
Male/Male
10/pack

- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing



Oxygen Bleed-In Adapter

Item# **RES019**
22mm with 6mm Rotating Stem
10/pack



RESPIRATORY DISPOSABLES Contract MS1087



Passive Heated Breathing Circuit

Item#

VC1114

For use with MR850 Humidifier

10/case



Non-Heated Breathing Circuit

For use with Respiration Trilogy

Item# **VC1001**



Active Adult, No Water Trap

22mm cuffs

30/case

OEM#

1073223 & 1073232



Passive Non-Heated Breathing Circuit

For use with Respiration Trilogy

Item# **VC1100**

With fixed exhalation port

Each

OEM# 1069467 & 1069468




Adult NIV Patient Circuit

Compatible with V60, No Filters

Item# **VC1160**

QTY?/case

Comparable to OEM# 1069210



RESPIRATORY DISPOSABLES Contract MS1087

Non-Heated Dual Limb Breathing Circuits

10/case

Item#	Size	Humidifier Limb	Water Trap
VC2001P	Pediatric	No	No
VC2001	Adult	No	No
VC2002P	Pediatric	Yes	No
VC2002	Adult	Yes	No
VC2003	Adult	Yes	Yes



Bacterial/Viral HME Filter

Item# **BF2200**
5/pack



Oxygen Tubing Connector

Item# **RES012**
Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Christmas Tree Connector

Item# **RES001**
Green without Swivel
10/pack

RES001-30 30/pack



Item# **RES002**
Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Circuit Connectors

Item# **RES014**
22mm OD x 22mm OD
With Lip
10/pack



Humidifier Tubing Adapter

Item# **RES016**
22mm ID x 22mm OD
with Lip
10/pack



RESPIRATORY DISPOSABLES Contract MS1087

Flex HME






Item# **RES027FLEX**
40/case



TV Range (mL)	300-1500
Moisture Output (mg H2O/L)	33 mg/L
Weight (g)	17
Resistance (cm H2O)	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (mL)	75
Connection (mm)	22M,15F


Heat Moisture Exchangers

30/case

Item#					
Item#	RES027 T-HME	RES027OP T-HME with Oxygen Port	RES027IL Inline HME — Large Volume	RES027NE Inline HME — Neonatal	RES027HV2 Foam HME — Adult/Child
TV Range (mL)	250-1000	250-1000	50-250	15-50	250-1000
Moisture Output (mg H2O/L)	30, Vt = 0.1L	29.50	31, Vt = 0.1L	30	29
Weight (g)	4.5	7	10	4.5	21
Resistance (cm H2O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
Dead Space (mL)	12.5	12	13	2.5	29
Connection (mm)	15F	15F 5.5mm Oxygen port	15M,22M/15F	15M,15F	22F/15M, 22M/15F

Adapter for T-HME

Item# **RES027A**
Each



Tracheostomy Mask

50/case

Item# **RES2130**
Adult Regular Size

RES2130L
Adult Large Size

RES2230
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O2 supply tubing

RESPIRATORY DISPOSABLES Contract MS1087

Tracheostomy Ties with Fastener Tabs

20/case


One-piece

• Easy on

Item#

RES029 Adult

RES029P Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

Two-piece

• Easy off

Item#

RES030 Adult

RES030P Pediatric



Color Coded Dual Channel Berman Oral Airway

50/pack

Item#	Size	Color
AOA40BC-50	40mm	Pink
AOA50BC-50	50mm	Blue
AOA60BC-50	60mm	Black
AOA70BC-50	70mm	White
AOA80BC-50	80mm	Green
AOA90BC-50	90mm	Yellow
AOA100BC-50	100mm	Red
AOA110BC-50	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

Oropharyngeal Guedel Airway

10/pack

Item#	Size	Color
AOP500GC	50mm	Blue
AOP601GC	60mm	Black
AOP702GC	70mm	White
AOP803GC	80mm	Green
AOP904GC	90mm	Yellow
AOP1005GC	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

RESPIRATORY DISPOSABLES Contract MS1087

Standard Finger Pulse Oximeter

Item# **RES5100**
Each



- 1 display mode
- Digital LED display
- SpO2, PR, and pulse bar

Advanced Finger Pulse Oximeter

Item# **RES5101**
Each



- 4 display modes
- Dual color OLED display
- SpO2, PR, pulse bar and waveform

RESPIRATORY DISPOSABLES Contract MS1087

Firesafe Thermal Shut-Off Valve

Item#
RES010
Bidirectional
10/case




- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire
- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients




FireSafe Nozzle

Item#
RES010FND-20
20/pack



FireSafe Tubing Connector

Item#
RES010TC 100/case
RES010TC-10 10/case



RESPIRATORY DISPOSABLES Contract MS1087

Volumetric Incentive Spirometer

12/case

Item# **RES055**
2500ml

RES056
5000ml

RES055P
Pediatric



SCAN FOR VIDEO



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise.
- Dual-sided, clear and legible identification of volumes.
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration.
- Compact and ergonomic design.

Aerosol Mask

50/case

Item# **RES2110**
Adult Size

Item# **RES2110F**
Adult Size
w/ Flex Tube

Item# **RES2210**
Pediatric Size



- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged

Three Ball Spirometer

Item# **RES057**
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient




- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

Aerosol Tubing

Item# **TUB100-6**
Corrugated Aerosol Tubing
100ft in 6in segments
Each



Item# **TUB100-6B**
Blue Corrugated Aerosol Tubing
100ft in 6in segments
Each



Valved T-Adapter, Aerosol, 22mm ID x 22mm OD

Item# **NEB9001T**
30/case



RESPIRATORY DISPOSABLES Contract MS1087

Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty
RES1100	None	Standard	50
RES1104	4 ft	Standard	50
RES1107N	7 ft	Standard	50
RES1107U	7 ft	Universal	50
RES1107NU	7 ft	Universal	50
RES1107	7 ft	Standard	50
RES1115U	15 ft	Universal	50
RES1115	15 ft	Standard	50
RES1125	25 ft	Standard	25
RES1135	35 ft	Standard	20
RES1150U	50 ft	Universal	20
RES1150Z	50 ft	Standard	20



Soft Cannula

Item#	Tubing Length	Connector	Qty
RES1100S	None	Standard	50
RES1107SU	7 ft	Universal	50
RES1107S	7 ft	Standard	50
RES1114S	14 ft	Standard	50
RES1114SU	14 ft	Universal	50
RES1115S	15 ft	Standard	50
RES1125S	25 ft	Standard	25
RES1125SU	25 ft	Universal	25

Flared Cannula

Item#	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

Adult Cannula with Ear Cushions

Standard Touch, Standard Tip

Item#	Tubing Length	Connector	Qty
RES1100EC	None	Standard	50
RES1104EC	None	Standard	50
RES1107ECU	7 ft	Universal	50
RES1107EC	7 ft	Standard	50
RES1114EC	14 ft	Standard	50
RES1115EC	15 ft	Standard	50
RES1125EC	25 ft	Standard	25

Cannula Ear Cushions

Item# **RES020**
50 Pairs/case

• Add cushions to existing cannula



RESPIRATORY DISPOSABLES Contract MS1087

High Flow Soft Cannula

25/case

Length	Item#
4ft	RES1104SHF
7ft	RES1107SHF
15ft	RES1115SHF
25ft	RES1125SHF




- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

7 ft Cannula With CO₂ Sampling Line

50/case

Item#	Size	Connector	Luer	Filter
RES8007	Adult	Standard	Male	No
RES8007F	Adult	Standard	Male	Yes
RES8007I	Infant	Standard	Male	No
RES8007IF	Infant	Standard	Male	Yes
RES8007IU	Infant	Universal	Male	No
RES8007P	Pediatric	Standard	Male	No
RES8007PF	Pediatric	Standard	Male	Yes
RES8007PU	Pediatric	Universal	Male	No
RES8007U	Adult	Universal	Male	No
RES8107	Adult	Standard	Female	No
RES8107F	Adult	Standard	Female	Yes
RES8107I	Infant	Standard	Female	No
RES8107IF	Infant	Standard	Female	Yes
RES8107P	Pediatric	Standard	Female	No
RES8107PF	Pediatric	Standard	Female	Yes




RESPIRATORY DISPOSABLES Contract MS1087

Secure Cannula Stickers

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No

Humidoflo High Flow Nasal Cannula

10/case

Item#	Size	F&P Adapter
HFN001AS	Small	No
HFN001AM	Medium	No
HFN001AL	Large	No
HFN001ASA	Small	Yes
HFN001AMA	Medium	Yes
HFN001ALA	Large	Yes

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respirationics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)

With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerable due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O₂ concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.



RESPIRATORY DISPOSABLES Contract MS1087

Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector
RES2100NT	Adult	None	Standard
RES2100U	Adult	7 ft	Universal
RES2100	Adult	7 ft	Standard
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203A	Adult	None	Standard
RES2203A7	Adult	7 ft	Standard
RES2203I	Infant	None	Standard
RES2203I7	Infant	7 ft	Standard
RES2203P7	Pediatric	7 ft	Standard

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



High Concentration Oxygen Mask


50/case

Item#	Size	Tubing Length
RES2202A	Adult	None
RES2202A7	Adult	7 ft
RES2202I	Infant	None
RES2202I7	Infant	7 ft
RES2202P	Pediatric	None
RES2202P7	Pediatric	7 ft

Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent
RES2102I	Infant	7 ft	Standard	No
RES2102IM	Infant	None	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No
RES2102S	Adult	7 ft	Standard	Yes
RES2102SM	Adult	None	Standard	Yes
RES2102U	Adult	7 ft	Universal	No
RES2102	Adult	7 ft	Standard	No



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

RESPIRATORY DISPOSABLES Contract MS1087

Oxygen Supply Tubing

Item#	Length	Connector	Color	Bore	Qty
RES3007SB	7 ft	Standard	None	Smooth	50
RES3007U	7 ft	Universal	None	Channeled	50
RES3007	7 ft	Standard	None	Channeled	50
RES3007G	7 ft	Standard	Green	Channeled	50
RES3015G	15 ft	Standard	Green	Channeled	50
RES3015	15 ft	Standard	None	Channeled	50
RES3025G	25 ft	Standard	Green	Channeled	25
RES3025U	25 ft	Universal	None	Channeled	25
RES3025	25 ft	Standard	None	Channeled	25
RES3030	30 ft	Standard	None	Channeled	20
RES3030G	30 ft	Standard	Green	Channeled	20
RES3040G	40 ft	Standard	Green	Channeled	20
RES3050G	50 ft	Standard	Green	Channeled	20
RES3050	50 ft	Standard	None	Channeled	20

Universal Connector



Adult Air Entrainment Mask

50/case

Item#	Size	Kit	Tubing Length	Connector	Flex Tube
RES2400A	Adult	Mask Only	7 ft	Standard	Yes
RES2400AK	Adult	Yes	7 ft	Standard	Yes
RES2400AKU	Adult	Yes	7 ft	Universal	Yes
RES2400AU	Adult	Mask Only	7 ft	Universal	Yes
RES2400P	Pediatric	Mask Only	7 ft	Standard	Yes
RES2400PK	Pediatric	Yes	7 ft	Standard	Yes
RES2400PKU	Pediatric	Yes	7 ft	Universal	Yes
RES2400PU	Pediatric	Mask Only	7 ft	Universal	Yes



Disposable Nebulizer Kits with Mouthpiece


50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- T-piece
- Mouthpiece
- Corrugated tube
- 7ft supply tube

- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal



Item#	Size	Accessories
RES091	Adult	All
RES091NA	Adult	No corrugated tube

Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**

Disposable
Single Patient
50/case

Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



Disposable Nebulizer Kits with Mask

50/case

All include 6cc jet nebulizer medicine cup

Accessory options

- Mask
- 7ft supply tube

Item#	Size	Accessories
RES090	Adult	With Mask
RES092	Pediatric	All



- Disposable single patient use for long term delivery procedures
- 6cc bottle is easy to assemble and disassemble

Medicine Cup

Item# **RES090A**
50/case

MDD Contract MS1097

Nebulizer Compressor

Item# **NEB100**
Each



Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

Features

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Peak Flow Meter

Item# **RES052**
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**
Pediatric Mouthpiece
for Peak Flow Meter
10/case

Valved T-Adapter

Item# **VEN001**
22F-18F-22M/15F
10/pack



SperaChamber

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Item#
RES195-10
10/case

RES195-50
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask
10/case

Item#	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)
RES195ML	Large (Adult)



Adult Mask

