

 belluscura®

# DISCOV-R™

## Continuous & Pulse Flow Portable Oxygen Concentrator

Item#

**RES70051001** DISCOV-R with 1 Battery (Ascent Package)

**RES70051002** DISCOV-R with 2 Batteries (Peak Package)



The world's  
smallest  
**continuous flow  
concentrator**

Airgonomically™ designed using Belluscura's patented ModulAir™ technology, the compact DISCOV-R weighs less than 6.5 lbs, producing more oxygen per pound than any other portable oxygen concentrator.

### Features:

- **4 continuous flow settings**
- **8 pulse flow settings**
- **Weights only 6.5 lbs**
- **FAA compliant for air travel**



*Weights less than  
6.5 lbs. with battery!*



#### Lightweight

Weights only 6.5 pounds (2.95kg) with battery



#### User-Friendly

User-replaceable oxygen enrichment cartridge and direct-charge batteries



#### Airgonomically™ Designed

Designed to be small and discreet for easy travel



#### Bluetooth Capability

Connect the DISCOV-R to your mobile device through the Nomad Biometrics app



#### Multiple Flow Settings

8 pulse flow oxygen settings + 4 continuous flow settings delivering up to 2000mL/min



#### All-In-One Solution

With both pulse dose and continuous flow settings, the DISCOV-R can be used as a home unit as well as for traveling



#### FAA Compliant

Take your device with you when you travel by air



## Ascent Package

---

- DISCOV-R™ Unit
- 8-Cell Battery
- Pair of User-Replaceable Cartridges
- Air Intake Filter
- AC Power Adapter
- DC Power Adapter
- Carry Bag with Handle & Shoulder Strap
- Backpack
- Cart
- Standard 3-Year Warranty

## Peak Package

---

- DISCOV-R™ Unit
- Two (2) 8-Cell Batteries
- Pair of User-Replaceable Cartridges
- Air Intake Filter
- AC Power Adapter
- DC Power Adapter
- Carry Bag with Handle & Shoulder Strap
- Backpack
- Cart
- Standard 3-Year Warranty

Pulse and Continuous Flow Portable Oxygen Concentrator. The DISCOV-R weighs 6.5 lbs\*, with a lightweight design that is perfect for patients looking to replace their heavy, metal oxygen tanks and live a more mobile, active lifestyle.

\*Weight with battery.

## Specs

---

**Dimensions** L x W x H: 9.9 in x 4.1 in x 7.9 in

---

**Weight** ≤ 6.5 lbs (2.95 kg)

---

**Auditory Output** < 39 dB(A) at Level 2

---

**Charge Duration** ≤ 4.5 Hrs

---

**Charge Time** ≤ 6 Hrs

---

**To order or for more information**

P 312.997.9980 E [sales@sunsethcs.com](mailto:sales@sunsethcs.com) W [SUNSETHCS.COM](http://SUNSETHCS.COM)

