








# Heat/Moisture Exchangers (HMEs)



	 <b>RES027</b> T-HME	 <b>RES027OP</b> T-HME with Oxygen Port	 <b>RES027IL</b> Inline HME — Large Volume	 <b>RES027UOP</b> HME with Oxygen Port	 <b>RES027NE</b> Inline HME — Neonatal	 <b>RES027FLEX</b> Flex HME
TV Range (ml)	250-1000	150-1250	50-250	50-1000	15-50	300-1500
Moisture Output (mg H <sub>2</sub> O/L)	30, Vt = 0.1L	36.3 mg/L	31, Vt = 0.1L	13.2 mg/L	30	33
Weight (g)	4.5	7.4	10	5	4.5	17
Resistance (cm H <sub>2</sub> O)	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.2 @ 30 LPM	0.4 @ 10 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	0.5 @ 10 LPM	0.5 @ 30 LPM 1.34 @ 60 LPM
Dead Space (ml)	12.5	7	13	17	2.5	75
Connection (mm)	15F	15F, 5.5mm Oxygen port	15M,22M/15F	15F	15M,15F	22M,15F