



Contracted Supplier



vizient
Awarded Supplier





Sales Representatives



Egle Sventnickaite
872-258-1654
esventnickaite@sunsethcs.com



Mason Moss
262-788-1862
mmoss@sunsethcs.com



Owen Ruth
312-533-2451
owenr@sunsethcs.com



Cindy Geighes
312-533-2478
cgeighes@sunsethcs.com




Hunter Baudino
262-788-1855
hbaudino@sunsethcs.com



Wyatt Kanemoto
262-788-8941
wkanemoto@sunsethcs.com

 *Director, Strategic Accounts*
Lisa Iannazzone RN, BSN
636-222-3134
lisa@sunsethcs.com

 *Manager, Product & Business Development*
Tom Munar
312-533-2476
tmunar@sunsethcs.com



 *Contract Manager, Strategic Accounts*
Gina Rosenberg
914-837-5215
grosenberg@sunsethcs.com

 *RRT Clinical Specialist*
Krystle Bodine BS, RRT
kbodine@sunsethcs.com

Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Returns/Distribution Center

2725 Alft Ln
Elgin IL 60124

West Returns/Distribution Center

6750 W Van Buren St
Phoenix AZ 85043

Content

NIV

	Page
NIV Masks	4
CPAP/NIV Masks	5
Ventilators	6

Oxygen Regulators

	Page
Microdial Flowmeter	7
Regulators	8

Respiratory Disposables

	Page
Spirometers	9
Oxygen Cannula and Accessories	9
Oxygen Masks	10
Air Entrainment Masks	11
Oxygen Supply Tubing	11
Breathing Circuits	12
Trach Ties	13
HMEs	13
Trach Masks	13
Adapters	13
SpO2 Sensors	14

Medication Delivery

	Page
Nebulizers	15-16
Nebulizer Kits	16
Aerosol Masks	16
Inhaler Chamber	16

Closed Suction Catheters

	Page
Standard Switch	17
Deluxe Switch	18

Laryngeal Masks

	Page
Silicone Disposable	19
Silicone Reusable	19
Flexible Silicone Disposable	19
Flexible Silicone Reusable	19
PVC Disposable	20
Flexible PVC Disposable	20
Silicone Disposable Dual Seal	20
PVC Disposable Dual Seal	20

Resuscitation Devices

	Page
Resuscitators	21
Breathing Masks	21



Contracted Supplier

Respiratory Disposables • Contract PP-NS-1741
Resuscitation Devices • Contract PP-NS-1762
Medication Delivery • Contract PP-NS-1753
Laryngeal Masks • Contract PP-OR-2217



vizient
Awarded Supplier

NIV • Contract MS1049
Respiratory Disposables • Contract MS1087
MDD • Contract MS1097

SleepNet MiniMe 2 Pediatric Nasal Mask V

10/pack

Vented Mask with Headgear

Size	Item#
Extra Small	CM60216
Small	CM60217
Medium	CM60218
Large	CM60221

Non-Vented Mask with Headgear

Size	Item#
Extra Small	CM60256
Small	CM60257
Medium	CM60258
Large	CM60261

novaplus⁺

Vented Mask with Headgear

Size	Item#
Extra Small	VCM055XS
Small	VCM055S
Medium	VCM055M
Large	VCM055L

Non-Vented Mask with Headgear

Size	Item#
Extra Small	VCM056XS
Small	VCM056S
Medium	VCM056M
Large	VCM056L

HCPCS#
MASK - A7034
HEADGEAR - A7035

Replacement Headgear

Size	Item#
Extra Small/Small	HG027S
Medium/Large	HG027L



Ray Pediatric Masks Comparable to Nonny masks V

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	CM061S
Medium	CM061M
Large	CM061L
Extra Large	CM061XL

Replacement Cushions

Size	Item#
Small	CUCM061S
Medium	CUCM061M
Large	CUCM061L
Extra Large	CUCM061XL

Headgear Only

Size	Item#
Small	HG061S
Large	HG061L

SCAN CODE to see
Ray Full Face Mask
on YouTube



Nasal Mask with Headgear, Tube, and Replaceable Cushions


Comes with	Item#
S&M Cushions	CM060-KITS
L&XL Cushions	CM060-KITL

Replacement Cushions

Size	Item#
Small	CUCM060S
Medium	CUCM060M
Large	CUCM060L
Extra Large	CUCM060XL

Headgear Only

Size	Item#
Small	HG060S
Large	HG060L




Available Size Guide



SCAN CODE to see
Ray Nasal Mask
on YouTube



MiniMe Pediatric Nasal Mask



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

Headgear Size	Item#	OEM#	
Small	CM004S	60213	V
Medium	CM004M	60214	V P
Large	CM004L	60215	V

HCPCS# A7034



Pediatric SleepWeaver® CPAP Mask V

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

Vivo 45LS Portable Life Support Ventilator



Item#
VNT230016

The Vivo 45LS is designed to maximize independence and mobility for adult and pediatric patients from 11 lbs/5 kg.

Packed with features like:

- Auto-EPAP*
- HFNT*
- AHI* plus built-in humidifier* with heated wire circuit*
- eSync
- Built in information iButton
- 24+ hour battery
- Can connect to the securely hosted Breas cloud solution, EveryWare

*Modification in accordance with FDA's guidance, Enforcement Policy for Ventilators and Accessories and Other Respiratory Devices During the Coronavirus Disease 2019 (COVID-19) Public Health Emergency.

Vivo 65 Dual Limb Life-Support Ventilator



Item#
VNT224516

The Vivo 65 is an advanced homecare ventilator designed to deliver comfortable and effective (life-support) treatment to adult and pediatric patients from 11 lbs/5 kg.

The Vivo 65 can be used for a wide variety of patients thanks to a comprehensive set of modes, circuits and accessories. The dual limb circuit allows for measurement of exhaled volumes providing additional security.

The extensive monitoring capabilities, including EveryWare, help obtain a better insight into the quality of ventilation. The Vivo 65 is an excellent choice for mechanical ventilation at home, in the hospital and in long term care facilities.

Microflow Microdial Flowmeter P

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160048			A
RES8160050	400 kPa (58 psi)	350 – 550 kPa	E



DISS Inlet • DISS Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160049	400 kPa (58 psi)	345 – 500 kPa	A
RES8160051	350 kPa (51 psi)	250 – 450 kPa	E

*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
Range A	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
Range E	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

Assured, repeatable accuracy to enable safe, effective weaning.

What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned



Mini Oxygen Regulator CGA 870 Style

Item# **RES046**

Pediatric • Barbed

0-4 LPM • **Each**

0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$, $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1, 1 $\frac{1}{2}$, 2, 3, 4 LPM

HCPCS# E1353

For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.

Barbed hose outlet
for easy connection to
the oxygen tubing



"T" Handle for easy
mounting — no tools
needed

Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty

Mini Oxygen Regulator CGA 540 Style

Item# **RES047**

Pediatric • DISS

0-4 LPM • **Each**

0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$, $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1, 1 $\frac{1}{2}$, 2, 3, 4 LPM

HCPCS# E1353

For use with M60 and larger cylinders: M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.

Threaded
DISS outlet



Volumetric Incentive Spirometer P

12/case

Item#
RES055P

Pediatric

- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design



Pediatric, Infant, and Neonate Cannula V

50/case

Item #	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No



Ear Cushions for Cannula V P

Item# **RES020**
Cushions Only
50 Pairs/case



- Add cushions to existing cannula



Secure Cannula Stickers V P

Item# **RES080**
25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



High Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length
RES2202I	Infant	None
RES2202I7	Infant	7 ft
RES2202P	Pediatric	None
RES2202P7	Pediatric	7 ft



Medium Concentration Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector
RES2203P	Pediatric	None	Standard
RES2100P	Pediatric	7 ft	Standard
RES2203I	Infant	None	Standard
RES2203I7	Infant	7 ft	Standard
RES2203P7	Pediatric	7 ft	Standard



- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather Oxygen Mask V

50/case

Item #	Size	Tubing Length	Connector	Safety Vent
RES2102I	Infant	7 ft	Standard	No
RES2102IM	Infant	None	Standard	No
RES2102P	Pediatric	7 ft	Standard	No
RES2102PS	Pediatric	7 ft	Standard	Yes
RES2102PSM	Pediatric	None	Standard	Yes
RES2102PSU	Pediatric	7ft	Universal	Yes
RES2102PU	Pediatric	7 ft	Universal	No

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



Air Entrainment Mask

With Flex Tube

50/case

Item#	Kit	Tubing Length	Connector
RES2400P	Mask Only	7 ft	Standard V
RES2400PK	Yes	7 ft	Standard V P
RES2400PKU	Yes	7 ft	Universal V
RES2400PU	Mask Only	7 ft	Universal V P



Oxygen Supply Tubing **V**

Item #	Length	Connector	Bore	Color	Qty
RES3007	7ft	Standard	Channelled	None	50
RES3007SB	7 ft	Standard	Smooth	None	50
RES3007U	7 ft	Universal	Channelled	None	50
RES3015	15 ft	Standard	Channelled	None	50
RES3015G	15 ft	Standard	Channelled	Green	50
RES3025	25 ft	Standard	Channelled	None	25
RES3025G	25 ft	Standard	Channelled	Green	25
RES3025U	25 ft	Universal	Channelled	None	25
RES3030	30 ft	Standard	Channelled	None	20
RES3030G	30 ft	Standard	Channelled	Green	20
RES3040G	40 ft	Standard	Channelled	Green	20
RES3050	50 ft	Standard	Channelled	None	20
RES3050G	50 ft	Standard	Channelled	Green	20

- 6-channel
- PHT/DEHP-free




Universal Connector

Non-Heated Breathing Circuits


For use with Respiration Trilogy
22mm cuffs
30/case

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
VC1001P	1073220 & 1073229	Active	Yes	Pediatric	15mm
VC1004P	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#
VC1001P
Active Pediatric
OEM#
1073220
1073229



Item#
VC1004P
Passive Pediatric
OEM#
1073224
1073233



To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Non-Heated Dual Limb Breathing Circuits V

30/case

REF	Humidifier Limb	Water Trap
VC2001P	No	No
VC2002P	Yes	No




Passive Heated Breathing Circuit P V

Item# **VC1003P**
**With Swivel
Exhalation Port**
For use with Philips
Respiration Trilogy
Each

Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port

OEM# 1076046,
1073227, 1073235




Tracheostomy Ties with Fastener Tabs **P**

20/case

One-piece

- Easy on

Item#
RES029P Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton

Two-piece

- Easy off

Item#
RES030P Pediatric



Heat Moisture Exchanger **V**

30/case

Item#
RES027NE-10
Inline HME — Neonatal



TV Range (mL)	15-50
Moisture Output (mg H₂O/L)	30
Weight (g)	4.5
Resistance (cm H₂O)	0.5 @ 10 LPM
Dead Space (mL)	2.5
Connection (mm)	15M,15F

Tracheostomy Mask with O₂ Connector **V**

50/case

Item#
RES2230
Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing

Oxygen Bleed-In Adapter **P**

Item# **RES019** 10/pack
22mm with 6mm
Rotating Stem



Flexible Connector **P V**

50/case

Item# **TRAA004P**
15mm F x 15mm M
Length 4.9-7.5cm
OEM# 3215



SpO2 Sensors


NEW

Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**

Similar to Nellcor MAX-N

20/box




Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**

1.5ft cord • Non-replaceable

20/box




Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**

3ft cord • Non-replaceable

20/box



Coban Tape Roll

Item# **RES7202CT**

Use for sensor life extension
or sensor cord stability


Each roll 1" x 177"

Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**

1.5ft cord

20/box



Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**

1.5ft cord

20/box



Compressor Nebulizers

Item#

NEB100 **V**

Each

Includes

- Compressor
- Disposable T-Piece Kit
- Foam Filters

NEB100-R

- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Sunny the Bear

NEB300BEAR

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

NEB300BEAR-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



Gordon the Fish

NEB300FISH

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Foam Filters (BF018 x5)

NEB300FISH-R

- Disposable Pediatric Mask (RES2210)
- Disposable T-Piece Kit (RES091)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)



Handheld Compressor Nebulizer **P**

Item#

NEB200

Each

Included

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (approximately 3 microns)
- No replacement filter required



Handheld Mesh Nebulizer **P**

Item#

NEB401

Each

With internal rechargeable battery

Included

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual

- Small particle size: MMAD 2.2µm*
- Efficient: Low Residual Liquid
- Convenient: Small and Portable
- Quiet: ≤50dB (at one meter)
- Modern: Rechargeable Lithium Battery



Handheld Mesh Nebulizer

Item#
NEB402

With internal
rechargeable battery

Included
Medication cup
USB cord
Mask
Mouthpiece
Portable bag
Manual

Particle Size
MMAD < 5µm

Nebulization Rate
0.25ml/min-0.9ml/min

Medication Cup Capacity
10ml (Max)

Device Weight
Approx. 1/4 lb

Device Size
Approx. 1 1/2 x 2 x 3 7/8 inches



Disposable Nebulizer Kit V

Item# **RES092**
Pediatric Size
50/case

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



Aerosol Mask V P

50/case

RES2210
Pediatric Size

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



SperaChamber V P

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Item#
RES195-10
10/case

RES195-50
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



SperaChamber with Mask 10/case

Item#	Size
RES195MS	Small (Infant)
RES195MM	Medium (Pediatric)

Mask Only

Item#	Size	Qty
RES195MASK-S	Small (Infant)	Each
RES195MASK-S-10	Small (Infant)	10/box
RES195MASK-M	Medium (Pediatric)	Each
RES195MASK-M-10	Medium (Pediatric)	10/box



Closed Ventilation Suction System

Deluxe Switch (One-handed operation)

Sterile and latex-free

Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port

Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.

Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.



Closed Ventilation Suction System

Standard Switch

Sterile and latex-free

24 Hour Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

72 Hour Neonatal/Pediatric

With Y-Adapter

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr

Silicone Disposable Laryngeal Masks P

Item #	Size
ALM1001	#1
ALM1015	#1.5
ALM1002	#2
ALM1025	#2.5
ALM1003	#3

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



Silicone Reusable Laryngeal Masks P

Item #	Size
ALM1001R	#1
ALM1015R	#1.5
ALM1002R	#2
ALM1025R	#2.5
ALM1003R	#3

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



Flexible Silicone Disposable Laryngeal Masks P

Item #	Size
ALM3001	#1
ALM3015	#1.5
ALM3002	#2
ALM3025	#2.5
ALM3003	#3



Flexible Silicone Reusable Laryngeal Masks P

Item #	Size
ALM3001R	#1
ALM3015R	#1.5
ALM3002R	#2
ALM3025R	#2.5
ALM3003R	#3



- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness

Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

PVC Disposable Laryngeal Masks P

Item #	Size
ALM2001	#1
ALM2015	#1.5
ALM2002	#2
ALM2025	#2.5
ALM2003	#3

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



Flexible PVC Disposable Laryngeal Masks P

Item #	Size
ALM4001	#1
ALM4015	#1.5
ALM4002	#2
ALM4025	#2.5
ALM4003	#3



Silicone Disposable Dual Seal Laryngeal Masks P

Item #	Size
ALM5015	#1.5
ALM5002	#2
ALM5025	#2.5
ALM5003	#3



PVC Disposable Dual Seal Laryngeal Masks P

Item #	Size
ALM6003	#3



Size	Patient Weight
#1	Up to 5 kg
#1.5	5-10 kg
#2	10-20 kg
#2.5	20-30 kg
#3	30-50 kg

Child Resuscitator P

12/case

Item #	Pop-Off	Mask
RES8110H	Yes	None
RES8111	Yes	#2
RES8112	Yes	#3



Infant Resuscitator P

12/case

Item #	Pop-Off	Mask
RES8120H	Yes	None
RES8121	Yes	#0
RES8122	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

Breathing Mask P

20/case

Item #	Size	Hooks
SBC01N	Neonatal (#0)	No
SBC01I	Infant (#1)	No
SBC01P	Pediatric (#2)	No
SBC02N	Neonatal (#0)	Yes
SBC02I	Infant (#1)	Yes
SBC02P	Pediatric (#2)	Yes



