


Heat/Moisture Exchangers (HMEs)



	 RES027 T-HME	 RES027OP T-HME with Oxygen Port	 RES027IL Inline HME with Gas Sampling Port	 RES027UOP HME with Oxygen Port	 RES027NE Inline HME — Neonatal	 RES027FLEX Flex HME
TV Range (ml)	150-1250	150-1250	70-600	50-1000	5-50	250-1500
Moisture Output (mg H ₂ O/L)	36.4 at VT = 500 ml	36.3	34.5 at VT = 250 ml 32.5 at VT = 500 ml	13.2	11 at VT = 250 ml	30
Weight (g)	6.6	7.4	10	5	1.9	15.5
Resistance (cm H ₂ O)	0.1 @ 30 LPM	0.2 @ 30 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM	2.3 @ 30 LPM	0.5 @ 30 LPM 1.4 @ 60 LPM
Dead Space (ml)	15	7	21	17	2.2	90
Connection (mm)	15F	15F, 5.5mm Oxygen port	22M/15F-15M	15F	15F-15M	22F, 15F