

Portable Oxygen Concentrator

Lightweight and compact, the X-PLOR® is ideal for patients who want to replace their heavy oxygen tanks without missing a moment. The X-PLOR is also Bluetooth®-capable and can be connected to your mobile device through Belluscura's newly developed Nomad Biometrics™ App. With a rechargeable battery, easy-to-use interface, and high oxygen output, supplementary oxygen has never been more approachable or affordable.



Item# **RES52051001**

X-PLOR

Item# **RES52051002**

X-PLOR with Extra Battery



Included

X-PLOR® Unit

8-Cell Battery

User-Replaceable Sieve Bed Cartridge

Air Intake Filter

AC Power Adaptor

DC Power Adaptor

Carry Bag With Handle & Shoulder Strap

Features

- Oxygen output 82%-92% purity
- User-Replaceable Sieve Bed Cartridge
- Airgonomically™ Designed to be lightweight and compact for easy travel
- Long-Lasting, Direct Charge Battery: Up to 4.5 hours of use when fully charged
- Pulse flow settings 1 through 5 delivering 215 to 900 ml/min
- No Breath Detected (NBD) mode is a battery-saving mode that runs at Auto-pulse15 BPM at the units current setting
- FAA Compliant



Pulse and No Breath Detection mode Portable Oxygen Concentrator. The X-PLOR weighs 3.75 lbs*, with a lightweight design that is perfect for patients looking to replace their heavy, metal oxygen tanks and live a more mobile, active lifestyle.

Belluscura's patented, award-winning ModulAir® technology integrated with their proprietary ultra-rapid processing and consumer-replaceable sieve technology allows them to make smaller, smarter, and more efficient oxygen concentrators, resulting in reduced costs and broader potential treatment options.

The Nomad Biometrics™ mobile app tracks your X-PLOR oxygen concentrator unit's oxygen purity, battery life, alarm history, and battery cycle. It is Bluetooth®-compatible with an easy-to-use interface. By downloading the Nomad Biometrics app, you can enjoy the experience of oxygen therapy on the go. Available in the iOS store and Google Play.

X-PLOR is FAA compliant, meaning it is just as safe and easy to travel with as it is reliable!

*Weight with battery.



Metrics

Dimensions 7.3 in x 2.9 in x 7.6 in (18.5 cm x 7.4 cm x 19.3 cm)

Weight ≤ 3.75 lb (1.7 kg)

Auditory Output < 39 dB(A) at Level 2

Charge Duration ≤ 4.5 Hrs

Charge Time ≤ 6 Hrs

Specifications

	Setting 1	Setting 2	Setting 3	Setting 4	Setting 5
Total Volume Per Minute (mL/min)	215	425	635	750	900
No Breath Detected Mode	15 BPM Auto pulse				

Accessories

X-PLOR® User-Replaceable Sieve Bed Cartridge
RES50073



X-PLOR® 8-Cell Battery
RES50188



X-PLOR® DC Power Adaptor
RES53012



X-PLOR® AC Power Adaptor
RES53011



X-PLOR® Carry Bag
RES52191



X-PLOR® Backpack
RES52291



To order or for more information

P 312.997.9980 E sales@sunsethcs.com W SUNSETHCS.COM

