




# Acute Care

## Product Catalog



**Contracted Supplier**



**vizient.**  
Awarded Supplier





## Sales Representatives



**Egle Sventnickaite**  
872-258-1654  
esventnickaite@sunsethcs.com



**Mason Moss**  
262-788-1862  
mmoss@sunsethcs.com



**Owen Ruth**  
312-533-2451  
owenr@sunsethcs.com




**Cindy Geighes**  
312-533-2478  
cgeighes@sunsethcs.com



**Hunter Baudino**  
262-788-1855  
hbaudino@sunsethcs.com



**Wyatt Kanemoto**  
262-788-8941  
wkanemoto@sunsethcs.com



**Director, Strategic Accounts**  
**Lisa Iannazzone** RN, BSN  
636-222-3134  
lisa@sunsethcs.com

**Manager, Product & Business Development**  
**Tom Munar**  
312-533-2476  
tmunar@sunsethcs.com

**Contract Manager, Strategic Accounts**  
**Gina Rosenberg**  
914-837-5215  
grosenberg@sunsethcs.com

**RRT Clinical Specialist**  
**Krystle Bodine** BS, RRT  
kbodine@sunsethcs.com

### Corporate Office

Phone 872-275-3210

Fax 312-997-9985

141 W Jackson Blvd Ste 1950  
Chicago IL 60604

### Central Returns/Distribution Center

2725 Alft Ln  
Elgin IL 60124

### West Returns/Distribution Center

6750 W Van Buren St  
Phoenix AZ 85043

# Content

## **NIV** **Page**

Adult NIV Masks	4-9
Pediatric NIV Masks	6-9
CPAP Tubing	10
Connectors and Ports	10
Headgear	11
Chinstraps	12

## **Microdial Flowmeter** **Page**

Microflow	13
-----------	----

## **Respiratory Disposables** **Page**

Vent Circuits	14-17
Lynx Venturi Oxygen Delivery	18
Oxygen Masks	19-20
Air Entrainment Masks	20
Oxygen Cannula and Accessories	21-23
Oxygen Supply and CPAP Tubing	24
Firesafe	25
Oxygen Connectors and Adapters	26
Oxygen Accessories	27
SpO2 Sensors	28
Bacteria Filters and Humidifiers	29-30
Oxygen Filters	30-31
Spirometers	32
Drainage Bags	32
Aerosol Tubing	32
Concentrator Filters	33-34
Berman Airways	35
Guedel and Nasopharyngeal Airways	36
HMEs and Accessories	37
Suction Catheter Kits	38
Suction Catheters	39
Trach Ties and Connectors	39
Suction Canisters and Liners	40
Suction Canister/Liner Kits	41
Suction Canister Brackets and Wall Mounts	42
Suction Tubing, Yankauers, and Bulk Roll	43
Solidifiers	44

## **Medication Delivery** **Page**

Inhaler Chambers	45
Peak Flow Meters	45
T-Adapters	45
Compressor Nebulizers	46
Handheld and Mini Nebulizers	46
Nebulizer Kits	47
Aerosol Masks	47

## **Closed Suction Catheters** **Page**

Standard Switch	48
Deluxe Switch	49-50

## **Resuscitation Devices** **Page**

Resuscitators	51
Manometers	52
PEEP Valves	52
Reservoir Bags	53
Pocket Resuscitator Set	53
Breathing Masks	53

## **Laryngeal Masks** **Page**

Silicone Disposable	54
Silicone Reusable	54
Flexible Silicone Disposable	54
Flexible Silicone Reusable	54
PVC Disposable	55
Flexible PVC Disposable	55

## **Anesthesia** **Page**

PeDIA	56
-------	----

## **Vizient**

NIV • Contract MS1049  
Respiratory Disposables • Contract MS1087  
MDD • Contract MS1097

## **Premier**

Respiratory Disposables • Contract PP-NS-1741  
Closed Suction Catheters • Contract PP-NS-1571  
Resuscitation Devices • Contract PP-NS-1762  
Medication Delivery • Contract PP-NS-1753  
Laryngeal Masks • Contract PP-OR-2217  
Medical Gas Pipeline Equipment • PP-FA-2138





## Full Face Mask

### With Standard Headgear

#### Standard Elbow, Non-Vented

Size	Item#
Small	<b>CM041S</b>
Medium	<b>CM041M</b>
Large	<b>CM041L</b>

#### Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	<b>CM041SEE</b>
Medium	<b>CM041MEE</b>
Large	<b>CM041LEE</b>

#### Single Limb Elbow, Vented, AAV

Size	Item#
Small	<b>CM041SEEV</b>
Medium	<b>CM041MEEV</b>
Large	<b>CM041LEEV</b>

**novaplus<sup>+</sup>**

#### Standard Elbow, Non-Vented

Size	Item#
Small	<b>VCM041S</b>
Medium	<b>VCM041M</b>
Large	<b>VCM041L</b>

#### Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	<b>VCM041SEE</b>
Medium	<b>VCM041MEE</b>
Large	<b>VCM041LEE</b>

#### Single Limb Elbow, Vented, AAV

Size	Item#
Small	<b>VCM041SEEV</b>
Medium	<b>VCM041MEEV</b>
Large	<b>VCM041LEEV</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM005S-10</b>
Medium	<b>CUCM005M-10</b>
Large	<b>CUCM005L-10</b>



### With Cap Style Headgear

#### Standard Elbow, Non-Vented

Size	Item#
Small	<b>CM041SC</b>
Medium	<b>CM041MC</b>
Large	<b>CM041LC</b>

#### Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	<b>CM041SEEC</b>
Medium	<b>CM041MEEC</b>
Large	<b>CM041LEEC</b>

#### Single Limb Elbow, Vented, AAV

Size	Item#
Small	<b>CM041SEEV</b>
Medium	<b>CM041MEEV</b>
Large	<b>CM041LEEV</b>

- Unobtrusive design offers a clear line of sight
- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone cushion increases comfort while maintaining optimal seal

## Sunset Deluxe Nasal CPAP Masks

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	<b>CM006S</b>
Medium	<b>CM006M</b>
Large	<b>CM006L</b>

**HG003** Headgear - Fits Most  
**HG003XL** Headgear - XL

**CMSIZE01**  
Deluxe Sizing Guide

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
CUSHION - A7032

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



### Mask without Headgear

Size	Item#
Small	<b>CM006SNH</b>
Medium	<b>CM006MNH</b>
Large	<b>CM006LNH</b>

### Replacement Cushions

Size	Item#
Small	<b>CUCM006S</b>
Medium	<b>CUCM006M</b>
Large	<b>CUCM006L</b>



## Nasal Mask V

### With Standard Headgear

Size	Item#
Small	CM008S
Medium	CM008M
Large	CM008L

novaplus+

Size	Item#
Small	VCM008S
Medium	VCM008M
Large	VCM008L



### Standard Elbow, Non-Vented

Size	Item#
Small	CM008SVNV-10
Medium	CM008MVNV-10
Large	CM008LVNV-10

- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

## Total Face Mask V

### Standard Elbow

Size	Item#
Small	CM040S
Medium	CM040M
Large	CM040L

### Single Limb Elbow, Non-Vented, AAV

Size	Item#
Small	CM040SEE
Medium	CM040MEE
Large	CM040LEE



## Sunset Deluxe Nasal Pillow CPAP Mask V

### Mask with Headgear and Small, Medium, and Large Pillows

Item# CM036

### Headgear Only HG036

HCPCS#  
MASK - A7034  
HEADGEAR - A7035  
PILLOWS - A7033

### Replacement Pillows

Size	Item#
Small	CM036SP
Medium	CM036MP
Large	CM036LP



- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

### Mask without Headgear CM036NH

## Sunset Clearsight Deluxe Nasal Mask

### Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM110S V
Medium	CM110M V
Large	CM110L V

### Replacement Cushions

Size	Item#
Small	CUCM110S V
Medium	CUCM110M V
Large	CUCM110L V

HCPCS#  
MASK - A7034  
HEADGEAR - A7035



*Clearsight design provides unobstructed view*

## Sol Full Face CPAP Mask

### Mask with Headgear and Replaceable Cushion


10/case

Size	Item#
Small	<b>CM205S</b>
Medium	<b>CM205M</b>
Large	<b>CM205L</b>

#### Replacement Cushions

Size	Item#
Small	<b>CUCM205S-10</b>
Medium	<b>CUCM205M-10</b>
Large	<b>CUCM205L-10</b>

Easy Pop Cushions



### A premium mask without the premium price

#### Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

#### Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

#### Quick release connector

Rapidly and easily separates or connects mask and the tube

#### Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

## Ray Pediatric Masks Comparable to Nonny masks

### Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#
Small	<b>CM061S</b>
Medium	<b>CM061M</b>
Large	<b>CM061L</b>
Extra Large	<b>CM061XL</b>

#### Replacement Cushions

Size	Item#
Small	<b>CUCM061S</b>
Medium	<b>CUCM061M</b>
Large	<b>CUCM061L</b>
Extra Large	<b>CUCM061XL</b>

#### Headgear Only

Size	Item#
Small	<b>HG061S</b>
Large	<b>HG061L</b>

SCAN CODE to see Ray Full Face Mask on YouTube



### Nasal Mask with Headgear, Tube, and Replaceable Cushions


Comes with	Item#
S&M Cushions	<b>CM060-KITS</b>
L&XL Cushions	<b>CM060-KITL</b>

#### Replacement Cushions

Size	Item#
Small	<b>CUCM060S</b>
Medium	<b>CUCM060M</b>
Large	<b>CUCM060L</b>
Extra Large	<b>CUCM060XL</b>

#### Headgear Only

Size	Item#
Small	<b>HG060S</b>
Large	<b>HG060L</b>



Available Size Guide



SCAN CODE to see Ray Nasal Mask on YouTube



## Sleepnet Innova

### Nasal Mask with Headgear

#### Vented

Size	Item#
Small/Medium	<b>CM50180</b>
Medium/Large	<b>CM50181</b>

**novaplus<sup>+</sup>** V

Size	Item#
Small/Medium	<b>VCM034SM</b>
Medium/Large	<b>VCM034ML</b>

#### Non-Vented

Size	Item#
Small/Medium	<b>CM50594</b>
Medium/Large	<b>CM50595</b>

**novaplus<sup>+</sup>** V

Size	Item#
Small/Medium	<b>VCM034NVSM</b>
Medium/Large	<b>VCM034NVML</b>

Innova® CPAP masks feature patented AIR®gel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA



## Sleepnet Ascend Mask Kit

### Full Face System

Item# **CM037F**

Each

Vented

#### Mask Only

Size	Item#
Small	<b>CM037S</b> V
Medium	<b>CM037M</b> V
Large	<b>CM037L</b> V

#### Replacement Cushions

Size	Item#
Small	<b>CM037CS</b> V
Medium	<b>CM037CM</b> V
Large	<b>CM037CL</b> V



### Nasal System

Item# **CM037N**

Each

Vented

#### Mask Only

Size	Item#
Small	<b>CM037NS</b> V
Medium	<b>CM037NM</b> V
Large	<b>CM037NL</b> V

#### Replacement Cushions

Size	Item#
Small	<b>CM037NCS</b>
Medium	<b>CM037NCM</b>
Large	<b>CM037NCL</b>



## Sleepnet Phantom®2 and iQ®2 V

### Phantom®2

Item# **CM050**

Replacement Cushion  
**CUCM050**

Sleepnet's patented AIR®gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR®gel® allow the mask to seal with minimal tension. By design AIR®gel® is twice as soft as competing gel products.



Replacement Frame/Tubing Assembly  
**CM051**

### iQ®2

Item# **CM049**

Replacement Cushion  
**CUCM049**

- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.



Ask us about our 30-day mask guarantee

- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.



## Veraseal 2 Disposable Gel Mask with Headgear V



### Standard Elbow, Non-Vented



Size	Item#
Small	<b>CM50933</b>
Medium	<b>CM50934</b>
Large	<b>CM50935</b>
Extra Large	<b>CM50936</b>



Size	Item#
Small	<b>VCM053S</b>
Medium	<b>VCM053M</b>
Large	<b>VCM053L</b>
Extra Large	<b>VCM053XL</b>



### Single Limb Elbow, Non-Vented, AAV



Size	Item#
Small	<b>CM50937</b>
Medium	<b>CM50938</b>
Large	<b>CM50939</b>
Extra Large	<b>CM50940</b>



Size	Item#
Small	<b>VCM054S</b>
Medium	<b>VCM054M</b>
Large	<b>VCM054L</b>
Extra Large	<b>VCM054XL</b>



### Single Limb Elbow, Vented, AAV



Size	Item#
Small	<b>CM50827</b>
Medium	<b>CM50828</b>
Large	<b>CM50829</b>
Extra Large	<b>CM50830</b>



Size	Item#
Small	<b>VCM052S</b>
Medium	<b>VCM052M</b>
Large	<b>VCM052L</b>
Extra Large	<b>VCM052XL</b>

- The Veraseal® 2 disposable full face mask brings AIR°gel® Technology to the acute care environment for a giant step forward in patient comfort. Three color coded elbow designs make choosing the right mask for CPAP, single, or dual limb circuit NIV systems simple. The Veraseal® 2 fits quickly and easily while enhancing patient comfort. Both vented (with leak ports) and non-vented (without leak ports) versions are available for CPAP, Bilevel, and ventilator therapy.
- The AIR°gel® cushion provides the utmost in gentle cushioning to protect from pressure points and irritation. AIR°gel® promotes advanced skin protection, superior performance, and unbeatable comfort, all in a lightweight mask.
- The active headgear connector moves as the patient moves for a reliable seal.
- Breathable headgear lets patients enjoy soft, lightweight stability they can count on for a comfortable fit.

## MiniMe 2 Pediatric Nasal Mask V



### Vented

Size	Item#
Extra Small	<b>CM60216</b>
Small	<b>CM60217</b>
Medium	<b>CM60218</b>
Large	<b>CM60221</b>

### Non-Vented

Size	Item#
Extra Small	<b>CM60256</b>
Small	<b>CM60257</b>
Medium	<b>CM60258</b>
Large	<b>CM60261</b>



### Vented

Size	Item#
Extra Small	<b>VCM055XS</b>
Small	<b>VCM055S</b>
Medium	<b>VCM055M</b>
Large	<b>VCM055L</b>

### Non-Vented

Size	Item#
Extra Small	<b>VCM056XS</b>
Small	<b>VCM056S</b>
Medium	<b>VCM056M</b>
Large	<b>VCM056L</b>



- Maximum comfort, minimum weight. Patented AIR°gel makes MiniMe masks some of the most gentle fitting masks on the market.
- The shapeable shell mask can be easily molded to fit individual needs.
- 100% designed and manufactured in the USA.

## Pediatric SleepWeaver® CPAP Mask **V**

Item#  
**CM101325**



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.


- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

## MiniMe Pediatric Nasal Mask



Headgear Size	Item#	OEM#	
Small	<b>CM004S</b>	60213	<b>V</b>
Medium	<b>CM004M</b>	60214	<b>V P</b>
Large	<b>CM004L</b>	60215	<b>V</b>

HCPCS# A7034



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

## SleepWeaver® Advance Nasal Mask **V**

Item#  
**CM017B**

Mask with Adjustable Headgear

- Easy to store and pack for travel

OEM# 100274  
 HCPCS# A7034



## NIV Tubing

Item#

**TUB0018** 18in**TUB02** 2ft**TUB004** 4ft**TUB06-10** 6ft V**TUB08** 8ft V P**TUB010** 10ft V P

19mm with 22mm cuffs

**10/case**

- Sunset's most popular item
- ISO standard 22mm cuffs
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

## Slim CPAP Tube V P

Item#

**TUB006SS**6ft x 15mm with  
22mm cuffs

HCPCS# A7037



### Highly durable

Strain relief cuff decreases  
separation and tearing

### Virtually no memory

Designed to be the most flexible  
tube available

### Minimal weight

One of the lightest tubes on the  
market at 2.6 ounces

### Great value

Affordable and sustainable


## Heated Tube for Respironics P

Item#

**TUB06HPR**Heated Tube for Philips  
Respironics Dreamstation  
Series and System One

OEM# HT15

HCPCS# A4604



- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off

## Standard Elbow, Non-Vented

Item# **CAP009****10/case**

## Single Limb Elbow, Non-Vented, AAV

Item# **CAP008****10/case**

## Swivel Connector Silencer

Item#

**VEN002** Each P**VEN002-10** 10/pack V Contract MS1087

## Single Use Exhalation Port V

Item# **VEN003****10/pack**



## Standard 4-Point Headgear **V**

Item#  
**HG003S** Small  
**HG003M** Medium  
**HG003L** Large



## Cap Style Headgear **V**

Item# **HG002**



## Universal Neoprene Chinstrap

Item#  
**CS001**  
 Universal Chinstrap  
 HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows



## Deluxe Chinstrap

Size	Item#
26in	<b>CS004S</b>
28in	<b>CS004L</b>

Compare to  
 Respironics Deluxe  
 OEM# 302425  
 HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort



## Comfort Chinstrap

Item# **CS005**  
 Compare to  
 ResMed Sullivan  
 OEM# 16015  
 HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use—attaches to itself at the top of the head
- Soft cloth chin cup adds comfort



## Premium Chinstrap

Item# **CS006**  
 Compare to  
 Respironics Premium  
 OEM# 1012911  
 HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.



## Ruby Style Adjustable Chinstrap




Size	Item#
Regular	<b>CS007-ADJ</b>
Extra Large	<b>CS007-ADJXL</b>

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



## Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	<b>CS007S</b>
22-23in	<b>CS007M</b>
23-24in	<b>CS007L</b>

Compare to Puresom Ruby

OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

## Halo Style Chinstrap V

Item# **CS025**

Compare to BreatheWear Halo Chinstrap

OEM# H3407

HCPCS# A7036



## Microflow Microdial Flowmeter P

2024 Innovation Celebration Winner



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160048</b>	345 – 500 kPa	<b>A</b>
<b>RES8160050</b>	350 – 550 kPa	<b>E</b>



DISS Inlet • DISS Outlet

Item#	Inlet Pressure Range	Flow Range*
<b>RES8160049</b>	345 – 500 kPa	<b>A</b>
<b>RES8160051</b>	250 – 450 kPa	<b>E</b>

\*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
<b>Range A</b>	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
<b>Range E</b>	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

## Assured, repeatable accuracy to enable safe, effective weaning.

### What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.

### Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned






## Non-Heated Breathing Circuits

For use with Respiration Trilogy  
22mm cuffs


**30/case**

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
<b>VC1001</b>	1073223 & 1073232	Active	No	Adult	22mm
<b>VC1001P</b>	1073220 & 1073229	Active	Yes	Pediatric	15mm
<b>VC1004</b>	1069443 & 1069444	Passive	Yes	Adult	22mm
<b>VC1004P</b>	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#  
**VC1001** V P  
Active Adult  
OEM#  
1073223  
1073232




Item#  
**VC1001P**  
Active Pediatric  
OEM#  
1073220  
1073229




To use active circuits with  
Respiration Trilogy, simply  
cut ends off pressure and  
exhalation lines.



Item#  
**VC1004** V P  
Passive Adult  
OEM#  
1069443  
1069444



Item#  
**VC1004P** V P  
Passive Pediatric  
OEM#  
1073224  
1073233



## Adult Passive Heated Circuit with Exhalation Port **V**


**30/case**

For use with Philips  
 Respironics Trilogy

Item#  
**VC1003**


**VC1003P**  
 Pediatric Size

Comparable to OEM# 1076046,  
 1073227, 1073235



### Includes:

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port



## Passive Circuit w/Exhalation Port and Bacterial/Viral Filter **V P**

With 22mm cuffs

Item#  
**VC1005**  
**Each**

OEM# 1147378, 1132340, 7615



## Adult Passive Circuit with Bacteria Filter **V**

Comparable to Breas 007936

Item#  
**VC1006**  
**Each**



## Adult Passive Circuit with Bacteria Filter and T-Adapter **V**

Comparable to Breas 007615


Item#  
**VC1007**  
**Each**



## Passive Non-Heated Breathing Circuits **V** **P**

### Exhalation Port



Exhalation Port	Item#	OEM#
Swivel	<b>VC1000</b>	1073228 & 1073236
Fixed	<b>VC1100</b>	1069467 & 1069468



## Passive Heated Breathing Circuit **V** **P**

Item#  
**VC1114**

For use with  
 MR850 Humidifier  
**10/case**



## Non-Heated Dual Limb Breathing Circuits **V** **P**

**30/case**

Item#	Size	Humidifier Limb	Water Trap
<b>VC2001P</b>	Pediatric	No	No
<b>VC2001</b>	Adult	No	No
<b>VC2002P</b>	Pediatric	Yes	No
<b>VC2002</b>	Adult	Yes	No
<b>VC2003</b>	Adult	Yes	Yes





## Adult NIV Patient Circuit **V** **P**



Item# **VC1160**

Compatible with V60, No Filters

**10/case**

Comparable to

OEM# 1069210



## Adult NIV Single Limb Patient Circuit with Main Flow Filter **V** **P**

Item# **VC1161**

Includes 1 BF2200 22mm Bacteria Disk Filter

**10/case**




## Adult NIV Single Limb Patient Circuit with Main Flow Filter and Exhalation Port Filter **V** **P**

Item# **VC1162**

Includes 2 BF2200 22mm Bacteria Disk Filters

**10/case**



## Lynx Venturi Oxygen Delivery V P


### Low Concentration

Item# **RES2400**

Each

Adapter Only **RES240LC**

- Design replaces multiple venturi adapters
- Connects to oxygen tubing and regulates flow
- Provides precision FiO<sub>2</sub> delivery from 24-60% with less than 2% variance during mild, moderate, and severe respiratory distress
- Excellent CO<sub>2</sub> washout



### Unprecedented simplicity

No training required. Intuitive design makes it easy for nurses to operate without assistance, freeing up RTs to work where most needed.

### Do no harm

Quicker de-escalation of O<sub>2</sub> levels and prevention of excessive CO<sub>2</sub> in the patient.

### Reduction in waste

Remains with patient throughout EMS and hospital journey, replacing multiple adapters which may get lost or contaminated during transition.

### High Concentration

Item# **RES2401**


Each

Adapter Only **RES241HC**

- Allows rapid de-escalation of FiO<sub>2</sub> in the EMS ambulance or ED to prevent CO<sub>2</sub> buildup.
- Replaces bulkier high concentration venturis and non-rebreather bag systems
- Accurate assessment of the severity of hypoxemia (PaO<sub>2</sub>/FiO<sub>2</sub>) allows early intervention and better outcomes
- Excellent CO<sub>2</sub> washout

Single flow (flush) with an adjustable window shutter delivers oxygen concentrations of 65% • 80% • 100%

Capabilities could include replacing non-rebreathers and/or partial rebreather products



### Unmatched precision

Provides specific levels of oxygen with less than 2% variance. Comparable products vary by as much as 30%.

## Simple O<sub>2</sub> Mask V P

With 7ft Safety Tubing  
Elongated

Item# **RES2100**  
Adult Size  
50/case




- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

## Medium Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	V	P
<b>RES2100NT</b>	Adult	None	Standard	V	P
<b>RES2100P</b>	Pediatric	7 ft	Standard	V	P
<b>RES2100U</b>	Adult	7 ft	Universal	V	P
<b>RES2203I</b>	Infant	None	Standard	V	
<b>RES2203P</b>	Pediatric	None	Standard	V	P
<b>RES2203A</b>	Adult	None	Standard	V	P
<b>RES2203A7</b>	Adult	7 ft	Standard	V	P
<b>RES2203I7</b>	Infant	7 ft	Standard	V	P
<b>RES2203P7</b>	Pediatric	7 ft	Standard	V	P

- 5-12 LPM
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal



## High Concentration Oxygen Mask

50/case

Item#	Size	Tubing Length	V	P
<b>RES2202A</b>	Adult	None	V	P
<b>RES2202A7</b>	Adult	7 ft	V	
<b>RES2202I</b>	Infant	None	V	P
<b>RES2202I7</b>	Infant	7 ft	V	P
<b>RES2202P</b>	Pediatric	None	V	
<b>RES2202P7</b>	Pediatric	7 ft	V	



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

## Non-Rebreather Oxygen Mask

50/case

Item#	Size	Tubing Length	Connector	Safety Vent	
RES2102	Adult	7 ft	Standard	No	V P
RES2102I	Infant	7 ft	Standard	No	V P
RES2102IM	Infant	7 ft	Standard	No	V
RES2102P	Pediatric	7 ft	Standard	No	V P
RES2102PS	Pediatric	7 ft	Standard	Yes	V P
RES2102PSM	Pediatric	None	Standard	Yes	V P
RES2102PSU	Pediatric	7ft	Universal	Yes	V P
RES2102PU	Pediatric	7 ft	Universal	No	V P
RES2102S	Adult	7 ft	Standard	Yes	V P
RES2102SM	Adult	None	Standard	Yes	V P
RES2102U	Adult	7 ft	Universal	No	V P

- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape



## Adult Air Entrainment Mask

With Flex Tube

50/case

Item#	Size	Kit	Tubing Length	Connector	
RES2400A	Adult	Mask Only	7 ft	Standard	V P
RES2400AK	Adult	Yes	7 ft	Standard	V P
RES2400AKU	Adult	Yes	7 ft	Universal	V P
RES2400AU	Adult	Mask Only	7 ft	Universal	V P
RES2400P	Pediatric	Mask Only	7 ft	Standard	V
RES2400PK	Pediatric	Yes	7 ft	Standard	V P
RES2400PKU	Pediatric	Yes	7 ft	Universal	V
RES2400PU	Pediatric	Mask Only	7 ft	Universal	V P





## Humidoflo High Flow Nasal Cannula **V** **P**

### 10/case

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respironics V60
- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)



Size	Item#
Small	<b>HFN001AS</b>
Medium	<b>HFN001AM</b>
Large	<b>HFN001AL</b>

### With F&P Adapter

Size	Item#
Small	<b>HFN001ASA</b>
Medium	<b>HFN001AMA</b>
Large	<b>HFN001ALA</b>

## Adult Oxygen Cannula

Item#	Tubing Length	Connector	Qty	
<b>RES1100</b>	None	Standard	50	<b>V</b> <b>P</b>
<b>RES1100N</b>	None	Standard	50	<b>V</b>
<b>RES1104</b>	4 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1107</b>	7 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1107N</b>	7 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1107NU</b>	7 ft	Universal	50	<b>V</b>
<b>RES1107U</b>	7 ft	Universal	50	<b>V</b> <b>P</b>
<b>RES1110</b>	10 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1114N</b>	14 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1114U</b>	14 ft	Universal	50	<b>V</b> <b>P</b>
<b>RES1115</b>	15 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1115U</b>	15 ft	Universal	50	<b>V</b> <b>P</b>
<b>RES1125</b>	25 ft	Standard	25	<b>V</b> <b>P</b>
<b>RES1135</b>	35 ft	Standard	20	<b>V</b>
<b>RES1150Z</b>	50 ft	Standard	20	<b>V</b>
<b>RES1150U</b>	50 ft	Universal	20	<b>V</b> <b>P</b>



- 1-6 LPM
- 6-channel supply tubing for safety

## Soft Cannula

Item#	Tubing Length	Connector	Qty	
<b>RES1100S</b>	None	Standard	50	<b>V</b> <b>P</b>
<b>RES1107S</b>	7 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1107SU</b>	7 ft	Universal	50	<b>V</b> <b>P</b>
<b>RES1114S</b>	14 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1114SU</b>	14 ft	Universal	50	<b>V</b> <b>P</b>
<b>RES1115S</b>	15 ft	Standard	50	<b>V</b> <b>P</b>
<b>RES1125S</b>	25 ft	Standard	25	<b>V</b> <b>P</b>
<b>RES1125SU</b>	25 ft	Universal	25	<b>V</b>



- 1-6 LPM
- 6-channel supply tubing for safety

## Pediatric, Infant, and Neonate Cannula

50/case

Item#	Size	Tubing Length	Connector	Touch	Ear Cushions
RES1100P	Pediatric	None	Standard	Standard	No
RES1107I	Infant	7 ft	Standard	Standard	No
RES1107NE	Neonate	7 ft	Standard	Standard	No
RES1107P	Pediatric	7 ft	Standard	Standard	No
RES1107PEC	Pediatric	7 ft	Standard	Standard	Yes
RES1107PECU	Pediatric	7 ft	Universal	Standard	Yes
RES1107SI	Infant	7 ft	Standard	Soft	No
RES1107SNE	Neonate	7 ft	Standard	Soft	No
RES1107SP	Pediatric	7 ft	Standard	Soft	No
RES1107SPU	Pediatric	7 ft	Universal	Soft	No
RES1114P	Pediatric	14 ft	Standard	Standard	No
RES1114PU	Pediatric	14 ft	Universal	Standard	No
RES1114SP	Pediatric	14 ft	Standard	Soft	No
RES1114SPU	Pediatric	14 ft	Universal	Soft	No

- 1-6 LPM
- 6-channel supply tubing for safety

## Adult Cannula with Ear Cushions

Item#	Tubing Length	Connector	Qty
RES1100EC	None	Standard	50
RES1104EC	4 ft	Standard	50
RES1107EC	7 ft	Standard	50
RES1107ECU	7 ft	Universal	50
RES1114EC	14 ft	Standard	50
RES1115EC	15 ft	Standard	50
RES1125EC	25 ft	Standard	25

## Flared Cannula V

Item#	Tubing Length	Connector	Qty
RES1100F	None	Standard	50
RES1107F	7 ft	Standard	50
RES1107FU	7 ft	Universal	50
RES1114F	14 ft	Standard	50
RES1114FU	14 ft	Universal	50
RES1125F	25 ft	Standard	25
RES1150F	50 ft	Standard	20

- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

## High Flow Soft Cannula V

25/case


Item#	Tubing Length
RES1104SHF	4ft
RES1107SHF	7ft
RES1115SHF	15ft
RES1125SHF	25ft



## 7 ft Cannula With CO<sub>2</sub> Sampling Line **V**

50/case

Item#	Size	Connector	Luer	Filter
<b>RES8007</b>	Adult	Standard	Male	No
<b>RES8007F</b>	Adult	Standard	Male	Yes
<b>RES8007I</b>	Infant	Standard	Male	No
<b>RES8007IF</b>	Infant	Standard	Male	Yes
<b>RES8007IU</b>	Infant	Universal	Male	No
<b>RES8007P</b>	Pediatric	Standard	Male	No
<b>RES8007PF</b>	Pediatric	Standard	Male	Yes
<b>RES8007PU</b>	Pediatric	Universal	Male	No
<b>RES8007U</b>	Adult	Universal	Male	No
<b>RES8107</b>	Adult	Standard	Female	No
<b>RES8107F</b>	Adult	Standard	Female	Yes
<b>RES8107I</b>	Infant	Standard	Female	No
<b>RES8107IF</b>	Infant	Standard	Female	Yes
<b>RES8107P</b>	Pediatric	Standard	Female	No
<b>RES8107PF</b>	Pediatric	Standard	Female	Yes




## Secure Cannula Stickers **V** **P**

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



## Ear Cushions for Cannula **V** **P**

Item# **RES020**

Cushions Only

50 Pairs/case



- Add cushions to existing cannula

## Oxygen Supply Tubing

Item#	Length	Connector	Bore	Color	Qty	
RES3007	7 ft	Standard	Standard	None	50	V P
RES3007G	7 ft	Standard	Standard	Green	50	V P
RES3007SB	7 ft	Standard	Smooth	Green	50	V
RES3007U	7 ft	Universal	Standard	None	50	V P
RES3015	15 ft	Standard	Standard	None	50	V P
RES3015G	15 ft	Standard	Standard	Green	50	V
RES3025	25 ft	Standard	Standard	None	25	V P
RES3025G	25 ft	Standard	Standard	Green	25	V P
RES3025U	25 ft	Universal	Standard	None	25	V P
RES3030	30 ft	Standard	Standard	None	20	V
RES3030G	30 ft	Standard	Standard	Green	20	V
RES3040G	40 ft	Standard	Standard	Green	20	V P
RES3050	50 ft	Standard	Standard	None	20	V
RES3050G	50 ft	Standard	Standard	Green	20	V P

- 6-channel
- PHT/DEHP-free



**Universal Connector**

## CPAP Tubing



Item#

**TUB008** 8ft

**TUB010** 10ft

19mm with 22mm cuffs

HCPSC# A7037



- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free





## Firesafe Thermal Shut-Off Valve V

Item#  
**RES010-10**  
Bidirectional  
10/case




- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire
- Smaller
- Laser etched
- Expiration date
- Human readable UDI
- Designed to last 5 years



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

### Fatalities per 100,000 patients



## FireSafe Nozzle V P

Item#  
**RES010FND**  
10/pack



## FireSafe Tubing Connector

Item#  
**RES010TC** 100/case V  
**RES010TC-10** 10/case



## Christmas Tree Connector

Item# **RES001** V P  
Green without Swivel  
10/pack



**RES001-30** P  
30/pack

Item# **RES002** V P  
Green with Swivel  
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

## Oxygen Tubing Connector — Non-Swivel V P

Item# **RES012**  
Fits 3/16in and 1/4in  
ID Tubing  
10/pack



- Barbed ends ensure secure, leak-free connections

## Tubing Reducer Connector V P

Item# **RES011**  
Fits 3/16in, 1/4in,  
5/16in, 3/8in, and 7/16in  
ID Tubing  
10/pack



- Tapered on both ends to connect multiple sizes of tubing, works with suction, oxygen, and catheter tubing

## Humidifier Tubing Adapter V P

Item# **RES016**  
22mm ID x  
22mm OD  
with Lip  
10/pack



## Circuit Connectors V P

Item# **RES013**  
22mm OD x  
15mm ID/  
22mm OD  
10/pack



Item# **RES014**  
22mm OD x  
22mm OD  
With Lip  
10/pack



## Oxygen Tubing Connector — Swivel

Item# **RES018** V P  
Male/Male  
10/pack



**RES018-25** P  
25/pack

## Oxygen Tubing Y-Connector V P

Item# **RES017**  
Fits 3/16in and 1/4in ID Tubing  
(All 3 sides are the same)  
10/pack



- Connects two oxygen sources for higher flows

## Oxygen Bleed-In Adapter V P

Item# **RES019**  
22mm with 6mm  
Rotating Stem  
10/pack



## Finger Pulse Oximeter V P

Item# **RES5102**  
Each



- Oversized 1.3" OLED display face
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- Includes lanyard and manual

## Wrenches P

Item# **RES031**  
With Chain  
10/pack



## Water Traps V P

Item# **RES037**  
Push On End Caps  
Each



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

## Oxygen Liter Meter V P

Item# **RES034**  
0-10 Liter/Minute  
Each




- Oxygen pen
- Works with tanks and concentrators

## Humidifier Connectors V P

Item# **HUM008**  
L-Connector 10in  
10/pk

**HUM008L**  
L-Connector 21in  
10/pk

**HUM008S**  
Straight Connector 9in  
10/pk




## SpO2 Sensors

**NEW**

### Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**


Similar to Nellcor MAX-N  
**20/box**



### Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**


1.5ft cord • Non-replaceable  
**20/box**



### Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**

3ft cord • Non-replaceable  
**20/box**



### Coban Tape Roll

Item# **RES7202CT**

Use for sensor life  
extension or sensor  
cord stability


**Each roll 1" x 177"**



### Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**

1.5ft cord  
**20/box**



### Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)


Item# **RES7204**

1.5ft cord  
**20/box**





## Suction Machine Inline Filters **V** **P**



**Item# BF1412**  
 Final Bacteria Filter •  
 1/4in to 1/2in Stepped  
 Barb on Both Sides  
**10/pack**



**Item# BF1312**  
 1/3in to 1/2in  
**1/pack**



**Item# BF3812**  
 3/8in to 1/2in  
**10/pack**



**Item# BF516TS**  
 5/16in Barbed with  
 Threaded NPT  
**10/pack**

## Oxygen Filters



**Item# BF316**  
 Final Bacteria Filter •  
 3/16 Stepped Barb and  
 1/4in Straight Barb on  
 Both Sides  
**10/pack** **V** **P**



**Item# BF516**  
 Final Bacteria Filter •  
 1/4in Barbed  
**10/pack**  
 OEM# P0750X030 **V** **P**



**Item# BF1438** **P**  
 1/4in to 3/8in  
**10/pack**

## HEPA Filter for 3B Stratus **V** **P**

**Item# BF6020**  
 HEPA Filter  
**Each**  
 OEM# STR1020



## 6in HEPA Capsule with Removable Sticker **V** **P**

**Item# BFH6S**  
 Threaded  
**Each**



## Disposable Bubble Humidifier 500 ml **P**

**Item# RES050**  
 Up to 15 LPM • 6 psi • 500 ml  
**50/case**  
 HCPCS# E0560  
 OEM# HUM001

- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents

- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient



**Extended Life Intake Bacteria Filter** V P

Item# **OF1025D**  
Each

**Foam Cabinet Filter** V P

Item# **OF9001**  
6 3/4" x 33/8" x 1/2"  
**10/pack**  
OEM# 2000489



For Invacare  
Mobilaire and  
HomeFill oxygen  
concentrators

**3B Stratus 5 Style Filter** P

Item# **OF6030**  
Foam Cabinet Filter  
**10/pack**  
OEM# STR1030

**Invacare Platinum 5 and 10 Style Filter** V P

Item# **OF9002**  
Foam Cabinet Filter  
9" x 27/8" x 1/2"  
**10/pack**  
OEM#1107412

**Extended Life Bacteria Filter with Muffler** V P

Item# **OF6011**  
22mm fittings  
Each  
Use if SolAiris was  
manufactured  
before Nov '03



For Drive/DeVilbiss  
SolAiris and  
525 Series oxygen  
concentrators

**DeVilbiss/Drive SolAiris III, 303 Series Style Filter** V P

Item# **OF5004**  
Foam Cabinet Filter  
33/4" x 33/4" x 1/2"  
**10/pack**  
OEM# 30302605

**Invacare Perfecto2 Style Filter** V P

Item# **OF9004**  
Foam Cabinet Filter  
67/8" x 33/4" x 1/2"  
**10/pack**  
OEM# 1143492



## Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units **V P**

Item# **BF2200**  
 22mm Bacteria  
 Disk Filter  
**5/pack**



Attach between the machine and tube for additional protection against bacteria and allergens

### Filter Specifications

Bacterial Filtration Efficiency: (BFE) 99.99+%  
 Viral Filtration Efficiency: (VFE) 99.9+%  
 NaCl Particulate Filtration Efficiency: 96+%  
 Flow Resistance: Approximately 8.8 mmH2O at 30 LPM  
 Internal Volume: 66.3 mL without fittings / 78.3 mL with fittings

### Connections:

- Patient End: 22 mm OD / 15 mm ID
- Machine End: 22 mm ID

## Plastic Intake Filter with Head **V P**

Item# **BF017**  
 For DeVilbiss Nebulizers  
**5/pack**  
 HCPCS# A7013



## Foam Nebulizer Filter **V P**

Item# **BF018**  
**10/pack**  
 HCPCS# A7013



## Plastic Intake Filter for Respirationics InnoSpire **V P**


Item# **BF026S**  
**10/pack**  
 HCPCS# A7013



- Mini
- MicroElite
- MiniElite

## EverFlo Intake Filter **V P**

Item# **BF1001**  
**Each**  
 Respirationics OEM# 1038831



## Bacterial/Viral HME Filter – Adult **V P**

Item# **BF104**



- Bacterial efficiency 99.9999%
- Viral efficiency 99.999%
- Humidification efficiency 37.5 mg/l H2O (500 ml tidal volume)
- Coupled with a hydrophobic bacterial viral filter membrane

## Volumetric Incentive Spirometer V P

12/case

Item#  
**RES055**  
2500ml

**RES056**  
5000ml

**RES055P**  
Pediatric



SCAN FOR VIDEO



- Chambers for ball and volumetric
- Tube
- Mouthpiece

The Volumetric Incentive Spirometer helps patients perform and monitor their own post-surgical breathing exercises. The patients are encouraged to breathe in from the device as slowly and as deeply as possible.

- Adjustable target indicator provides a volumetric guideline for adequate lung exercise
- Dual-sided, clear and legible identification of volumes
- Inspiratory indicator encourages patients to maintain an appropriate rate inspiration
- Compact and ergonomic design

## Three Ball Spirometer

Item#  
**RES057-10** P  
10/pack

**RES057** V P  
12/case

- Improves the patient's lung expansion by training their inhaled lung volume
- Three color coded balls inside the chamber offer incentive for the patient



- Three color-code balls in three chambers
- Chamber to lung capacity: 600cc/sec, 900cc/sec, and 1200cc/sec
- Compact design and break-resistant plastic
- Each set includes incentive spirometer, disposable mouth piece and flexible tube

## Drainage Bag with Y-Adapter V P

Item# **RES071**  
750cc  
50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

## Drainage Bag with T-Adapter V P

Item# **RES070**  
750cc  
50/case




## Aerosol Tubing

Item#  
**TUB100-6**  
Corrugated Aerosol Tubing  
100ft in 6in segments  
Each



**TUB100-6B**  
Blue Corrugated Aerosol Tubing  
100ft in 6in segments  
Each





### ResMed S9 and AirSense 10 Series Filters **V** **P**

**Item# CF2107**

Disposable Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107-1** 1/pack OEM# 36850

**CF2107-2** 2/pack OEM# 36851

**CF2107-50** 50/pack OEM# 36853

**Item# CF2107HA**

Hypoallergenic Filter

**3/pack** • 2 1/8" x 1 3/8"

HCPCS# A7038

**CF2107HA-1** 1/pack OEM# 36855

**CF2107HA-2** 2/pack OEM# 36856

**CF2107HA-50** 50/pack

OEM# 36858



### ResMed AirSense 11 Series Filters **P**

**Item# CF2111**

Disposable Filter

**3/pack**

HCPCS# A7038

**CF2111-1** 1/pack OEM# 39300

**CF2111-2** 2/pack OEM# 39301

**CF2111-50** 50/pack OEM# 39303

**Item# CF2111HA-3**

Hypoallergenic Filter

**3/pack**

1 1/4" x 1 3/16"

HCPCS# A7038

**CF2111HA-1** 1/pack

**CF2111HA-2** 2/pack



### 3B Luna G3 Style Filters **V** **P**

**Item# CF8006-6**

Disposable Filter

**6/pack**

HCPCS# A7038

**CF8006-1** 1/pack

**CF8006-2** 2/pack

**Item# CF8006F-1**

Reusable Filter

**1/pack**

HCPCS# A7039



### 3B Luna and Luna II Style Filters **P**

**Item# CF8005**

Disposable Filter

**12/pack** • 2" x 1 3/16"

HCPCS# A7038

**Item# CF8005F-1**

Reusable Filter

**1/pack** • 2" x 1 3/16"

HCPCS# A7039

**CF8005F-2** 2/pack



### 3B RESmart Style Filter **P**

**Item#**

**CF8004**

Disposable Filter

**3/pack** • 3 5/16" x 1 5/8"

HCPCS# A7038



## Respironics DreamStation 2 Series Filters V P

### Item# **CF1012**

Disposable Filter

**6/pack**

OEM# 1142831

HCPCS# A7038

**CF1012-1** 1/pack OEM# 1142828

**CF1012-2** 2/pack OEM# 1142829

**CF1012-3** 3/pack

**CF1012-50** 50/pack



### Item# **CF1012F**

Reusable Filter

**2/pack**

OEM# 1142687

HCPCS# A7039

**CF1012F-1** 1/pack



## Respironics DreamStation Series Filters V P

### Item# **CF1010**

Disposable Filter

**6/pack**

OEM# 1122519

HCPCS# A7038

**CF1010-1** 1/pack

**CF1010-2** 2/pack



### Item# **CF1010F**

Reusable Filter

**2/pack**

OEM# 112446 (1-pack)

HCPCS# A7039

**CF1010F-1** 1/pack

**CF1010F-6** 6/pack



## 90-Degree CPAP Tubing Elbow V P

Item# **CAP1027**  
Each

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90-degree angle, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

Comparable to OEM# 37394



## Color Coded Dual Channel Berman Oral Airway V P

50/pack

Item#	Size	Color
<b>AOA40BC-50</b>	40mm	Pink
<b>AOA50BC-50</b>	50mm	Blue
<b>AOA60BC-50</b>	60mm	Black
<b>AOA70BC-50</b>	70mm	White
<b>AOA80BC-50</b>	80mm	Green
<b>AOA90BC-50</b>	90mm	Yellow
<b>AOA100BC-50</b>	100mm	Red
<b>AOA110BC-50</b>	110mm	Orange



- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## White Dual Channel Berman Oral Airway V P

10/pack

Item#	Size
<b>AOA40BD</b>	40mm
<b>AOA60BD</b>	60mm
<b>AOA80BD</b>	80mm
<b>AOA90BD</b>	90mm
<b>AOA100BD</b>	100mm



## Oropharyngeal Guedel Airway V P

10/pack

Item#	Size	Color
<b>AOP500GC</b>	50mm	Blue
<b>AOP601GC</b>	60mm	Black
<b>AOP702GC</b>	70mm	White
<b>AOP803GC</b>	80mm	Green
<b>AOP904GC</b>	90mm	Yellow
<b>AOP1005GC</b>	100mm	Red



- Designed to overcome associated risks of loose or detached bite blocks
- Smooth surface & rounded tip for improved patient comfort
- 100% Polyethylene (PE)
- Latex free

## Nasopharyngeal Airway V

PVC

10/pack

Item#	Size	Length
<b>ANP1472PNE</b>	14fr	72mm
<b>ANP1685PNE</b>	16fr	85mm
<b>ANP18105PNE</b>	18fr	105mm
<b>ANP20115PNE</b>	20fr	115mm
<b>ANP22125PNE</b>	22fr	125mm
<b>ANP24130PNE</b>	24fr	130mm
<b>ANP26140PNE</b>	26fr	140mm
<b>ANP28155PNE</b>	28fr	155mm
<b>ANP30165PNE</b>	30fr	165mm
<b>ANP32170PNE</b>	32fr	170mm
<b>ANP34175PNE</b>	34fr	175mm



## Flex HME **V** **P**






Item# **RES027FLEX**  
 40/case



<b>TV Range (mL)</b>	300-1500
<b>Moisture Output (mg H2O/L)</b>	33 mg/L
<b>Weight (g)</b>	17
<b>Resistance (cm H2O)</b>	0.5 @ 30 LPM 1.34 @ 60 LPM
<b>Dead Space (mL)</b>	75
<b>Connection (mm)</b>	22M,15F

## Heat Moisture Exchangers **V** **P**

30/case

Item#					
<b>Item#</b>	<b>RES027</b> T-HME	<b>RES027OP</b> T-HME with Oxygen Port	<b>RES027IL</b> Inline HME — Large Volume	<b>RES027NE</b> Inline HME — Neonatal	<b>RES027HV2</b> Foam HME — Adult/Child
<b>TV Range (mL)</b>	250-1000	150-1250	50-250	15-50	250-1000
<b>Moisture Output (mg H2O/L)</b>	30, Vt = 0.1L	36.3 mg/L	31, Vt = 0.1L	30	29
<b>Weight (g)</b>	4.5	7.4	10	4.5	21
<b>Resistance (cm H2O)</b>	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.2 @ 30 LPM	0.4 @ 10 LPM	0.5 @ 10 LPM	0.52 @ 30 LPM 1.32 @ 60 LPM 2.54 @ 90 LPM
<b>Dead Space (mL)</b>	12.5	7	13	2.5	29
<b>Connection (mm)</b>	15F	15F 5.5mm O2 port	15M,22M/15F	15M,15F	22F/15M, 22M/15F

## Tracheostomy Mask with O<sub>2</sub> Connector **V** **P**

50/case

Item#  
**RES2130**  
 Adult Regular Size

**RES2130L**  
 Adult Large Size

**RES2230**  
 Pediatric Size



- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O<sub>2</sub> supply tubing

## Adapter for T-HME

Item#  
**RES027A** **V** **P**  
 Each  
**RES027A-30** 30/pack **P**



## Flexible Connector Length 4.9-7.5cm

50/case

Item#  
**TRAA004** **V** **P**  
 Adult — 22/15mm x 15mm  
**TRAA004P**  
 Pediatric — 15mm F x 15mm M





## Suction Catheter Kits P

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK305</b>
6fr	Y	DeLee	Pediatric	<b>TRATK306</b>
8fr	Y	DeLee	Pediatric	<b>TRATK308</b>
8fr	Y	Whistle	Adult	<b>TRATK108</b>
8fr	Chimney	Whistle	Adult	<b>TRATK108CV</b>
10fr	Y	Whistle	Adult	<b>TRATK110</b>
10fr	Chimney	Whistle	Adult	<b>TRATK110CV</b>
12fr	Y	Whistle	Adult	<b>TRATK112</b>
12fr	Chimney	Whistle	Adult	<b>TRATK112CV</b>
14fr	Y	Whistle	Adult	<b>TRATK114</b>
14fr	Chimney	Whistle	Adult	<b>TRATK114CV</b>
16fr	Y	Whistle	Adult	<b>TRATK116</b>
16fr	Chimney	Whistle	Adult	<b>TRATK116CV</b>



- Sterile
- Includes Catheter and Gloves

## Suction Catheter Kits with Cup P

50/case

Size	Valve Type	Tip	Intended Use	Item#
5fr	Y	DeLee	Pediatric	<b>TRATK405</b>
6fr	Y	DeLee	Pediatric	<b>TRATK406</b>
8fr	Y	DeLee	Pediatric	<b>TRATK408</b>
8fr	Y	Whistle	Adult	<b>TRATK208</b>
8fr	Chimney	Whistle	Adult	<b>TRATK208CV</b>
10fr	Y	Whistle	Adult	<b>TRATK210</b>
10fr	Chimney	Whistle	Adult	<b>TRATK210CV</b>
12fr	Y	Whistle	Adult	<b>TRATK212</b>
12fr	Chimney	Whistle	Adult	<b>TRATK212CV</b>
14fr	Y	Whistle	Adult	<b>TRATK214</b>
14fr	Chimney	Whistle	Adult	<b>TRATK214CV</b>
16fr	Y	Whistle	Adult	<b>TRATK216</b>
16fr	Chimney	Whistle	Adult	<b>TRATK216CV</b>



- Sterile
- Includes Catheter and Gloves

## Suction Catheters P

Size	Valve Type	Tip	Intended Use	Packaged	Item#	Qty
5fr	Y	DeLee	Pediatric	Coiled	<b>TRASC305C</b>	50
5fr	Y	DeLee	Pediatric	Straight	<b>TRASC305S</b>	100
6fr	Y	DeLee	Pediatric	Coiled	<b>TRASC306C</b>	50
6fr	Y	DeLee	Pediatric	Straight	<b>TRASC306S</b>	100
8fr	Y	DeLee	Pediatric	Coiled	<b>TRASC308C</b>	50
8fr	Y	DeLee	Pediatric	Straight	<b>TRASC308S</b>	100
10fr	Y	Whistle	Adult	Coiled	<b>TRASC110C</b>	50
10fr	Y	Whistle	Adult	Straight	<b>TRASC110S</b>	100
12fr	Y	Whistle	Adult	Coiled	<b>TRASC112C</b>	50
12fr	Y	Whistle	Adult	Straight	<b>TRASC112S</b>	100
14fr	Y	Whistle	Adult	Coiled	<b>TRASC114C</b>	50
14fr	Y	Whistle	Adult	Straight	<b>TRASC114S</b>	100
16fr	Y	Whistle	Adult	Coiled	<b>TRASC116C</b>	50
16fr	Y	Whistle	Adult	Straight	<b>TRASC116S</b>	100



- Sterile
- Single patient use
- Latex-free

## Tracheostomy Ties with Fastener Tabs V P

20/case


One-piece

- Easy on

Item#

**RES029** Adult

**RES029P** Pediatric



- Designed to provide secure stabilization of the tracheostomy tube
- Made from soft foam and cotton


Two-piece

- Easy off

Item#

**RES030** Adult

**RES030P** Pediatric



## Tracheostomy Swivel Elbow Connectors V P

Item# **RES021**

With suction port

15mm ID x 15mm OD

10/pack



## Disposable Rigid Suction Canister

Item#	Size	Direct to Regulator	Qty
<b>RES023S-24</b>	800cc	No	24
<b>RES025AW1200R</b>	1200cc	Yes	40
<b>RES025AW1200-40</b>	1200cc	No	40
<b>RES025AW2000</b>	2000cc	No	40
<b>RES025AW3000</b>	3000cc	No	40



- Prevents backflow and premature shut off.
- Includes 20-degree patient port for 1200cc, or straight port for 800c, to prevent tubing from kinking.
- Durability assurance for collecting and retaining aspirated fluids.
- Disposable.

## Flexible Disposable Liner

50/case

Item#	Size	Includes Solidifier
<b>LIN1000FD</b>	1000cc	No
<b>LIN1500FD</b>	1500cc	No
<b>LIN3000FD</b>	3000cc	No
<b>LIN1000FDS</b>	1000cc	Yes
<b>LIN1500FDS</b>	1500cc	Yes
<b>LIN3000FDS</b>	3000cc	Yes



- Part of a complete fluid management system.
- These soft, plastic liners fit into the reusable cannisters, creating less plastic waste and greater material utilization.
- Not made with natural rubber latex.

## Reusable Outer Canister for Disposable Liner with Built In On/Off Valve

10/case

### For Flexible Disposable Liner

Item#	Size
<b>LIN1000FDC</b>	1000cc
<b>LIN1500FDC</b>	1500cc
<b>LIN3000FDC</b>	3000cc

### For Semi-Rigid Disposable Liner

Item#	Size
<b>LIN1000SRDC</b>	1000cc
<b>LIN1500SRDC</b>	1500cc
<b>LIN3000SRDC</b>	3000cc



- Reusable outer cannisters provide a safe, convenient, and economical solution for fluid collection.
- Clear plastic cannisters are designed to hold the disposable liners.
- Reusable polycarbonate canister/bracket assembly has pre-attached tubing for connection to the liner.

### Semi-Rigid Disposable Liner

Item#	Size
<b>LIN1000SRD</b>	1000cc
<b>LIN1500SRD</b>	1500cc
<b>LIN3000SRD</b>	3000cc

## Canister/Liner Kits

Item#	Size and Type	Tubing Size and Length	Qty
<b>LIN1000FD618</b>	1000cc Flexible Liner	6mm x 1.8m	50
<b>LIN1000FD627</b>	1000cc Flexible Liner	6mm x 2.7m	50
<b>LIN1500FD618</b>	1500cc Flexible Liner	6mm x 1.8m	50
<b>LIN1500FD627</b>	1500cc Flexible Liner	6mm x 2.7m	50
<b>LIN3000FD618</b>	3000cc Flexible Liner	6mm x 1.8m	50
<b>LIN1000SRD618</b>	1000cc Semi-Rigid Liner Kit	6mm x 1.8m	72
<b>LIN1000SRD627</b>	1000cc Semi-Rigid Liner Kit	6mm x 2.7m	72
<b>LIN1500SRD518</b>	1500cc Semi-Rigid Liner Kit	5mm x 1.8m	72
<b>LIN1500SRD618</b>	1500cc Semi-Rigid Liner Kit	6mm x 1.8m	72
<b>LIN1500SRD627</b>	1500cc Semi-Rigid Liner Kit	6mm x 2.7m	72
<b>LIN3000SRD618</b>	3000cc Semi-Rigid Liner Kit	6mm x 1.8m	36
<b>RES023S-T</b>	800cc Disposable Suction Canister Kit	6mm x 1.8m	45
<b>RES025AW1200-T518</b>	1200cc Disposable Suction Canister Kit	5mm x 1.8m	30
<b>RES025AW1200-T618</b>	1200cc Disposable Suction Canister Kit	6mm x 1.8m	30
<b>RES025AW1200-T518V</b>	1200cc Disposable Suction Canister Kit	5mm x 1.8m 5mm x .46m Vacuum Tubing	30




**Four-Canister Roll Stand**

Item# **RES9100**  
30"  
Each


**Wall Plate**

Item# **RES9101**  
10/case

**Disposable Suction Canister Ring Bracket**

20/case

Item#	Size
<b>RES9102</b>	1200cc
<b>RES9103</b>	2000-3000cc

**Disposable Suction Canister Ring Bracket**

20/case

Item#	Size
<b>RES9102T</b>	1200cc
<b>RES9103T</b>	2000cc

**Disposable Suction Canister Ring Bracket Extended for Regulator with DISS**

20/case

Item#	Size
<b>RES9102D</b>	1200cc
<b>RES9103D</b>	2000-3000cc

**Disposable Suction Canister Ring Bracket with On/Off Valve**

20/case

Item#	Size
<b>RES9102V</b>	1200cc
<b>RES9103V</b>	2000-3000cc

**Semi-Rigid Suction Canister Ring Bracket Close to Wall**

Each

Item#	Size
<b>RES9104</b>	1000cc
<b>RES9105</b>	1500cc

**Semi-Rigid Suction Canister Ring Bracket Extended For Regulator**

Item# **RES9105E**  
1500cc  
Each





## Suction Tubing — Single Use Sterile

50/case


Item#	Diameter	Length
<b>RES025-31618</b>	5mm	46cm (18in)
<b>RES025-31672</b>	5mm	1.8m (6ft)
<b>RES025-316120</b>	5mm	3.0m (10ft)
<b>RES025-316144</b>	5mm	3.7m (12ft)
<b>RES025-316240</b>	5mm	6.1m (20ft)
<b>RES025-14072</b>	6mm	1.8m (6ft)
<b>RES025-14120</b>	6mm	3.0m (10ft)
<b>RES025-14144</b>	6mm	3.7m (12ft)
<b>RES025-14240</b>	6mm	6.1m (20ft)
<b>RES025-93272</b>	7mm	1.8m (6ft)
<b>RES025-932120</b>	7mm	3.0m (10ft)
<b>RES025-932144</b>	7mm	3.7m (12ft)
<b>RES025-932240</b>	7mm	6.1m (20ft)
<b>RES025-38120</b>	Large Bore 9.5mm	3.0m (10ft)



## Yankauer — Single Use Sterile

100/case

Item#	Tip	Vented
<b>TRAA002S</b>	Bulbous	No
<b>TRAA001S</b>	Bulbous	Yes
<b>TRAA004S</b>	Straight	No
<b>TRAA003S</b>	Straight	Yes





- These lightweight yankauers are designed to provide optimal control and visibility while helping to minimize tissue trauma.

## Bulk Roll — Non-Sterile

4/case

Item#	Diameter	Length
<b>RES025B-530</b>	5mm	30m
<b>RES025B-630</b>	6mm	30m
<b>RES025B-730</b>	7mm	30m

- Clear bulk tubing in rolls that can be cut into the exact length needed.
- Cut-to-size female connectors are located at 6ft intervals.
- Disposable accessories make it easy to add auxiliary items to your system.
- Flex Advantage liners, Guardian systems, Medi-Vac accessories are interchangeable to help you select the most effective setup for your needs.
- Allows you to adapt your system for various setups.
- Used to connect vacuum source to canisters and manifolds.



## Solidifier — Dissolvable Sachet

100/case

Item#	Size
<b>RES8001</b>	150cc
<b>RES8004</b>	450cc
<b>RES8015</b>	1500cc



### Economical

Pre-measured rupturing pouches eliminate the tearing and spillage associated with conventional solidifiers for medical waste.

### Particulate-Free

Self-activating technology reduces particulates and bacterial exposure, and eliminates harmful dust.

### Reduce Exposure

Helps to minimize the risk of cross contamination of infectious waste to medical staff.

## Solidifier with Treatment — Dissolvable Sachet

100/case

Item#	Size
<b>RES8002T</b>	250cc
<b>RES8005T</b>	500cc
<b>RES8007T</b>	750cc
<b>RES8015T</b>	1500cc



### Contain spills quickly with this astonishingly absorbent liquid solidifier packet

Super absorbent packets quickly absorb up to 400 times their own weight in aqueous solutions. These drop-in liquid absorbers offer the ideal solution for shipping, storing, and handling of aqueous solutions that must be contained quickly and effectively in the event of a leak.

### Smart Design

Specially formulated pouch bursts open in the presence of a wide range of fluids to quickly immobilize them. Solidified material is easily contained by an outer vessel, shipping container, or other such device, resulting in safer handling.

### Ideal for Medical Use

The most common applications for the use include medical shipping containers and bio-specimen transportation bags used for vials, test kits, and lab trays.

## SperaChamber **V** **P**

Anti-Static Valved Inhaler Chamber  
for Metered Dose Inhaler

Item# **RES195-10**  
10/case

Item# **RES195-50**  
50/case

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free



Adult Mask



### SperaChamber with Mask 10/case

Item#	Size
<b>RES195MS</b>	Small (Infant)
<b>RES195MM</b>	Medium (Pediatric)
<b>RES195ML</b>	Large (Adult)

### Mask Only

Item#	Size	Qty
<b>RES195MASK-S</b>	Small (Infant)	Each
<b>RES195MASK-S-10</b>	Small (Infant)	10/box
<b>RES195MASK-M</b>	Medium (Pediatric)	Each
<b>RES195MASK-M-10</b>	Medium (Pediatric)	10/box
<b>RES195MASK-L</b>	Large (Adult)	Each
<b>RES195MASK-L-10</b>	Large (Adult)	10/box

## Peak Flow Meter **V** **P**

Item# **RES052**  
10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD** **V** **P**  
Pediatric Mouthpiece  
for Peak Flow Meter  
10/case

## Valved T-Adapter, Aerosol **V** **P**

Item# **NEB9001T**  
22mm ID x 22mm OD  
10/case



## Valved T-Adapter — Adult **V** **P**

Item# **VEN001**  
22F-18F-22M/15F  
Each



## Compressor Nebulizers V P

Item#

**NEB100**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Foam Filters (BF018 x5)

**NEB100-R**

- Disposable Nebulizer Kit with Angled Mouthpiece (RES096)
- Reusable T-Piece Kit (RES093)
- Foam Filters (BF018 x5)

SCAN FOR VIDEO

**Sunny the Bear**Item# **NEB300BEAR****Gordon the Fish**Item# **NEB300FISH**

**NEB100-BAG** V P  
Carry Bag

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

## Handheld Compressor Nebulizer V P

Item#

**NEB200****Included:**

- Compressor
- Medication cup
- AC power adapter
- Air hose
- Elbow
- Mouthpiece
- Mask
- 4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required



## Mini Mesh Nebulizer with Rechargeable Battery V P

Item#

**NEB401****Included:**

- Main unit
- Medicine cup
- Disposable mouthpiece
- Disposable adult mask
- Disposable pediatric mask
- USB cord
- Carrying bag
- Instruction manual



**Small particle size**  
MMAD 2.2µm\*

**Efficient**

Low Residual Liquid

**Convenient**

Small and Portable

**Quiet**  
≤50dB (at one meter)

**Modern**

Rechargeable Lithium Battery

## Disposable Nebulizer Kit

Item#

**RES090** **V** **P**

Adult Size

**50/case**

- 6cc jet nebulizer bottle
- Mask
- 7ft supply tube

**RES090A** **V** **P**

No Accessories

**50/case**

**RES090NM** **P**

- 6cc jet nebulizer bottle
- 7ft supply tube

**50/case**

**RES092** **V** **P**

Pediatric Size

**50/case**

- Pediatric aerosol mask
- 6cc jet nebulizer bottle
- 7ft supply tube



## Disposable Nebulizer Kit with T-Piece **V** **P**

Item#

**RES091**

Disposable  
Single Patient

**50/case**

**RES091NA**

Disposable  
Single Patient,  
No Aerosol

**50/case**

HCPCS# A7003

- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

### Includes:

- 6cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube



## Reusable Nebulizer Kit with T-Piece **P**

Item# **RES093**

Reusable — Single Patient

**10/case**

HCPCS# A7005

### Includes:

- 10cc jet nebulizer bottle
- T-piece with mouthpiece and corrugated reservoir tube
- 7ft supply tube
- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°



## Disposable Nebulizer Kit with Angled Mouthpiece **V** **P**

Item# **RES096**

Disposable  
Single Patient

**50/case**

HCPCS# A7003

### Includes:

- 6cc jet nebulizer bottle
- Angled mouthpiece
- 7ft supply tube



## Aerosol Mask

**50/case**

Item#

**RES2110** **V** **P**

Adult Size

**RES2110F** **P**

Adult Size  
w/ Flex Tube

**RES2210** **V** **P**

Pediatric Size

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged





# Closed Ventilation Suction System

## Standard Switch

Sterile and latex-free

### 24 Hour Adult

#### Endotracheal 530mm With Swivel Connector

Item#	Description
VEN72453010M	10Fr, MDI Port
VEN72453012M	12Fr, MDI Port
VEN72453014M	14Fr, MDI Port
VEN72453016M	16Fr, MDI Port

#### Endotracheal 530mm With T-Piece

Item#	Description
VEN82453010M	10Fr, MDI Port
VEN82453012M	12Fr, MDI Port
VEN82453014M	14Fr, MDI Port
VEN82453016M	16Fr, MDI Port

#### Tracheostomy With Swivel Connector

Item#	Description
VEN72430510M	305mm, 10Fr, MDI Port
VEN72430514M	305mm, 14Fr, MDI Port
VEN72430516M	305mm, 16Fr, MDI Port

#### Tracheostomy With T-Piece

Item#	Description
VEN82430510M	305mm, 10Fr, MDI Port
VEN82430512	305mm, 12Fr
VEN82430512M	305mm, 12Fr, MDI Port
VEN82430514	305mm, 14Fr
VEN82430514M	305mm, 14Fr, MDI Port
VEN82430516	305mm, 16Fr
VEN82430516M	305mm, 16Fr, MDI Port
VEN82440512M	405mm, 12Fr, MDI Port
VEN82440514M	405mm, 14Fr, MDI Port
VEN82440516M	405mm, 16Fr, MDI Port

### 72 Hour Adult

#### Endotracheal 530mm With Swivel Connector

Item#	Description
VEN77253010M	10Fr, MDI Port
VEN77253012M	12Fr, MDI Port
VEN77253014M	14Fr, MDI Port
VEN77253016M	16Fr, MDI Port

#### Tracheostomy 305mm With Swivel Connector

Item#	Description
VEN77230510M	10Fr, MDI Port
VEN77230512M	12Fr, MDI Port
VEN77230514M	14Fr, MDI Port
VEN77230516M	16Fr, MDI Port

### 24 Hour Neonatal/Pediatric

#### With Y-Adapter

Item#	Description
VEN66000500	305mm, 5Fr
VEN66000500M	305mm, 5Fr, MDI Port
VEN66000600	305mm, 6Fr
VEN66000600M	305mm, 6Fr, MDI Port
VEN66000700	305mm, 7Fr
VEN66000700M	305mm, 7Fr, MDI Port
VEN66000800	305mm, 8Fr
VEN66000800M	305mm, 8Fr, MDI Port
VEN66001000	405mm, 10Fr
VEN66001000M	405mm, 10Fr, MDI Port

#### With Swivel Connector

Item#	Description
VEN72430505	305mm, 5Fr
VEN72430506	305mm, 6Fr
VEN72430507	305mm, 7Fr
VEN72430508	305mm, 8Fr

### 72 Hour Neonatal/Pediatric

#### With Y-Adapter

Item#	Description
VEN77230505	305mm, 5Fr
VEN77230506	305mm, 6Fr
VEN77230507	305mm, 7Fr
VEN77230508	305mm, 8Fr



## Closed Ventilation Suction System

# Deluxe Switch (one-handed operation) P

Sterile and latex-free

### Suction Control Switch

Ergonomic suction control switch with a built-in security lock to visually determine if the device is in the OFF/ON position. Offers a quick and easy, single-handed operation design, which allows you to keep your eye on the patient and do no harm. One-handed operation allows the free hand to support the ET or Trach device to ensure less trauma to the site.

### Irrigation Port

Offers a one-way valve to prevent backflow leakage and easy access for cleaning the closed suction system.

### MDI Port

Provides access for drug delivery as needed. Large and easy to grip port allows for rapid medication introduction.

### Catheter Tip

Atraumatic rounded catheter tip for minimal tracheal wall injury.


### Isolation Valve\*

Seals the irrigation chamber to provide a better cleaning of the catheter. Allows for continuous ventilation and prevents backflow of saline during cleaning to ensure patient comfort and safety.

\*Only available with 72 Hour catheter

### Black Ring Indicator

Ensures that the catheter is properly withdrawn from the ET or Trach device. Black ring must be visible in the plastic sheath after suctioning the patient.



# Deluxe Switch (One-handed operation) P Sterile and latex-free

## 24 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN160100A	10Fr
VEN160100AM	10Fr, MDI Port
VEN160120	12Fr
VEN160120M	12Fr, MDI Port
VEN160140	14Fr
VEN160140M	14Fr, MDI Port
VEN160160	16Fr
VEN160160M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN161100A	10Fr
VEN161100AM	10Fr, MDI Port
VEN161120	12Fr
VEN161140	14Fr
VEN161160	16Fr
VEN161160M	16Fr, MDI Port

## 72 Hour Adult

### Endotracheal 540mm

Item#	Description
VEN365100A	10Fr
VEN365100AM	10Fr, MDI Port
VEN365120	12Fr
VEN365120M	12Fr, MDI Port
VEN365140	14Fr
VEN365140M	14Fr, MDI Port
VEN365160	16Fr
VEN365160M	16Fr, MDI Port

### Tracheostomy 305mm

Item#	Description
VEN366100A	10Fr
VEN366100AM	10Fr, MDI Port
VEN366120	12Fr
VEN366120M	12Fr, MDI Port
VEN366140	14Fr
VEN366140M	14Fr, MDI Port
VEN366160	16Fr
VEN366160M	16Fr, MDI Port

## Neonatal/Pediatric

Item#	Description
VEN26000500	305mm, 5Fr
VEN26000500M	305mm, 5Fr, MDI Port
VEN26000600	305mm, 6Fr
VEN26000600M	305mm, 6Fr, MDI Port
VEN26000700	305mm, 7Fr
VEN26000700M	305mm, 7Fr, MDI Port
VEN26000800	305mm, 8Fr
VEN26000800M	305mm, 8Fr, MDI Port
VEN26001000	405mm, 10Fr
VEN26001000M	405mm, 10Fr, MDI Port



## Adult Resuscitator P

12/case

Item#	Pop-Off	Mask
<b>RES8100</b>	No	None
<b>RES8100H</b>	Yes	None
<b>RES8101</b>	Yes	#4
<b>RES8102</b>	Yes	#5



## Child Resuscitator P

12/case

Item#	Pop-Off	Mask
<b>RES8110H</b>	Yes	None
<b>RES8111</b>	Yes	#2
<b>RES8112</b>	Yes	#3



## Infant Resuscitator P

12/case

Item#	Pop-Off	Mask
<b>RES8120H</b>	Yes	None
<b>RES8121</b>	Yes	#0
<b>RES8122</b>	Yes	#1



- Single patient use
- Set includes PVC reservoir bag, PVC mask, and disposable oxygen tubing
- 100% Medical grade PVC material
- Latex-free components
- All-in-one intake valve connects to oxygen reservoir bag directly
- Diverter is able to attach PEEP valve

### Disposable Pressure Manometer P

Item# **RES8201**  
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

### Disposable Pressure Manometer, L-Shaped P

Item# **RES8201L**  
25/case



- Pressure range 0-60 cm H2O
- Indicator strip
- All models have the pressure range emboss printed on the housing
- Some users like to also have a color strip / indicator

### Reusable PEEP Valve, Orange P

2-10cm H2O  
10/case

Item#  
**RES8203R**  
15mm ID



**RES8203RL**  
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

### Reusable PEEP Valve, Blue P

5-20cm H2O  
10/case

Item#  
**RES8204B**  
15mm ID



**RES8204BL**  
30mm ID

- Available in orange color in 2-10cm H2O and blue color in 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Materials: Polycarbonate, silicone, stainless steel for spring

### Disposable PEEP Valve, Orange P

2.5-10cm H2O  
10/case

Item#  
**RES8205R**  
15mm ID



**RES8205RL**  
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel

### Disposable PEEP Valve, Blue P

5-20cm H2O  
10/case

Item#  
**RES8205B**  
15mm ID



**RES8205BL**  
30mm ID

- Available in 2.5-10cm H2O and 5-20cm H2O types with choice of 30mm or 15mm ID
- Adjustable from 0-20cm water pressure
- Single patient use
- Materials: Polycarbonate, silicone, stainless steel



**Reservoir Bag, PVC** P**25/case**

Item#	Size
<b>RES8300P</b>	0.5L
<b>RES8301P</b>	1L
<b>RES8302P</b>	2L

**Reservoir Bag, Silicone** P**25/case**

Item#	Size
<b>RES8300S</b>	0.5L
<b>RES8301S</b>	1L
<b>RES8302S</b>	2L



Suitable for reusable resuscitator.

**CPR Pocket Resuscitator Set** PItem# **RES8200****12/case**

- High transparency for better visualization
- Non-reusable, single use only
- 100% Medical level PVC material
- Latex-free components
- Economical choice
- Blow molding

**Breathing Mask** P**20/case**

Item#	Size	Hooks
<b>SBC01N</b>	Neonatal (#0)	No
<b>SBC01I</b>	Infant (#1)	No
<b>SBC01P</b>	Pediatric (#2)	No
<b>SBC01S</b>	Adult Small (#3)	No
<b>SBC01M</b>	Adult Medium (#4)	No
<b>SBC01L</b>	Adult Large (#5)	No
<b>SBC01XL</b>	Adult Extra Large (#6)	No
<b>SBC02N</b>	Neonatal (#0)	Yes
<b>SBC02I</b>	Infant (#1)	Yes
<b>SBC02P</b>	Pediatric (#2)	Yes
<b>SBC02S</b>	Adult Small (#3)	Yes
<b>SBC02M</b>	Adult Medium (#4)	Yes
<b>SBC02L</b>	Adult Large (#5)	Yes
<b>SBC02XL</b>	Adult Extra Large (#6)	Yes



## Silicone Disposable Laryngeal Masks P

Item#	Size
<b>ALM1001</b>	#1
<b>ALM1015</b>	#1.5
<b>ALM1002</b>	#2
<b>ALM1025</b>	#2.5
<b>ALM1003</b>	#3
<b>ALM1004</b>	#4
<b>ALM1005</b>	#5

- Medical grade silicone for super compatibility and patient comfort
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters
- Available with seven sizes for all patient use



## Silicone Reusable Laryngeal Masks P

Item#	Size
<b>ALM1001R</b>	#1
<b>ALM1015R</b>	#1.5
<b>ALM1002R</b>	#2
<b>ALM1025R</b>	#2.5
<b>ALM1003R</b>	#3
<b>ALM1004R</b>	#4
<b>ALM1005R</b>	#5

- Reusable Silicone Laryngeal Masks are designed to be sterilized and used up to 40 times
- Durable high temperature resistant material suitable for autoclave sterilization
- Serial no. is printed on each connector for easy tracking
- Additional user card for record of device usage times and patient information is supplied in each pouch



## Flexible Silicone Disposable Laryngeal Masks P

Item#	Size
<b>ALM3001</b>	#1
<b>ALM3015</b>	#1.5
<b>ALM3002</b>	#2
<b>ALM3025</b>	#2.5
<b>ALM3003</b>	#3
<b>ALM3004</b>	#4
<b>ALM3005</b>	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



## Flexible Silicone Reusable Laryngeal Masks P

Item#	Size
<b>ALM3001R</b>	#1
<b>ALM3015R</b>	#1.5
<b>ALM3002R</b>	#2
<b>ALM3025R</b>	#2.5
<b>ALM3003R</b>	#3
<b>ALM3004R</b>	#4
<b>ALM3005R</b>	#5

- Flexible Laryngeal Mask is specially designed for use in ophthalmic, ENT, dental, oro-maxillary, and other head and neck surgeries
- Wire reinforced tube reduces risk of kinking and airway tube occlusion
- Flexible tube allows for placement away from surgical field and improves surgical access while maintaining good airtightness



## PVC Disposable Laryngeal Masks P

Item#	Size
<b>ALM2001</b>	#1
<b>ALM2015</b>	#1.5
<b>ALM2002</b>	#2
<b>ALM2025</b>	#2.5
<b>ALM2003</b>	#3
<b>ALM2004</b>	#4
<b>ALM2005</b>	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



## Flexible PVC Disposable Laryngeal Masks P

Item#	Size
<b>ALM4001</b>	#1
<b>ALM4015</b>	#1.5
<b>ALM4002</b>	#2
<b>ALM4025</b>	#2.5
<b>ALM4003</b>	#3
<b>ALM4004</b>	#4
<b>ALM4005</b>	#5

- Medical grade PVC with good compatibility
- No epiglottis bars for easy insertion of fiberoptic scope, introducers and suction catheters



### Size Patient Weight and Age

#1	Up to 5 kg	Pediatric
#1.5	5-10 kg	Pediatric
#2	10-20 kg	Pediatric
#2.5	20-30 kg	Pediatric
#3	30-50 kg	Pediatric
#4	50-70 kg	Adult
#5	70-100 kg	Adult

## PeDIA

### Pediatric Device for Induction of Anesthesia

Item#	Size
<b>ANEP001B</b>	Blue
<b>ANEP001Y</b>	Yellow
<b>ANEP001P</b>	Purple

*PeDIA resolves the problems of traditional mask induction while fostering a sense of trust and patient satisfaction.*

For children, the process of undergoing surgery is already unsettling, and the process of mask-based induction is often the most traumatizing part. With the introduction of the revolutionary Pediatric Device for Induction of Anesthesia device, or “PeDIA” for short, that stress and anxiety can be replaced with something colorful and fun.

As the world’s first mask-free system and method for inhalation induction, the PeDIA provides a solution to this problem that anesthesia providers have been researching for decades. It replaces the mask with a colorful balloon, which actually allows the child to put themselves to sleep rather than being forced to submit to an unfamiliar mask.







