



Homecare Product Catalog

Policies

Ordering Options

Three Easy Ways To Order

Phone

(877) 5-SUNSET (877-578-6738)

Our customer service professionals are here to assist you
M-F 8:00 AM – 5:00 PM Central Time.

Fax

312-997-9985 — Accepted 24 hours a day.

Email

sales@sunsethcs.com

*Contact us for more information
on fulfillment and compliance
solutions that best fit your needs.*

Shipping

Ask your rep about shipping to various locations and
our drop shipping partners.

UPS

Whenever possible, orders ship via UPS. If another
carrier is desired, please contact your sales rep.

All orders are packed in the smallest boxes available
and combined with other products on the order to
reduce shipping costs.

All products ship prepaid and are added to the invoice.

All products ship from Elgin, IL or Phoenix, AZ.
Please refer to UPS.com for transit times.

Damaged Merchandise

Call us within 10 days. Keep original cartons for
UPS inspection.

We proudly partner with



Payment Methods

Payment Portal and ACH

Visit www.sunsethcs.com/paymentportal to pay online with ACH
or credit card.

Net 30-Day Payment Terms

For qualified companies. A 1.5% monthly finance charge will be
assessed on all past due invoices.

COD

For non-qualified companies.

Credit Card

We accept all major credit cards. A processing fee may apply.
After 10 days of invoice date, a surcharge of 2.5% will be applied.



Customer Service

Returns

Contact us within 7 days of receiving your shipment.

Call us for a return authorization number. There is no
restocking charge on orders returned within 20 days of invoice.
Orders returned after 20 days may be assessed a restocking
fee of 15%. No returns after 90 days.
(877) 5-SUNSET (877-578-6738)

Central Distribution Center

Sunset Healthcare Solutions
2725 Alft Ln
Elgin, IL 60124

Western Distribution Center

Sunset Healthcare Solutions
6750 W Van Buren St Ste 100
Phoenix, AZ 85043

Problems

Please call us immediately if you have any complaints (it is
helpful to ask for the customer service representative with
whom you placed your order).

Or email us at: customerservice@sunsethcs.com

Questions?

Our experienced customer service staff will be glad to provide
you with information. Just give us a call.

Pricing

All prices are in U.S. Dollars.

Taxes are added where applicable.

We make every effort to hold pricing throughout the year.
However, prices are subject to change without notice.

Content

CPAP

	Page
Masks and Replacement Parts	4-12
Sunset Healthcare Solutions	4-6
Sleepnet	7-9
Circadiance	10-11
InnoMed	12
Airway Management	12
Bleep	13
Sleep Apnea Appliances	14-15
Tubing	16-17
Filters	18-26
Respironics	18-19
ResMed	20-21
Fisher & Paykel	22
React Health/3B	22
AEIOMed	23
Beyond	23
Breas/GE	23
CareFusion/Viasys	24
Covidien Nellcor Puritan Bennett	24
DeVilbiss	25
Evo Medical	25
Healthdyne	26
Invacare	26
Löwenstein	26
Resvent	26
ProBasics	26
Humidifier Chambers	27
Chinstraps	28-29
Headgear	29
Cash Items	30-37
Cleaners and Accessories	30-31
Battery/Power Inverter	32
Transcend Micro™ Auto CPAP	33
Mask Liners and Soft Headgear Covers	34
Self Care	35
Connectors, Adapters, and Inline Filter	35
Pillows and Storage	36
Tubing Wraps	37

Oxygen

	Page
Pulse Oximeter	38
Oxygen Sensors	39
Firesafe	40
Concentrators	41-42
Masks	43
Cannula and Cannula Accessories	44-45
Tubing	45
Microdial Flowmeter	46
Regulators and Conserver	47
Cylinder Carriers	49-50
Connectors	51
Accessories and Cleaners	51-53
Analyzers	54
Celon Wet Bands	54
Bacteria and HEPA Filters	55
Concentrator Cabinet Filters	56-62
Caire/AirSep/SeQual	56
CHAD	57
Covidien Nellcor Puritan Bennett	57
Dynarex	57
OxyGo	57
DeVilbiss/Drive	58
Invacare	59
Rhythm P2	59
Rhythm HC, Jumao, Attendant, Lifestyle, Medacure	60
Nidek	60
React Health/3B	60
Respironics	61
Healthdyne	62
Replacement Parts	63-64

Other

	Page
Disinfectants	65
Aerosol	66
Tabletop Nebulizers	67
Handheld Nebulizers	68
Nebulizer Kits	69
Nebulizer and Anesthesia Filters	70
Suction	71-73
Closed Suction System	74-75
Tracheostomy	76
HMEs	77
Ventilator Filters, Masks, and Connectors	78
Ventilator Circuits	79-81
Anesthesia Masks	80
StrongArm Mobility and Comfort	82-83
Bath Safety	83
Maternity	84-85
Wound Care	86-89
Equipment Covers	90-91

Sunset Healthcare Solutions

Sol Full Face CPAP Mask

Mask with Headgear and Replaceable Cushion


Size	Item#
Small	CM205S
Medium	CM205M
Large	CM205L

Replacement Cushions

Size	Item#
Small	CUCM205S
Medium	CUCM205M
Large	CUCM205L

HG205 Replacement Headgear - Fits Most

Ask us about
our 30-day
mask guarantee



Sol without Headgear

Size	Item#
Small	CM205SNH
Medium	CM205MNH
Large	CM205LNH

HCPCS#
MASK - A7030
HEADGEAR - A7035
CUSHION - A7031



Available to add
to your formulary
in SleepGlad

*A premium mask without
the premium price*

Easy pop cushions

Cushions attach or detach from mask in an easy pop; no need to painstakingly align pieces

Quiet vent

Honeycomb design reduces airflow erosion intensity for gentle and quiet flow

Quick release connector

Rapidly and easily separates or connects mask and the tube

Soft and light headgear

Skin-friendly fabric, light and stable design helps patients comfortably adapt to treatment

Sol comes in attractive retail packaging



Ray Pediatric CPAP/NIV Masks

• Comparable to Nonny masks

Mask with Headgear, Tube, and Replaceable Cushion

Size	Item#	OEM#
Small	CM061S	AG-PEDFF-S
Medium	CM061M	AG-PEDFF-M
Large	CM061L	AG-PEDFF-L
Extra Large	CM061XL	AG-PEDFF-XL

Replacement Cushions

Size	Item#	OEM#
Small	CUCM061S	AG-PEDFF-CS
Medium	CUCM061M	AG-PEDFF-CM
Large	CUCM061L	AG-PEDFF-CL
Extra Large	CUCM061XL	AG-PEDFF-CXL

Replacement Headgear

Size	Item#	OEM#
Sm/Med	HG061S	AG-PEDFF-HGS
Lg/XL	HG061L	AG-PEDFF-HGL

HCPCS#
MASK - A7030 • HEADGEAR - A7035 • CUSHION - A7031



Nasal Mask with Headgear, Tube, and Replaceable Cushions

Size	Item#	OEM#
S&M Cushions	CM060-KITS	AG-PEDKIT-S
L&XL Cushions	CM060-KITL	AG-PEDKIT-L

Replacement Cushions

Size	Item#	OEM#
Small	CUCM060S	AG-PEDKIT-CS
Medium	CUCM060M	AG-PEDKIT-CM
Large	CUCM060L	AG-PEDKIT-CL
Extra Large	CUCM060XL	AG-PEDKIT-CXL

Replacement Headgear

Size	Item#	OEM#
Small	HG060S	AG-PEDKIT-HGS
Large	HG060L	AG-PEDKIT-HGL

HCPCS#
MASK - A7034 • HEADGEAR - A7035 • CUSHION - A7032



Sunset Healthcare Solutions

Sunset Deluxe Full Face CPAP Mask

Mask with Flex Tube, Headgear, and Replaceable Cushion

Size	Item#
Small	CM005S
Medium	CM005M
Large	CM005L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

CMSIZE01
 Deluxe Sizing Guide

HCPCS#
 MASK - A7030
 HEADGEAR - A7035
 CUSHION - A7031

Unobtrusive design
 offers a clear
 line of sight



- 360° dual-swivel elbow port allows for free movement while connected to the tubing
- Quick-release clips make it easy to remove the mask without resetting the headgear
- Ultra-lightweight silicone replaceable cushion increases comfort while maintaining optimal seal

Replacement Cushions

Size	Item#
Small	CUCM005S
Medium	CUCM005M
Large	CUCM005L

Mask without Headgear

Size	Item#
Small	CM005SNH
Medium	CM005MNH
Large	CM005LNH

Sunset Deluxe Nasal CPAP Masks

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM006S
Medium	CM006M
Large	CM006L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

CMSIZE01
 Deluxe Sizing Guide

HCPCS#
 MASK - A7034
 HEADGEAR - A7035
 CUSHION - A7032

Same features as the full face with 360° dual-swivel elbow port, quick release clips, and ultra-lightweight replaceable cushions.



Mask with Headgear, Adjustable Forehead T, and Replaceable Cushion

Size	Item#
Small	CM106S
Medium	CM106M
Large	CM106L

Replacement Cushions

Size	Item#
Small	CUCM006S
Medium	CUCM006M
Large	CUCM006L



Mask without Headgear

Size	Item#
Small	CM006SNH
Medium	CM006MNH
Large	CM006LNH

Sunset Healthcare Solutions

Sunset Deluxe Nasal Pillow CPAP Mask

Mask with Headgear and Small, Medium, and Large Pillows

Item# **CM036**

Headgear Only
HG036

HCPCS#
MASK - A7034
HEADGEAR - A7035
PILLOWS - A7033

Replacement Pillows

Size	Item#
Small	CM036SP
Medium	CM036MP
Large	CM036LP

Mask without Headgear
CM036NH



- Easy to fit, assemble, and maintain
- Three different pillow sizes to fit more people more easily
- 360° rotational elbow with flexible silicone tube
- Quiet vent design

Sunset Clearsight Deluxe Nasal Mask

Mask with Headgear and Replaceable Cushion

Size	Item#
Small	CM110S
Medium	CM110M
Large	CM110L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

Replacement Cushions

Size	Item#
Small	CUCM110S
Medium	CUCM110M
Large	CUCM110L

HCPCS#
MASK - A7034
HEADGEAR - A7035



Clearsight design provides unobstructed view

Mask Series	Mask Type	Replaceable Cushion	Flex Tube Included	Headgear Included
CM005	Full Face	✓	✓	✓
CM006	Nasal - Forehead Support	✓	✓	✓
CM106	Nasal - Adjustable	✓		✓
CM007	Full Face			✓
CM008	Nasal			✓
CM110	Nasal	✓		✓
CM036	Nasal Pillow	✓	✓	✓

Sunset Classic CPAP Masks

Full Face Mask with Headgear

Size	Item#
Small	CM007S
Medium	CM007M
Large	CM007L
Extra Large	CM007XL

HG003 Headgear - Fits Most
HG003XL Headgear - XL

HCPCS#
MASK - A7030
HEADGEAR - A7035



- Conventional forehead pad removed to create a clear line of sight
- Minimalist design and ultra-lightweight silicone
- Quick-release clips make it easy to remove
- 360° dual-swivel elbow port allows for free movement
- One-piece mask and cushion design is durable and easy

Nasal Mask with Headgear

Size	Item#
Medium	CM008M
Large	CM008L

HG003 Headgear - Fits Most
HG003XL Headgear - XL

HCPCS#
MASK - A7034
HEADGEAR - A7035



- Great value for cash patients or those with lower insurance reimbursement
- Headgear clips allow for quick removal of mask — no further adjustments required after initial fitting
- One-piece mask and cushion design is durable and easy

Sleepnet

Phantom®2 and iQ®2

Phantom®2

Item#
CM050

Replacement
Cushion
CUCM050



iQ®2

Item#
CM049

Replacement
Cushion
CUCM049



*Ask us about
our 30-day
mask guarantee*

Replacement Frame/Tubing Assembly
CM051

- Sleepnet's patented AIR°gel® technology provides maximum comfort while minimizing skin irritation and pressure points on the nasal bridge. Unlike many other gel masks, the unique qualities of AIR°gel allow the mask to seal with minimal tension. By design AIR°gel is twice as soft as competing gel products.
- Sleepnet has the only mask that you can shape to each individual face. Flex, Fit, Forget!
- Cushions may be removed and replaced without purchasing a new mask.
- Magnetic Clips — Easy on and off!
- 360° Swivel — Less tube drag for active sleepers.
- One size cushion — Larger fitting range.



Sleepnet

Innova

Nasal Mask with Headgear

Size	Item#
Small/Medium	CM034SM
Medium/Large	CM034ML

HCPCS#
MASK - A7034
HEADGEAR - A7035



Full Face Mask with Headgear

Size	Item#
Small	CM035S
Medium	CM035M
Large	CM035L

HCPCS#
MASK - A7030
HEADGEAR - A7035



Innova® CPAP masks feature patented AIR^ogel with Advanced Cushion Technology and Active Headgear Connector™. Together, they provide an optimal fit, maximum comfort, and freedom.

- The Touchless Spacebar's unique design eliminates any contact with the patient's forehead
- The quick-release headgear is lightweight, breathable, and slides off easily with one hand
- Made in the USA

MiniMe Pediatric Nasal Mask



MiniMe 2 Pediatric Nasal Mask

Mask with Headgear

Size	Item#	OEM#
Extra Small	CM027XS	60216
Small	CM027S	60217
Medium	CM027M	60218
Large	CM027L	60221

Replacement Headgear

Size	Item#
Extra Small/Small	HG027S
Medium/Large	HG027L

HCPCS#
MASK - A7034
HEADGEAR - A7035

- Non-vented version available



Mask with Headgear


Headgear Size	Item#	OEM#
Small	CM004S	60213
Medium	CM004M	60214
Large	CM004L	60215

Replacement Headgear

Size	Item#	OEM#
Small	HG010S	55241
Medium	HG010M	55242
Large	HG010L	55243

HCPCS#
MASK - A7034
HEADGEAR - A7035

- Non-vented version available



MASK SIZE	HEAD CIRCUMFERENCE
SMALL	14-16in
MEDIUM	17-18in
LARGE	19-21in

Measure around head above eyebrows for head circumference

Sleepnet

Ascend

Nasal Mask with Headgear and All Cushions

Item# **CM037N**


Replacement Cushions

Size	Item#
Small	CM037NCS
Medium	CM037NCM
Large	CM037NCL

Nasal Mask with Headgear and 1 Cushion

Size	Item#
Small	CM037NS
Medium	CM037NM
Large	CM037NL

HCPCS#
MASK - A7034
HEADGEAR - A7035
CUSHION - A7032



- One frame with three sizes of replaceable cushions: S, M, L
- Unique flexible spacebar design provides superior stability, holding a seal without contacting the skin
- Sleepnet's patented AIR°gel minimizes pressure points and skin irritation while offering maximum comfort with minimum weight
- Lightweight, breathable quick-release headgear slides off easily with one hand
- 360-degree rotating elbow for optimal tube positioning and easier mobility
- Active Headgear Connector moves as the patient moves for a reliable seal
- Made in the USA

Full Face Mask with Headgear and All Cushions

Item# **CM037F**

Replacement Cushions

Size	Item#
Small	CM037CS
Medium	CM037CM
Large	CM037CL

Full Face Mask with Headgear and 1 Cushion

Size	Item#
Small	CM037S
Medium	CM037M
Large	CM037L

HCPCS#
MASK - A7030
HEADGEAR - A7035
CUSHION - A7031



iQ Nasal Mask

Item#
CM001

iQ Nasal Mask with Headgear

OEM# TMS-889

HCPCS#
MASK - A7034
HEADGEAR - A7035



Mojo Full Face Mask

Mask with Headgear

Size	Item#	OEM#
Small	CM003SH	TMS-835
Medium	CM003MH	TMS-836
Large	CM003LH	TMS-837

Replacement Headgear

HG003
OEM# TMS-89

HCPCS#
MASK - A7030
HEADGEAR - A7035



Mojo 2 Full Face Mask

Mask with Headgear

Size	Item#	OEM#
Small	CM50817	50817
Medium	CM50818	50818
Large	CM50819	50819
Extra Large	CM50820	50820

- Maximum comfort, minimum weight. Patented AIR°gel minimizes pressure points and skin irritation
- Forehead adjustment screw removes pressure from the bridge of the nose for personalized fit
- The shapeable shell of each mask can be easily molded to fit individual needs: bend it, mold it, shape it!
- Lightweight, breathable headgear disconnects from mask easily with one hand thanks to the quick-release magnetic headgear clips
- Oxygen port allows easy delivery of supplemental oxygen
- 360° rotating elbow for optimal tube positioning and easier mobility



Circadiance

SleepWeaver® Advance Nasal Mask



- One size fits most (Small Tan available for smaller faces)
- Includes Advance headgear that can be adjusted to fit a broader range of patients

Mask with Adjustable Headgear

Mask Color	Item#	OEM#
Blue	CM017B	100274
Camo	CM017C	100281
Leopard	CM017L	100285
Pink	CM017P	100277
Tan	CM017T	100332

CM017T-S
Small Tan
OEM# 101224

HCPCS#
MASK - A7034
HEADGEAR - A7035

HG017
Headgear Only
OEM# 100435

CM017TNH
Tan — No Headgear
OEM# 100336

- Easy to store and pack for travel

SleepWeaver® Élan Nasal Mask

Mask Kit

Mask Color	Item#	OEM#
Tan	CM022T-KIT	100563
Blue	CM022B-KIT	100782

Kits Include:

- Mask
- Small Cushion
- Regular Cushion
- Large Cushion
- Headgear
- Feather Weight Tube®
- Overhead Tube Tether Strap
- Quick Start Illustration Guide
- Instructional DVD

HCPCS#
MASK - A7034
HEADGEAR - A7035
TUBE - A7037

Replacement Cushions

Size and Color	Item#	OEM#
Small Blue	CUCM022B-S	100743
Regular Blue	CUCM022B-R	100654
Large Blue	CUCM022B-L	100664
Small Tan	CUCM022T-S	100742
Regular Tan	CUCM022T-R	100653
Large Tan	CUCM022T-L	100663
Small Leopard	CUCM022L-S	100746
Regular Leopard	CUCM022L-R	100657
Large Leopard	CUCM022L-L	100667

TUB0018F
Replacement
Feather Weight Tube®
OEM# 100477

HG022
Headgear Only
OEM# 100630

Mask with Headgear only

Size and Color	Item#	OEM#
Small Blue	CM022B-S	100728
Regular Blue	CM022B-R	100519
Large Blue	CM022B-L	100482
Small Tan	CM022T-S	100727
Regular Tan	CM022T-R	100518
Large Tan	CM022T-L	100481



All SleepWeaver Masks

- Made of soft cloth
- Allow sleeping on side, stomach, and even wearing glasses
- Easy to store and pack for travel
- Hand or machine washable gentle cycle, mild detergent, hang to dry



- Detachable 90° elbow with 360° rotation allows for over-the-head hose placement

Circadiance

SleepWeaver® Anew™ Full Face Mask

Blue Mask with Headgear

Size	Item#	OEM#
Small	CM026B-S	100955
Regular	CM026B-R	100960
Large	CM026B-L	100965

Headgear

Size	Item#	OEM#
Small	HGCM026-S	101145
Reg/Lg	HGCM026-R	101084

TUB0018F

Feather Weight Tube®
OEM# 100477

HPCPS#
MASK - A7027
HEADGEAR - A7035
TUBE - A7037



- Soft material wicks away moisture and leaves no pressure marks on the face
- Loose fit provides ultimate comfort and eliminates strap marks

- Open line of sight accommodates eyeglasses and is ideal for people who like to fall asleep while reading or watching TV
- 90° elbow swivel connector rotates 360° allowing users to sleep in more natural positions; easy to remove and clean

Pediatric SleepWeaver® CPAP Mask

Item#
CM101325



The launch of the SleepWeaver Advance Pediatric Soft Cloth CPAP Mask marks the development of a revolutionary interface for CPAP use in pediatric patients. This new configuration extends the range of SleepWeaver products to include pediatric users with delicate facial features.

- Same proven cloth technology enhanced to accommodate children 2 through 7 years of age

SleepWeaver® 3D®

Item#
CM027


Blue Mask without Headgear

Headgear	Item#
With Large Headgear	CM027-LB
With Regular Headgear	CM027-RB

100/packs

Headgear	Item#
With Large Headgear	CM027-L-100
With Regular Headgear	CM027-R-100

- New wing-shaped stay routes the headgear under the ears, ensuring comfort and a great fit
- Permanent 360° swivel provides reliable performance in any position, including over the headboard
- Breathable cloth design eliminates fogging that is common in hard plastic masks - nothing breathes like cloth, keeping skin dry and comfortable
- For those who want to read or watch TV while in bed, SleepWeaver 3D can be worn with glasses
- Users prefer a soft cloth against the face over rubber, silicone and other gels. Cloth masks eliminate the need for a cushion or nasal pad



Innomed

Aloha® Nasal Pillows System

Item#

CM014

All Sizes Kit

OEM# ALO100

Replacement Headgear

HG014

Aloha Headgear

OEM# ALO130

Replacement Pillows

Size	Item#	OEM#
Small	CM014SP	ALO121
Medium	CM014MP	ALO123
Large	CM014LP	ALO125

Frame

Size	Item#	OEM#
Small	CM014FS	ALO131
Medium	CM014FM	ALO133
Large	CM014FL	ALO135

HCPCS#
MASK - A7034
HEADGEAR - A7035
PILLOWS - A7033

- One of the highest consumer-rated CPAP masks on the market
- Limited headgear reduces points of contact on the face and prevents interference with the eyes

- Removable micro-fleece headgear wraps are included for comfort
- Pillow angle can be adjusted for a more customized and natural fit



Airway Management

Tap Pap

No headgear needed!

Item#

CM024NH

Mask Kit with Small, Medium, and Large Pillows

OEM# PAP-NP1-000

Replacement Pillows

Size	Item#	OEM#
Small	CUCM024S	PAP-NP1-101
Medium	CUCM024M	PAP-NP1-102
Large	CUCM024L	PAP-NP1-103

CM024MP Replacement Mouthpiece

OEM# PAP-NP1-201

HCPCS#
MASK - A7034
PILLOWS - A7033
HEADGEAR - A7035



- Mouthpiece molds to the upper teeth to provide a custom and secure fit even if the bottom teeth are not in contact with the mouthpiece
- Unique nut and bolt joint allows the nasal pillows to be set to the desired angle
- Provides a stable fit regardless of the sleep position

- Patients report not being able to feel the mask on their face
- 3 sizes of nasal pillows included (Small, Medium, and Large)
- 90-day compliance guarantee for new users of the mask, 30-day guarantee for experienced users

Bleep

Bleep Eclipse Starter Kit

Eclipse Frame and Two 32/Packs of Halos

Item#
CAP100589

The BLEEP Eclipse™ Starter Kit includes the headgear-free Eclipse™ mask frame and two 32-packs of Halos™ (32-day total supply) for immediate use. This clinically-proven NO LEAKS PAP solution is the only one of its kind on the market. No headgear, no mask lines and no broken hair. The comfortable Eclipse mask frame includes a diffuser and flexible tube and connects to Halos in lieu of headgear. Halos affix to the nasal area without inserting into the nostrils, significantly reducing irritation and claustrophobia, driving patient comfort and health, and preventing face marks. Halos are hypoallergenic, with gentle, surgical grade adhesive, and free from BPA, corn or latex. Halos utilize innovative MagSeal™ technology for rapid connection and disconnection. BLEEP sleep solutions are covered by Medicare and insurance, and the only CPAP solutions fully made in the USA.

HPCS# A7034



Bleep Halos for Eclipse Mask Frame

32/Pack

Item#
CAP100590

The Eclipse™ Halos™ 32-pack provides a 16-day supply of disposable nasal interfaces for use with the Eclipse mask frame (sold separately). Halos affix to the nasal area without inserting into the nostrils, significantly reducing irritation and claustrophobia, driving patient comfort and health, and preventing face marks. This clinically-proven NO LEAKS PAP solution is the only one of its kind on the market. Halos are hypoallergenic, with gentle, surgical grade adhesive, and free from BPA, corn or latex. Halos utilize innovative MagSeal™ technology for rapid connection and disconnection. Apply Halos on dry skin after cleansing the nasal area with Dream Preps™ or a mild astringent such as witch hazel or rubbing alcohol. BLEEP sleep solutions are the only CPAP solutions fully made in the USA.

HPCS# A7033




Sleep Apnea Appliances

Smart Nora Contact-Free, Hassle-Free Snoring Solution

Item#
CAP8000
Each

NEW

Turns any pillow into an anti-snore pillow



Easy to set up and use

What's Included

- 1. Pebble** This goes on the patient's nightstand and listens for early signs of snoring.
- 2. Felt box** The felt box functions as a travel bag. It also houses the silent base. Box can be stored under the bed.
- 3. Expander** The expander goes under the pillow between the pillow and the pillowcase.
- 4. Universal charger** Smart Nora comes with a universal charger and interchangeable plug-heads.
- 5. USB cable for pebble**



Hush Strips®

Item#
CAP1025
Box of 32 Strips

NEW

Hush Strips® are designed to improve sleep quality by reducing mouth breathing and snoring during sleep

Hush Strips® sleep strips are made from medical grade nonwoven fabric and gentle easy-release adhesive, which moves with the skin flexibly, making them easy to reposition and remove.


Easy to apply residue-free material

The simple, non-invasive and easy to use anti-snoring strips are made of medical-grade hypoallergenic material which doesn't leave residue and is safe for both kids and adults.



How to use

1. Wash and dry lips
2. Remove strip from the liner
3. Close your mouth and turn lips slightly inward
4. Apply strip on lips and press to seal



Airway Management

myTAP



Item# **CM039**



Precision-fit, immediate, simple and affordable solution for snoring and mild to moderate sleep apnea.

- Offers immediate treatment
- Can be fitted in 15 minutes in a sleep lab, physician's office or DME
- Clinically proven in over 40 independent peer reviewed studies to treat snoring and obstructive sleep apnea (OSA)
- Patients control and remold trays with Precision-Fit, enabling overall comfort, compliance and success: 90% of patients achieve ideal fit upon first try, no bite registration needed
- Thin, sleek trays create a comfortable, low-profile appliance
- Simple design makes it an affordable and convenient option
- Includes an AM Aligner to help offset potential bite changes
- Includes a mouth shield accessory to prevent dry mouth or excess saliva. A secondary benefit is the promotion of nasal breathing, the best natural way to breathe

AirAvant

Bongo Rx EPAP

Item#

CAPBNG500

Bongo Rx all sizes starter kit

- All four sizes: small, medium, large, and extra large
- Travel case
- Drying stand
- Headgear
- No machine, tube, electricity, or batteries needed
- Use at home, traveling, camping... use anytime, anywhere
- Small enough to fit in a shirt pocket
- Reusable and easy to clean with simply soap and water
- Designed, molded, and assembled in the USA
- Clinically tested to demonstrate safety and effectiveness
- Supplemental lifestyle option for your compliant patients
- Patient friendly option for your non-compliant patients



The Bongo Rx is a simple device with soft nasal seals and is available in various sizes (S, M, L, and XL). The Bongo Rx is an EPAP (expiratory positive airway pressure) device and works as follows:

- During inhalation, small valves open to allow normal breathing through the device.
- During exhalation, the valves close and direct exhaled air through small vents to generate EPAP, which keeps the airway open.

Annual Replenishment Pack

Includes 4 pairs, travel case, drying stand, and headgear

Size	Item#
Small	CAPBNG502
Medium	CAPBNG503
Large	CAPBNG504
Extra Large	CAPBNG505



Single Replenishment Unit

Includes 1 pair and headgear

Size	Item#
Small	CAPBNG512
Medium	CAPBNG513
Large	CAPBNG514
Extra Large	CAPBNG515



CPAP Tubing

Durable CPAP Tubing — Gray



Item#

TUB06

6ft

19mm with 22mm cuffs

HCPCS# A7037


- Sunset's most popular item
- Compatible with nearly all CPAP, BiPAP and BiLevel systems for every manufacturer on the market
- Smooth interior and a corrugated exterior for excellent flexibility and unrestricted interior flow characteristics
- Injection molded cuffs to reduce tearing
- Low memory to avoid kinks
- BPA-free and latex-free

All Lengths

TUB0018	18in
TUB002	2ft
TUB004	4ft
TUB06	6ft
TUB006.5	6ft 6in — Compatible with ResMed tube wraps
TUB007	7ft
TUB008	8ft
TUB010	10ft

Compatible with Sunset's new Quick Release Connector!

Item#

CAP1016

Slim CPAP Tube

Item#

TUB006SS

6ft x 15mm with 22mm cuffs

HCPCS# A7037



Highly durable

Strain relief cuff decreases separation and tearing

Virtually no memory

Designed to be the most flexible tube available

Minimal weight

One of the lightest tubes on the market at 2.6 ounces

Great value


Affordable and sustainable

Pressure Port Tubing

Item# **TUB006PP**

6ft with Pressure Sensor Lines 22mm Cuffs

HCPCS# A7037



- Built-in sensor line for use with the Puritan Bennett: GoodKnight 420S, 420E and 425 and KnightStar 330

Heated CPAP Tubing

Heated Tube for Respirationics

Item#


TUB06HPR

Heated Tube for
Philips Respironics
Dreamstation Series
and System One

OEM# HT15

HCPCS# A4604

- 510(k)
- UL Recognized
- Strain relief cuffs
- Auto safety shut-off



TUBHC900ICON208

Heated Tube for Fisher & Paykel ICON

OEM# 900ICON208

Hybernite Heated CPAP Tube 6ft

Item#

TUB006ROC *MOQ applies*

- Copper wire embedded in the tubing wall generates uniform heat and maintains the air temperature inside the tube to prevent rainout
- Stand-alone kit can be used in the sleep lab or at home with any CPAP or BiPAP machine with humidifier
- Includes heated circuit tube and independent power supply

Compatible with machines that use standard CPAP Tubing.
Not compatible with machines that use slim 15mm tubing
such as newer Respironics and ResMed machines.



TUB006HBT

Hybernite Replacement Heated CPAP Tube 6ft

Respironics

DreamStation 2



Item# CF1012

Disposable Filter

6/pack

OEM# 1142831

HCPCS# A7038

CF1012-1 1/pack OEM# 1142828

CF1012-2 2/pack OEM# 1142829

CF1012-3 3/pack

CF1012-50 50/pack

Item# CF1012F

Reusable Filter

2/pack

OEM# 1142687

HCPCS# A7039

CF1012F-1 1/pack

DreamStation



Item# CF1010

Disposable Filter

6/pack

OEM# 1122519

HCPCS# A7038

CF1010-1 1/pack

CF1010-2 2/pack

Item# CF1010F

Reusable Filter

2/pack

OEM# 112446 (1-pack)

HCPCS# A7039

CF1010F-1 1/pack

CF1010F-6 6/pack

M-Series



Item# CF1007

Ultrafine Filter - Tabbed

6/pack • 1³/₄" x 7⁷/₈"

HCPCS# A7038

CF1007-1 1/pack

CF1007-2 2/pack

Item# CF1007F

Foam Filter

2/pack • 1⁵/₈" x 7⁷/₈"

OEM# 1029330

HCPCS# A7039

CF1007F-1 1/pack

PR System One



Item# CF1009NT

Ultrafine Filter - No Tab

6/pack • 1³/₄" x 7⁷/₈"

OEM# 1029331

HCPCS# A7038

CF1009NT-1 1/pack

CF1009NT-2 2/pack

Item# CF1009

Ultrafine Filter - Tabbed

6/pack • 1³/₄" x 7⁷/₈"

HCPCS# A7038

CF1009-1 1/pack

CF1009-2 2/pack

Item# CF1009F

Foam Filter

2/pack • 1⁵/₈" x 7⁷/₈"

OEM# 1029330

HCPCS# A7039

CF1009F-1 1/pack

*All filters can be ordered in any pack quantity.
Contact your sales rep for details.*

Respironics

REMstar Pro, REMstar Plus



Item# CF1006
Ultrafine Filter
2/pack • 3 $\frac{1}{2}$ " x 1 $\frac{1}{2}$ "
OEM# 1005945
HCPCS# A7038

Item# CF1006F
Foam Filter
2/pack • 3 $\frac{5}{8}$ " x 1 $\frac{5}{8}$ "
OEM# 1005964
HCPCS# A7039

BiPAP Pro, Plus, Duet LX, Synchrony



Item# CF1005
Ultrafine Filter
6/pack • 4 $\frac{3}{4}$ " x 3 $\frac{1}{4}$ "
OEM# 622219
HCPCS# A7038

Item# CF1005F
Foam Filter
2/pack • 4 $\frac{3}{4}$ " x 2 $\frac{3}{8}$ "
OEM# 622220
HCPCS# A7039

Aria, Virtuoso, Duet, Quartet



Item# CF1004
Ultrafine Filter
6/pack • 3 $\frac{1}{2}$ " x 2 $\frac{3}{16}$ "
OEM# 532311
HCPCS# A7038

Item# CF1004F
Foam Filter
2/pack • 3 $\frac{1}{4}$ " x 2 $\frac{1}{8}$ "
OEM# 532247
HCPCS# A7039

REMstar



Item# CF1003
Ultrafine Filter
6/pack • 5 $\frac{1}{2}$ " x 4 $\frac{3}{4}$ "
OEM# 362522
HCPCS# A7038

Item# CF1003F
Foam Filter
1/pack • 5 $\frac{3}{4}$ " x 5 $\frac{1}{4}$ "
OEM# 362521
HCPCS# A7039

Solo, Aria LX, Virtuoso LX, REMstar LX



Item# CF1002
Ultrafine Filter
6/pack • 3 $\frac{3}{8}$ " x 2 $\frac{3}{4}$ "
OEM# 622017
HCPCS# A7038

Item# CF1002F
Foam Filter
2/pack • 3 $\frac{5}{16}$ " x 2 $\frac{5}{16}$ "
OEM# 622018
HCPCS# A7039

SleepEasy/BiPAP



Item# CF1001
Ultrafine Filter
6/pack • 7 $\frac{3}{4}$ " x 4 $\frac{7}{8}$ "
OEM# 302064
HCPCS# A7038

ResMed

AirSense 11



Item#

CF2111

Disposable Filter

3/pack

HCPCS# A7038

CF2111-1 1/pack OEM# 39300

CF2111-2 2/pack OEM# 39301

CF2111-50 50/pack OEM# 39303

Item#

CF2111HA-3

Hypoallergenic Filter

3/pack

1 1/4" x 1 3/16"

HCPCS# A7038

CF2111HA-1 1/pack

CF2111HA-2 2/pack

S9 and AirSense 10 Series



Item#

CF2107

Disposable Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107-1 1/pack OEM# 36850

CF2107-2 2/pack OEM# 36851

CF2107-50 50/pack OEM# 36853



Item#

CF2107HA

Hypoallergenic Filter

3/pack • 2 1/8" x 1 3/8"

HCPCS# A7038

CF2107HA-1 1/pack OEM# 36855

CF2107HA-2 2/pack OEM# 36856

CF2107HA-50 50/pack
OEM# 36858

AirMini



Item#

CF2108

Disposable Filter

2/pack • 1 1/4" x 1 3/16"

OEM# 38815

HCPCS# A7038

Adapt SV



Item#

CF2106

Disposable Filter

3/pack • 1 1/4" x 1 1/2"

OEM# 26926

HCPCS# A7038

Item#

CF2106HA

Hypoallergenic Filter


3/pack • 1 1/4" x 1 1/2"

OEM# 26927

HCPCS# A7038

ResMed

S8 Series



Item#
CF2105
Disposable Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33916
HCPCS# A7038

CF2105-1 1/pack OEM# 33915
CF2105-2 2/pack OEM# 33973
CF2105-50 50/pack OEM# 33918



Item#
CF2105HA
Hypoallergenic Filter
3/pack • 1³/₈" x 1³/₈"
OEM# 33920
HCPCS# A7038

CF2105HA-1 1/pack OEM# 33919
CF2105HA-2 2/pack OEM# 33975
CF2105HA-50 50/pack OEM# 33922

S7 Series



Item#
CF2103
Disposable Filter
3/pack • 1¹/₄" x 1¹/₂"
OEM# 30907
HCPCS# A7038

CF2103-1 1/pack OEM# 30919
CF2103-2 2/pack
CF2103-50 50/pack OEM# 30908



Item#
CF2103HA
Hypoallergenic Filter
3/pack • 1¹/₄" x 1¹/₂"
HCPCS# A7038

CF2103HA-1 1/pack OEM# 30933

S6 Series



Item#
CF2102
Disposable Filter
3/pack • 5³/₄" x 5⁵/₈"
OEM# 21935
HCPCS# A7038

CF2102-50 50/pack
OEM# 21936

Fisher & Paykel

SleepStyle Series



Item#
CF4004-2
Disposable Filter
2/pack
HCPCS# A7038



ICON Series



Item#
CF4003
Disposable Filter
4/pack
2 1/16" x 5/8"
HCPCS# A7038



CF4003-1 1/pack
OEM# 900ICON503

CF4003-6 6/pack
OEM# 900ICON212

HC230 Series, HC600 Series



Item#
CF4002
Disposable Filter
4/pack • 2 3/4" x 7/8"
OEM# 900HC240
HCPCS# A7038
CF4002-2 2/pack



HC200, HC220, HC221



Item#
CF4001
Disposable Filter
4/pack • 5" x 3/4"
OEM# 900HC222
HCPCS# A7038
CF4001-2 2/pack



React Health/3B

Luna G3



Item# **CF8006-6**
Disposable Filter
6/pack
HCPCS# A7038
CF8006-1 1/pack
CF8006-2 2/pack



Item# **CF8006F-1**
Reusable Filter
1/pack
HCPCS# A7039

Luna and Luna II



Item#
CF8005
Disposable Filter
12/pack • 2" x 1 3/16"
HCPCS# A7038



Item#
CF8005F-1
Reusable Filter
1/pack • 2" x 1 3/16"
HCPCS# A7039
CF8005F-2 2/pack



RESmart




Item#
CF8004
Disposable Filter
3/pack • 3 5/16" x 1 5/8"
HCPCS# A7038



AEIOMed

Everest Aura, Everest III and Integrated



Item# **CF8001**
Foam Filter
2/pack • 2" x 1⁹/₁₆"
OEM# 100446
HCPCS# A7039

Beyond

ResPlus



Item# **CF9005F-2**
Reusable Filter
2/pack
HCPCS# A7039
CF9005F-1 1/pack

Breas/GE

iSleep 20 Series



Item# **CF7003**
Ultrafine Filter
2/pack • 2¹/₄" x 7⁷/₈"
OEM# 004153
HCPCS# A7038

Item# **CF7003F**
Foam Filter
Each • 2³/₄" x 1"
OEM# 004154
HCPCS# A7039

PV10



Item# **CF7002**
Ultrafine Filter
2/pack • 2³/₄" x 7⁷/₈"
OEM# 2B01975
HCPCS# A7038

Item# **CF7002F**
Foam Filter
Each • 2⁵/₈" x 7⁷/₈"
OEM# 2B01972
HCPCS# A7039

PV1



Item# **CF7001**
Ultrafine Filter
2/pack • 2¹/₂" x 1¹/₄"
HCPCS# A7038



Item# **CF7001F**
Foam Filter
3/pack • 2¹/₂" x 1¹/₄"
HCPCS# A7039

CareFusion/Viasys

PureSom



Item#
CF9003
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item#
CF9003F
Foam Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7510
HCPCS# A7039

Orion and Pegasus



Item#
CF8002
Ultrafine Filter
6/pack • 2⁷/₈" x 1³/₄"
OEM# 774688
HCPCS# A7038

Item#
CF8002F
Foam Filter
Each • 2⁷/₈" x 1³/₄"
OEM# 774689
HCPCS# A7039

Covidien Nellcor Puritan Bennett

GoodKnight 314/318



Item#
CF6005F
Foam Filter
2/pack • 6¹/₄" x 1"
OEM# S-290609-00
HCPCS# A7039

GoodKnight 420/425



Item#
CF6003
Ultrafine Filter
2/pack • 1³/₁₆" x 3¹/₁₆"
OEM# M-413950-04
HCPCS# A7038

Item#
CF6003F
Foam Filter
Each • 1³/₁₆" x 3¹/₁₆"
OEM# M-413950-02
HCPCS# A7039

GoodKnight 418 A/G/P/S



Item#
CF6002
Ultrafine Filter
4/pack • 5³/₁₆" x 3³/₄"
OEM# Y-413560-01
HCPCS# A7038

GoodKnight 418



Item#
CF6001
Ultrafine Filter
4/pack • 2¹/₄" x 1²/₅"
OEM# S-645030-00
HCPCS# A7038

Item#
CF6001F
Foam Filter
6/pack • 2⁷/₁₆" x 1¹/₂"
OEM# S-645029-00
HCPCS# A7039

DeVilbiss

IntelliPAP



Item#
CF5004
Ultrafine Filter
4/pack • 1³/₄" x 1¹/₄"
OEM# DV51D-603
HCPCS# A7038

Item#
CF5004F
Foam Filter
4/pack • 1⁵/₈" x 1³/₁₆"
OEM# DV51D-602
HCPCS# A7039

Horizon LT



Item#
CF5003
Ultrafine Filter
4/pack • 4¹/₄" x 1⁵/₁₆"
OEM# 8000D-603
HCPCS# A7038

Item#
CF5003F
Foam Filter
4/pack • 4¹/₄" x 1⁵/₁₆"
OEM# 8000D-602
HCPCS# A7039

Evo Medical

RemRest



Item#
CF8003F
Foam Filter
3/pack • 2³/₁₆" x 1¹/₁₆"
OEM# 21117
HCPCS# A7039

Comfort PAP



Item#
CF9003
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item#
CF9003F
Foam Filter
3/pack • 1³/₈" x 1¹/₄"
OEM# 7510
HCPCS# A7039



Enhance your branding!

FREE PRIVATE LABELING

Sunset's free service allows you to put your logo and contact information on our products.

Increase brand awareness and allow patients to easily contact you.

Work with our marketing team to customize the layout with your logo, phone number, website address, and reorder information.

Available on all CPAP tubes and filters.
Minimum orders apply.
Certain restrictions may apply.
Contact your sales rep for details.

Healthdyne

Quest



Item# CF3003
Ultrafine Filter
with Curves
5/pack • 4¹/₂" x 4³/₈"
OEM# 7302
HCPCS# A7038



Item# CF3003F
Foam Filter
with Curves
5/pack • 4¹/₂" x 4¹/₄"
OEM# 7301
HCPCS# A7039

Tranquility



Item# CF3001
Ultrafine Filter
5/pack • 5¹/₈" x 1⁷/₈"
OEM# 7170
HCPCS# A7038

Löwenstein

prisma



Item# VF008-12
Disposable Filter
12/pack
OEM# WM29652



Item# VF008F-2
Reusable Filter
2/pack
OEM# WM29928

Resvent

iBreeze



Ultrafine Filter HCPCS# A7038
CF9901-1 1/pack
CF9901-2 2/pack
CF9901-6 6/pack
CF9901-50 50/pack

Invacare

Polaris LT



Item# CF9002
Ultrafine Filter
3/pack • 5" x 2"
OEM# 1082428
HCPCS# A7038

Item# CF9002F
Foam Filter
3/pack • 5" x 2"
OEM# 1082995
HCPCS# A7039


Polaris EX



Item# CF9001
Ultrafine Filter
3/pack • 2³/₁₆" x 1"
HCPCS# A7038

Item# CF9001F
Foam Filter
3/pack • 2¹/₄" x 1¹¹/₁₆"
OEM# 1124211
HCPCS# A7039

Polaris TR



Item# CF8001
Foam Filter
2/pack • 2" x 1⁹/₁₆"
OEM# 100446
HCPCS# A7039

ProBasics

Zzz-PAP



Item# CF9003
Ultrafine Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7511
HCPCS# A7038

Item# CF9003F
Foam Filter
3/pack • 1³/₁₆" x 1⁷/₁₆"
OEM# 7510
HCPCS# A7039

Respironics

Humidifier Chamber for PR System One



Item# **HUM1136502**
Humidifier Chamber for PR System
One 60 Series
Each HCPCS# A7046

Humidifier Chamber for DreamStation



Item# **HUM1122520**
Humidifier Chamber
for DreamStation
Each HCPCS# A7046

Fisher & Paykel

Humidifier Chambers for HC100 Series



Item# **HUMHC300**
Extended Life Humidifier Chamber
for the Fisher & Paykel HC100 Series
Each OEM# HC 300



Item# **HUMHC325S**
Replacement Chamber for the
Fisher & Paykel HC100 Series
Each OEM# HC 325S
HCPCS# A7046

Compatible with:

- HC100
- HC105
- HC150



Humidifier Chambers for HC200 Series



Item# **HUMHC355**
Extended Life Humidifier Chamber for
the Fisher & Paykel HC200 Series
Each OEM# HC 355



Item# **HUMHC385**
Replacement Chamber for the Fisher
& Paykel SleepStyle HC200 Series
Each OEM# HC 385S
HCPCS# A7046



Compatible with:

• HC200	• HC231	• HC236
• HC201	• HC232	• HC238
• HC220	• HC233	• HC242
• HC221	• HC234	• HC244
		• HC254

Humidifier Chamber for HC600 Series



Item# **HUMHC365**
Replacement Chamber for the
Fisher & Paykel SleepStyle 600 Series
Each OEM# HC 365S • HCPCS# A7046

Compatible with:

- HC604
- HC608



Humidifier Chamber for ICON Series



Item# **HUMHC900ICON200**
Replacement Chamber for the
Fisher & Paykel ICON Series
Each OEM# 900ICON200 • HCPCS# A7046

Compatible with:

- ICON Premo
- ICON Novo
- ICON Auto



Sunset Chinstraps

Universal Neoprene Chinstrap



Item# **CS001**

Universal Chinstrap

HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to itself at the top of the head
- Can be used with any nasal mask or pillows

Adjustable Neoprene Chinstrap



Size	Item#
6in	CS003S
7in	CS003R
8in	CS003M
9in	CS003L
10in	CS003XL

Compare to Puresom
Adjustable Ruby

OEM# TMS-09 ADJ
HCPCS# A7036

- Latex-free
- Thin adjustable strap attaches at the top of the head. Thicker bottom strap attaches to itself at the back of the head.
- Great for hard-to-fit patients

For fitting, measure along one side of the head, as indicated by the dashed line above.

Comfort Chinstrap



Item# **CS005**

Compare to
ResMed Sullivan

OEM# 16015
HCPCS# A7036

- Latex-free
- One size fits most
- Easy to use —attaches to itself at the top of the head
- Soft cloth chin cup adds comfort

Basic Chinstrap



Item# **CS002**

Compare to
Respironics Basic

OEM# 302175
HCPCS# A7036

- Latex-free
- One size fits most
- Attaches to the top strap of headgear
- Cups the chin for comfort
- Works best with any nasal mask

Deluxe Chinstrap



Size	Item#
26in	CS004S
28in	CS004L

Compare to
Respironics Deluxe

OEM# 302425
HCPCS# A7036

- Latex-free
- Attaches to itself at the top of the head, two smaller support straps keep the system in place
- Extra-wide strap distributes pressure to add comfort

Premium Chinstrap



Item# **CS006**


Compare to
Respironics Premium

OEM# 1012911
HCPCS# A7036

- One size fits most
- Attaches to itself at the top of the head to keep the chin lifted. Smaller support strap attaches to a split section of the chinstrap and wraps around the head for added stability.

Sunset Chinstraps and Headgear

Ruby Style Chinstrap



For fitting, measure completely around the head, as indicated by the dashed line above.

Size	Item#
21.5-22.5in	CS007S
22-23in	CS007M
23-24in	CS007L

Compare to Puresom Ruby
OEM# TMS-09 • HCPCS# A7036

- Latex-free
- One piece stretches to slide over the top of the head — no need to Velcro
- Can be used with any nasal mask or pillows

Halo Style Chinstrap



Item# **CS025**

Compare to BreatheWear
Halo Chinstrap

OEM# H3407
HCPCS# A7036

Ruby Style Adjustable Chinstrap



Size	Item#
Regular	CS007-ADJ
Extra Large	CS007-ADJXL

HCPCS# A7036

- Latex-free
- Fully adjustable
- Super soft materials for extra comfort



Universal Mesh Headgear



Item# **HG001**

Compare to Respiroics
SimpleStrap Headgear

OEM# 302328
HCPCS# A7035

- Latex-free
- One size fits most
- Can be used with most 4-point nasal or full face masks

Universal Headgear



Size	Length	Item#
Regular	22.5in	HG003
Extra Large	24.5in	HG003XL

HCPCS# A7035

- Latex-free
- Can be used with most 4-point nasal or full face masks

Cleaners

Sunset Cleaning Wipes — Flowpack

Item# **CAP1003S**

Flowpack of 64 wipes
12/case

- Ingredients are safe and gentle for daily cleaning of all types of CPAP and oxygen masks and equipment



Private Label CPAP Wipes

Order the flowpack wipes with your logo and contact info on a lid sticker, with a minimum quantity commitment. Talk to your sales rep for more info.



Citrus II Cleaning Wipes

Item# **CAP1003**

Canister of 62 wipes
12/case



Cleaning Wipes — Single Packs

Item# **CAP1005**

Citrus II Box of 12 single-pack wipes
12/case

- Convenient individually wrapped wipes are great for patient samples or travel use



CleanSmart CPAP Disinfectant Spray

Item# **CAP1018-L**

16oz • **6/case**

Item# **CAP1018-S**

20z • **9/case**

- Kills 99.9% of viruses & bacteria
- Included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2
- Active ingredient Hypochlorous Acid does not leave a chemical residue
- EPA Registration Number 89896-2
- Formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment



Citrus II Mask Cleaning Spray

Item# **CAP1001**

1.5oz • **24/case**

Item# **CAP1002**

8oz • **12/case**


- All-natural cleanser for daily mask cleaning
- Available in 1.5oz travel size or 8oz non-aerosol spray cans



Filter and Check Valve Kit for SoClean 2 CPAP Cleaning Device

Item# **CAP1007S-KIT**
Each

- With filter, check valve, retail packaging and instruction manual
- Replace every 6 months



CPAP Cash Items

CPAP Tube Cleaning Brush

Item#
CAP1013
Standard 6ft
Each

- Brushes all the way through a 6ft CPAP tube without leaving debris in the middle
- Stainless steel wire and soft nylon bristles



Handheld CPAP Pressure Meter

Item#
CAP003
Each

- Measures flow 0-18 cm/H2O
- Connects directly to the CPAP tubing without an adapter
- Lightweight device used to verify the pressure of a CPAP unit



SleepRes V-Com™

Item#
CAP1026
Each

NEW

Helping patients adjust to CPAP on their own

Adding V-Com to a patient's therapy is like Training Wheels for CPAP™. This simple device can mean the difference between your patient giving up on PAP or continuing their lifesaving treatment.

Within minutes, a respiratory therapist can show a patient what the PAP pressure feels like with and without the V-Com device inline. This feedback mechanism allows them to better understand how their airflow should feel.

As you know, comfort equals adherence and adherence equals better treatment outcomes.



The V-Com Difference

Improve early PAP tolerance with V-Com by providing your patients with initial comfort that doesn't compromise their treatment.

Traditionally, when a patient struggles with the constant flow of air pressure, pressure settings are reduced.

Reducing pressure settings to increase a patient's initial tolerance improves comfort, but at a cost.

Decreasing pressure settings also means decreasing expiratory positive airway pressure (EPAP) which is principally responsible for therapy.

The V-Com effect is flow-dependent with its maximum comforting effect during the high inspiratory flows. As a result, V-Com has minimal effect on expiratory positive airway pressure.

90-Degree CPAP Tubing Elbow

Item#
CAP1027
Each

NEW

The Sunset 90-Degree CPAP Tubing Elbow Adapter connects standard 22mm or slim 15mm CPAP tubes to most CPAP machines at a 90°, relieving stress on the tubing connection and creating a more convenient pathway to the bed.

OEM# 37394



CPAP Battery and Power Inverter/Charger

Portable Outlet UPS CPAP Battery

Item#

CAP10082

Includes battery, wall charger, and travel bag

- Autoswitch technology: When power goes out, your CPAP stays on
- Three-prong AC outlet
- Two USB ports
- Digital display
- Convenient travel bag
- Meets TSA/FAA battery guidelines
- 1-year warranty
- Rated for over 500 charge cycles
- Recharges in less than 4 hours
- May be used for other devices powered by AC or USB
- AC and USB ports may be used simultaneously



Specs

- 2.875 lbs
- 5.375" wide x 9" long x 1.5" tall
- Rated Capacity: 159Wh/14.4Ah
- Rated Voltage: 14.5V
- Charge Input: 16.8V/2A
- USB Output: 5V/4.2A
- AC Output: 110v AC 60hz/100w
- AC Adapter Input: 100V-240V/50-60Hz

DC to AC Power Inverter, Car Charger, Lighter Plug-in

Item# **CAPDC2AC150**

Recharge your Portable Outlet CPAP Battery using a car's 12V lighter outlet

- Can also power most CPAPs and other small electronics
- 150 Watts
- 3-prong outlet
- 2 USB ports



CPAP Cash Items

Transcend Micro™ Auto CPAP

Item#
CAPSOM506001
Each

NEW

- **Small:** About the size of a baseball
- **Lightweight:** Weighs less than a mobile phone
- **Works with any CPAP mask:** Use your existing mask
- **Quiet operation:** The Transcend WhisperSoft™ muffler is included

- Smallest and lightest travel CPAP device
- Compatible with any CPAP mask
- Unplug and sleep anywhere with more power options
- Adjustable GentleRise and AirRelief technology
- Stay dry, all the time
- AirMist HME – waterless humidification

Included:

- Transcend Micro CPAP
- Transcend Micro AirFlex™ 6-foot hose
- Transcend WhisperSoft™ Muffler
- AC power adapter
- 2 Year manufacturer's warranty



Transcend Micro™ Accessories

Item#

CAPSOM505001	PowerAway® Battery
CAPSOM505035	Plug Adapters (AU/EU/UK)
CAPSOM505010	Power Supply
CAPSOM505030	AirMist Starter Pack (2 filters and 1 holder)
CAPSOM505009	AirMist Pack (4 filters and 1 holder)
CAPSOM505014	AirMist 4 Pack
CAPSOM505015	AirMist 8 Pack
CAPSOM505016	AirMist 12 Pack
CAPSOM505008	Travel Bag
CAPSOM505017	PureFresh Filters (2 Pack)
CAPSOM505033	WhisperSoft Micro Muffler and AirFlex 6ft Slim Hose
CAPSOM505007	Filter Cover



Item# **CAPSOM505001** PowerAway® Battery

Extended Run Time

- **Long-Lasting Power:** Offers 2-3 nights or 17.7 hours of run time on a single charge at a pressure setting of 9cm, ensuring restful sleep without frequent recharging.
- **Quick Charge Time:** Recharges in less than 4 hours, so it's ready to go when you are.

Reliability and Safety

- **FAA-Approved:** Cleared for in-flight use, making it an essential for travelers
- **Backup Power:** Provides a permanent backup power source during outages, ensuring uninterrupted sleep therapy.
- **Long Life Span:** Capable of 1,000+ full charge/discharge cycles for long-term reliability.

Item#
CAPSOM505042 Essentials Pack

Included:

- NEW Micro WhisperSoft Muffler and Hose
- (2) Hose Clips
- AirMist Pack (4 filters and 1 holder)
- Travel Power Supply Plug Adapters
- Micro filters (2 pack)
- Transcend Micro Travel Pouch



Item#
CAPSOM505043 Resupply Pack

Included:

- NEW Micro WhisperSoft Muffler & Hose
- AirMist Pack (4 filters and 1 holder)
- Micro filters (2 pack)



CPAP Cash Items

Snugz CPAP Mask Liners

Item#
CAP7001


Machine-Washable Liners
for Most Full Face Masks
2/pack

CAP7002

Machine-Washable Liners
for Most Nasal Masks
2/pack

- 2 per pack, 90-day supply
- Machine washable
- Fits MOST nasal or full face CPAP masks (Not recommended for ResMed Mirage Quattro or Respirationics Amara View and Dreamwear)
- Blend of spandex, moisture-wicking polyester and lycra fibers
- Latex-free


Snugz is a machine washable, one-size-fits-most liner for CPAP masks that is designed to help reduce noisy leaks, prevent redness and irritation, and feel great against the skin. The liner helps to grip the skin and keep the mask in place around the nose and mouth, reducing the need for a chinstrap. The improved seal is especially beneficial for high-pressure/high-flow PAP patients and can reduce nightly sleep disturbances and AHI.



Soft Headgear Strap Covers

Item#
CAP2004

- Reduces strap marks
- No bulky and uncomfortable velcro
- Reduces strap irritation on sensitive skin
- Thick padding for extra comfort
- No need for clips – just slip wraps around headgear straps
- Easy and helpful for patients with arthritis in hands or fingers
- Can be machine washed and dried
- Wraps securely grip headgear and slide to adjust
- One size fits all
- Made from 230g Coral Velvet



RemZzzs CPAP Mask Liners

Item#
CAP3001(S,M,L)

Liners for Full Face ResMed
& Respirationics

OEM#	K1-FS	SMALL
	K2-FM	MEDIUM
	K3-FL	LARGE

CAP3002(S,M,L)

Liners for Full Face Fisher
& Paykel, Hans Rudolph,
and Sunset

OEM#	K4-FS	SMALL
	K5-FM	MEDIUM
	K6-FL	LARGE

CAP3003(S,M,L)

Liners for Nasal
Respirationics

OEM#	K7-NS	SMALL
	K8-NM	MEDIUM
	K9-NL	LARGE

CAP3004(S,M,L)

Liners for Nasal ResMed
and Fisher & Paykel

OEM#	K10-NS	SMALL
	K11-NM	MEDIUM
	K12-NL	LARGE

CAP3005L

Liners for Respirationics
FitLife Mask

OEM#	6A-FLK	SMALL AND LARGE
	6B-FLK	EXTRA LARGE

CAP3006


Liners for Nasal Pillows

OEM# 7A-NPK

CAP3007

Liners for Respirationics
DreamWear and Amara View

OEM# 6C-AVK



30-Day Supply

Boxes sold individually.
Full face and nasal mask
RemZzzs include 30 facial
liners and 2 forehead liners.

- Increases compliance by adding a comfortable and absorbent pad to CPAP mask
- Reduces air leaks and protects skin from red marks



CPAP Cash Items

CPAP Moisture Therapy

Item# **CAP6100L**
1oz Tube • 18/case

- Formulated to prevent chafing, irritation, dryness or cracking from using a CPAP mask or other delivery device
- Apply to the face at the beginning of treatment and as needed
- Petroleum-free with aloe vera, vitamins A & E, and argan oil



Pur-Sleep CPAP Aromatherapy

Item# **CAP5001**
VaporClear Aromatherapy Starter Kit
OEM# VAPOR

Kit Includes:

- 1 aluminum diffusion stand
- 4 reusable diffusion pads
- 1.5ml bottle of Clear (peppermint, lime)
- 1.5ml bottle of Calm (eucalyptus, lavender, rosemary)



30ml Aromatic Refills

Scent/OEM#	Item#
Black Cherry	CAP5004-BCH
Calm	CAP5004-CLM
Clear	CAP5004-CLR
Dark Chocolate	CAP5004-DCH
Deep	CAP5004-DEE
French Vanilla	CAP5004-FRV
Fresh	CAP5004-FSH
Peace	CAP5004-PCE
Spice	CAP5004-SPC
Strawberry	CAP5004-STR

Item# **CAP5002**
Reusable Diffusion Pads
10/pack • OEM# DP10


Item# **CAP5003**
Replacement Aluminum Diffusion Tray
Each • OEM# TRAY



- Aromatherapy products specifically for use with CPAP
- Set up the diffuser right below the CPAP filter intake, add a drop of essential oil to the pad and the fragrance is drawn into the CPAP machine

Quick Tube Release Connector

Item# **CAP1016**
Each



Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item# **BF2200**
22mm
5/pack
Attach between the machine and tube for additional protection against bacteria and allergens



Circuit Connectors

Item# **RES013**
22mm OD x
15mm ID/
22mm OD
10/pack



Item# **RES014**
22mm OD x
22mm OD
With Lip
10/pack



Humidifier Tubing Adapter

Item# **RES016**
22mm ID x 22mm OD • with Lip
10/pack



Oxygen Bleed-In Adapter

Item# **RES019**
22mm with 6mm
Rotating Stem 10/pack



RES019-CAP
With Stem Cap 10/pack



Swivel Exhalation Port

Item# **VEN002**
For CPAP and
Ventilator Circuits



CPAP Cash Items

CPAP Pillows

CPAPMax Pillow

Item# **CAP4003**

Each OEM# 15-501R-1-920



CPAPMax Replacement Cover
CAP4003RC

Each OEM# 1-630-101-1-933

BEFORE

Traditional pillows interfere with mask performance and can cause a stiff neck, shifting, leaks, and facial discomfort.



AFTER

The CPAPMax Pillow offers superb mask accommodation, excellent support, and unparalleled comfort.



Cool 3D Mesh Air-Flow
Reduces heat and perspiration



Fiberfill
Soft microfiber pillow top with hypoallergenic fiberfill

10-in-1 Flip Pillow

Item# **CAP4006**

Each



The most versatile wedge and bed pillow available!

10 support positions with 1 pillow! Comfort and align your neck, shoulder and lumbar region when used as back support cushion. Use it as a knee or leg pillow wedge or even in an office chair for added support. Filled with a soft hypo-allergenic fiber-fill.

10-in-1 Flip Pillowcase Cover

CAP4006CN Navy

CAP4006CW White

Each




CoolPAP Pillow

Item# **CAP4007**

Each

Works with any style and brand of mask



CoolPAP Pillow Cover

CAP4007C

Each



CPAP Bedside Storage Bag

Item#

CAP1015

- Well-made and stylish
- Designed to hold CPAP supplies such as a mask and tube
- Weighted lid allows bag to hang off top of nightstand or tabletop



CPAP Cash Items

Deluxe Tubing Wrap

Item#

CAP20026ft Nylon
Tubing Wrap**Each**Compare to
ResMed
Tubing Wrap

- 6ft nylon wrap with fleece lining maximizes heat retention and reduces noise
- Zipper closure for ease of use



Comfort Tubing Wrap

Item#

CAP20016ft Velour
Tubing Wrap**Each**Compare to
Respironics
Tubing Wrap

- 6ft plush velour material adds comfort while preventing rainout
- Zipper closure for ease of use

**Both Tubing Wraps
come in attractive
retail packaging**

CPAP HCPCS and Replacement Schedule

HCPCS	Item	Recommended Replacement Schedule
A4604	Heated Circuit Tubing	1 per 3 months
A7027	Combination Oral/Nasal Mask	1 per 3 months
A7028	Replacement Oral Cushion for Combination Oral/Nasal Mask	2 per month
A7029	Replacement Nasal Pillows for Combination Oral/Nasal Mask — Pair	2 per month
A7030	Full Face Mask	1 per 3 months
A7031	Replacement Full Face Mask Cushion or Liner	1 per month
A7032	Replacement Nasal Cushion or Liner	2 per month
A7033	Replacement Pillows — Pair	2 per month
A7034	Nasal or Pillows Mask	1 per 3 months
A7035	Headgear	1 per 6 months
A7036	Chinstrap	1 per 6 months
A7037	Tubing	1 per 3 months
A7038	Filter — Disposable	2 per month
A7039	Filter — Nondisposable	1 per 6 months
A7044	Oral Interface	
A7045	Replacement Exhalation Port	
A7046	Replacement Water Chamber	1 per 6 months
A9279	Monitoring Feature/Device	
E0470	BiLevel Device without Backup Rate	
E0471	BiLevel Device with Backup Rate	
E0561	Humidifier — Non-Heated	
E0562	Humidifier — Heated	
E0601	CPAP Device	

Due to the rapid changes in reimbursement, this information is for reference only. Sunset Healthcare Solutions cannot guarantee its comprehensiveness or accuracy. Sunset urges you to contact a reimbursement specialist for up-to-date information and advice.

Oxygen Products

Oxygen	Page
Pulse Oximeter	38
Oxygen Sensors	39
Firesafe	40
Concentrators	41-42
Masks	43
Cannula and Cannula Accessories	44-45
Tubing	45
Microdial Flowmeter	46
Regulators and Conserver	47
Cylinder Carriers	49-50
Connectors	51
Accessories and Cleaners	51-53
Analyzers	54
Celon Wet Bands	54
Bacteria and HEPA Filters	55
Concentrator Cabinet Filters	56-62
Caire/AirSep/SeQual	56
CHAD	57
Covidien Nellcor Puritan Bennett	57
Dynarex	57
OxyGo	57
DeVilbiss/Drive	58
Invacare	59
Rhythm P2	59
Rhythm HC, Jumao, Attendant, Lifestyle, Medacure	60
Nidek	60
React Health/3B	60
Respironics	61
Healthdyne	62
Replacement Parts	63-64

Pulse Oximeter

Finger Pulse Oximeter

Item# **RES5102**

- Oversized 1.3" OLED display face is the largest of any finger pulse oximeter
- Features large numeric font
- Five grades of brightness adjustment
- Displays SpO2, PR and PI parameters and pulse wave graph
- Offers buzzer, demo and visual alarm functions
- Easy to navigate menu sets alarm limits or switches buzzer on/off
- CE, ISO and FDA approved
- Includes lanyard and manual




SpO2 Sensors

NEW


Disposable SpO2 Sensor Neonatal/Adult

Item# **RES7100**
Similar to Nellcor MAX-N
20/box



Neo-M-1.5 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7200**
1.5ft cord • Non-replaceable
20/box



Neo-M-3 Disposable SpO2 Sensor Neonatal/Adult (<3 kg or >40 kg)

Item# **RES7201**
3ft cord • Non-replaceable
20/box




Coban Tape Roll

Item# **RES7202CT**
Use for sensor life
extension or sensor
cord stability
Each roll 1" x 177"



Infant-M-1.5 Disposable SpO2 Sensor Infant (3-20 kg)

Item# **RES7203**
1.5ft cord
20/box



Pediatric-M-1.5 Disposable SpO2 Sensor Pediatric (10-50 kg)

Item# **RES7204**
1.5ft cord
20/box



Firesafe

Firesafe Thermal Shut-Off Valve


Item#
RES010
Each

Bidirectional



- Now smaller
- Laser etched expiration date
- Human readable UDI
- Designed to last 5 years

- A thermal fuse designed to stop the flow of gas in the event that the downstream oxygen tubing is ignited
- Helps protect against oxygen tubing fires
- Helps keep patients safe no matter what method of oxygen delivery they use
- Offers a cost-effective solution to lower the risk of serious patient injury in an oxygen fire



Home oxygen fires represent a much higher risk in the United States than the United Kingdom, where the fitting of thermal fuses for oxygen is mandatory.

Fatalities per 100,000 patients



* 1 death was reported among 73% of the patient population (2013-17)

FireSafe Kits

Item#
RES010KIT
2 Firesafe Valves
1 Tubing Connector

Item#
RES010CANNULAKIT
2 Firesafe Valves
1 Tubing Connector
1 Adult Soft Cannula

Item#
RES010VAKIT
2 Firesafe Cannula Valves
1 Tubing Connector
1 Cannula with Female Port
1 7ft Oxygen Supply Tube



NEW

FireSafe Tubing Connector

Item#
RES010TC
100/case

RES010TC-10
10/case



FireSafe Nozzle

Item#
RES010FND-20
20/pack



Belluscura Oxygen Concentrators

X-PLOR® Portable Oxygen Concentrator

Item#

RES52051001**X-PLOR****NEW****RES52051002****X-PLOR with Extra Battery**

Lightweight and compact, the X-PLOR® is ideal for patients who want to replace their heavy oxygen tanks without missing a moment. The X-PLOR is also Bluetooth®-capable and can be connected to your mobile device through Belluscura's newly developed Nomad Biometrics™ App. With a rechargeable battery, easy-to-use interface, and high oxygen output, supplementary oxygen has never been more approachable or affordable.



Accessories

Item#

RES50073 User-Replaceable Sieve Bed Cartridge**RES50188** 8-Cell Battery**RES53012** DC Power Adaptor**RES53011** AC Power Adaptor**RES52191** Carry Bag**RES52291** Backpack

What's Included

- X-PLOR® Unit
- 8-Cell Battery
- User-Replaceable Sieve Bed Cartridge
- Air Intake Filter
- AC Power Adaptor
- DC Power Adaptor
- Carry Bag With Handle & Shoulder Strap

DISCOV-R™ Continuous & Pulse Flow Portable Oxygen Concentrator

NEW

Item#

RES70051004**DISCOV-R with 1 Battery (without Cart and Backpack)****RES70051005****DISCOV-R with 2 Batteries (without Cart and Backpack)****RES70051001****DISCOV-R with 1 Battery (with Cart and Backpack)****RES70051002****DISCOV-R with 2 Batteries (with Cart and Backpack)**

Airgonomically™ designed using Belluscura's patented ModulAir™ technology, the compact DISCOV-R weighs less than 6.5 lbs, producing more oxygen per pound than any other portable oxygen concentrator.

- 4 continuous flow settings (0.5, 1.0, 1.5, 2.0 LPM)
- 8 pulse flow settings
- Weighs only 6.5 lbs
- FAA compliant for air travel

Replacement Parts

Item#

RES73011 AC Power Adapter**RES73021** DC Power Adapter

What's Included

- DISCOV-R™ Unit
- 8-Cell Battery OR Two (2) 8-Cell Batteries
- Pair of User-Replaceable Cartridges
- Air Intake Filter
- AC Power Adapter
- DC Power Adapter
- Carry Bag with Handle & Shoulder Strap
- Standard 3-Year Warranty



Oxygen Concentrators

Drive DeVilbiss Oxygen Concentrators

Item#
RES525DS
5 Liter
OEM# 525DS

Item#
RES1025DS
10 Liter
OEM# 1025DS



- Built upon a legacy of reliable, field-proven concentrator models, the DeVilbiss 5 Liter and 10 Liter Oxygen Concentrators were designed with enhanced quality, durability and simplicity in mind
- Patented DeVilbiss Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator
- Exclusive DeVilbiss OSD (Oxygen Sensing Device), standard on every unit, ensures patient safety and reliability for longer service intervals
- Readily accessible patient controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms Audible and visual high/low pressure, low flow, low oxygen, power fail, Oxygen Sensing Device

Drive DeVilbiss iG02 Portable Oxygen Concentrator

Item#
RES125D
iG02

Item#
RES125D-XB
iG02 with
Extra Battery



- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market

Covers for Concentrators



25×15x30"
1.5mil
250 bags/roll

Item#	Color
EC251530C	Clear
EC251530T	Tan Tint

Oxygen Masks

Simple O2 Mask

With 7ft Safety Tubing

Item# **RES2100**

Adult Size

50/case

Item# **RES2200V**

Pediatric Size

50/case

HCPSC# A4620



- 5-12 LPM
- Medium concentration
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Partial Rebreather / Non-Rebreather O2 Mask with 7ft Tubing

Oxygen High Concentration
Mask with 7ft Safety Tubing
and Reservoir Bag

Item# **RES2101V**

Adult Size

50/case

OEM# 8120-7-50

HCPSC# A4620



- 5-15 LPM
- Safety vent
- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal

Non-Rebreather O2 Mask

With 7ft Safety Tubing
& Reservoir Bag

Item# **RES2102**

Adult Size

50/case

Item# **RES2202V**

Pediatric Size

50/case

HCPSC# A4620



- 10+ LPM
- Low-resistance check valves prevent rebreathing and allow exhaled gas to escape

Oxygen Cannula

Cannula with Tubing

Adult Cannula with Tubing

Length	Material	Item#	Qty
4ft	Standard	RES1104	50/case
4ft	Soft	RES1104S	50/case
7ft	Standard	RES1107	50/case
7ft	Soft	RES1107S	50/case
15ft	Standard	RES1115	50/case
25ft	Standard	RES1125	25/case
25ft	Soft	RES1125S	25/case



Item#

RES1207V

Pediatric Cannula
with 7ft Supply Tube
50/case

**RES1307V**

Infant Cannula
with 7ft Supply Tube
50/case



HCPCS# A4615


- Clean, fresh look
- 1-6 LPM
- Curved nasal prongs for added comfort
- 6-channel safety tubing designed to prevent kinks
- Soft, lightweight material is memory free for long-term comfort

Cannula without Tubing

Adult Cannula without Tubing

Port	Material	Item#	Qty
Male	Standard	RES1100	50/case
Male	Soft	RES1100S	50/case
Female	Soft	RES1100FS	50/case

HCPCS# A4615



Dual-Lumen (Demand) Cannula

Item# **RES4147V**

7ft

25/case

HCPCS# A4615



- 1-6 LPM
- Designed for use with dual port oxygen conservation delivery systems
- Dual tubing set provides independent tubes for oxygen delivery and sensing of patient breathing
- Smoothbore tubing

High Flow Soft Cannula

Length	Item#	Qty
4ft	RES1104SHF	25/case
7ft	RES1107SHF	25/case
15ft	RES1115SHF	25/case
25ft	RES1125SHF	25/case

HCPCS# A4615



- Now softer
- 6-15 LPM
- Maintains saturation for high-flow patients without a mask
- 6-channel supply tubing for safety

Oxygen Cannula

Humidoflo High Flow Nasal Cannula

10/case

Size	Item#
Small	HFN001AS
Medium	HFN001AM
Large	HFN001AL

With F&P Adapter

Size	Item#
Small	HFN001ASA
Medium	HFN001AMA
Large	HFN001ALA

HFN001-FPA F&P Adapter



With Nasal High Flow (NHF) therapy, patients benefit from warmed, humidified oxygen at a flow rate greater than those seen with conventional oxygen delivery devices, and more than peak inspiratory airflow. NHF treatment is better tolerated due to the heated and humidified air. The NHF apparatus is also smaller and more comfortable on the face than an NIV mask. Patients can continue eating, drinking and talking during treatment, leading to higher compliance. Delivery respiratory gases range from 2 – 60 L/min and O₂ concentrations range from 21 – 100%. NHF is intended for patients in homes and long-term care facilities and may help prevent intubation.

- Tube Side Quick Switch connects tube to either side of the patient
- Kink-free prongs and flexible, low-noise tube
- Compatible with any ventilator with a standard 22mm ISO connector, including Respirationics V60

- Suitable alternative to Philips AC611
- Can be used with F&P 850 and 950 with HFN001-FPA adapter
- Comparable to F&P OptiFlow cannula (OPT942, OPT944, OPT946)

See page 77 for
LUIA Ventilator




Secure Cannula Stickers

Item# **RES080**

25 pairs/case

- Holds cannula in place
- Easy to use
- Flexible, sheer fabric
- Designed to let skin breathe and peel off without leaving residue
- Latex-free
- Individually wrapped pairs



Comfort Cannula with Ear Cushions

Item# **RES1107EC**

7ft Cannula with
Ear Cushions

25/case

HCPCS# A4615

Item# **RES020**

Cushions Only
50 Pairs/case



- Add cushions to existing cannula

Oxygen Tubing

Oxygen Supply Tube



Length	Item#	Qty
7ft	RES3007	50/case
15ft	RES3015	50/case
25ft	RES3025	25/case
40ft	RES3040	20/case
50ft	RES3050	20/case

HCPCS# A4616

- Crush resistant with low memory to prevent kinking
- Great durability and comfort level

Oxygen Supply Tube — Green




Length	Item#	Qty
7ft	RES3007G	50/case
15ft	RES3015G	50/case
25ft	RES3025G	25/case
40ft	RES3040G	20/case
50ft	RES3050G	20/case

HCPCS# A4616

- Latex-free
- 1-6 LPM
- 6 channels

Microdial Flowmeter

Microflow Microdial Flowmeter



DISS Inlet • Barbed Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160048			A
RES8160050	400 kPa (58 psi)	350 – 550 kPa	E



DISS Inlet • DISS Outlet

Item#	Inlet Pressure	Inlet Pressure Range	Flow Range*
RES8160049	400 kPa (58 psi)	345 – 500 kPa	A
RES8160051	350 kPa (51 psi)	250 – 450 kPa	E

*Volumetric flow rates in standard liters per minute

Dial Position	Off	1	2	3	4	5	6	7	8	9	10	11
Range A	0	0.02	0.03	0.05	0.08	0.12	0.2	0.3	0.5	0.75	1	3
Range E	0	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09	0.1	1

Assured, repeatable accuracy to enable safe, effective weaning.

What is Microflow used for?

Microflow is primarily used to supply supplementary oxygen to infants with bronchopulmonary dysplasia (BPD).

BPD is common in very low birth weight and pre-term infants that have needed mechanical ventilation or have been treated with high oxygen concentrations due to respiratory distress syndrome (RDS)/chronic lung disease (CLD).

The risk of RDS/CLD is inversely related to birth weight and gestational age.

Infants may need very low supplementary oxygen for many months to compensate for underdeveloped lungs.



Microflow benefits

- Enables baby to receive oxygen therapy at home, putting families together sooner and increasing parent satisfaction
- Less time in hospital, lower cost of care
- Accurate, repeatable steps
- Micro step changes in flow
- Enables low flow rates to achieve target saturation
- Enables oxygen flow to be weaned

Oxygen Regulators and Conserver

Mini Oxygen Regulator CGA 870 Style



For use with E size or smaller cylinders: A, B, C, D, JD, E, M2, M4, M6, ML6, M7, M9, M15, M22, or M24.



Item# RES042
Adult • Barbed
0-8 LPM • **Each**
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

Item# RES043
Adult • Barbed
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

Item# RES046
Pediatric • Barbed
0-4 LPM • **Each**
0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$, $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1, $1\frac{1}{2}$, 2, 3, 4 LPM

HPCPS# E1353



Pediatric models have a green dial cap

Mini Oxygen Regulator CGA 540 Style



For use with M60 and larger cylinders: M60, M90, M, MM, M122, M250, H, J, T or other cylinders with CGA 540 style valve.

Item# RES044
Adult • DISS
0-8 LPM • **Each**
0, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8 LPM

Item# RES045
Adult • DISS
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

Item# RES047
Pediatric • DISS
0-4 LPM • **Each**
0, $\frac{1}{32}$, $\frac{1}{16}$, $\frac{1}{8}$, $\frac{1}{4}$, $\frac{1}{2}$, $\frac{3}{4}$, 1, $1\frac{1}{2}$, 2, 3, 4 LPM

HPCPS# E1353



Both Regulators Feature

- Pressure gauge reads to 3000 PSI
- Easy-click flow selector knob
- 2-year warranty



Private label laser etching now available!

Ask your rep for details.

Oxygen Regulator CGA 540 Style Long Stem



Item# RES045L
Adult • DISS
0-15 LPM • **Each**
0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 LPM

Bonsai Pneumatic Conserver

Item# RES049
Single Lumen • Barbed Outlet
Each



- Delivers a uniform oxygen pulse in increments: 1, 2, 3, 4, 5, and 6
- Two continuous flow (CF) settings preset at 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14-40 breaths per minute
- Senses breath/delivers oxygen in first half of inspiratory cycle
- Operating pressure range: 200-3000 PSI
- Light weight and quiet operation
- Easy to operate single selector knob
- Compact design fits easily in virtually any oxygen carry bag
- Limited 2-year warranty

Bags for Regulators and Conserver



Item# EC1012C
10×12"
1.5mil, Clear
1000 bags/box

Sample Oxygen Cylinder Specifications*

Use the handy chart below to determine which accessories to purchase for your oxygen cylinder.

* This information is intended to be used as a guide. Dimensions may vary by manufacturer. Sunset Healthcare Solutions does not sell oxygen tanks.



Older Name		A	B			C	D	JD	E		
Newer Name	M2	M4	M6	ML6	M7	M9	M15	M22	M24	M60	M/MM/M122
Diameter (in.)	2.5	3.2	3.2	4.3	4.3	4.3	4.3	5.3	4.3	7.3	8
Height (in.)	5.3	8.5	11.5	7.6	9.1	11	16.5	16.5	25.5	23	36
Empty Weight (lb)	0.7	1.6	2.2	2.8	3.3	3.7	5.3	8	7.9	22.3	39.5
Capacity (L)	42	113	164	165	198	255	425	640	680	1738	3455
Transport Method	CSBM6, CSBM6BP	CSBM6, CSBM6BP	CSBM6, CSBM6BP, CSBM6Flex, CSBM6-Carrier	CSBML6M7, CSBML6-Hardshell	CSBML6M7	CSBM9, CSBM6BP, CSBM6Flex, RES038B	CSBD, CSBD8P, CSBD1Roller, CSBD2Roller, RES038B, RES060		CSBE, CSBD1Roller, CSBD2Roller, RES038B, RES060	Not Portable	Not Portable
Regulator Type	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES042, RES043	RES044, RES045	RES044, RES045

Oxygen HCPCS

HCPCS	Item
A4615	Cannula, nasal
A4616	Oxygen tubing, per foot
A4620	Oxygen mask, variable concentration
E0424	Stationary compressed gaseous oxygen system, rental
E0431	Portable gaseous oxygen system, rental
E0443	Portable oxygen contents, gaseous (for use only with portable gaseous systems when no stationary gas or liquid system is used), 1 month supply = 1 unit
E0445	Oximeter device for measuring blood oxygen levels non-invasively
E0560	Humidifier, durable for supplemental humidification during IPPB treatment or oxygen delivery
E1353	Oxygen supply regulator
E1354	Oxygen supply wheeled cart
E1355	Oxygen supply stand/rack
E1390	Oxygen concentrator, single delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1391	Oxygen concentrator, dual delivery port, capable of delivering 85% or greater oxygen concentration at the prescribed flow rate
E1392	Portable oxygen concentrator, rental
K0738	Portable gaseous oxygen system, rental; home compressor used to fill portable oxygen cylinders; includes portable containers, regulators, flow meter, humidifier, cannula or mask and tubing

Cylinder Bags

Comfort Shoulder Bags

Cylinder	Item#
ML6/M7	CSBML6M7
C/M9	CSBM9
B/M6	CSBM6
D	CSBD

- Comfort flex shoulder strap
- Double zipper for easy access
- Machine washable



CASH ITEM



CSBML6M7

CSBM9

CSBM6

CSBD

Flex O2 Bag

Item#
CSBM6Flex
Reconfigurable Bag
for M6 and M9
Cylinders
Each

- Adjusts to 5 carrying modes



CASH ITEM



Backpack

Cylinder	Item#
B/M6/C/M9	CSBM6BP
D	CSBDBP

- Fully adjustable and comfortable backpack straps
- Machine washable
- Converts to either backpack or shoulder bag



CASH ITEM



Wheelchair/Scooter Carrier

Item# **CSBE**
For E Cylinders
Each

- Multi-part stabilizer straps for easy chair attachment
- Machine washable



CASH ITEM



Cylinder Carriers

Single Tank Rolling Bag/Backpack

Item#
CSBD1Roller
For D/E Cylinders
Each



Double Tank Rolling/Shoulder Bag

Item#
CSBD2Roller
For D/E Cylinders
Each

- May also be used with shoulder strap



ML6 Hardshell Tote

Item#
CSBML6-Hardshell
For Homefill Cylinders
Each



Cylinder Rack

Item# **RES038B**
For D, E, or M9 Cylinders —
Lightweight, Durable
Plastic Composite Rack
Each



Oxygen Cylinder Carrier

Item#
CSBM6-Carrier
Soft Polyester Tote
for M6 Cylinders
Each

- Tote handles
- Padded shoulder strap
- Fastening strap for mobility
- Holds up to 8 M6 tanks



D/E Cart

Item#
RES060
4/case

- Accommodates one D or E tank
- Sturdy yet lightweight (4.18 LBS welded metal construction)
- Adjustable height handle with PVC grip
- Large EVA Foam wheels have just the right amount of resilience for good traction and maneuverability
- Elegant, grey satin finish
- Bottom fork rest with rubber caps



Oxygen Connectors and Accessories

Christmas Tree Connector

Item# **RES001**
Green without Swivel
10/pack



Item# **RES002**
Green with Swivel
10/pack



- Tapered and barbed nut and nipple assembly allows secure tubing connections
- For standard DISS threaded oxygen outlets
- Improved threading is easier to connect to regulators or flow meters

Oxygen Tubing Y-Connector

Item# **RES017**
Fits 3/16in and 1/4in ID Tubing
(All 3 sides are the same)
10/pack



- Connects 2 oxygen sources for higher flows

Oxygen In Use Signs

Item# **RES003**
6in x 9in
With String
50/pack



- Double-sided (English and Spanish)

Washers with Rubber Ring

Item# **RES035**
Brass
25/pack



Item# **RES036**
Aluminum
25/pack



Oxygen Accessory Kit

Item# **RES1000**
RES002 Christmas Tree Connector - Swivel
RES012 Tubing Connector - Non-Swivel
RES017 Tubing Y-Connector
RES018 Tubing Connector - Swivel
RES032 Wrench
RES035 Brass Washer

Each



CUSTOMIZABLE KIT



- Order all common oxygen accessories in one easy-to-handle kit
- Contents can be customized to fit your needs

Oxygen Tubing Connector — Non-Swivel

Item# **RES012**
Fits 3/16in and 1/4in
ID Tubing
10/pack



- Barbed ends ensure secure, leak-free connections

Oxygen Tubing Connector — Swivel

Item# **RES018**
Male/Male **10/pack**

RES018-25 25/pack

RES018KIT-1
Female/Male **Each**



- Enhance patient mobility and comfort by minimizing twisting and tangling of oxygen supply tubing

Oxygen Accessories

O2 Starter Kit

Item# RES02KIT

RES1107 Cannula with 7ft Tube x2
 RES3015 15ft Supply Tube x2
 RES002 Connector — Swivel
 RES003 Oxygen In Use Sign
 RES012 Tubing Connector — Non-Swivel
 RES017 Tubing Y-Connector
 RES018 Tubing Connector — Swivel
 RES032 Wrench
 RES035 Brass Washer



CUSTOMIZABLE KIT



Wrenches and Chains

Item# RES031

With Chain **10/pack**

Item# RES032

Without Chain **10/pack**



- For use with post valve oxygen cylinders: A, B, C, D, JD or E tanks
- Made of tough, stress-resistant nylon

Item# RES031M

Metal Wrench with Chain for Post Valve Oxygen Cylinders and CGA 540 Style Regulators **5/pack**



Item# RES031-CHAIN

12in Ball Chain with Coupling Connector **50/pack**

- Add chain to an oxygen valve wrench



NEW

Item# RES031-BW

Cylinder Wrench with Two-Ring Bungee and Two Yoke Washers **10/pack**



Disposable Bubble Humidifier 500 ml

Item#

RES050

Up to 15 LPM • 6 psi • 500 ml

50/case

HCPCS# E0560

OEM# HUM001



- Multi-threaded lid ensures a secure fit, preventing leaks and service calls
- Audible back pressure alarm activates at 6 psi
- No aluminum fitting in stem eliminates aluminum shreds often found in lower quality models
- Soft plastic diffuser reduces clogging due to mineral build-up
- Translucent jar provides improved visibility of jar contents
- Humidifier bottle lid is colored green/yellow, making it easy to tell patients on the phone how to change the bottle
- Bottle fill line is very easy to see
- Finger grips at the top of the bottle make the lid easier to screw on/off for the patient

Humidifier Connectors

Item# HUM008

L-Connector 12in • **10/pk**

Item# HUM008L

L-Connector 21in • **10/pk**

Item# HUM008S

Straight Connector 9in • **10/pk**



Water Traps

Item# RES037

Push On End Caps • **Each**

Item# RES037T

Threaded End Caps • **Each**



- Keeps water out of cannula
- Connects between cannula and supply tubing
- Removable end caps allows drainage of accumulated water
- 9cc water capacity in horizontal position

Accessories and Cleaners

Homefill Tubing Connector


Item# **RES039**
10/pack

- Compatible with Invacare, Drive, and Respironics homefill machines
- Durable, coiled tubing




Homefill Racks


Item# **RES039-RACKB**
For Invacare Perfecto2
4/case



RES039-RACKW
For Invacare Platinum
4/case



RES039-RACKD
For DeVilbiss 525/1025 Series
4/case



Easily attaches onto concentrator and is designed to conveniently hold Homefill and transfill oxygen systems.

Oxygen Liter Meter

Item# **RES034**
0-10 Liter/Minute
Each



- Oxygen pen
- Works with tanks and concentrators

RoEzIt Dermal Care

Item# **RES6104L**
1oz Tube • 18/case

- Apply before beginning oxygen therapy and as needed during treatment
- Moisturizes the nose, lip, and over the ear where friction may cause discomfort from cannula
- Petroleum-free with aloe vera and vitamins A and E



CASH ITEM

Germicidal Deodorizing Cleaner

Item# **RES6100**
22oz spray bottles
12/case

- Ready-to-use spray cleans, disinfects and deodorizes
- Formulated for use on hard, non-porous surfaces
- Bactericidal, Tuberculocidal, Fungicidal, Mildewstat
- Contains no abrasives, bleach, phosphates or phosphorous
- Effective against SARS-CoV-2 (COVID 19 virus) on hard non-porous surfaces



Germicidal Deodorizing Cleaner — Gallon

Item# **RES6102**
1 Gallon

- No-scrub, hospital strength germicidal cleaner
- Cleans, deodorizes, and disinfects surfaces, equipment and non-critical instruments



Oxygen Accessories

Smart Check Analyzer

Item#
RES7003
Each

NEW

The ultimate oxygen + flow + pressure analyzer.

The Smart Check is an all-in-one oxygen analyzer for testing continuous and pulse (POC) concentrators. Measure Oxygen, flow, pressure and pulse volume with a single device. Use the Smart Check to also test CPAP, BIPAP and nebulizers.

The Smart Check is the only device that also offers a built-in flashlight, backlight, and integrated tube management system.

With its fast temperature response, SmartCheck quickly adapts to temperature changes, ensuring accurate results in less time.

Offering all these features at a very competitive price, SmartCheck delivers exceptional performance at an incredible value.



Features:

- LED Display
- Calibration Check Button
- Sample Inlet
- Flashlight
- Sample Outlet
- Battery Cover

Maxtec UltraMax O2 Oxygen Analyzer

Item#
RES7000
Each

NEW

The UltraMaxO2 helps oxygen concentrator service technicians looking to save on costs and time when checking patient O2 concentrators.

With integrated oxygen, flow, and outlet pressure monitoring in a single device, this handheld unit is easy to operate, store, and transport and does not require a traditional electrochemical sensor, which reduces the overall maintenance and minimizes the cost of ownership.



- Reduced cost of ownership with the ultrasonic sensor
- No in-field calibration required
- Quick set-up and readouts
- Ability to check outlet pressure

Celon Wet Bands

NEW

Create a perfect seal to protect oxygen tank valves

**Easy to apply, safe and effective
Biodegradable and sustainable material**

Item#
WB66120
66 x 120mm
for CGA 540 Valves
Bucket of 1,500



Item#
WB30540
30.5 x 40mm
for CGA 870 Valves
Bucket of 11,500



Item#
WB30525
30.5 x 25mm
for CGA 870 Valves
Bucket of 18,000



- Eliminates dust, dirt & dampness
- Hand applied
- Heat-free solution
- Creates a secure and wrinkle-free seal

- Professional finish
- Environmentally friendly
- Biodegradable

Oxygen Bacteria And HEPA Filters

6in HEPA Capsule with Removable Sticker

Item# **BFH6S**
Threaded
Each



EverFlo Intake Filter

Item# **BF1001**
Each
Respironics OEM# 1038831



Extended Life Bacteria Filter

Item# **BF5000**
Nonthreaded
Each



Extended Life Intake Bacteria Filter

Item# **BF1103**
Each



Final Bacteria Filter

Item# **BF516**
1/4in Barbed on Both Sides
10/pack



Respironics SimplyGo Style Inlet Filter

Item# **BF8000**
Each




HEPA Filter for React Health/3B Stratus

Item# **BF6020**
HEPA Filter
Each
OEM# STR1020




Intake Bacteria Filter

Item# **BF2222**
22mm fittings
5/pack



Extended Life Bacteria Filter with Muffler

Item# **OF6011**
Each



Extended Life Intake Bacteria Filter

Item# **OF1025D**
Each



Intake Bacteria Filter

Item# **BF7000**
Jumao JMC5A, Attendant 8X899,
Lifestyle LM5A, Medacure OC5L,
Rhythm LM5CA, Drive DECO2 5L
Concentrators
Each



Final Bacteria Filter — 1/4in to 3/8in

Item# **BF1438**
1/4in to 3/8in
Stepped Barb
on Both Sides
10/pack



Final Bacteria Filter — 3/16in

Item# **BF316**
3/16 Stepped Barb and 1/4in
Straight Barb on Both Sides
10/pack



Millennium M10 Style HEPA Filter

Item# **BF1002**
Long Life Filter
Each



Caire/AirSep/SeQual

SeQual Integra



Item# **OF11005**
Foam Cabinet Filter
10/pack
75/8" x 35/8"-5 x 3/4"
OEM# 2185



Item# **BF5000**
Extended Life
Bacteria Filter
Each
OEM# 2033



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack

AirSep Newlife, QuietLife, Intensity



Item# **OF8002**
Foam Cabinet Filter
10/pack
51/4" x 47/8" x 7/8"
OEM# F1-50-08-05



Item# **BF1438**
Final Bacteria Filter
1/4in to 3/8in
Barbed
10/pack
OEM# F1-30-03



Item# **OF8001**
Felt Filter
10/pack
17/8" OD x 3/4" ID x 1"
OEM# M1 161-1



KIT


Item# **OF8002SP**
Newlife Kit
2 - OF8002
1 - BF1438
1 - OF8001
Each

CAIRE Companion 5



Item# **BF1103**
Extended Life Intake
Bacteria Filter
Each
OEM# 2185

SeQual Eclipse



Item# **OF11006**
Foam Cabinet Filter
10/pack
21/2" x 21/2" x 3/4"

CHAD

Total 02



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
Healthdyne OEM# 100-6631-00
Respironics OEM# H620




Item# **OF3003**
Foam Intake Filter
for use with BF2222
25/8" OD x 3/4" ID x 23/8"
10/pack
OEM# H613

Dynarex

Resp-02



Item# **OF6013-10**
Foam Cabinet Filter
5.75" x 3.625"
10/pack



Item# **OF6012**
Air Intake Filter
Each
OEM# RPOC-14AIF

Covidien Nellcor Puritan Bennett

590 and 592A



Item# **OF6001**
Foam Cabinet Filter
63/4" x 41/4" x 3/8"
10/pack
OEM# 492672



Item# **BF5000**
Extended Life
Bacteria Filter
Each



KIT

Item# **OF6001SP**
590 and 592 Kit
2 - OF6001
1 - BF5000
1 - BF316
Each



Item# **BF316**
Final Bacteria Filter
3/16 Stepped Barb and
1/4in Straight Barb
10/pack
OEM# 492141

OxyGo

OxyHome



Item# **OF8006-10**
Foam Cabinet Filter
63/4" x 41/4" x 3/8"
10/pack
OEM# 1170-5401

DeVilbiss/Drive

SolAiris III, 303 Series



Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 30302605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# BF5000
Extended Life Bacteria Filter
Each

SolAiris V, 505, 515 Series



Item# OF5005SP
SolAiris V Kit
2 - OF5005
1 - OF6011
1 - BF1438
Each



Item# OF5005
Foam Cabinet Filter
5 1/2" x 35/8" x 5/8"
10/pack
OEM# 50502694



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# BF5000
Extended Life Bacteria Filter
Each



Item# OF6011
Extended Life Bacteria Filter—22mm fittings
Each
Use if SolAiris was manufactured before Nov '03

1025 Series



Item# OF1025D
Extended Life Bacteria Filter
Each



Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 303DZ-605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651

525 DS Series




Item# OF5004
Foam Cabinet Filter
33/4" x 33/4" x 1/2"
10/pack
OEM# 303DZ-605



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651



Item# RES515LF607
Low Output Flow Meter
Each
OEM# 515LF-607



Item# OF6011
Extended Life Bacteria Filter — 22mm fittings
Each

All MC Units



Item# OF5001
Foam Cabinet Filter
10/pack
6 1/8" x 3 1/2" x 1/4"
OEM# 444-514



Item# BF5000
Extended Life Bacteria Filter • **Each**



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# MC29D-651

Invacare

Mobilaire 3 and 5



Item# **OF9001SP**
Mobilaire Kit
2 - OF9001
1 - BFH6S
1 - BF516
Each



Item# **OF9001**
Foam Cabinet Filter
6³/₄" x 3³/₈" x 1¹/₂"
10/pack
OEM# 2000489



Item# **BFH6S**
HEPA Filter
6in Capsule
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack
OEM# P0750X030

Platinum 5 and 10



Item# **OF9002SP**
Platinum Kit
2 - OF9002
1 - BFH6S
1 - BF516
Each



Item# **OF9002**
Foam Cabinet Filter 9" x 2⁷/₈" x 1¹/₂"
10/pack
OEM# 1107412



Item# **BFH6S**
HEPA Filter 6in Capsule with Sticker
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter 1/4" Barbed
10/pack
OEM# P0750X030



Item# **HUM008S**
Straight Humidifier Connector
w/ 9in Tube
10/pack

Perfecto2



Item# **OF9004SP**
Perfecto2 Kit
2 - OF9004
1 - BFH6S
1 - BF516
Each



Item# **OF9004**
Foam Cabinet Filter
6⁷/₈" x 3³/₄" x 1¹/₂"
10/pack
OEM# 1143492



Item# **BFH6S**
HEPA Filter 6in Capsule
Each
OEM# 2000806



Item# **BF516**
Final Bacteria Filter
1/4in Barbed
10/pack
OEM# P0750X030

HomeFill



Item# **OF9001**
Foam Cabinet Filter
6³/₄" x 3³/₈" x 1¹/₂"
10/pack
OEM# 2000489

Item# **RP008**
2000psi Pressure
Switch
Each
OEM# 1089467

Rhythm P2

Rhythm P2, P2-E6, and P2-E7



Item# **OF7002-5**
Cotton Intake Filter
5/pack
OEM# P2FC-1

Rhythm HC, Jumao, Attendant, Lifestyle, Medacure

Rhythm LM5CA, Jumao JMC5A, Attendant 8X899, Lifestyle LM5A, Medacure OC5L

• Rhythm LM5CA



Item# OF7001-10
Foam Cabinet Filter
10/pack
OEM# LM209069

• Jumao JMC5A, Attendant 8X899, Lifestyle LM5A, Medacure OC5L, Drive DEC02



Item# OF7000
Foam Cabinet Filter
10/pack



Item# BF7000
Intake Filter
Each
OEM# LM209066



Item# BF516
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# LM209068

React Health/3B

Stratus 5



Item# OF6030
Foam Cabinet Filter
10/pack
OEM# STR1030



Item# BF6020
HEPA Filter
Each
OEM# STR1020



Item# BF516
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# P0750X030

Rhythm HC

Rhythm HC LM5BA



Item# OF8000
Foam Cabinet Filter
10/pack
OEM# LM5BA-ITF



Item# OF6011
Extended Life
Bacteria Filter —
22mm fittings
Each
OEM# LM5BA-ILF



Item# BF516
Final Bacteria Filter,
1/4in Barbed
10/pack
OEM# LM5BA-BF

Nidek

Nuvo



Item# OF11003
Foam Cabinet Filter
9 1/2" x 4" x 1/4"
10/pack
OEM# 9250-1025



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in Barbed
10/pack
OEM# 7631-1053

Item# OF11003SP
Nuvo Kit
2 - OF11003
1 - BF1438
1 - BF5000
Each



Item# BF5000
Extended Life
Bacteria Filter
Each
OEM# 9250-1028

Mark IV, V, V+



Item# BF1438
Final Bacteria Filter
1/4in to 3/8in
Barbed
10/pack
OEM# 7631-1053

Respironics

EverGo



Item#
OF1005SP
EverGo Kit
2 - OF1005
1 - BF316
Each



Item# **OF1005**
Foam Cabinet Filter
2³/₁₆" x 11¹/₁₆" x 3³/₈"
10/pack
OEM# 900-150



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

SimplyGo and SimplyGo Mini



Item# **BF8000**
Inlet Filter for SimplyGo
Each
OEM# 1095248



Item# **RES1008**
Compressor Filter for
SimplyGo Mini
Each

Millennium



Item#
OF1001SP
Millennium Kit
2 - OF1001
1 - BF2222
1 - OF3003
1 - BF316
Each



Item# **OF1001**
Foam Cabinet Filter
8¹/₄" x 4" x 3³/₈"
(3/16" notches)
10/pack
OEM# H619



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
OEM# H620



Item# **OF3003**
Foam Intake Filter
for use with BF2222
2⁵/₈" OD x 3³/₄" ID x 2³/₈"
10/pack
OEM# H618



Item# **BF316**
Final Bacteria Filter
3/16 Stepped Barb,
1/4in Straight Barb
10/pack
OEM# H621



Item# **OF1002**
Felt Filter
6/pack



Item# **BF1002**
Millennium M10 Style
Long Life Filter
Each

EverFlo



Item#
OF1003SP
EverFlo Kit
1 - BF1001
1 - BF316
Each



Item# **OF1006SP**
Cabinet Filter Kit
(No Cover)
2 Filters, 1 Casing
OEM# 1055625



Item# **OF1006**
Cabinet Filter
5¹¹/₁₆" x 4¹/₂" x 7⁷/₁₆"
10/pk



Item# **OF1006-CASE**
Foam Cabinet Filter
Casing
10/pack



Item# **BF1001**
EverFlo Intake Filter
Each
OEM# 1038831

Item#
HUM008
Humidifier
L-Connector
With 12in
Tube
10/pk



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

Healthdyne

Alliance



Item# **OF3002**
Foam Cabinet Filter
8¹/₄" x 3³/₄" x 1¹/₂"
10/pack
OEM# 100-0078-10



Item# **BF2222**
Intake Bacteria Filter
22mm fittings
5/pack
OEM# 100-6631-00



Item# **OF3003**
Foam Intake Filter
for use with BF2222
2⁵/₈" OD x 3³/₄" ID x 2³/₈"
10/pack
OEM# H618



Item# **BF316**
Final Bacteria Filter
3/16in Stepped Barb
and 1/4in Straight Barb
10/pack
OEM# H621

Repair Parts

NEW

Pilot Valve + 1 O-Ring for Invacare

Item# RP001

For Invacare Platinum (IRC10LX02, IRC5LX02), Perfecto2 (IRC5P02) and Perfecto2V (IRC5P02V)

Each

OEM# 2003052



O-Ring for Invacare Pilot Valve

Item# RP002

25/pack

OEM# 1085112




PE Valve for Invacare

Item# RP003

For Invacare Platinum (IRC10LX02, IRC5LX02), Perfecto2 (IRC5P02) and Perfecto2V (IRC5P02V)

Each

OEM# 1143578



On/Off Switch

Item# RP004

For Invacare Platinum and PerfectO2, Drive/DeVilbiss, Lifestyle, Direct Supply

Each

OEM# 1139950



Power Cord for Invacare

Item# RP009

For Invacare Perfecto2, Platinum 5 & 10

Each

OEM# 1108764



Flow Meter for EverFlo 0-5L

Item# RP012

Each

OEM# 0209184



Hour Meter

Item# RP005

For Invacare Platinum and PerfectO2, Lifestyle, Direct Supply (Meter Only - No Clip)

Each

OEM# 2002939



O2 Conc Hose Barb Fitting 1/8" Conc Port

Item# RP006

For Invacare and EverFlo Post 2013

Each

OEM# 1104151



O2 Conc Fitting Receptacle Port

Item# RP007

For Invacare HomeFill

Each

OEM# 1104152



2000psi Pressure Switch

Item# RP008

For Invacare HomeFill

Each

OEM# 1089467




4-Way Valve Rebuild Gasket Kit

Item# RP010

For Invacare Platinum 10, Platinum 5, PerfectO2V, and PerfectO2

Each



Concentrator Capacitors

Each

Item#	Size	Compatibility	OEM#
RP013	10MFD	PerfectO2	113125002
RP014	15MFD	Platinum Series	2003019
RP015	22MFD	PerfectO2V	1191077
RP016	30MFD	Platinum 10	1128740



Repair Parts

NEW

Compressor Rebuild Kit

Item# **RP018**

For Invacare Perfecto2V and
Drive/DeVilbiss 525/1025, 303DZ
Each

OEM# 1191916



Pressure Regulator

Item# **RP019**

For Invacare Platinum 10 / Perfecto2V
Each

OEM# MC29D-612



Flow Meter

Item# **RP020**

For Drive/DeVilbiss 525/1025
Each

OEM# 505DZ-607



Cooling Fan

Item# **RP021**

For Drive/DeVilbiss 525/1025
Each

OEM# 525DD-634



Exhaust Muffler

Item# **RP022**

For Drive/DeVilbiss 525/1025
Each

OEM# 515A-705



Circuit Breaker

Item# **RP023**

For Drive/DeVilbiss 525/1025
Each

OEM# 525D-613



Oxygen Outlet Port

Item# **RP024**

For Drive/DeVilbiss 525/1025
5/case

OEM# 525DD-606



O2 Quick Connect L-Shape Adapter

Item# **RP025**

For Respironics Trilogy
2/pack



Cylinder Dust Cap

Item# **RP026**

For Invacare HomeFill
10/pack



Rubber Dust Cap

Item# **RP027**

For Invacare HomeFill
10/pack



O-Ring for Homefill Tank

Item# **RP028**

For Invacare HomeFill
10/pack



Other

Other

Page

Disinfectants	65
Aerosol	66
Tabletop Nebulizers	67
Handheld Nebulizers	68
Nebulizer Kits	69
Nebulizer and Anesthesia Filters	70
Suction	71-73
Closed Suction Systems	74-75
Tracheostomy	76
HMEs	77
Ventilator Filters, Masks, and Connectors	78
Ventilator Circuits	79-81
Anesthesia Masks	80
StrongArm Mobility and Comfort	82-83
Bath Safety	83
Maternity	84-85
Wound Care	86-89
Equipment Covers	90-91

Other HCPCS

A6402	Gauze - sterile, pad size 16 sq. in. or less, without adhesive border, each dressing
A4483	Heat and Moisture Exchanger — Disposable
A4626	Tracheostomy Cleaning Brush
A4628	Oropharyngeal Suction Catheter (Yankauer)
A4629	Tracheostomy Care Kit - for established tracheostomy
A7000	Suction Canister — Disposable
A7002	Tubing for Suction Pump
A7003	Nebulizer Set - Small Volume—Disposable
A7005	Nebulizer Set - Small Volume—Nondisposable

A7013	Nebulizer Compressor Filter— Disposable
A7015	Aerosol Mask used with Nebulizer
A7525	Tracheostomy Mask
E0570	Nebulizer with Compressor
E0600	Respiratory Suction Pump - home model, portable or stationary, electronic
E2601	General Use Wheelchair Seat Cushion - width less than 22 inches, any depth
K1005	Disposable collection and storage bag for breast milk, any size, any type, each

Disinfectants

Simpli Germicidal Disposable Wipes

Item#

RES6108W

Canister of 160 wipes

NEW



- Bactericidal, Tuberculocidal, Fungicidal, Virucidal*
- Premoistened wipes containing a quaternary/alcohol-based solution.
- Made in the USA
- Kills *SARS-CoV-2, Covid 19, Mycobacterium bovis – BCG (Tuberculosis—TB) in 1 Minute
- Kills Candida Auris in 2 Minutes

CleanSmart Disinfectant

Item#

CAP1019-C

23oz Spray Bottle

6/case

- CleanSmart is EPA listed as a broad-spectrum and hospital disinfectant
- EPA Registration Number 89896-2
- CleanSmart is formulated for disinfecting surfaces of external medical equipment and external surfaces of respiratory equipment
- CleanSmart is included on the EPA's List N: Disinfectants for Use Against SARS-CoV-2, the virus that causes COVID-19
- Active ingredient Hypochlorous Acid does not leave a chemical residue



Aerosol

Aerosol Mask

Item# **RES2110**

Adult Size

50/case

Item# **RES2210**

Pediatric Size

50/case

HCPCS# A7015

- Clear, soft vinyl for patient comfort
- Adjustable nose clip for secure seal
- Connector accepts 22mm aerosol tubing
- Individually packaged



Peak Flow Meter

Item# **RES052**

10/case

50-800 L/min

- High performance
- Quality controlled
- Color areas indicator
- Slip design
- One-way valve
- Latex-free
- Phthalate-free



Item# **RES052MPD**

Pediatric Mouthpiece
for Peak Flow Meter

10/case



Drainage Bag with T-Adapter

Item# **RES070**

750cc • 50/case



Drainage Bag with Y-Adapter

Item# **RES071**

750cc • 50/case



- Designed to ensure patients receive mist rather than liquid accumulated in the tubing
- Includes hanger and bottom drainage port
- Works with 22mm corrugated aerosol tubing

SperaChamber

Anti-Static Valved Inhaler Chamber
for Metered Dose Inhaler

Chamber Only

Item#

RES195-10 10/pack

RES195-50 50/pack

Chamber with Mask

10/pack

Item#

RES195MS Small (Infant)

RES195MM Medium (Pediatric)

RES195ML Large (Adult)

- Inhalation valve
- Universal adapter
- Spiral vents
- Seal with exhalation safety valve
- Latex-free and DEHP-free
- Dog mask appeals to kids



Mask Only

Size	Item#
Small	RES195S
Medium	RES195M
Large	RES195L



Tabletop Compressor Nebulizers

Compressor Nebulizers

Item#

NEB100

- Disposable Nebulizer Kit
- Foam Filters

NEB100-R

- Disposable Nebulizer Kit
- Reusable Nebulizer Kit
- Foam Filters x 5



NEB100-BAG
Carry Bag

- Rust-resistant
- Quiet (below 60dBA)
- 6ml medication capacity
- Average nebulization rate above 0.2ml/min
- Particle size: respirable fraction 0.5 to 5um 0.5um 80%
- Compressor pressure range 30-45 PSI
- Operation pressure range 12-14 PSI
- Liter flow range 8-10 LPM

Sunny the Bear



NEB300BEAR

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

NEB300BEAR-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



NEW

NEB300-BAG Bag for Nebulizer Compressor

To include bag in configuration, just add "B" to Item Number:

NEB300FISH-B
NEB300FISH-RB
NEB300BEAR-B
NEB300BEAR-RB

Gordon the Fish



NEB300FISH

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Foam Filters x 5

NEB300FISH-R

- Disposable Pediatric Mask
- Disposable T-Piece Kit
- Reusable T-Piece Kit
- Foam Filters x 5



Covers for Nebulizers



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Handheld Nebulizers

Handheld Mesh Nebulizer

Item#
NEB402

With internal rechargeable battery

Included

Medication cup
USB cord
Mask
Mouthpiece
Portable bag
Manual



Particle Size MMAD < 5µm
Nebulization Rate 0.25ml/min-0.9ml/min
Medication Cup Capacity 10ml (Max)
Device Weight Approx. 1/4 lb
Device Size Approx. 1 1/2 x 2 x 3 7/8 inches

Handheld Compressor Nebulizer

Item#
NEB200

Included

Compressor
Medication cup
AC power adapter
Air hose
Elbow
Mouthpiece
Mask
4 AA batteries

- Runs on 4 AA batteries or AC power adapter
- Portable (small and weighs only 1/2 pound)
- Quiet (60dBA)
- Fine particle size (Approximately 3 microns)
- No replacement filter required



Nebulizer Kits

Disposable Nebulizer Kit with T-Piece

Item# **RES091**
Disposable
Single Patient
50/case
HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - T-piece with mouthpiece and corrugated reservoir tube
 - 7ft supply tube



- Green, convex cone is easy to see and remove for simple cleaning and reassembly
- Easy to read medication jar provides a leak-proof seal

Disposable Nebulizer Kit with Angled Mouthpiece

Item# **RES096**
Disposable
Single Patient
50/case
HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - Angled mouthpiece
 - 7ft supply tube



Reusable Nebulizer Kit with T-Piece

Item# **RES093**
Reusable
Single Patient
10/case
HCPCS# A7005

- Includes:**
- 10cc jet nebulizer bottle
 - T-piece with mouthpiece and corrugated reservoir tube
 - 7ft supply tube



- Latex-free
- Dishwasher safe
- High capacity (10cc)
- CC volume markings
- Very low residual volume (0.7cc)
- 1/4 turn cap provides leakproof seal
- Works upright or at angles up to 45°

Disposable Nebulizer Kit with Mask

Item# **RES090**
Adult Size
50/case

Item# **RES092**
Pediatric Size
50/case



HCPCS# A7003

- Includes:**
- 6cc jet nebulizer bottle
 - Mask
 - 7ft supply tube

Disposable single patient use for long term delivery procedures

6cc bottle is easy to assemble and disassemble



Nebulizer Filters

Plastic Intake Filter with Head

Item# **BF016**

10/pack

HCPCS# A7013



Plastic Intake Filter with Head

Item# **BF017**

5/pack

HCPCS# A7013



Foam Filter

Item# **BF018**

10/pack

HCPCS# A7013



Intake Filter for 50 PSI Compressor

Item#

BF020

3/pack

NEW

Compatible with CHAD Drive
50 PSI Compressor
OEM# 18450-FILTER



Plastic Intake Filter for Respironics InnoSpire

Item# **BF026S**

10/pack

HCPCS# A7013



- Mini
- MicroElite
- MiniElite

Nebulizer Filter Compatibility

Caire

Medi-Mist	BF018
-----------	-------

CHAD

Drive 50psi	BF020
-------------	-------

Drive/DeVilbiss

PulmoNeb	BF017
PulmoNeb LT	BF017
Sunrise Compact	BF018
Pulmo-Aide	BF017
Pulmo-Aide Compact	BF017
PulmoMate	BF017
Power Neb Ultra 18080	BF017
CHAD 50psi	BF020

Evo

AeroMax	BF018
SportsNeb	BF018
Super SportsNeb	BF018

Invacare

Passport	BF018
Envoy	BF018
Envoy Jr.	BF018
Select	BF018

John Bunn

Nebulite EV	BF018
LT2	BF018
LX	BF018

Medline

SportsMist	BF018
------------	-------

Respironics

Inspiration 626 & 929	BF016
MisterNeb	BF016
Inspiration Elite	BF016
MiniElite	BF018
InnoSpire	BF026S

Salter

AIRE Plus	BF018
-----------	-------

Schuco

2000	BF018
3000	BF018

Sunset

NEB100 Series	BF018
NEB300 Series	BF018

Suction

Suction Machine

Item#

SU100DC

With Internal
Rechargeable
Battery

SU100AC

Without Battery

SU100DPA

Includes DC Mobile
Power Adapter
and Internal
Rechargeable Battery



Comes with

Machine
Carry Bag
800cc Suction Canister
AC Adapter & Power Cord
6ft Suction Tubing
3.5in Suction Tube Connector
4.5in Suction Tube Connector
BF1438S Bacteria Filter

- Effective removal of respiratory fluids
- Durable yet lightweight
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty

Accessories and replacement parts

SU100DC-PA	DC Mobile Power Adapter (for car power outlet)
SU100DC-Adapt	Replacement Power Adapter (without cord)
SU100DC-Cord	Replacement Power Cord (without adapter)
SU100DC-BAG	Carrying/Storage Bag
SU100DC-BATTERY	Internal rechargeable battery replacement

Suction Kit for SU100DC

Item# **RES026S-SS**

RES023AW Suction Canister
RES024 Suction Tubing Connector 3.5in
RES024M Suction Tubing Connector 4.5in
BF1438S Bacteria Filter
RES025 Suction Tubing 1/4in x 6ft

HCPCS#
CANISTER - A7000
TUBING - A7002



Suction Canister

Item#

RES023AW

800cc • Single Use
48/case

HCPCS# A7000

- Float valve shutoff



Suction Filter Component Kit for SU100DC

Item# **RES024-MKIT**

RES024 Suction Tubing Connector 3.5in
RES024M Suction Tubing Connector 4.5in
BF1438S Bacteria Filter



RES026S-SS and RES024-MKIT are compatible with:

Drive/DeVilbiss

Vacu-Aide® QSU Quiet Suction Unit	7314P-D
	7314P-D-EXF
	7314D-D
Vacu-Aide® Suction Unit	7305D-D
	7305P-D-EXF
	7305D-D-EXF
	7305P-D

Sunset

Suction Machine (Portable Unit)	SU100DC
Suction Machine (Stationary Unit)	SU100AC

Covers for Suction Machines



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Suction

Suction Kit

Item# **RES026S-MULTI**

RES023AW Suction Canister 800cc
 RES024XL Suction Tube Connector 16in
 RES025 Suction Tubing 6ft x 1/4in
 BF516TS Bacteria Filter
 RES024S Suction Tube Connector 1.5in

NEW



Suction Component Kit

Item# **RES024-MULTI**

RES024S Connector Tube 1.5in
 RES024XL Connector Tube 16in
 BF516TS Bacteria Filter

NEW



RES026S-MULTI and RES024-MULTI are compatible with:

Precision Medical

Easy Vac	PM60
Power Vac	PM61
Easy Go	PM66 series PM65 PM65E PM65HG

Compass/Roscoe

Roscoe Portable Suction	50004 (AC model) 50006 (DC model)
Cliq Single Cylinder Aspirator	ASP-ROS

Medline

Heavy Duty	HCS7000
------------	---------

Invacare

Aspirator durable suction unit	IRC1135 IRC1136
--------------------------------	-----------------

Drive/DeVilbiss

Heavy-Duty Suction Machine	18600
VacuMax™ Go Suction Machine	18615


Allied Healthcare

Schuco Portable Aspirator with Sturdy Legs	S430A
Schuco Portable Suction Unit with Metal Base	S130A
Schuco Portable Aspirator with Molded Plastic Base	S330A
Gomco OptiVac® AC/DC Portable Suction Unit & 305 Aspirator	G180 G180CE G180CE-JAR G180-JAR 01-22-3005 01-22-0305
Gomco® Tabletop Aspirator Model 300	01-22-0300
Gomco® General Use Aspirators Models 3020 & 3040	01-62-3020 01-62-3040
Gomco® 405 Aspirator	01-12-0405 01-12-4005 01-23-0405
Gomco® Stand-Mounted Models 4032 & 4042	01-52-4032 01-52-4042
Gomco Table Top Gastric Aspirator 270	01-22-0270

Suction

Suction Machine Inline Filters

Item# **BF1412**
Final Bacteria Filter
1/4in to 1/2in
Stepped Barb on
Both Sides
10/pack



Item# **BF1438S**
1/4in to 3/8in
10/pack



Item# **BF3812**
3/8in to 1/2in
10/pack



Item# **BF1312**
1/3in to 1/2in
1/pack



Item# **BF516TS**
5/16in Barbed with
Threaded NPT
10/pack



Suction Tubing

1/4in

Length	Item#	Qty
18in	RES025S	50/case
6ft	RES025	10/case



3/16in

Length	Item#	Qty
18in	RES025-31618	100/case
20in	RES025-31620	100/case
6ft	RES025-31672	50/case

Suction Tube Connectors

10/pack

Size	Item#
1.5in	RES024S
3.5in	RES024
4.5in	RES024M
6in	RES024L
10in	RES024L10
16in	RES024XL



Rubber Suction Elbow Connector Drive DeVilbiss Style

Item#
RES028
6/pack



Suction Elbow Connector Precision Medical Style

Item#
RES033
6/pack



Bulb Tip Yankauers

Item# **TRAA001S**
Vented • **50/case**


Item# **TRAA002S**
Non-Vented • **50/case**

HPCPS# A4628



- Sterile
- Single patient use
- Not made with natural latex rubber

Standard Closed Ventilation Suction System



Standard Closed Ventilation Suction System

24 Hour Adult

Endotracheal 540mm

Item#	Description
VEN46001000A	10Fr
VEN46001000AM	10Fr, MDI Port
VEN46001200	12Fr
VEN46001200M	12Fr, MDI Port
VEN46001400	14Fr
VEN46001400M	14Fr, MDI Port
VEN46001600	16Fr
VEN46001600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN46011000A	10Fr
VEN46011000AM	10Fr, MDI Port
VEN46011200	12Fr
VEN46011200M	12Fr, MDI Port
VEN46011400	14Fr
VEN46011400M	14Fr, MDI Port
VEN46011600	16Fr
VEN46011600M	16Fr, MDI Port

72 Hour Adult

Endotracheal 540mm

Item#	Description
VEN56051000A	10Fr
VEN56051000AM	10Fr, MDI Port
VEN56051200	12Fr
VEN56051200M	12Fr, MDI Port
VEN56051400	14Fr
VEN56051400M	14Fr, MDI Port
VEN56051600	16Fr
VEN56051600M	16Fr, MDI Port

Tracheostomy 305mm

Item#	Description
VEN56061000A	10Fr
VEN56061000AM	10Fr, MDI Port
VEN56061200	12Fr
VEN56061200M	12Fr, MDI Port
VEN56061400	14Fr
VEN56061400M	14Fr, MDI Port
VEN56061600	16Fr
VEN56061600M	16Fr, MDI Port

Neonatal/Pediatric

NEW

With DSE Connection

Item#	Description
VENCSS6-305	6Fr, 305mm
VENCSS8-203	8Fr, 203mm
VENCSS8-355	8Fr, 355mm
VENCSS10-203	10Fr, 203mm
VENCSS10-405	10Fr, 405mm



Tracheostomy

Tracheostomy Mask

Item#
RES2130
Adult Size
50/case
HCPCS# A7525


- Latex-free
- 360° swivel in 2 places allows access from either side of the patient
- Comes with removable adapter to fit with aerosol or O₂ supply tubing



RES2230
Pediatric Size
50/case

Adapter for T-HME

Item# **RES027A**
Each



Flexible Trach Adapter

Item#
TRAA003
22mm Cuff
Each



OEM# 1073902, 1147898

NEW

Flexible Connector

50/case

NEW

Item#
TRAA004 Adult
22/15mm x 15mm
Length 4.9-7.5cm
OEM# 3222



TRAA004P Pediatric
15mm F x 15mm M
Length 4.9-7.5cm
OEM# 3215



Bacterial/Viral HME Filter

Item#
BF104
Each

NEW

Bacterial Filtration Efficiency (BFE)
99.99+%

Viral Filtration Efficiency (VFE)
99.99+%

Flow Resistance 0.09 (30L/min kpa)

Dead Space 38ml

Weight 30.6g

Tidal Volume 150-1500ml

Moisture Output 35.1 mg/L



Tracheostomy Ties with Fastener Tabs

20/case

One-piece • Easy on

Item#
RES029P-11
1 in x 11.5 in
RES029-17
1 in x 17 in



NEW

Two-piece • Easy off

Item#
RES030P-10
0.78 in x 10 in
RES030P-14
0.78 in x 14 in
RES030P-18
0.78 in x 18 in
RES030-20
1 in x 20 in



Sunset's trach ties deliver value and comfort, with soft straps and light velcro closures for customization. Secure, but flexible for patient movement, the ties help prevent sliding, skin irritation and accidental extubation. 60% Polyester, 20% Cotton, 20% Nylon

HMEs

Heat Moisture Exchangers

30/case*

*RES027FLEX 40/case

*RES027UOP 50/case



NEW

	RES027FLEX Flex HME — Foam	RES027 T-HME	RES027IL Inline HME with Gas Sampling Port	RES027NE Inline HME — Neonatal	RES027OP T-HME with O2 Port	RES027UOP HME with O2 Port OEM# 1190
TV Range (mL)	250-1500	250-1000	70-600	15-50	150-1250	50-1000
Moisture Output (mg H₂O/L)	30 mg/L	30, Vt = 0.1L	34.5 at VT = 250 ml 32.5 at VT = 500 ml	30 mg/L	36.3 mg/L	13.2 mg/L
Weight (g)	15.5	4.5	10	4.5	7.4	5
Resistance (cm H₂O)	0.5 @ 30 LPM 1.4 @ 60 LPM	0.32 @ 30 LPM 0.81 @ 60 LPM 1.62 @ 90 LPM	0.8 @ 30 LPM 2.1 @ 60 LPM	0.5 @ 10 LPM	0.2 @ 30 LPM	0.4 @ 30 LPM 0.9 @ 60 LPM 1.5 @ 90 LPM
Dead Space (mL)	90	12.5	21	2.5	7	17
Connection (mm)	22F ,15F	15F	22M/15F-15M	15M,15F	15F 5.5mm O2 port	15F

Ventilator

Movair / Löwenstein LUISA®

Item# **VNT31390**
Each



LUISA with High Flow Therapy provides long term care with easy serviceability. Routine patient maintenance consisting of only filter changes and a device wipe down provides uncomplicated upkeep.

With simple, routine maintenance and no adapters required between circuits, the LUISA provides clinicians and patients' peace of mind with an all-encompassing system in a compact, home ventilator with the added benefit of High Flow Therapy.

- High Flow Therapy: all circuit systems
- Connectivity: customized patient download reports
- Uncomplicated: no adapters, intuitive operation
- Battery run time: up to 18 hours
- Monitoring: FiO₂ and SpO₂ measurement
- Safety: alarm displayed in two languages
- Display: 10" angled touchscreen
- Functionality: 100 - 3,000 mL Vt without additional flow sensor
- Responsive: adjustable target volume speed
- Lightweight: ~ 8 lbs.
- Positioning: vertical or horizontal; Circuit orientation options (left or right side)

Includes:

- Exhalation Module
- Leakage Circuit
- Power Supply/Power Cord
- O₂ Nozzles
- Fine and Course Filters
- Protection Bag
- Accessories Bag
- USB Drive

Accessories and Replacement Parts

VNT31404	Exhalation Module - Disposable
VNT31569	Power Supply
VNT27027	Power Cord
VNT30669	Oxygen Connection Nozzle
VF29652	Set of 12 Pollen/Fine Filters - 12/PK
VF29928	Set of 2 Air/Coarse Dust Filters - 2/PK
VNT31417	Protective Bag
VNT31440	Accessories Bag
VNT31414	USB Flash Drive

Ventilator Filters, Masks, And Connectors

NIV Mask with Headgear — Non-vented

Each

Item# **CM038M** Medium • **CM038L** Large



Circuit Connectors

Item# **RES013**

15mm ID x
22mm OD
10/pack



Item# **RES014**

22mm OD x
22mm OD
10/pack



Reusable Filter for Respiration Trilogy

Item# **CF1009F**

2/pack

OEM# 1029330



Filter for ResMed Astral Ventilator

Item# **VF003-4**

4/pack

OEM# 27939

NEW



Ventilator Felt Filter

Item# **VF006**

For use with NPB LP6 & LP10

10/pack

OEM# L-002917-000



Filters for Löwenstein LUISA

NEW

Item# **VF008-12**

Disposable filter
12/pack

OEM# WM29652



Item# **VF008F-2**

Reusable Filter
2/pack

OEM# WM29928



Breas 45LS Style Air Inlet Filter

NEW

Item# **VF010-5**

Disposable filter
5/pack

OEM# 007103



Item# **VF010F-5**

Reusable Filter
5/pack

OEM# 007104



Breas Vivo 45 Style Cooling Air Filter

NEW

Item# **VF011F-5**

Reusable filter
5/pack

OEM# 007105



Inline Bacterial Viral Filter for CPAP/BIPAP/Vent Units

Item#

BF2200

22mm

5/pack

Attach between the machine and tube for additional protection against bacteria and allergens



Valved T-Adapter — Adult

Item# **VEN001**

22F-18F-22M/15F

Each



Swivel Exhalation Port (For Ventilator Circuits & CPAP)

Item# **VEN002**

Each



Fixed Exhalation Port

Item# **VEN003**

For Ventilator Circuits and CPAP

22 mm OD x 15 mm ID

10/pack

OEM# 312149



Air Inlet Vent Filter for Philips Respironics V60

Item# **VF002**

For use with Philips Respironics V60

5/pack



Reusable Filter for Respiration Trilogy EV300

For use with EV300

Item#

Qty

VF007-1

1/pack

VF007-10

10/pack

OEM# 1134591



NEW

Proximal Inline Filter for Newport HT50

Item# **VF009-5**

5/pack

NEW



Ventilator Circuits without Humidifier Chamber

Passive Heated Breathing Circuits

Item#

VC1114

For use with MR850 Humidifier

10/case

OEM# RT114

**VC1003****With Swivel Exhalation Port**

For use with Philips Respironics Trilogy

Each**Includes:**

- 1 VC1114 Passive Heated Breathing Circuit
- 1 VEN002 Swivel Exhalation Port

OEM# 1076046,
1073227, 1073235

Adult NIV Patient Circuit

Item#

VC1160Compatible with V60,
No Filters**10/case**

OEM# 1069210



Passive Non-Heated Breathing Circuits

For use with Respironics Trilogy
Each**Exhalation Port****Item#****OEM#**

Swivel

VC1000

1073228 & 1073236

Fixed

VC1100

1069467 & 1069468



Passive Circuit w/Exhalation Port and Bacterial/Viral Filter

With 22mm cuffs

Each**Size****Tube Diameter****Item#**

Adult

19mm

VC1005

Pediatric

15mm

VC1005P

OEM# 1147378, 1132340, 7615



Adult Passive Circuit with Bacteria Filter

Item#

VC1006

Comparable to Breas 007936

Each

Adult Passive Circuit with Bacteria Filter and Fixed Exhalation Port

Item#

VC1007

Comparable to Breas 007615

Each

Ventilator Circuits without Humidifier Chamber

Cough Assist Kits

All kits come with 1 BF2200 Bacterial/Viral Filter and 1 RES013 Circuit Connector

Item#

CAKIT001

With 1 TUB06 6ft Tube

Each

OEM# 1090835



Item#

CAKIT002

With 1 TUB06 6ft Tube and
1 RES019 Oxygen Enrichment Adapter

Each



CAKIT003

With 1 TUB06 6ft Tube
and 1 TRAA003 Flexible
Trach Adapter

Each



NEW

ABM Respiratory Care BiWaze Cough System



Item# **OF8003-3**

Air Inlet Foam Filter
3/pack

Kits With Tube and SBC01 Series Breathing Mask

With 1 TUB06 6ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	CAKIT004I	1090830
Toddler	CAKIT004T	1090831
Adult Small	CAKIT004S	1090832
Adult Medium	CAKIT004M	1090833
Adult Large	CAKIT004L	1090834



With 1 TUB04 4ft Tube

Mask Size	Item#	Comparable to OEM#
Infant	CAKIT005I	325-9237
Toddler	CAKIT005T	325-9236
Adult Small	CAKIT005S	325-9235
Adult Medium	CAKIT005M	325-9234
Adult Large	CAKIT005L	325-9217



CAKIT006

with 1 CAMP001 Mouthpiece
and 1 TUB06 6ft Tube

Each

OEM# 1090836



NEW

Anesthesia Masks

Inflatable Air Cushion Breathing Mask with Valve - No Hooks

20/case

Mask Size	Item#
Infant	SBC01I
Toddler	SBC01T
Adult Small	SBC01S
Adult Medium	SBC01M
Adult Large	SBC01L



- Ergonomic cushion design provides an optimal seal with minimal pressure
- Clear mask permits easy patient observation
- Available in a variety of sizes for the best possible fit
- Inflation port for optimized seal


Ventilator Circuits without Humidifier Chamber

Non-Heated Breathing Circuits


For use with Respiration Trilogy
22mm cuffs
30/case

Item#	OEM#	Mode	Water Trap	Size	Tube Diameter
VC1001	1073223 & 1073232	Active	No	Adult	22mm
VC1001P	1073220 & 1073229	Active	Yes	Pediatric	15mm
VC1004	1069443 & 1069444	Passive	Yes	Adult	22mm
VC1004P	1073224 & 1073233	Passive	Yes	Pediatric	15mm


Item#
VC1001
Active Adult
OEM#
1073223
1073232




Item#
VC1001P
Active Pediatric
OEM#
1073220
1073229




To use active circuits with Respiration Trilogy, simply cut ends off pressure and exhalation lines.



Item#
VC1004
Passive Adult
OEM#
1069443
1069444



Item#
VC1004P
Passive Pediatric
OEM#
1073224
1073233



StrongArm®

StrongArm® Comfort Cane

12/case

Item#	Color
BM001	Black
BM001-BL	Blue
BM001-PR	Purple

**Color pack
(one of each)**

BM001-3MC



- Ergonomic handgrip supports forearm and relieves pressure from wrist
- Offset bend centers grip over tip for maximum stability when cane is at the user's side
- Lightweight and strong
- 12 height adjustments
- Comfort Fit Handle
- Supports up to 500 lbs
- For left and right arm
- Self-standing base: Larger 'footprint' greatly improves traction and overall stability
- Safely supports user walking on grass, gravel, sand, carpet, hard surfaces, and stairs



**Now comes in
retail packaging**

StrongArm® Grabber Reacher Tool

Item#
BM002
12/case

Your ultimate companion for effortless grabbing and reaching! Engineered with precision and designed for practicality, this innovative tool is perfect for those hard-to-reach items and provides extra support for heavier items. Grab with Ease and Confidence in Each Squeeze.



- Transfers weight for improved support
- Rotating jaw for vertical or horizontal use
- Integrated hook and magnetic retriever
- Durable, lightweight aluminum construction
- Folds for easy storage and portability
- 32" of added reach

StrongArm® Clear View Universal Walker Tray

12/case

Item#	Color
BM003-BL	Blue
BM003-BK	Grey
BM003-PL	Purple



Unobstructed view for safer mobility

The Clear View Walker Tray is designed to enhance the user's ability to see obstacles ahead and reduce the risk of falls. Its transparency helps maintain the user's line of sight, providing added safety and confidence while moving. Compatible with the leading brand walkers.

Transparent

See where you are walking. Avoid tripping hazards.

Universal Fit

Tray works with the leading brand walkers.

Convenience

Built-in cup holder and storage bag.

StrongArm®

StrongArm® Universal Clip-On Light

Item#
BM004
4/box

NEW



Featuring a super-bright LED light, this clip-on accessory significantly enhances visibility, providing added safety during evening strolls or walking throughout the house. Effortlessly clip onto canes, walkers, rollators, crutches, grabber reachers, knee scooters, wheelchairs, and more, ensuring you always have reliable illumination on any mobility device.



StrongArm® LED Magnifying Glass

Item#
BM005
10/box

NEW



- Handheld 4x LED Magnifier
- 6 LED Lights
- Compact, Lightweight, and Portable
- Lens Diameter: 2.5"
- Weighs Only 1.7 oz
- Uses 2 AAA Batteries (Included)

Bath Safety

PreserveTech™ 360° Swivel Bath Chair

Item#
BSE001

Compact and customizable, the PreserveTech™ Swivel Bath Chair provides a safe transfer in and out of the tub for those with limited mobility. The seat rotates a full 360° in either direction and automatically locks at 90° intervals to help prevent skin irritation and friction wounds during transfers. The seat height can be adjusted from 15"-22", and the padded armrests can be raised out of the way.



Splash Defense™ Transfer Bench

Item#
BSE002

Splash Defense™ Transfer Bench is the first of its kind to allow the shower curtain liner to remain inside the bath tub while a user is bathing. This unique product enhancement is essential to help prevent water from spilling out onto the bathroom floor, creating a danger zone for users and caretakers. The Splash Defense™ Transfer Bench provides security for the caretaker, while giving users back their independence and allowing them to bathe safely and privately again.



PreserveTech™ 360° Universal Raised Toilet Seat

Item#
BSE003

Introducing the PreserveTech™ Universal Raised Toilet Seat. With this truly universal raised toilet seat users can set three different seat heights (2, 4 or 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space. This item also includes a lid with a recessed grip for effortless opening.



Maternity

Motif Twist Double Electric Breast Pump

Item#
MAT2300T1

Motif Twist is a convenient, portable option for every mom. This pump is lightweight and simple to use with a quiet motor. The Twist is battery-powered with AA batteries, making it a great option for moms on-the-go. Simply twist the knob to select different suction and massage settings and use the green LED light timer to track pumping sessions. Includes a storage bag and bottle stands.

- Expression Mode: Works to express milk from the breast as efficiently as possible
- Closed System: Backflow prevention keeps the motor clean, preventing contamination
- Lightweight: Weighs 0.7 lb. — perfect for portability
- Massage Mode: Helps stimulate let-down by mimicking a baby's nursing pattern
- Dual Power Options: Allows you to pump stationary or on battery power with AA batteries
- Quiet Motor: 45 dB hum making this a discreet option for moms



Motif Luna Double Electric Breast Pump

Item#
MATAAA0013

Sophisticated, strong, and soothing, the Motif Luna is Motif's most powerful breast pump. With a battery life of 2+ hours between charges, the Luna was built for modern motherhood. The Luna is available in both battery and non-battery versions. Lightweight and portable, the Luna weighs less than 2 pounds. With a suction strength of 280 mmHg, the Luna is proven to provide More Milk in Less Time®.

- Very quiet
- Customizable pumping sessions
- Compact size and light weight
- Less guilt, pain, & anxiety
- Convenient nightlight & digital screen

NEW



Maternity

Motif Milk Storage Bags

Designed to make breastmilk storage as simple as possible, these bags feature an easy-pour spout for a no-mess experience. Easily transfer milk from the bottle to the bag, seal with the leak-proof double zipper lock, label, and store in your fridge or freezer. No spills. No hassle. Made with your busy day in mind. Comes in a box of 40, 90, 100 or 300 count.

HCPCS: K1005

Item#

MATAAC000806

100/pack

Other pack counts are available.
Ask your rep for details.

NEW



- Easy-pour spout design
- BPA-free
- Leak-proof with double zipper lock
- Stores up to 8 oz. of milk
- Write-on label
- Freezer safe
- Lays flat for compact freezer storage
- Self stand for refrigerator storage

Motif BiliTouch™ Phototherapy Blanket

Item#

MAT2300EAB

Motif BiliTouch™ Phototherapy Blanket
with Hard Case
Each

Motif BiliTouch™ Disposable Covers

MAT2300EAB-DC

25/case

Motif Phototherapy Eye Shield

25/box

Size	Item#
Small	MAT2300ES-S
Large	MAT2300ES-L

The Motif BiliTouch™ phototherapy blanket is a biliblanket designed by healthcare professionals to treat newborn infants diagnosed with jaundice (hyperbilirubinemia). This portable, convenient phototherapy treatment device allows your patients to continue to hold, feed, and care for their baby while reducing high bilirubin levels at home.



- Convenient, ultra-portable phototherapy treatment with battery-powered option
- Patient-centered, lightweight design facilitates bonding with newborn infant
- LED light source, no fiber optic cables for easy positioning

Wound Care

ZeniFOAM

Polyurethane Foam
Wound Dressing
10/box

ZeniFOAM is an absorbent, non-adherent sterile polyurethane foam dressing used to absorb moderate to heavy wound drainage to protect sensitive and fragile skin. ZeniFOAM is conformable with a colored waterproof backing to prevent strike through drainage and provide protection from external contaminants.



Item#	Size
WC10022	2" x 2"
WC10045	4" x 5"
WC10066	6" x 6"
WC10088	8" x 8"

ZeniFOAM GENTLE

Polyurethane Foam Wound
Dressing with Silicone
Adhesive, No Border
10/box

ZeniFOAM GENTLE is a sterile, silicone-faced, polyurethane foam dressing without an adhesive border designed to absorb moderate to heavy wound drainage. ZeniFOAM GENTLE provides maximum absorption while protecting sensitive or fragile skin and reduction of painful dressing changes.



Item#	Size
WC20022	2" x 2"
WC20044	4" x 4"
WC20066	6" x 6"

ZeniFOAM GENTLE BORDER

Polyurethane Foam Wound
Dressing with Silicone
Adhesive and Border
10/box

ZeniFOAM GENTLE BORDER is a sterile, silicone-faced, polyurethane foam dressing with a silicone adhesive border designed to absorb moderate to heavy wound drainage. ZeniFOAM GENTLE BORDER is a multi-layered dressing that provides maximum absorption while protecting sensitive or fragile skin and reducing painful dressing changes.



Item#	Size
WC30000H	Heel
WC30022	2" x 2"
WC30033	3" x 3"
WC30044	4" x 4"
WC30066	6" x 6"
WC30077S	Sacral 7" x 7"
WC30099S	Sacral 9" x 9"

ZeniPAD PLUS

Multi-Layer Composite Dressing
10/box

ZeniPAD PLUS is a sterile, multi-layered, waterproof composite wound dressing used for absorption on minimal to moderately draining wounds. It has multiple layers comprised of a non-adherent contact layer, a non-woven pad for absorption, a breathable polyurethane barrier film and a conformable adhesive fabric tape.



Item#	Size
WC40022	2" x 2"
WC40044	4" x 5"
WC40066	6" x 6"
WC400410	4" x 10"

Wound Care

ZeniPLAST

Hydrocolloid Wound Dressing 10/box

ZeniPLAST is a sterile, hydrocolloid wound dressing used to help maintain an optimal moist wound environment conducive to wound healing. ZeniPLAST interacts with wound exudate to form a soft gel which helps to protect and prevent damage to newly formed tissue during dressing removal.



Item#	Size
WC50022	2" x 2"
WC50044	4" x 4"
WC50066	6" x 6"

ZeniFIBER

Calcium Alginate Wound Dressing 10/box

ZeniFIBER is a sterile, non-woven pad or rope made from calcium alginate fibers. Upon contact with wound exudate, ZeniFIBER becomes a fibrous gel that helps to maintain a moist wound environment that is conducive to wound healing.




Item#	Size
WC60022	2" x 2"
WC60045	4" x 5"
WC60066-Rope	6" Rope
WC60066	6" x 6"
WC60088	8" x 8"

ZeniCONTACT

Silicone Wound Contact Dressing 10/box

ZeniCONTACT is a sterile, open mesh, soft silicone wound contact layer that allows exudate to pass through to the secondary dressing. The silicone is gentle and non-adherent to the wound while preventing secondary dressings from adhering to fragile tissue. ZeniCONTACT dressings can be used as topical ointment carriers or in conjunction with negative pressure wound therapy.




Item#	Size
WC70047	4" x 7"

ZeniABSORB

Super Absorbent Dressings 10/box

ZeniABSORB is a sterile dressing for use on heavily exuding wounds, consistent of an oxygen permeable, waterproof nonwoven layer, a super absorbent pad layer and hydro-penetrating nonwoven layer. The layered construction facilitates fluid management to provide the optimal moist wound environment and locks in exudate for maximum absorption.



Item#	Size
WC80045	4" x 5"
WC80069	6" x 9"

Wound Care

ZeniCOL

Bovine Collagen Dressing and Powder 10/box

ZeniCOL is a high quality, sterile, non-hydrolyzed pure collagen dressing composed from 100% native bovine origin. Collagen is used to promote angiogenesis, decrease MMP levels in chronic wounds, stimulate collagen deposition and formation of granulation tissue in the wound.



Item#	Size
WC90022	2" x 2"
WC90044	4" x 4"
WC9001GP	Powder - 1 gram

ZeniCOL Ag

Collagen Matrix Dressing with Silver 10/box

ZeniCOL Ag Collagen Dressing with Silver is an advanced wound care dressing made of animal derived collagen, sodium alginate, carboxymethylcellulose (CMC), ethylenediaminetetraacetic acid (EDTA) and silver chloride. The topically applied wound dressing transforms into a soft gel sheet when in contact with wound exudate to maintain a moist wound environment that aids in the formation of granulation tissue and epithelialization and removes zinc to inhibit the activity of MMPs. The silver chloride in the dressing is an effective antimicrobial intended to prevent colonization of the dressing.



Item#	Size
WC90022AG	2" x 2"
WC90044AG	4" x 4"

ZeniFOAM Ag

Foam Wound Dressings with Silver 10/box

ZeniFOAM Ag is a highly absorbent, non-adherent, conformable foam dressing without an adhesive border used for the management of moderately to heavily draining wounds. Each dressing contains 0.35-0.40mg/cm² ionic silver as the active component to provide an optimal moist wound healing environment combined with effective antibacterial activity and exudate management.



Item#	Size
WC10022AG	2" x 2"
WC10045AG	4" x 5"
WC10066AG	6" x 6"
WC10088AG	8" x 8"

ZeniFOAM Ag GENTLE

Foam Dressings with Silver 10/box

ZeniFOAM Ag GENTLE is a sterile, silver impregnated, silicone-faced, polyurethane foam dressing without an adhesive border designed to absorb moderate to heavy wound drainage. ZeniFOAM Ag GENTLE provides maximum absorption while protecting sensitive or fragile skin and reducing painful dressing changes. ZeniFOAM Ag GENTLE contains 0.25-0.35mg/cm² of silver ions and has been shown to have antibacterial effectiveness for up to 7 days, as demonstrated in vitro.



Item#	Size
WC20044AG	4" x 4"

Wound Care

ZeniFOAM Ag GENTLE BORDER

Silver, Silicone Foam Dressings with Adhesive Border
10/box

ZeniFOAM Ag GENTLE BORDER is a sterile, silicone-faced, polyurethane foam dressing with a silicone adhesive border designed to absorb moderate to heavy wound drainage. ZeniFOAM Ag GENTLE BORDER contains 0.25-0.35mg/cm³ of silver ions and has been shown to have antibacterial effectiveness for up to 7 days, as demonstrated in vitro. ZeniFOAM Ag GENTLE BORDER is a multi-layered dressing that provides maximum absorption while protecting sensitive or fragile skin and reducing painful dressing changes.



Item#	Size
WC30044AG	4" x 4"
WC30066AG	6" x 6"

ZeniFIBER Ag

Calcium Alginate with Silver
10/box

ZeniFIBER Ag is a silver alginate dressing composed of natural fibers from seaweed combined with silver to provide exudate control and reduction of bacterial growth in the dressing while maintaining a moist environment and reducing wound odor.



Item#	Size
WC60012AGR	12" Rope
WC60022AG	2" x 2"
WC60045AG	4" x 5"
WC60066AG	6" x 6"
WC60088AG	8" x 8"

Equipment Covers

Bags for Regulators, Conservers, and Small Items



Item#
EC1012C
10×12"
1.5mil, Clear
1000 bags/box

Covers for Nebulizers and Suction Machines



20×24"
2mil
500 bags/roll

Item#	Color
EC2024C	Clear
EC2024T	Tan Tint

Respiratory Set Up Bags with Draw Tape



Item#
EC1216RB
12×16"
1mil
500 bags/box

Covers for Concentrators



25×15×30"
1.5mil
250 bags/roll

Item#	Color
EC251530C	Clear
EC251530T	Tan Tint

Covers for CPAP/Vents and Walkers



Item#
EC2430C
24×30"
1.5mil, Clear
500 bags/roll

Equipment Covers

Covers for Wheelchairs, Walkers, and Commodes




Item#	Thickness	Size	Color	Quantity
EC302035C	1.5mil	30x20x35"	Clear	200 bags/roll
EC302035T	1.5mil	30x20x35"	Tan Tint	200 bags/roll
EC5045	1mil	50x45"	Clear	250 bags/roll

- Protect equipment and materials from dust and moisture
- Packed on convenient rolls
- Low Density Polyethylene
- Packaged in dispenser carton

Mattress Covers



Item#	Vented	Thickness	Size	Color	Quantity
EC39990C	Yes	1.5mil	39x9x90"	Clear	100 bags/box
EC7252C	No	1mil	72x52"	Clear	200 bags/roll
EC7252T	No	1mil	72x52"	Tan Tint	100 bags/roll



Connect with us



Your sales representatives



Sean O'Donnell

Direct 312-637-6558
Email sodonnell@sunsethcs.com



Madison Gwilliam

Direct 312-533-2459
Email mgwilliam@sunsethcs.com



Evan Berger

Direct 262-788-8948
Email eberger@sunsethcs.com



John Bisceglie

Direct 312-533-2462
Email jbisceglie@sunsethcs.com



Jon Kmiec

Direct 312-533-2482
Email jkmiec@sunsethcs.com



Phil Cosens

Direct 312-533-2470
Email pcosens@sunsethcs.com



Brian Lipson

Direct 312-533-2463
Email blipson@sunsethcs.com



Jonny Deutsch

Direct 872-259-3691
Email jdeutsch@sunsethcs.com



Michelle Chojnacki

Direct 312-533-2486
Email mchojnacki@sunsethcs.com



Chris Crawford

Direct 872-258-1650
Email ccrawford@sunsethcs.com



Pavle Pantovic

Direct 872-282-0751
Email ppantovic@sunsethcs.com



Danielle Dubowski

Direct 312-533-2479
Email ddubowski@sunsethcs.com



Phone 877-578-6738
Fax 312-997-9985
Email sales@sunsethcs.com
Visit www.sunsethcs.com

Corporate Headquarters
141 W Jackson Blvd Ste 1950
Chicago IL 60604

Central Distribution Center
2725 Alft Ln
Elgin IL 60124

Western Distribution Center
6750 W Van Buren St Ste 100
Phoenix AZ 85043

National Accounts



Nikki Drake

Direct 312-533-2469
Email ndrake@sunsethcs.com



Gabrielle Murphy

Direct 312-638-8680
Email gmurphy@sunsethcs.com

International Accounts



Bryan Vigna

Direct 312-533-2455
Email bvigna@sunsethcs.com